Gebrauchsanweisung
Operating Instructions
Manual de instrucciones
Mode d'emploi
Istruzioni d'uso

STERILE EO

CE 0123

Rx only

Medicon eG · Gänsäcker 15
 78532 Tuttlingen · Germany
 Tel.: +49 (0) 74 62 / 20 09-0
 Fax: +49 (0) 74 62 / 20 09-50
 E-Mail: sales@medicon.de
 Internet: www.medicon.de

DE

MEDICON AURIKULAR REKO FIXIERDRAHT

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EN

MEDICON RECONSTRUCTION FIXING WIRE

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

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- General notes
- Intended use
- Indications
- Contraindications
- Possible side effects and complications
- Single use product
- Application and handling
- Pre-operative and post-operative aspects
- Sterile product
- Disposal
- Liability
- Signs and symbols

1. GENERAL NOTES

The MEDICON reconstruction fixing wire with double reinforcement is made of implant steel DIN ISO 5832-1, the needles are made of stainless steel DIN EN ISO 7153-1. The materials are biocompatible, corrosion-resistant and non-toxic in a biological environment. The surface is chemically passive, the material is antimagnetic.

The MEDICON reconstruction fixing wire with double reinforcement is not reusable and is supplied in sterile state. Please observe the following instructions. These will ensure accurate and reliable functioning of the product.

2. INTENDED USE

The MEDICON reconstruction fixing wire is used for reconstruction of the auricle by surgeons with sufficient experience in the required surgical techniques in hospitals and clinical practices. The MEDICON reconstruction fixing wire is a single-use product and may not be reused.

3. INDICATION

The MEDICON reconstruction fixing wire is used for reconstruction of the auricle using autogenic rib cartilage.

4. CONTRAINDICATIONS

Allergic reactions to the chrome and nickel constituents in the implant steel have not been observed in clinical practice but cannot be fully excluded.

5. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

In rare cases, the MEDICON reconstruction fixing wire may extrude through the skin if strong and long-term mechanical stress acts on the reconstructed auricle. This can be caused by wearing a helmet or a close-fitting cap or hat. In addition to the undesirable effects and complications already mentioned, the surgical procedure can lead to other problems such as infections, pain etc. that cannot be attributed to the implant.

6. SINGLE USE PRODUCT

The MEDICON reconstruction fixing wire is intended for single use on one patient only and may not be reused.

7. APPLICATION AND HANDLING

Bent or damaged needles must be discarded as they can break. To avoid damage, the reconstruction fixing wire should not be held on the connection points of needle and thread. Hold the needle only on the middle third. Only suitable needle holders with smooth mouth should be used, e.g. MEDICON needle holder REF 10.09.61 or REF 10.09.63.

8. PRE-OPERATIVE AND POST-OPERATIVE ASPECTS

Preoperative

Before the use of the product, the surgeon should thoroughly discuss the desired operation result with the patient. Particular attention must be given to post-operative behaviour and a potentially required follow-up care.

Post-operative

The patient must be instructed to immediately inform the surgeon about unusual in situ changes. A LOT number is indicated on the label of the packaging. To ensure 100% traceability, we recommend noting down the LOT number in the patient records.

9. STERILE PRODUCT

The MEDICON reconstruction fixing wire is a sterile product for single use. Products from damaged or broken sterile packaging may not be used, as sterility cannot be guaranteed. The product may not be reesterilised. The product may not be used after the end of the expiry date.

10. DISPOSAL

Needles as well as wire leftovers must be collected in solid and break-resistant disposable containers that must be tightly closed and stored, transported and disposed safe from unauthorised access.

11. LIABILITY

In the event of discrepancies between the english and the german version of these instructions for use only the German version shall be applicable. Only the latest revision of the instructions for use applies. Due to constant technical development the contents of these MEDICON instructions for use are updated regularly. Please use the MEDICON extranet www.medicon.de/extranet to ensure that you are using the current version. The current version date is printed on each edition of the instructions for use. MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use, non-observance of post-operative instructions and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

12. SIGNS AND SYMBOLS

- | | |
|--|---|
| | Follow instructions |
| | CE label according to directive 93/42/EEC |
| | Warning |
| | Single use product |
| | Manufacturer |
| | Order number |
| | Batch number |
| | Not reesterilizable |
| | Sterilized with ethylene oxide |
| | Do not use if packaging is damaged |
| | Keep away from sunlight |
| | Store in a dry place |
| | Limit humidity |
| | Limit temperature |

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