Hämorrhoiden-Ligator
Nr. 50.93.78

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- Anwendung
- Demontage / Montage
- Reinigung
- Sterilisation
- Lagerung
- Montageplan

Anwendung

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Hemorrhoid Ligator
No. 50.93.78

Contents

- Use
- Dismantling / assembly
- Cleaning
- Sterilization
- Storage
- Assembly plan

Use

This instrument is used solely for applying ligature rubber rings (Art. No. 50.93.75) and should also only be used for the intended purpose.

Dismantling

- See assembly plan -
Undo milled screw **D** in the direction of arrow a and remove slide **B**. Unscrew suction head **A** in the direction of arrow a.

Assembly

- See assembly plan -
Screw on suction head **A** in the direction of arrow b. Plug slide **B** onto suction head **A** and slide rod **C** into clamping device **E**. Turn milled screw **D** in the direction of arrow b and tighten.

Cleaning

Preparation for decontamination

The instruments must be placed on instrument carriers suitable for machine-based cleaning processes. The instrument carriers (e.g. strainer baskets) must be designed so that subsequent cleaning in the ultrasonic device or in the cleaning and disinfection device is not impaired by areas inaccessible to ultrasound or water.

Cleaning

Suitable alkaline cleaners or those with a neutral pH can be used. Selection of the cleaning agent is based on the properties of the instruments as well as on national guidelines and recommendations. Elevated chloride concentrations in the water can cause pitting or stress corrosion cracking on the instruments. Such corrosion can be minimized by using demineralized water or alkaline cleaning agents.

Thermal disinfection

Using demineralized water, thermal disinfection is carried out at temperatures from 80-95°C, with a corresponding soaking time according to the A0 concept, EN ISO 15883-1.

Drying

Sufficient drying must be ensured using the cleaning and disinfection device.

Cleaning agents

The parameters specified by the manufacturer of the cleaning agent for the concentration, temperature and soaking time must be complied with and it must be possible to check automatic dosing systems.

Maintenance, inspection and testing

Following cleaning/disinfection, the instruments must be macroscopically clean, i.e. free of visible residue. This is checked by visual inspection. Critical areas such as joints and jaw serrations require especially careful checking. Inadequately cleaned or clogged instruments must be cleaned again and must then be rinsed sufficiently. Before functional testing, instruments with movable parts must have cooled down and been lubricated with instrument maintenance oil. Instruments such as clamps and needle holders which have a detent should only be closed up to the first detent (danger of stress cracking). Defective instruments (showing hairline cracks, distortion or wear) must be

replaced because their functionality can no longer be guaranteed or can no longer be guaranteed sufficiently.

Sterilization

Only cleaned and disinfected instruments may be sterilized.



Only the sterilization methods listed may be used for sterilization.

- Fractionated vacuum method with sufficient product drying
- Autoclaving in accordance with DIN EN 13060 or DIN EN 13485 and DIN EN 554, as well as other processes validated in accordance with ISO 14937:2000 and USP 23.
- Sterilization time and temperature:
MEDICON eG recommends as follows:
At least 5 minutes holding time at 134°C

Storage

Instruments can corrode due to the storage conditions. To prevent this, instruments must be stored dry and protected from dust. To prevent condensation from forming, major temperature fluctuations should be avoided. Chemicals must not be stored with the instruments. The storage period depends on the type of packaging.

Assembly plan

Medicon eG, as the company marketing these products, does not assume any liability for direct or subsequent damage resulting from improper use or handling, in particular as a result of failure to observe the instructions for use or through incorrect care or maintenance.