
**Gebrauchsanweisung
Operating Instructions
Manual de instrucciones
Mode d'emploi
Istruzioni per l'uso**

Medicon
...DieBLAUEN.

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Medicon ... DieBLAUEN.

Dismantable rongeurs with silicone handle

Flat slider
Ergonomic design and balance
One-piece dismantling technology

Carefully read the operating instructions prior to clinical application and keep them safe and at hand. The instructions contained therein must be followed.

Rx only

Medicon eG · Gänsäcker 15 · D-78532 Tuttlingen · Germany
Tel.: +49 (0) 74 62 / 20 09-0 · Fax: +49 (0) 74 62 / 20 09-50
E-mail: sales@medicon.de · Internet: www.medicon.de

TABLE OF CONTENTS

1.	General notes	10
2.	Purpose and intended use	10
3.	Possible adverse effects and complications	10
4.	Application and handling	10-12
	4.1 Dismantling instruction	
	4.2 Assembly instructions	
	4.3 Handling of the 3-step control element	
5.	Decontamination, cleaning and sterilisation	12-15
	5.1 Reprocessing instructions	
6.	Liability	16
7.	Description of symbols and icons	16

1. GENERAL NOTES

These instruments are reusable. They are supplied non-sterile and must be cleaned, disinfected and sterilised prior to their first use and any further use. Before each use, the instruments must be visually examined, checked for damage, and suitably maintained. Please observe the following notes. This will ensure that the instrument functions faultlessly and reliably.

2. PURPOSE AND INTENDED USE

The demountable rongeurs "Medicon ...Die Blauen" are used exclusively for removing bone, cartilage and tissue from the skull and spinal column. The rongeurs may be used only by surgeons with sufficient training and experience in neurosurgery, orthopaedics or traumatology in hospitals and practices.

No materials foreign to the body, such as implants, wires or instruments, may be grasped.

MEDICON offers two versions of the demountable rongeurs "Medicon ...Die Blauen". The "Classic" and "Comfort" lines are each available with wide, narrow or extra-narrow cutting plate. The rongeurs of the "Comfort" line are additionally equipped with a grip adjustment that permits selecting from among three different levels of grip width and jaw opening ③.

Rongeurs with wide cutting plate are suitable for removal of hard bone and tough tissue. The rongeurs with wide cutting plate feature a silver push-button on the "Classic" and in addition a silver control element on the "Comfort" model.

The rongeurs with narrow cutting plate are labelled "for delicate use only" and may be used only for removal of thin or porous bone structures and soft tissues. The rongeurs with narrow cutting plate feature a golden push-button on the "Classic" and in addition a golden control element on the "Comfort" model.

The rongeurs with extra-narrow cutting plate are marked with the orange movable silicone grip part ② and the lettering "for cervical use only". They may be used only on the cervical spine for removal of thin or porous bone structures and soft tissue. The rongeurs with extra-narrow cutting plate feature a black push-button on the "Classic" and in addition a black control element on the "Comfort" model.

When selecting the cutting width, the relation to the quantity and structure of the material to be cut must always be taken into account.

3. POSSIBLE ADVERSE EFFECTS AND COMPLICATIONS

Possible adverse effects and complications linked to the demountable rongeurs "Medicon ...Die Blauen" include:

- Any use contrary to the intended use may harm the patient and/or damage the instruments. Safe use cannot be guaranteed with damaged instruments
- Nerve, tissue or vessel damage, infections or other unintended injuries during the surgical procedure

In addition to the adverse effects and complications already mentioned, the surgical procedure can lead to other problems such as nerve damage, infections, pain etc. that cannot be attributed to the instruments.

4. APPLICATION AND HANDLING

Even when used normally, the instruments are subject to wear and mechanical stresses, but even more so if used too forcefully. In order to prevent failure or mechanical damage to the instruments during surgery, they must be checked before use to ensure

that they are mechanically intact, that there are no deformations, and that the parts are fully functioning. Damaged instruments must not be used and are to be replaced.

- As an option, the rongeurs are available with a 3-step control element (see figure section 4.3). This allows three settings for the grip width/jaw opening ⑥.
- Rongeurs may be used only for their intended purpose.
- All the rongeurs are equipped with a fully integrated ejector ⑦ and an external guide ⑧.
- If any cartilage or bone residues still remain in the jaw opening ⑨, remove them before closing the rongeur. Do not use any hard objects to do this, because this might damage the blades.
- Rongeurs may not be used for leverage or turning. This would damage or deform the rongeur.
- The amount of material that can be punched is limited. Under no circumstances may more bone or cartilage be removed than fits into the jaw cavities of the shaft and main component.
- Conduct a function test and visually inspect the rongeurs before each use. Damaged rongeurs must not be used under any circumstances.
- When the spring is unhooked, the grips of the rongeurs must not be actuated under any circumstances, as otherwise the silicone could be damaged.

4.1 DISMANTLING INSTRUCTIONS


For rongeurs with adjustable grip width / jaw opening ⑥, the 3-step control element must be placed into the initial position before dismantling (see section 4.3).


The rongeurs must be dismantled for cleaning and disinfection.

Press the "push" button ①.

1. Pull the movable grip part ② downwards.
2. Push the slide ③ back up to the stop.


Fold the slide ③ open.


4.2 ASSEMBLY INSTRUCTIONS


For rongeurs with adjustable grip width/jaw opening ⑥, the 3-step control element must be placed into the initial position before assembly (see section 4.3).

1. Introduce the slide ③ into the external guide ⑧ of the main component ④.
2. Push the slide ③ forward until the two marking arrows ⑨ face each other.

- Push the movable grip ② upwards into the main component ④ until it snaps into position.
Do not hold the slide ③ tight!


Do not use handles while they are dismantled!

4.3 HANDLING OF THE 3-STEP CONTROL ELEMENT

By pressing the 3-step control element, grip width and jaw opening can be adjusted in three steps.

5. DECONTAMINATION, CLEANING AND STERILISATION

The following prion-specific precautions (v)CJD for instrument processing are indicated.

For persons:

- who suffer from (v)CJD or who are suspected of suffering from it [possible or clinically probable (v)CJD]
- who are related to a confirmed or suspected CJD patient (or a patient who died of CJD)
- who are receiving growth hormones, cornea or dura mater transplants
- with an unexplained, rapidly progressive CNS disease with and without dementia without concrete suspicion of CJD

If it is not possible to use single-use items, the invasive procedure must be planned in advance such that, among the instruments that are to be used, those are identified that can be fed into safe processing after careful consideration. All other instruments must be disposed of and incinerated after use. Instruments that can be safely processed owing to their structural and material properties are to be subjected to pre-cleaning taking into account the provisions of occupational safety.

Immerse the instruments in NaOH or NaOCl in an appropriate tub for 1 hour; the tub itself has to be adequately disinfected or disposed of after use. The resources used for the pre-cleaning (e.g. brushes) must be disposed of and incinerated after use. Following pre-cleaning and rinsing in water, clean the instrument under isolation in a machine using a validated prion-effective cleaning process (with proven cleaning power, alkaline pH values >10, application time >10 minutes, at a temperature of 55 °C and prion-inactivating cleaning agent, e.g. according to the ANSM list) and disinfect it with an appropriate, non-fixating thermochemical disinfection process. Owing to the possible impairment of the subsequent sterilisation, thermal drying should be avoided. After that, the WD should be put on an empty cycle. Steam sterilisation at 134 °C with a hold time of 18 minutes.

Alternatively:

In case of an unclear neurological diagnosis or suspicion of vCJD, the potentially contaminated instruments should be stored as evidence after isolated pre-cleaning and rinsing until the diagnosis is conclusively clarified. If (v)CJD is ruled out, routine processing according to these reprocessing instructions can take place. In case of a confirmed or inconclusive diagnosis (CJD/vCJD), the instruments must be disposed of and incinerated in a manner that precludes injury (e.g. processed with prion-specific measures, as far as acceptable in connection with CJD).

Storage as evidence:

Store the instruments dry on trays in an alkaline-resistant and steam-sterilisable container (e.g. V4A steel material number 1.4401) with a tightly sealing lid or in a suitable disposable container. The container must be secured with seals and must be stored with permanent and clear markings (under quarantine).

In case of an unidentifiable risk of (v)CJD:

Preventive routine processing: clean the instruments using a validated prion-effective automated cleaning process (with proven cleaning performance, alkaline pH values >10, contact time >10 minutes, at a temperature of 55 °C and prion-inactivating cleaning agent, e.g. according to the ANSM list) and disinfect them with an appropriate, non-fixating thermochemical disinfection process. Followed by steam sterilisation at 134 °C with a hold time of 18 minutes.

The respective national statutory regulations and guidelines regarding processing, disposal, and incineration are to be complied with.

Patient safety and occupational safety:

Process the surgical instruments intended for use in the central nervous system only by automated means. Instruments made of stainless steels may not be stored in physiological saline (NaCl) solution, as prolonged exposure will lead to corrosion problems such as pitting and stress corrosion. Only cleaned and disinfected instruments may be sterilised.

Limitations on reprocessing:

Immediately sort out instruments that are damaged, e.g. showing loss of function, wear, corrosion, cracks, deformations, loose parts, and fractures.

5.1 REPROCESSING INSTRUCTIONS

Place of use:

- Immediately after use, remove any gross contamination from the instruments with a disposable cloth or paper
- Do not use any fixing agents or hot water (>40 °C / 104 °F), as this would result in fixation of residues and might affect the success of the cleaning
- Feed the instruments into reprocessing immediately
- Dry disposal is to be preferred

Preparation for decontamination:

The rongeurs must be put into the washer-disinfector device (WD) in disassembled state, on instrument carriers suitable for machine cleaning, so that they can be rinsed. The instrument carriers (e.g. sieve trays) must be such that subsequent cleaning in an ultrasonic device or in the washing and disinfection device is not impaired by "rinsing shadows" or "sound shadows" that may render any areas inaccessible to ultrasound or rinsing.

Pre-cleaning:

- Immerse the instruments into cold water for 5 minutes
- Brush the instruments (using plastic brushes) under running cold water until all visible dirt has been removed
- Rinse the interior cavities, threads and drilled holes with the water gun for 10 seconds, and then brush them again

Automated cleaning:

Use appropriate pH-neutral or alkaline cleaning agents. The type of cleaning agent depends on the national regulations and recommendations. The washer-disinfector (WD) must meet the requirements of DIN EN ISO 15883-1.

- 1. Prewashing: 1 minute with softened cold water, without additives
- Emptying
- 2. Pre-rinse: 3 minutes with softened cold water, without additives
- Emptying
- Cleaning: with softened water, heating to >55 °C / 131 °F and 5 minutes washing/cleaning, adding of the cleaning agent at 45 °C / 113 °F, alkaline cleaning agent, dosage 0.5 %
- Emptying
- Neutralisation: 3 minutes with warm water (>40 °C) including the addition of neutraliser, dosage 1 ml/l
- Emptying
- Final rinse: 2 minutes with warm, fully deionised water (>40 °C) (with additives)
- Emptying

**Disinfection:**

Thermal disinfection, A_0 value 3000:

- Fully demineralised water, thermal disinfection is carried out at temperatures $>80\text{ }^{\circ}\text{C}$ / $176\text{ }^{\circ}\text{F}$ with a corresponding application time according to the A_0 concept of standard DIN EN ISO 15883 and guideline DGKH, DGSV and AKI (e.g. A_0 3000 = $90\text{ }^{\circ}\text{C}$ / $194\text{ }^{\circ}\text{F}$ and 5 minutes application time). The operator is responsible for the A_0 value to be implemented.

Drying:

Ensure sufficient drying in the WD. The instruments must be removed from the WD immediately after completion of the cleaning and disinfection programme. If necessary, use of compressed air for drying is recommended due to its good and rapid effect (RKI recommendation).

Manual cleaning/disinfection:

As a matter of principle, the detergents and disinfectants to be used must be suited for the manual cleaning and disinfection of instruments, and they must also be mutually compatible. The effectiveness of the disinfectant must be proven. When choosing the disinfectant and procedure, refer to the relevant lists and recommendations of the Robert Koch Institute (RKI) and the German Society for Hygiene and Microbiology (DGHM).

Pre-cleaning:

- Immerse the instruments into cold water for 5 minutes
- Brush the instruments (using plastic brushes) under running cold water until all visible dirt has been removed
- Rinse the interior cavities, threads and drilled holes with the water gun for 10 seconds, and then brush them again
- Remove the instruments and rinse them with cold, fully desalinated water

Cleaning / disinfection:

- Place the instruments into a bath with an approved detergent and disinfectant
- The instruments must be completely immersed into the solution
- The contact times, temperatures and concentrations specified by the manufacturer of the cleaning or disinfectant agent must be absolutely complied with
- Remove the instruments and rinse them with cold demineralised water for two minutes
- Repeat the cleaning process if any residual contamination can still be seen on the instrument
- If necessary, perform cleaning and disinfection in an ultrasonic bath with sonication

Each day, freshly prepared working solutions must be used. In cases of extreme contamination, the working solution must be changed more frequently. A high level of contamination in the ultrasonic bath will impair the cleaning effect and promote corrosion. The cleaning solution must be regularly replaced, depending on the conditions of use. The decisive criterion is optically recognisable contamination. In any case, the bath must be changed frequently, at least once a day. Observe the valid national guidelines.

Drying:

Manual drying must be carried out using compressed air and a lint-free cloth. Use of compressed air for drying is recommended due to its good and quick effect (RKI recommendation).


Maintenance, control and inspection:

Following cleaning/disinfection, the instruments must be macroscopically clean, i.e. free from visible dirt and residues. The inspection is performed visually. Critical areas such as the sliding surfaces, guide, plug-in connection, dismantling mechanism, and spring component need to be checked particularly carefully. Insufficiently cleaned rongeurs must be cleaned again and then rinsed and dried sufficiently. Before sterilisation, the rongeurs must be assembled according to the assembly instructions. The sliding surfaces, guide, plug-in connection and dismantling mechanism must be treated with steam-permeable instrument spraying oil (e.g. Medicon instrument care oil).

REF 46.00.40 Medicon instrument care oil

After assembly, repeatedly open and close the grip parts to spread the oil evenly on the sliding surfaces, guide, plug-in connection and dismantling mechanism, in order to avoid metal-to-metal friction and related functional and fretting corrosion problems. The instruments must be cooled before functional testing, and oiled with an instrument spraying oil (e.g. Medicon instrument care oil).

REF 46.00.40 Medicon instrument care oil

Defective instruments (with hairline cracks, deformation or wear) must be replaced, as they cannot fulfil their function, or cannot fulfil it safely enough. Similarly, any rusty instruments must be removed, as they can cause corrosion on intact instruments through the spread of extraneous rust.

Packaging:

The instruments are to be placed in an appropriate sterile barrier system. The sterile barrier system must meet the following requirements:

- DIN EN 868
- DIN EN ISO 11607
- suitable for steam sterilisation (steam permeability)
- sufficient temperature resistance up to 138 °C / 281 °F

Sterilisation accessories and sterilisation wrap must be adapted to the package content and the sterilisation procedure applied.

Sterilisation:

For sterilisation, the following sterilisation procedure is to be used, taking into account all applicable national requirements:

- Fractionated vacuum procedure, fractionated three times and with sufficient product drying time
- Steam steriliser validated according to DIN EN 13060 or DIN EN 285 and according to DIN EN ISO 17665-1
- Sterilisation time and temperature: At least 5 minutes hold time at 134 °C

Reaching an SAL (Sterility Assurance Level) of 10⁻⁶ is indispensable.

Storage:

Reprocessed sterile instruments are to be stored in dry, dust-protected, low-germ, dark and cool rooms, free of vermin, in an appropriate reusable sterilisation container. To prevent condensation, extreme temperature fluctuations should be avoided during storage. No chemicals may be stored together with instruments. Walls, floors and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must have a ground clearance of at least 30 cm.

The permissible on-site storage period depends on the type of sterile barrier system used and on the storage conditions. The permissible storage period is to be specified by the operator.

Additional information on reprocessing:

A validated automated cleaning and disinfection procedure is always to be preferred to manual cleaning because of the greater safety of the process. Thorough cleaning also helps to maintain value and is a prerequisite for successful sterilisation. Please take note of the following points during automated processing:

- For effective automated processing, load the wire mesh baskets in such a way that they can easily be mechanically rinsed. Mesh baskets must not be overloaded
- Avoid "rinsing shadows" cast by large instruments
- To avoid any damage, the mechanical sensitivity of the instruments must be taken into account when they are placed down or stored

The times and temperatures stated in these reprocessing instructions are minimum requirements which must be adhered to. Should any downward deviation be needed for procedural reasons, this must be validated by the operator. Exceeding the specified times and temperatures is basically possible. However, this will lead to increased stress on the material, which may cause the instruments to wear out prematurely. We accept no responsibility for use of other sterilisation procedures.

Information on validating the processing procedure:

Validation is carried out with the following devices, materials and chemicals:

Cleaning and disinfection equipment:	Type Miele Desinfektor G 7735 CD and 7836 CD, trolley for surgical instruments
Alkaline cleaning agent:	neodisher® FA · Dr. Weigert GmbH & Co. KG
Neutral detergent:	Endozime · Fa. Ruhof (enzymatic) Decondex 23 Neutrazym · Borer Switzerland
Neutraliser:	neodisher® Z · Dr. Weigert GmbH & Co. KG
Water gun:	Selecta
Cleaning brushes:	Plastic or nylon bristles
Ultrasonic bath:	Sonorex
Sterilisers:	MMM Vaculab 969 S 3000, MMM Selectomat S 3000 Stiefenhofer KS 666-2ED, H+P Varioclav 400E

Minimising the risk of transmission of (v)CJD:

 Based on national statutory guidelines and regulations, we recommend using the alkaline cleaning agent neodisher® SeptoClean instead of neodisher® FA for the routine automated cleaning procedure in case of indicated prion-specific precautions, as it has been shown to have prion-inactivating properties (see ANSM list).

Note:

The person processing the instrument bears the responsibility that the processing actually performed with the equipment, materials and staff at the reprocessing facility achieves the desired results. Validation and routine monitoring of the process are normally required. If the devices, materials and chemicals described above are not available, it is up to the persons processing the instrument to have their procedures appropriately validated.

Observe the notes and regulations of the relevant national guidelines and standards. Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilised prior to dispatch. MEDICON eG reserves the right to make changes to these instructions due to new findings.

Disposal:

Observe the applicable national laws when disposing of the rongeurs and their accessories!

6. LIABILITY

In the event of any discrepancies between the non-German and the German version of these operating instructions, only the German version is authoritative. Only the latest revision of the operating instructions applies. Due to constant technical development, the contents of these MEDICON operating instructions are updated regularly. Please use the

MEDICON extranet: www.medicon.de/extranet

to ensure that you are using the current version. The version date of each edition of the operating instructions is stamped on it. MEDICON eG assumes no liability for any damage caused by improper use, incorrect postoperative handling, care or maintenance, or non-compliance with the restrictions on use and other guidelines in the operating instructions. MEDICON eG assumes no liability for defects in case of changes or repairs to the product without prior written consent from MEDICON eG either, nor in case of repairs not carried out by workshops authorised by MEDICON eG or by the MEDICON Repair Service.

7. DESCRIPTION OF SYMBOLS AND ICONS

	Observe the operating instructions
	Please note
	CE label acc. to Directive 93/42/EEC
	Manufacturer
Rx only	Subject to medical prescription (US federal law)
	Non-sterile
	Batch code
	Item number

PLEASE NOTE: According to US federal law, in the USA this product may be bought only by a physician or hospital or upon prescription!

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© Copyright 01/2018, MEDICON eG, Tuttlingen, Germany