Gebrauchsanweisung
Operating instructions
Manual de instrucciones
Mode d'emploi
Istruzioni per l'uso

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... DieBLAUEN.2

Table of Contents

1.	General information	10
2.	Purpose and proper use	10
3.	Possible side effects and complications	10
4.	Use and handling 10 –	11
	4.1 Dismantling instruction	11
	4.2 Assembly instructions	11
5.	Use, Cleaning, Disinfection, Sterilization 12 –	15
6.	Note	15
7.	Liability	15
8.	Signs and symbols	16

Rx only

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions in this document have to be followed.

Medicon eG · Gänsäcker 15 · D-78532 Tuttlingen · Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50 · E-Mail: sales@medicon.de · www.medicon.de

1. General information

The dismantable rongeurs "MEDICON ...Die Blauen.2" are not supplied sterile and must be cleaned, disinfected and sterilized before every use. Please observe the following instructions. This will guarantee that the rongeurs function faultlessly and reliably. The dismantable rongeurs with silicone handle feature two blue silicone handle sections.

2. Purpose and proper use

The dismantable rongeurs "MEDICON ...Die Blauen.2" are used exclusively for gripping and holding bone, cartilage and tissue. The rongeur must only be used by surgeons with sufficient training and experience in hospitals and practices. The rongeurs are reusable after application of appropriate treatment procedures.

No materials foreign to the body, such as implants, wires or instruments, may be grasped.

3. Possible side effects and complications

The possible undesirable effects and complications linked to the dismantable rongeurs, include:

- Using the rongeurs for purposes other than those for which it was designed may damage its components. It may not be safe to use it if any of the components become damaged.
- Nerve, tissue or vascular damage, infections or other undesirable injuries during the surgical procedure.

In addition to the undesirable effects and complications already mentioned, the surgical procedure could lead to other problems such as nerve damage, infections, pain etc. that cannot be attributed to the instruments.

4. Use and handling


- The dismantable Rongeurs may only be used for their intended purpose.
- All demountable rongeurs are equipped with a jaw (4), shaft (3), slider (2), dismantling mechanism (1), and blue silicone ring handles.
- Rongeurs may not be used for leverage or turning. This will damage or deform the rongeur.
- The amount of material that can be punched is limited. Under no circumstances may more bones or cartilage be removed than fit in the jaw cavities of the shaft and main component.
- Overloading must be avoided.

The dismantable rongeurs "MEDICON ...Die Blauen.2" are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force. To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional.


4.1 Dismantling instruction


The rongeurs must be dismantled prior to cleaning and disinfection.

Opening of the disassembly mechanism ① by 90° rotation to the left or to the right.


Move silicone ring handle ② in the direction of the arrow.


Open slider ③

4.2 Assembly instructions


The rongeurs must be assembled prior to sterilization.

Connect slider ② to shaft ①. By appropriate positioning of the silicone ring handle, the slider latches ③.

Move silicone ring handle ② in the direction of the arrow.

Move disassembly mechanism ① back to the starting position.

In the initial position, the recess of the disassembly mechanism points upwards.

5. Use, Cleaning, Disinfection, Sterilization

The following vCJD prion-specific protective measure is indicated when processing instruments.

In the event of diagnosis of a definite or probable case of vCJD.

If it is not possible to use disposable products, the instrument used, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

If prion contamination is suspected:

In the event of suspected prion contamination, incineration of the instrument is recommended according to the vCJD task force final report.

If vCJD is excluded:

Continue to use after instrument processing is completed. Otherwise, the instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

In the case of a non-identifiable vCJD illness.

Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions – e.g. mechanical alkaline cleaning combined with steam sterilization.

If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (CNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134°C.

Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress-corrosion cracking.

Only cleaned and disinfected instruments may be sterilized.

Reprocessing restrictions

Frequent reprocessing has few effects on these instruments. The end of their working life is normally determined by wear and damage through use. If there are any reprocessing restrictions, this is noted in the instructions for use of the specific product.

Reprocessing Instructions

Location of use

- Directly after use, gross contamination should be removed from the instruments with a disposable cloth/paper.
- Fixing agents or hot water (>40°C) must not be used as this leads to fixation of residues and can affect cleaning outcome.
- The instruments should be conveyed promptly to processing.
- Dry disposal preferred.

Preparation for decontamination

The rongeurs must be put in the washing and disinfection device so that they can be rinsed on appropriate suitable instrument carriers. The instrument carriers (e.g. sieve trays) must be such that subsequent cleaning in an ultrasonic device or in the washing and disinfection device is not impaired by any areas that are inaccessible to ultrasonics or rinsing.

Pre-cleaning

- Place the instruments in cold water for 5 minutes;
- Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
- Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again;

**Cleaning: machine**

The cleaning and disinfection device (CDD) must meet DIN EN ISO 15883-1 requirements.

- Pre-rinse 1: 1 minute with cold demineralised water, without additive;
- Emptying;
- Pre-rinse 2: 3 minutes with cold demineralised water, without additive;
- Emptying;
- Cleaning: with demineralised water, heat to 55°C and wash/clean for 10 minutes, add cleansing agent at 45°C, alkaline cleansing agent, strength 0.5%;
- Emptying.
- Neutralization: 3 minutes with warm water (>40°C) with addition of neutralizer, strength 1ml/l;
- Emptying;
- Final rinse: 2 minutes with warm deionised water (>40°C) (without any additive);
- Emptying

Disinfection

Thermal Disinfection A0 value 3000:

Completely desalinated water, the thermal disinfection is carried out at temperatures > 80 °C and corresponding application time according to the A0 concept, DIN EN ISO 15883 (e.g. A0 3000 = 90°C and 5 minutes application time). The operator is responsible for the implemented A0 value.

Drying

Complete drying must be ensured by the CDD. The instruments must be removed from the CDD promptly once the cleaning and disinfection program has finished. If necessary, compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Cleaning/disinfection: manual

The cleansing agents and disinfectants employed must be suitable for manual cleaning and disinfection of instruments and must be compatible with each other. The disinfectant must be of tested efficacy. When selecting the disinfectant agent and method, the relevant lists and recommendations of the Robert Koch Institute (RKI) and the German Society for Hygiene and Microbiology (DGHM) should be noted.

Pre-cleaning

- Place the instruments in cold water for 5 minutes;
- Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
- Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again;
- Thorough rinsing with demineralised water

Cleaning/disinfection

- Place the instruments in a bath with a tested cleansing and disinfectant agent;
- The instruments must be completely covered with the solution;
- The application times, temperatures, and concentration stated by the manufacturer of the cleansing/ disinfectant agent must always be observed;
- Remove the instruments and rinse for 2 minutes with cold deionised water;
- Repeat the cleansing process if visible contamination is still present on the instrument;
- If necessary, clean and disinfect the instruments in an ultrasonic bath.

Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed sooner. A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

Drying

Dry manually with compressed air and a lint free cloth. Compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Maintenance, inspection and testing

After cleaning/disinfection the instruments must be macroscopically clean, i.e., free from visible residues and soiling. Inspection is performed visually. Critical areas such as the sliding surfaces, hinges and dismantling mechanism need to be checked particularly carefully.

Insufficiently cleaned instruments must be cleaned again and then adequately rinsed and dried. Before sterilization, the rongeurs must be assembled according to the assembly instructions. The sliding surfaces, hinges and dismantling mechanism must be treated with steam-permeable instrument spraying oil (e.g. Medicon [REF] 46.00.40). After assembly, repeatedly open and close the grip parts to spread the oil evenly on the sliding surfaces, hinges and dismantling mechanism, in order to avoid metal-to-metal friction and related functional and fretting corrosion problems.

Before functional testing, the sliding surfaces, hinges and dismantling mechanism must be lubricated with instrument maintenance oil Medicon [REF] 46.00.40.

Defective instruments (hairline cracks, deformation or wear) must be replaced as they no longer fulfil their function or do so without adequate safety. Corroded instruments must also be replaced as they may corrode intact instruments by extraneous rust.

Packing

The instruments must be placed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868,
- DIN EN ISO 11607,
- suitable for steam sterilization (vapor permeable),
- adequate temperature resistance up to 138°C.

Sterilization equipment and sterilization wrapping must match both the wrap contents and the employed sterilization method.

Sterilization

Taking into account the respective national regulations, the following method must be employed for sterilization:

- Vacuum autoclave with triple vacuum and adequate drying of the products (Vacuum minimal 15 Minutes)
- Steam sterilizer complying with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: **at least 5 minutes hold time at 134°C.**

It is essential to attain a Sterility Assurance Level of 10^{-6} .

Storage

Reprocessed sterile instruments must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments.

The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority.

Further information on reprocessing

Validated machine cleaning and disinfection is always preferred over manual cleaning because of the greater certainty of the method. Good cleaning helps to preserve value and is a precondition of successful sterilization. During machine processing, the following points should be noted:

- Correct loading of the trays for rinsing is a precondition for effective machine processing. Trays must not be overloaded.
- Rinsing shadows due to large instruments must be avoided.
- The instruments must be placed or stored based on their susceptibility to mechanical damage in order to prevent them from becoming damaged.

Operating instructions

The times and temperatures specified in these reprocessing instructions are minimum requirements and must not be less than those stated here. If they are to be reduced for technical reasons, this must be validated by the operating authority. Exceeding the stated times and temperatures is always possible but leads to increased stress on the material, which may result in premature ageing of the instruments. The use of other sterilization methods is outside our responsibility.

Information on validation of the processing

Validation was performed with the following equipment, materials, and chemicals:

Cleaning and disinfection device:	Type Miele Disinfector G 7735 CD and 7836 CD compartmented cart for surgical instruments
Cleaning agent alkaline:	neodisher® FA, Dr. Weigert GmbH & Co. KG
Cleaning agent enzymatic:	Endozime, Ruhof (enzymatic) Decondex 23 Neutrazym, Borer, Switzerland neodisher® Z, Dr. Weigert GmbH & Co. KG
Neutraliser:	Selecta
Water pistol:	Plastic/nylon brushes
Cleaning brush:	Sonorex
Ultrasonic bath:	MMM Vaculab 969 S 3000, MMM Selectomat S 3000, Stiefenhofer KS 666-2ED, H+P Varioclav 400E
Steriliser:	

6. Note

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method.

If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device.

Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch.

Medicon eG reserves the right to modify these instructions whenever new information is obtained.

Removal

Observe the legal regulations or the applicable national guidelines for disposal of the Medicon rongeurs, its components and accessories!!

7. Liability

MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines in the instructions for use.

The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

CAUTION:

For USA, federal law restricts this device to purchase by or on behalf of physician or hospital.

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