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
please read these instructions carefully before clinical application and keep them safe and close to hand. the instructions must be followed.

The MEDICON retractor X-ray transparent spinal spreader system is manufactured from PEEK and may be used only by surgeons of proper training and experience in spine surgery (neurosurgeon, orthopaedic surgeon, etc.).

table of contents

- General notes
- Packaging contents
- Purpose and proper use/indications
- Contraindications
- Possible side effects and complications
- MR Instructions
- Instructions for assembly and disassembly
- Decontamination, cleaning and sterilization
- Repair
- Proper disposal
- Exclusion of liability
- Signs and symbol

1. GENERAL NOTES

The retractor X-ray transparent spinal spreader system is supplied unsterile and must be disinfected, cleaned and sterilised prior to use. Please observe the following instructions. This will guarantee that the implant functions faultlessly and reliably.

2. PACKAGING CONTENTS

The packaging contents of the retractor Spinal Spreader System (57.62.05) include:

- Spindle (905.72.00)
- Stabiliser with right-handed thread (905.72.01) (R)
- Stabiliser with left-handed thread (905.72.02) (L)
- Valve support (x4 57.62.10)
- U-ring spanner (57.62.15)
- Hexagon socket screwdriver (57.62.16)

Valves to be used in combination (not included): Only PEEK valves from MEDICON eG may be used with this spreader system. (Art.-Nr. 57.59.63 – 57.59.92).

3. PURPOSE AND PROPER USE/INDICATIONS

The retractor Spinal Spreader System is used for spreading and keeping open surgical access during lumbar, dorsal and cervical spine surgery as well as with spreading the surgical access during surgical exposures of longer sections of the spine requiring permanent repositioning of adjacent soft tissues. In particular, it is used in spinal canal stenosis, resection procedures of extensive tumours or stabilisation surgery in which intraoperative imaging is used. Due to the complete avoidance of metallic materials the spreader system allows artefact-free intraoperative imaging.

4. CONTRAINDICATIONS

Applications contrary to the intended use of the product.

5. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

With excessive loads due to pulling and pushing forces during the spreading action, irritation, local circulatory disorders in the skin and damage to the soft tissue and muscles may occur. This may cause tissue necrosis, impaired wound healing and infections. Furthermore, an overloading of the spreader system may lead to breakage of the instrument or individual components.

WARNINGS AND PRECAUTIONS

So that there is no overloading of the spreader system or tensions being applied to soft tissues, it is the responsibility of the surgeon to extend the surgical wound if necessary, in order to bring about a corresponding degree of relief.

6. MR INSTRUCTIONS Instruments

The use of medical devices in the vicinity of an MR poses a hazard. Individual medical devices must not be located near the equipment during the application of these procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MR environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MR devices.

7. INSTRUCTIONS FOR ASSEMBLY AND DISASSEMBLY

The stabilisers (2 905.72.01 and 3 905.72.02) must be screwed onto the spindle 1 (905.72.00) at the respective ends corresponding to their labelling (L or R) (Images 1 + 2). Please note that both the stabilisers and the spindle have different threads (left-/right-handed threads). Prior to insertion into the wound opening the PEEK-valves must be inserted into the valve support 4 (57.62.10) from below at a 90° angle to the direction of insertion (Image 3). In doing so, the base of the valve bolt is countersunk into the recess at the bottom of the valve support. Following insertion, the valves can be pivoted as required (Image 4). After placing the valves incl. valve holders in the wound opening, the stabilisers mounted at the spindle can be inserted into the valve holders by turning the spindle. The wound spreading is achieved by turning the spindle using the U-ring spanner 5 (57.62.15) at the left or right of the gripping surface of the spindle (Image 5) or using the hexagon socket screwdriver 6 (57.62.16) at the respective ends of the spindle (Image 6).

The spreader system is released in the opposite direction of the wound spreading. The removal of the valves from the valve holder is achieved by another 90° turn. Disassembly of the system is performed in the reverse order.

8. DECONTAMINATION, CLEANING AND STERILIZATION The following prion-specific protective measure (vCJD) is indicated when preparing the instrument.

In the event of a diagnosis of definite or probable vCJD:

If the use of disposable products is not possible, the instruments used (which had become contaminated or whose contamination cannot be ruled out) must be disposed of as waste for incineration.

In suspected cases:

If prion contamination is suspected, it is recommended that any medical instruments be incinerated in accordance with the final report by the vCJD task force.

If vCJD infection is ruled out:

Re-use after final treatment. Otherwise the instruments used (that had become contaminated or where contamination cannot be ruled out) must be disposed of as waste for incineration.

In the event that vCJD infection is not apparent.

Even if there is no known existence of a prion disease, we should, in preparation, set up two procedures with at least partial effectiveness against prion contamination, for example, an alkaline mechanical cleaning process combined with steam sterilisation.

If an alkaline mechanical cleaning process is not used, or another cleaning process with proven effectiveness against prion contamination is deployed, and medical devices that come into contact with risk tissue (CNS, eyes, lymphatic tissue) are used, the Robert Koch Institute (IRK) recommends an extended sterilisation time of 18 minutes at 134°C.

Instruments made of stainless steel must not be placed in physiological saline solution (nacl solution) as prolonged contact leads to corrosion such as pitting and stress-corrosion cracking. Only cleaned and disinfected implants may be sterilized.

Limitations on Reprocessing

Frequent reprocessing has little effect on these instruments. Their service life normally comes to an end as a result of wear and damage by use.

Location of use

- Directly after use, gross contamination should be removed from the instruments with a disposable cloth/paper.
- Fixing agents or hot water (>40°C) must not be used as this leads to fixation of

residues and can affect cleaning outcome.

- The instruments should be conveyed promptly to processing.
- Dry disposal preferred.

Preparation for decontamination

The Instrument has to be disassembled for reprocessing (see Chapter 7. Instructions for assembly and disassembly). The instrument must be placed on appropriate instrument carriers so that they can be rinsed in the washing and disinfection device. The instrument carriers (e.g. sieve trays) must be such that subsequent cleaning in an ultrasonic device or in the washing and disinfection device is not impaired by any areas that are inaccessible to ultrasonics or rinsing.

Cleaning: automatically

The cleaning and disinfection device must conform to the requirements of DIN EN ISO 15883-1.

- Pre-rinse 1: 1 minute with cold demineralised water with no additives;
- Drain;
- Pre-rinse 2: 3 minutes with cold demineralised water with no additives;
- Drain;
- Cleaning: use demineralised water, heat to >55°C, wash and clean for 10 minutes, add cleaning agent at 45°C, alkaline cleaning agent, dosage 0.5% (pH 11.4 – 11.9);
- Drain;
- Neutralization: 3 minutes with warm deionised water (>40°C) with neutralizer, dosage 1 ml/l;
- Drain;
- Final rinse: 2 minutes with warm deionised water (>40°C) (no additives);
- Drain

Disinfection

Thermal Disinfection A0 value 3000: Completely desalinated water, the thermal disinfection is carried out at temperatures > 80 °C and corresponding application time according to the A0 concept, DIN EN ISO 15883 and guideline DGKH, DGSV and AKI (e.g. A0 3000 = 90°C and 5 minutes application time). The operator is responsible for the implemented A0 value."

Drying

Sufficient drying time must be provided by the cleaning and disinfection device. The instruments must be removed from the cleaning-disinfection device immediately after the cleaning and disinfection program has finished. If necessary, the use of compressed air for drying is recommended because of its effective and quick effect (RKI recommendation).

Cleaning/disinfection: manual

The cleaning agents and disinfectants employed must be suitable for manual cleaning and disinfection of instruments and must be compatible with each other. The disinfectant must be of tested efficacy. When selecting the disinfectant agent and method, the relevant lists and recommendations of the Robert Koch Institute (RKI) and the German Society for Hygiene and Microbiology (DGHM) should be noted.

Pre-cleaning

- Place the instruments in cold water for 5 minutes;
- Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
- Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again;
- Thorough rinsing with demineralised water

Cleaning/disinfection

- Place the instruments in a bath with a tested cleansing and disinfectant agent;
- The instruments must be completely covered with the solution;
- The application times, temperatures, and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed;
- Remove the instruments and rinse for 2 minutes with cold deionised water;
- Repeat the cleansing process if visible contamination is still present on the instrument;
- If necessary, clean and disinfect the instruments in an ultrasonic bath.

Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed sooner. A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

Drying

Dry manually with compressed air and a lint free cloth. Compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Maintenance, control and inspection

The instruments must be macroscopically clean, i.e. free of visible residues, after cleaning and disinfection. The instruments must be visually inspected in this respect. Inadequately cleaned parts must be cleaned, thoroughly rinsed and dried again. Defective instruments (hairline cracks, deformation or wear) must be replaced, as they cannot fulfil their function, or cannot fulfil it safely enough. Similarly, any rusty instruments must be removed as they can cause corrosion through the spread of rust.

Packing

The instruments must be packed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868 and DIN EN ISO 11607
- suitable for steam sterilization (steam permeable)
- adequate temperature resistance to 138°C

Sterilization accessories and sterilization packaging must be appropriate for the content of the package and for the sterilization process.

Sterilization

Use the following sterilization procedures subject to the applicable national requirements:

- Fractionated vacuum method with triple fractionation and sufficient product drying.
- Steam sterilization in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: **at least 5 minutes holding time at 134°C.**

It is essential that a SAI (Sterility Assurance level) of 10⁻⁶ is reached.

Storage

Reprocessed sterile instruments must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments.

The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority.

Further information on reprocessing

A validated machine cleaning and disinfection procedure is always to be preferred to manual cleaning because of the greater safety of the process. Thorough cleaning also helps to maintain value and is a prerequisite for successful sterilisation.

Please take notice of the following points regarding mechanical processing:

- Load the wire mesh baskets in such a way that they can easily be mechanically rinsed. Mesh baskets must not be overloaded.
- Avoid „rinsing shadows“ caused by large instruments.
- The mechanical sensitivity of the instruments must be taken into account when they are placed down or stored to avoid any damage occurring.

The times and temperatures stated in these reprocessing instructions are minimum requirements, they must be reached. Should any downward deviation be needed for procedural reasons, this must be validated by the operator. Exceeding the specified times and temperatures is basically possible. However, this will lead to increased stress on the material, which can cause the instruments to wear out prematurely.

We accept no responsibility for the use of other sterilisation procedures.

Information on validating the reprocessing procedure

Validation is carried out with the following apparatus, materials and chemicals:

Rx only

Cleaning and disinfection equipment: Type Miele Desinfektor G 7735 CD
Trolley for surgical instruments
Alkaline detergent: neodisher® Fa. Dr. Weigert GmbH & Co. KG
Neutralizing detergent: neodisher® Z, Dr. Weigert GmbH & Co. KG
Water gun: Selecta
Cleaning brushes: bristles of plastic or nylon
Ultrasonic bath: Sonorex
Sterilisers: MMM Vaculab 969 S 3000,
MMM Selectomat S 3000,
Stiefenhofer KS 666-2ED, H+P Varioclav 400E

12. SIGNS AND SYMBOLS

Symbol	Meaning
	Manufacturer
	Warning
	Follow instructions
	MR unsafe
	CE label accord. To directive 93/42/EEC

Note

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method. If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device. Medicon eG reserves the right to modify these instructions whenever new information is obtained. Please refer to the instructions and requirements of the relevant national regulations and standards as well as the instructions for use included with the medical device, if applicable.

9. REPAIR

If the above-described apparatus, materials and chemicals are not available, it is up to the persons preparing the instrument to have their own procedures validated. Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch.


10. PROPER DISPOSAL

To avoid the risk of infection to third parties, the instruments must be cleaned and sterilised before disposal. In addition, the instruments must be disposed of in specific and appropriately labelled containers, in order to avoid any danger of third parties suffering injury from cuts.

11. EXCLUSION OF LIABILITY

In the event of discrepancies between the English and the German version of these instructions for use only the German version shall be applicable. The information in these instructions for use applies only to the product with which the instructions for use are supplied. The instructions for use for the current product in each case can be downloaded on the Medicon extranet at www.medicon.de/extranet or requested directly from MEDICON eG. MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

CAUTION:
for USA, federal law restricts this device to purchase by or on behalf of physician or hospital.

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