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English

Rev.: 10-12-2025

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Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

Medicon implants are manufactured from PEEK or titanium alloy. Both materials are biocompatible, corrosion-proof and non-toxic in the biological environment. The surface is chemically passive and the material is antimagnetic.

#### Table of contents

1. General notes
2. Purpose and proper use/Indications
3. Contraindications
4. Possible side effects and complications
5. Single use product
6. MR Note
7. Implant selection
8. Application/handling of Implants
9. Pre-operative and post-operative aspects
10. Instruments
11. Decontamination, cleaning and sterilization
12. Liability
13. Signs and symbols

#### 1. GENERAL NOTES

The Medicago system thoracic is supplied unsterile and must be disinfected, cleaned and sterilised prior to use. Please observe the following instructions. This will guarantee that the implant functions faultlessly and reliably.

#### 2. PURPOSE AND PROPER USE / INDICATIONS

The Medicago system thoracic is used for the reconstruction of the anterior spine to ensure the correct alignment and mechanical stability of one or more spine segments (in patients with complete skeletal growth) until the bone segments have healed. It allows both minimally invasive and conventional reconstruction of vertebral height when treating vertebral disk-related acute and chronic instability, stenosis or deformity, and to correct previous unsuccessful arthrodesis (pseudarthrosis). Mechanical fixation of the neighboring vertebral bodies with plate or rod systems is always recommended. The Medicago system thoracic efficiently replaces the support function of a bone graft and has to be filled with autologous bone material or substitute bone material, respectively. By filing the placeholder with small bone chips that are generated during the surgical procedure, a follow-up procedure may not be necessary in some cases.

#### 3. CONTRAINDICATIONS

The Medicago system thoracic may not be used in cases of:

- Patients that are unable to follow the instructions for postoperative care. This may be caused e.g. by psychological/mental or neurological problems of the patient
- Patients with unstable physical and/or mental health conditions
- Patients with inadequate bone tissue, or bone tissue of insufficient quality, and patients with circulatory disorders or latent infections
- Extreme material sensitivity, i.e. foreign-body reaction of the patient. For this case, appropriate tests must be carried out under all circumstances prior to implantation (even when there is only a suspicion of such condition)
- Acute infections

#### 4. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

- In many cases, possible complications are attributable not to the implants as such, but to incorrect patient selection, lack of practice or inadequate positioning of the implants.
  - Inadequate healing.
  - Delayed or insufficient fracture healing leading to possible implant breakage.
  - Breakage, bending or migration of the implant
  - Loosening of the implant
  - Material/foreign-body sensitivity of the patient with allergic reactions
  - Early or late infections
  - Insufficient or delayed bone stiffening
  - Complaints, pains, abnormal sensations or palpability of implants.
  - Degenerative changes or instability in segments that are adjacent to a stiffened spine
- Patients who put unnecessary strain on the implant due to illnesses such as alcoholism or drug abuse during the bone healing phase increase the risks of a sinking of the implant with secondary nerve constriction, malpositions or pseudarthrosis
- Insufficient bone formation, osteolysis, osteomyelitis, osteoporosis, inhibited revascularization or infections that can result in loosening, bending, rupture or breakage of the implants.
- Bone necrosis, osteoporosis and bone resorption can lead to pseudarthrosis

Implants may break if they are subjected to increased stress due to delayed or non-consolidated stiffening

- The use of different materials (e. g. when combined with plate fixation systems) can lead to corrosion
- If bone fusion does not occur, the implant will definitely lose its function. Apart from the side effects and complications mentioned above, the surgical intervention can also lead to problems such as nerve injuries, infections, pain, etc., which are not directly attributable to the implant.

When selecting patients, the following factors have been observed that promote possible adverse events and complications:

- Patient weight: The excess weight or obesity of a patient can increase the load on the implant in such a way that a failure is more likely.
- Profession or activity of the patient: Heavy lifting, strong muscle strain, great rotation bending of the neck, lateral flexion of the head or manual labour should be avoided in private or professional activities until to the bones have healed completely. Even after the patient has healed completely he or she may not be able to successfully perform the activities listed above any longer.
- Senility, mental illness, alcoholism or drug abuse: These circumstances may contribute to the patient ignoring the necessary restrictions or precautions that the implant calls for. This can lead to a failure of the implant or to other complications.
- Certain degenerative diseases
- In some cases a degenerative disease can have progressed so far at the time of implantation that the expected useful life of the pedicle screw and rod system is significantly decreased. In such cases orthopaedic aids can only delay degeneration or achieve a temporary stiltstand.
- Sensitivity to foreign bodies: Preoperative testing may not entirely exclude hypersensitivity or allergy. Such complications can occur even if the implant has already been in the body for some time.
- Smoking: In smokers, a higher rate of pseudarthrosis was found after surgical procedures where bone implants were used. The smoking-in-duced, progressive degeneration of adjacent segments can lead to a later clinical failure (accompanied by intermittent pain), even when an initially successful hardening of the bone has occurred and when clinical improvement has been demonstrated.

#### 5. SINGLE USE PRODUCT

- Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use. The physician should evaluate the risks of an explantation, as he bears responsibility. He or she must also ensure that explanted implants are disposed of according to the standard.

#### 6. MR Note

- **Instruments**  
The use of medical devices in the vicinity of an MR poses a hazard. Individual medical devices must not be located near the equipment during the application of these procedures.

#### Implants

The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MR environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MR devices.

#### 7. IMPLANT SELECTION

The surgeon is responsible for the correct selection of implants.

The following aspects are critical for the correct selection of implants:

- indication,
- body weight of the patient,
- general health condition, age and degree of activity and bone quality of the patient.

Incorrect implant selection can lead to premature implant loss and cause loosening, deformations and fractures of the implants. Failure to use the proper component to maintain adequate blood supply and provide rigid fixation may result in loosening, bending or fracturing of the product and/or bone. Before being implanted, all implants must be checked for function, deformation and mechanical integrity. Damaged implants must be replaced, as they can lead to a reduction in or even loss of function. Damaged implants may not be used and must be replaced. Medicon implants and instruments have been designed for being used together. Using implants or instruments from other manufacturers in conjunction with Medicon products entails unpredictable risks since the products will not be matched to each other. To avoid such risks, only Medicon products that are intended for combination with each other may be used together. This may also cause danger to the patient, the surgeon or third parties.

In cases of failed, delayed or insufficient fracture union, the incidence of implant

bending, implant failure or implant breakage cannot be excluded. Therefore, immobilization of the fracture site must be ensured until there is solid bone fusion. The incessant load changes, to which the implant will be subjected, can lead to fatigue fractures. It cannot be excluded that implants will crack, loosen, corrode, migrate through tissue and cause pain. The packaging label carries a "Lot No.". We recommend copying this "Lot No." to the patient records as this number allows tracing the production history of the implant back to the raw material.

#### 8. APPLICATION / HANDLING OF IMPLANTS

The Medicago system thoracic consists of intervertebral implants of different sizes for the thoracic spine. All components are available in titanium and PEEK. Please refer to the labelling on the packaging. Thanks to the design and the corresponding implant holder, a minimally-invasive implantation is possible in addition to the open, conventional implantation method. The safety and effectiveness of the Medicago system thoracic for spinal disk-related diseases and stenoses or deformations that require stiffening has been proven. For mechanical immobilization up to bone fusion, there should always be an additional internal fixation of the bridged segment. The implants serve to stabilize the spine until vertebral fusion has occurred. Make sure that the implants are positioned correctly with the aid of imaging procedures. Due to the proximity of vascular and neural structures as well as internal organs to the implantation location, the use of an intervertebral disc spacer poses the theoretical risk of a severe bleeding, neural damage or the erosion of internal organs. Such severe complications occur if vessels or internal organs erode, if they are damaged during the procedure/implantation or damaged due to migration after the implantation. The correct positioning of the implant prevents sinking of the cage. Implants can be damaged if the exerted forces are too high, due to which the implant may lose its function. Carving, scratching or bending of the implant during the operation can cause premature failure of the implant. As soon as the bone has healed, this danger is eliminated. Removal of the material is not required.

#### 9. PRE-OPERATIVE AND POST-OPERATIVE ASPECTS

Before the use of the product, the surgeon should thoroughly discuss the desired operation result with the patient. Particular attention must be given to post-operative aspects and a potentially required follow-up care. The patient must be informed that due to the limited stability the strain on the implant caused by body weight should be avoided. He or she must also be warned that non-observance of post-operative instructions may lead to a sinking of the implant, which may then require further surgical procedures to remove the implant. The patient must also be aware that a loss of function of the implant may still occur, even if he or she observes these restrictions. The patient must be instructed to inform the surgeon immediately if they notice unusual changes in the situs. The patient must be monitored carefully if there is a manifest change in the fixation region. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing. Post-operative follow-up care is extremely important. Further risks after an explantation should be avoided by adequate post-operative follow-up care. Until the formation of bone stiffening can be verified on the basis of an x-ray image, the patient should be encouraged to avoid excessive mechanical strain. In cases of a reduced bone quality, an additional external immobilisation may be necessary (e. g. by means of an orthosis).

#### 10. INSTRUMENTS

- Only the instruments that we specially develop and offer for plate systems may be used with the respective systems
- To avoid risks, the Medicago system thoracic may not be inserted using third-party instruments.
- The instruments intended for application of the system are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force.
- To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional.

#### 11. DECONTAMINATION, CLEANING AND STERILIZATION

The following prion-specific protective measure (vCJD) is indicated when preparing the instrument.

In the event of a diagnosis of definite or probable vCJD:

If the use of disposable products is not possible, the instruments used (which had become contaminated or whose contamination cannot be ruled out) must be disposed of as waste for incineration.

#### In suspected cases:

If prion contamination is suspected, it is recommended that any medical instruments be incinerated in accordance with the final report by the nvCJD task force.

#### If vCJD infection is ruled out:

Reuse after final treatment. Otherwise the instruments used (that had become contaminated or where contamination cannot be ruled out) must be disposed of as waste for incineration.

In the event that vCJD infection is not apparent.

Even if there is no known existence of a prion disease, we should, in preparation, set up two procedures with at least partial effectiveness against prion contamination, for example, an alkaline mechanical cleaning process combined with steam sterilisation.

If an alkaline mechanical cleaning process is not used, or another cleaning process with proven effectiveness against prion contamination is deployed, and medical devices that come into contact with risk tissue (CNS, eyes, lymphatic tissue) are used, the Robert Koch Institute (IRK) recommends an extended sterilisation time of 18 minutes at 134°C.

**Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress-corrosion cracking.**

**Cleaning solutions containing hydrogen peroxide and/or rinsing solutions with high alkalinity may result in colour changes. This can result in loss of the coding function.**  
**Only cleaned and disinfected implants may be sterilized.**

#### Limitations on reprocessing

Frequent reprocessing has few effects on these instruments. The end of their working life is normally determined by wear and damage through use.

#### Place of use

Implants with contaminations caused by blood and/or secretions must not be placed back into the implant storage tray. These must be disposed of.

- Directly after use, gross contamination should be removed from the instruments with a disposable cloth/paper.
- Fixing agents or hot water (>40°C) must not be used as this leads to fixation of residues and can affect cleaning outcome.
- The instruments should be conveyed promptly to processing.
- Dry disposal preferred.

#### Preparation for decontamination

Instruments must be removed from the storage tray for cleaning and disinfection. Articulated instruments must be processed in open position.

For rinsing, the instruments and the implant storage tray must be placed in instrument holders suitable for machine cleaning. The instrument holders (e.g. wire trays) must allow subsequent cleaning in an ultrasonic bath or in the cleaning and disinfection device (CDD) where they are not impeded by acoustic or rinsing shadows.

#### Pre-cleaning

- Place implant storage tray for 15 minutes in ultrasonic bath at 40°C with 0.5% alkaline cleaner (Neodisher Mediclean forte) and treat with ultrasonic waves.
- Remove implant storage tray and rinse thoroughly with cold water.
- Place the instruments in cold water for 5 minutes;
- Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
- Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again;
- Place instruments for 15 minutes in ultrasonic bath at 40°C with 0.5% alkaline cleaner and treat with ultrasonic waves;
- Remove instruments and rinse with cold water

A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

#### Cleaning: machine

The cleaning and disinfection device must conform to the requirements of DIN EN ISO 15883-1.

- Prerise 1: 1 minute with demineralized cold water with no additives;
- Drain;
- Prerise 2: 3 minutes with demineralized cold water with no additives;
- Drain;
- Cleaning: use demineralized water, heat to 55°C and wash and clean for 10 minutes, add cleaning agent at 45°C, alkaline cleaning agent, dosage 0.5%;
- Drain;
- Neutralization: 3 minutes with warm water (>40°C) with neutralizer, dosage 1 ml/l;
- Drain
- Final rinse: 2 minutes with warm demineralized water (>40°C) (no additives).
- Drain

#### Disinfection

Thermal disinfection A0 W-value 3000: Fully desalinated water, thermal disinfection is carried out at temperatures >80 °C with a corresponding application time according to A0 concept of standard DIN EN ISO 15883-1 and guideline DGHK, DGSV and AKI (e.g.

English

CE 0123

Rx only

Rev.: 10-12-2025

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A0 3000 = 90 °C and 5 minutes application time). The operator is responsible for the A0 value to be implemented.

**Drying**

Sufficient drying time must be provided by the cleaning and disinfection device. The instruments must be removed from the cleaning and disinfection device immediately after completion of the cleaning and disinfection programme. If necessary, the use of compressed air for drying is recommended because of its effective and quick effect (RKI recommendation).

**Cleaning/disinfection: manual**

The cleansing agents and disinfectants employed must be suitable for manual cleaning and disinfection of instruments and must be compatible with each other. The disinfectant must be of tested efficacy. When selecting the disinfectant agent and method, the relevant lists and recommendations of the Robert Koch Institute (RKI) and the German Society for Hygiene and Microbiology (DGHM) should be noted.

**Pre-cleaning**

- Place the instruments in cold water for 5 minutes;
- Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
- Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again;
- Place instruments for 15 minutes in ultrasonic bath at 40°C with 0.5% alkaline cleaner and treat with ultrasonic waves;
- Remove instruments and rinse with cold water

A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day.

**Cleaning/disinfection**

- Place the instruments in a bath with a tested cleansing and disinfectant agent;
- The instruments must be completely covered with the solution;
- The application times, temperatures, and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed;
- Remove the instruments and rinse for 2 minutes with cold deionised water;
- Repeat the cleansing process if visible contamination is still present on the instrument;
- If necessary, clean and disinfect the instruments in an ultrasonic bath. Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed sooner. National guidelines must be observed.

**Drying**

Dry manually with compressed air and a lint free cloth. Compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

**Maintenance, control and inspection**

After cleaning/disinfection the instruments must be macroscopically clean, i.e., free from visible residues and soiling. Inspection is performed visually. Insufficiently clean instruments must be cleaned again and then adequately rinsed and dried. Before functional testing, instruments with movable parts must have cooled down and been lubricated with instrument maintenance oil (**Medicon REF. 46.00.40**). Instruments such as clamps and needle holders, which have a detent, should only be closed up to the first detent (risk of stress cracking).

Defective instruments (hairline cracks, deformation or wear) must be replaced as they no longer fulfil their function or do so without adequate safety. Corroded instruments must also be replaced as they may corrode intact instruments by extraneous rust.

**Packing**

The implants must be packed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868
- DIN EN ISO 11607
- suitable for steam sterilization (steam permeable)
- adequate temperature resistance to 138°C

Sterilization accessories and sterilization packaging must be appropriate for the content of the package and for the sterilization process.

**Sterilization**

Use the following sterilization procedures subject to the applicable national requirements:

- Fractionated vacuum method with triple fractionation and sufficient product drying
- Steam sterilization in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.

- Sterilization time and temperature: at least 5 minutes holding time at 134°C  
**It is essential that a SAL (Sterility Assurance Level) of 10<sup>-6</sup> is reached.**

**Storage**

Reprocessed sterile instruments must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments. The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority.

**Further informations on reprocessing**

Validated machine cleaning and disinfection is always preferred over manual cleaning because of the greater certainty of the method. Good cleaning helps to preserve value and is a precondition of successful sterilization. During machine processing, the following points should be noted:

- Correct loading of the trays for rinsing is a precondition for effective machine processing. Trays must not be overloaded.
- Rinsing shadows due to large instruments must be avoided.
- The instruments must be placed or stored based on their susceptibility to mechanical damage in order to prevent them from becoming damaged.

The times and temperatures specified in these reprocessing instructions are minimum requirements and must not be less than those stated here. If they are to be reduced for technical reasons, this must be validated by the operating authority. Exceeding the stated times and temperatures is always possible but leads to increased stress on the material, which may result in premature ageing of the instruments. The use of other sterilization methods is outside our responsibility.

**Information on validating the reprocessing procedure**

Validation was performed with the following equipment, materials, and chemicals:

Cleaning and disinfection device: Type Miele Disinfector G 7735 CD and 7836 CD  
 compartmented cart for surgical instruments  
 Cleaning agent alkaline: neodisher® FA, Dr. Weigert GmbH & Co. KG  
 Cleaning agent enzymatic: Endozime, Ruhof (enzymatic)  
 Decondex 23 Neutrazym, Borer, Switzerland  
 Neutraliser: neodisher® Z, Dr. Weigert GmbH & Co. KG  
 Water pistol: Selecta  
 Cleaning brush: Plastic/nylon brushes  
 Ultrasonic bath: Sonorex  
 Steriliser: MMM Vaculab 969 S 3000, MMM Selectomat S 3000, Stiefenhofer KS 666-2ED, H+P Varioclav 400E

**Note**

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method. If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device. Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

**12. LIABILITY**

In the event of discrepancies between the english and the german version of these instructions for use only the German version shall be applicable. Only the latest revision of the instructions for use applies. Due to constant technical development the contents of these MEDICON instructions for use are updated regularly.

The current version date is printed on each edition of the instructions for use. MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use, non-observance of post-operative instructions and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

**13. SIGNS AND SYMBOLS**

Manufacturer

Batch number

Article number

Warning

Follow instructions

MR unsafe

CE label accord. To directive 93/42/EEC

Non reusable

**Rx only**

Prescription use only

Please contact MEDICON eG if you have further questions concerning the implants or instruments.

**CAUTION:**  
 For USA, federal law restricts this device to purchase by or on behalf of physician or hospital.

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