**Sagittale Säge nach SACHSE**  
**Nr. 68.24.02**

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**Sagittal Saw according to SACHSE  
No. 68.24.02**

For operating with motors with a coupling  
according to DIN 13940/ISO-DIS 3964

**Contents**

- Use
- Fitting the saw blades
- Cleaning
- Sterilization
- Storage
- Ambient conditions
- Service
- Warranty terms

**Use**

Sagittal saws are suitable for use in the following areas:

- Oral and maxillofacial surgery
- Implantology
- ENT surgery
- Plastic surgery
- Hand surgery
- Neurosurgery

The sagittal saw according to Sachse is a state-of-the-art micro saw. This saw can be operated with our micromotor system Servotronic EC 100 / 200 or with other motors which have an INTRA coupling (DIN 13940/ISO-DIS 3964). The saw also has an integrated irrigation tube which should always be connected. Saw blades are available in four sizes with various shaft lengths.

This product is re-usable and is supplied unsterile. Before using for the first time, the saw must be cleaned and sterilized according to the following instructions.

**This product should only be used by trained staff with specialist knowledge. The manufacturer's information and operating instructions must be observed for the micromotor system in use.**

**Fitting saw blades  
Inserting the saw blade**

**The motor should never be operated while inserting or changing the saw blade. As a precaution, we recommend disconnecting the saw handpiece from the motor.**

1. The cap nut of the clamping unit must be loosened before inserting.
2. When inserting the saw blade, observe the marking of the handpiece.
3. Tighten the cap nut carefully with the open-end wrench provided, until the saw blade sits tightly in the fixture.

**Only operate the drive motor clockwise, with a maximum speed of 20,000 rpm. In order to avoid any heat necrosis of the bone, we recommend continuous flushing throughout the entire use. Replace any worn saw blades in good time. We recommend that a new blade is always inserted before use.**

**The best result is achieved at motor speeds of at least 15,000 rpm. The saw has to be activated before it is put on the bone.**

**Cleaning**

Preparation for decontamination

The saw handpieces must be placed on an instrument carrier suitable for machine-based cleaning processes. The instrument carrier (e.g. strainer basket) must be designed so that subsequent cleaning in cleaning or disinfection equipment is not impaired by areas inaccessible to water. Before cleaning, remove the blade and cover.

**Cleaning**

Suitable alkaline cleaners or those with a neutral pH can be used. Selection of the cleaning agent is based on the properties of the saw handpieces as well as on national guidelines and recommendations. Elevated chloride concentrations in the water could cause pitting or stress corrosion cracking on the saw handpieces. Such corrosion can be minimized by using demineralized water or alkaline cleaning agents.

**Never clean the saw handpieces in an ultrasonic bath. This would have serious consequences on the functioning of the saw handpieces as well as on the safety of both user and patient.**

**Thermal disinfection**

Using demineralized water, thermal disinfection is carried out at temperatures from 80-95°C, with a corresponding soaking time according to the A0 concept, EN ISO 15883-1.

**Drying**

Sufficient drying must be ensured using the cleaning and disinfection device.

**Cleaning agents**

The parameters specified by the manufacturer of the cleaning agent for the concentration, temperature and soaking time must be complied with and it must be possible to check automatic dosing systems.

**Maintenance, inspection and testing**

Following cleaning/disinfection, the saw handpieces must be macroscopically clean, i.e. free of visible residue. This is checked by visual inspection. Critical areas such as joints and screw fixtures require especially careful checking. Inadequately cleaned or clogged saw handpieces must be cleaned again and must then be rinsed sufficiently.

Before the functional test, saw handpieces must have cooled down and have been lubricated with spray oil no. 63.69.95.

Defective saw handpieces (showing hairline cracks, distortion or wear) must be replaced because their functionality can no longer be guaranteed or can no longer be guaranteed sufficiently.

**Lubrication**

Lubricate before each sterilization or at least twice daily using spray oil no. 63.69.95.

1. Remove tools from the collet before lubricating.
2. Spray the inside and outside of the saw handpiece in order to remove the residue.
3. Cleaning and disinfection should be carried out using a clean damp cloth or else disinfectant cloths. We do not recommend any disinfectant products which contain acetone, chlorine or bleach.
4. Insert the spray can with spray adapter into the back of the saw handpiece and hold the aerosol can button depressed for approx. 1 second.

**Sterilization**

Only cleaned and disinfected saw handpieces may be sterilized.

Only the sterilization methods listed may be used for sterilization.

- Fractionated vacuum method with sufficient product drying
- Autoclaving in accordance with DIN EN 13060 or DIN EN 13485 and DIN EN 554, as well as other processes validated in accordance with ISO 14937:2000 and USP 23.
- Sterilization time and temperature:  
MEDICON eG recommends as follows:  
At least 5 minutes holding time at 134°C

After cleaning, disinfection and sterilization, insert a saw blade and leave the handpiece running for 10 to 15 seconds at a low speed in order to distribute the lubricant and remove any excess.

**Stains or discoloration may be caused by inadequate cleaning before sterilization or else by mineral deposits in the water used in autoclaving.**

**Make sure that the sterilizer and the water used are clean. After each sterilization cycle, remove the product from the sterilizer immediately in order to reduce the risk of corrosion.**

**Storage**

Saw handpieces can corrode due to the storage conditions. To prevent this, saw handpieces must be stored dry and protected from dust. To prevent condensation from forming, major temperature fluctuations should be avoided. Chemicals must not be stored with saw handpieces. The storage period depends on the type of packaging.

**Ambient conditions**

Temperature range – 40°C to 70°C, relative humidity between 10% and 100 %, atmospheric pressure between 50 kPa and 106 kPa.

**Service**

Never dismantle the saw handpieces. For all checks and repairs, we recommend contacting the supplier or MEDICON eG directly. MEDICON eG advises the user to have saw handpieces checked or overhauled at least once a year.

**Warranty terms**

MEDICON eG guarantees that the product is free of functional, material and processing defects. This warranty for the device extends for 12 months from the date of purchase.

In the event of justified complaints, MEDICON eG is responsible for repairing the product or replacing it free of charge. Any other claims of whatsoever kind are excluded, particularly any claims for damages.

MEDICON eG is not liable for defects, injuries and the consequences of the following:

- excessive wear and tear
- incorrect use
- disregarding the operating, assembly and care instructions
- unusual chemical, electrical and electrolytic effects
- defective connections for air, water and electricity

The warranty does not extend to fiber optic conductors or any plastic parts. The warranty becomes void if the damage and consequential damage occurs as a result of improper intervention or modifications to the product by any unauthorized third party. Warranty claims may only be asserted if a copy of the invoice or delivery note is presented with the product, clearly indicating the date of purchasing the product, its reference number and serial number.

**Medicon eG, as the company marketing these products, assumes no liability for direct or subsequent damage resulting from improper use or handling, in particular as a result of failure to observe the instructions for use or through incorrect care or maintenance.**

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