Winkelschraubendreher Nr. 68.30.01

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Contra-angle Screwdriver No. 68.30.01

Contents

- Use
- Inserting the screwdriver inserts
- Cleaning
- Sterilization
- Storage
- Ambient conditions
- Service
- Warranty terms

Use

For plate internal fixation in poorly accessible mandibular and midfacial regions, instead of transbuccal/transfacial soft tissue perforation directly through the intraoral, supra/intraorbital or coronal approach to the fracture. The contra-angle screwdriver is a product which has been produced with the utmost care. The contra-angle screwdriver can be operated with our micromotor system Servotronic EC 100 / 200 or with other motors which have an INTRA coupling (DIN 13940/ISO-DIS 3964).


Use the electromotor system solely for predrilling screw holes. For screwing in the screws, use the manual drive 68.30.05. We recommend using a maximum speed of 15,000 rpm for predrilling the screw holes.

This product is re-usable and is supplied unsterile. Before using for the first time, the contra-angle screwdriver must be cleaned and sterilized according to the following instructions.

This product should only be used by trained staff with specialist knowledge. The manufacturer's information and operating instructions must be observed for the contra-angle screwdriver in use.

Inserting the screwdriver inserts

Insert the right shape of screw insert into the holding device of the contra-angle screwdriver and press it in hard with your finger. The screw insert can also be pressed against an even surface. To remove the screw insert, proceed as for insertion, but from the rear.

Cleaning

Preparation for decontamination

The contra-angle screwdriver must be placed on an instrument carrier suitable for machine-based cleaning processes. The instrument carrier (e.g. strainer basket) must be designed so that subsequent cleaning in the cleaning and disinfection device is not impaired by areas inaccessible to water. Before cleaning, remove the blade and cover.

Cleaning

Suitable alkaline cleaners or those with a neutral pH can be used. Selection of the cleaning agent is based on the properties of the contra-angle screwdriver as well as on national guidelines and recommendations. Elevated chloride concentrations in the water can cause pitting or stress corrosion cracking on the contra-angle screwdriver. Such corrosion can be minimized by using demineralized water or alkaline cleaning agents.

Never clean the contra-angle screwdriver in an ultrasonic bath. This would have serious consequences on the functioning of the contra-angle screwdriver as well as on the safety of both user and patient.

Thermal disinfection

Using demineralized water, thermal disinfection is carried out at temperatures from 80-95°C, with a corresponding soaking time according to the A0 concept, EN ISO 15883-1.

Drying

Sufficient drying must be ensured using the cleaning and disinfection device.

Cleaning agents

The parameters specified by the manufacturer of the cleaning agent for the concentration, temperature and soaking time must be complied with and it must be possible to check automatic dosing systems.

Maintenance, inspection and testing

Following cleaning/disinfection, the contra-angle screwdriver must be macroscopically clean, i.e. free of visible residue. This is checked by visual inspection.

Critical areas such as joints and screw fixtures require especially careful checking. Inadequately cleaned or clogged instruments must be cleaned again and must then be rinsed sufficiently.

Before the functional test, the contra-angle screwdriver must have been allowed to cool down and have been lubricated with spray oil no. 63.69.95.

Defective contra-angle screwdrivers (showing hairline cracks, distortion or wear) must be replaced because their functionality can no longer be guaranteed or can no longer be guaranteed sufficiently.

Lubrication

Lubricate before each sterilization or at least twice daily using spray oil no. 63.69.95.

1. Spray the contra-angle screwdriver inside and outside with our spray oil on a lint-free cloth in order to remove any residue.
2. Cleaning and disinfection should be carried out using a clean damp cloth or else disinfectant cloths. We do not recommend any disinfectant products which contain acetone, chlorine or bleach.
3. Insert the spray can with spray adapter into the back of the saw handpiece and hold the aerosol can button depressed for approx. 1 second.
4. Place the contra-angle screwdriver on its back on a lint-free cloth to allow any excess oil to drip off.

Sterilization

Only cleaned and disinfected attachments may be sterilized.

Only the sterilization methods listed may be used for sterilization.

- Fractionated vacuum method with sufficient product drying
- Autoclaving in accordance with DIN EN 13060 or DIN EN 13485 and DIN EN 554, as well as other processes validated in accordance with ISO 14937:2000 and USP 23.
- Sterilization time and temperature:
MEDICON eG recommends as follows:
At least 5 minutes holding time at 134°C

After cleaning, disinfection and sterilization, leave the contra-angle screwdriver running for 10 to 15 seconds at a low speed in order to distribute the lubricant and remove any excess.

Stains or discoloration may be caused by inadequate cleaning before sterilization or else by mineral deposits in the water used in autoclaving.

Make sure that the sterilizer and the water used are clean. After each sterilization cycle, remove the product from the sterilizer immediately in order to reduce the risk of corrosion.

Storage

The contra-angle screwdriver can corrode as a result of storage conditions. To prevent this, it must be stored dry and protected from dust. To prevent condensation from forming, major temperature fluctuations should be avoided. Chemicals must not be stored with the contra-angle screwdriver. The storage period depends on the type of packaging.

Ambient conditions

Temperature range – 40°C to 70°C, relative humidity between 10% and 100 %, atmospheric pressure between 50 kPa and 106 kPa.

Service

Never dismantle the contra-angle screwdriver. For all checks and repairs, we recommend contacting the supplier or MEDICON eG directly. MEDICON eG advises the user to have the contra-angle screwdriver checked or overhauled at least once a year.

Warranty terms

MEDICON eG guarantees that the product is free of functional, material and processing defects. This warranty for the device extends for 12 months from the date of purchase.

In the event of justified complaints, MEDICON eG is responsible for repairing the product or replacing it free of charge. Any other claims of whatsoever kind are excluded, particularly any claims for damages.

MEDICON eG is not liable for defects, injuries and the consequences of the following:

- excessive wear and tear
- incorrect use
- disregarding the operating, assembly and care instructions
- unusual chemical, electrical and electrolytic effects
- defective connections for air, water and electricity

The warranty does not extend to fiber optic conductors or any plastic parts. The warranty becomes void if the damage and consequential damage occurs as a result of improper intervention or modifications to the product by unauthorized third parties. Warranty claims may only be asserted if a copy of the invoice or delivery note is presented with the product, clearly indicating the date of purchasing the product, its reference number and serial number.

Medicon eG, as the company marketing these products, does not assume any liability for direct or subsequent damage resulting from improper use or handling, in particular as a result of failure to observe the instructions for use or through incorrect care or maintenance.

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