Deutsch

CE 0123

Rx only

Rev.: 04-12-2025

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Rev.: 04-12-2025

Instructions for Use

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

Table of contents

1. General information
2. Intended use
3. Indications
4. Contraindications
5. Possible side effects and complications
6. Single use product
7. MR Notes
8. Application and handling
9. Pre- and postoperative behaviour
10. Decontamination, cleaning and sterilization
11. Liability
12. Signs and symbols

1. GENERAL INFORMATION

Medicon implants are manufactured from pure titanium or Titanium alloy. Both materials are biocompatible, corrosionproof and nontoxic in the biological environment.

The surface is chemically passive and the material is antimagnetic. The micro and mini implants are supplied unsterile and must be disinfected, cleaned and sterilised prior to use.

Please observe the following instructions for use/preparation. This will guarantee that the implant functions faultlessly and reliably.

2. INTENDED USE

The micro and mini implants that consist of plates and screws are intended for application in case of operations on the bony skull structures. The implants may be used only by surgeons with sufficient training and experience in oral, orthodontic and maxillofacial surgery in clinics and practices. The implants have been designed for supporting the normal bone healing process, not for replacing normal body structures or bearing the body weight in cases of insufficient bone healing. MEDICON implants are disposable products and may be used only together with the specified MEDICON combination and application instruments.

3. INDICATIONS

- fractures
- osteotomies, including orthognathic indications
- revisions in cases where other treatment methods or devices proved unsuccessful
- fixation and stabilisation of bone transplants for the alveolar ridge
- temporary fixation of suitable non-resorbable bone replacement materials for the alveolar ridge. There are no known interactions between non-resorbable bone replacement materials and our products.

4. CONTRAINDICATIONS

- Patients who are unable to follow the instructions for postoperative care. Such inability may be caused by psychological/mental or neuro-logical conditions of the patient.
- Patients with unstable physical and/or mental health conditions.
- Patients with insufficient or poor-quality bone tissue, circulatory disorders or latent infections.
- Oversensitivity to materials and foreign-body reactions. Preimplantation tests must be carried out even if such oversensitivity or reactions are only suspected.
- Acute infections.

5. POSSIBLE SIDE EFFECT AND COMPLICATIONS

- Insufficient adaptation of the implants can lead to bone union delay or failure.
- Complaints, pains, abnormal sensations or palpability of implants.

- Material/foreign-body sensitivity of the patient with allergic reactions.
- Stronger connective-tissue reaction in the fracture region or in the vicinity of the implant.
- Delayed or insufficient fracture healing leading to possible implant breakage.
- Inadequate healing.
- Insufficient bone formation, osteolysis, osteoporosis, reduced revascularization or infection, which can lead to loosening, bending, breakage or rupture of the implant.
- Breakage, bending, migration or loosening of the implant.
- Reduced bone density due to "stress shielding".
- Loosening of the implant as a result of insufficient screw tightening.
- Osseous necrosis, osteoporosis, inhibited revascularization, bone-resorption, and poor bone formation, which can cause premature loss of fixation leading to non-union.
- Extreme or repeated bending of the plates in the same place can result in plate breakage.
- Unstable comminuted fractures can lead to an increased connective tissue reaction in the fracture region.

This surgical procedure may cause not only the above-mentioned side effects and complications but also problems such as injuries to nerves, infections, pain etc., which are not necessarily caused by the implant.

6. SINGLE USE PRODUCT

Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist.

Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

7. MR NOTES

Instruments

The use of medical devices in the vicinity of an MR poses a hazard. Individual medical devices must not be located near the equipment during the application of these procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MR environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MRI devices.

8. APPLICATION AND HANDLING

Selection of implants

The surgeon is responsible for the correct selection of implants.

The following aspects are critical for the correct selection of implants:

- bone defect to be treated / intended bone repositioning,
- body weight of the patient,
- general health condition, degree of activity of the patient.

Incorrect implant selection can lead to premature implant loss and cause loosening, deformations and fractures of the implants. Before being implanted, all implants must be checked for function, deformation and mechanical integrity.

Scratches and other damage can reduce the stability of the product and result in early fatigue of the implant. If there is no complete recovery, or if recovery is delayed or inadequate, the possibility that the implant will bend, break or become loosened cannot be ruled out. The incessant load changes, to which the implant will be subjected, can lead to fatigue fractures. It cannot be excluded that implants will break, loosen, corrode, migrate through tissue and cause pain.

Medicon implants and instruments have been designed for being used together. Using implants or instruments from other manufacturers in conjunction with Medicon products entails unpredictable risks since the products will not be matched to each other. Patients, users and/or third parties can be endangered as a result. The surgeon is responsible for deciding whether the implants shall remain in the body after complete healing of the fracture or whether they shall be removed.

Application and handling of plates

The plates need to be adapted to the natural or planned bone contours as precisely as possible. Only use the instruments provided by the manufacturer for the purpose of adapting the plates.

During the plate bending procedure the plates are subjected to cold stresses, under which the titanium hardens and its ductility is decreased. Therefore it is absolutely necessary to achieve the required plate shape in as few bending steps as possible. Excessive bending can lead to postoperative plate breakage.

The success of an operation depends among other things on how the instruments are handled. The plates must not be bent carefully and cautiously, avoiding repeated and excessive bending. The plates must not be bent back and forth more than 2 – 3 times. Damage leads to a reduction in the product's strength and to premature fatigue.

Overaggressive application of the bending instruments can result in visible damage to the implant. When such damage has occurred, a new plate must be selected and adapted more cautiously. Deformed plate holes for the screws entail not only an increased risk of breakage at the holes; they also compromise the precise seating of the screw head in the plate. Therefore the plates must be bent cautiously. Adapted plates must be inspected for notches, deformed screw holes and other mechanical damage before they are applied to the patient.

Always select the plates that provide the closest fit to the specific situation and indication. If a plate needs to be shortened, only use the instruments supplied for this purpose. When shortening a plate, care must be taken that the portion of the plate that is cut off is not flung towards the patient or towards third parties during the cutting procedure.

The cut edges of a shortened plate must be deburred to prevent tissue injuries or irritations. Bending templates must not be implanted under any circumstances.

Application and handling of screws

Selftapping screws do not require precutting of threads. To apply the drill holes, only use the drills designated by Medicon for the respective screw so that the correct drill hole diameter for the screw is obtained. Only use drill bits with sharp cutting edges! During drilling, there is the risk of heat damage to the bone. Therefore, only work at low drill speed and always ensure sufficient cooling during the drilling procedure. Excessive application of force during the drilling procedure can lead to drill breakage, endangering the patient, the surgeon and third parties. Use the emergency screws intended for this purpose if, due to a failed drilling effort or a worn drill hole, the original screw cannot be anchored firmly in the bone. To ensure that the screwdriver blade is firmly seated in the screw head, axial pressure must be applied when the screwdriver is inserted into the screw head. This ensures the correct axial alignment of screw and screwdriver, and prevents slipping of the screwdriver blade and damage to the screw head.

Application of excessive force when tightening the screws can lead to intra-operative screw breakage. Wear of the screwdriver blade will impair the firm connection between blade and screw. Whenever such wear is detected, the screwdriver blade must be replaced by a new blade. Once the implantation has been completed, every screw and plate connection must be checked for stability. If necessary retighten the screws. Prior to implant removal, all screw heads must be cleaned carefully, using a scalpel or another suitable instrument, so that the screwdriver can be optimally seated in the screw head.

Self-drilling screws

It is normally not necessary to drill a pilot hole for self-drilling screws. Still, in dense and hard bones, e.g. at the calvarium or in the mandible, it may be necessary to prepare drill holes. Self-drilling screws are unsuitable for small bone fragments as the axial pressure on them can lead to fragment shifts or screw breakage.

Application and handling of instruments

The instruments intended for application of the system are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force. To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional. Use only accessories specifically approved by Medicon in order to avoid risks in connection with the compatibility of accessory parts with the products.

9. PRE- AND POSTOPERATIVE BEHAVIOUR

Preoperative

Prior to applying the products, the surgeon must thoroughly discuss with the patient the intended outcome of the operation. In such discussions, special attention must be paid to postoperative care measures that may become necessary.

Postoperative

The patient must be instructed to immediately inform the treating physician about unusual in situ changes. The treating physician has to consider possible implications such as implant

failure and discuss required measures for further healing with the patient. The surgeon is responsible for deciding whether the implants shall remain in the body after complete healing or whether they shall be removed.

For his decision, the surgeon has to weigh risks and advantages of leaving the implants in the body. It cannot be excluded that implants will break, loosen, corrode, migrate through tissue and cause pain. In cases of failed, delayed or insufficient bone union, the incidence of implant bending, implant failure or implant breakage cannot be excluded.

Therefore, immobilization of the fracture site must be ensured until there is solid bone fusion. The incessant load changes, to which the implant will be subjected, can lead to fatigue fractures. The treating physician takes decisions regarding the type, duration and intensity of physical activities after the surgery.

The packaging label carries a "Lot No.". To ensure 100% traceability, we recommend noting down the LOT number in the patient records.

Disposal

Observe the applicable national guidelines when disposing of the MEDICON system and accessories.

10. DECONTAMINATION, CLEANING AND STERILIZATION

The following vCJD prion-specific protective measure is indicated when processing instruments.

In the event of diagnosis of a definite or probable case of vCJD.

If it is not possible to use disposable products, the instrument used, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

If prion contamination is suspected:

In the event of suspected prion contamination, incineration of the instrument is recommended according to the vCJD task force final report.

If vCJD is excluded

Continue to use after instrument processing is completed. Otherwise, the instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

In the case of a non-identifiable vCJD illness.

Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions – e.g. mechanical alkaline cleaning combined with steam sterilization.

If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (CNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134°C.

Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress-corrosion cracking. Only cleaned and disinfected implants may be sterilized.

Reprocessing restrictions

Frequent reprocessing has a slight impact on implants

Location of use

Implants that have been contaminated by blood and/or secretions must not be placed back in the implant storage tray. They must be disposed of.

Preparation for decontamination

The implant storage tray must be placed in instrument holders suitable for machine cleaning. The instrument holders (e.g. wire trays) must allow subsequent cleaning in an ultrasonic bath or in the cleaning and disinfection device (CDD) where they are not impeded by acoustic or rinsing shadows.

Pre-cleaning

For the protection against crosscontamination during manual cleaning in an ultrasonic bath, a fresh and unused cleaning agent must always be used.

- Place implant storage tray for 15 minutes in ultrasonic bath at 40°C with 0.5% alkaline cleaner (Neodisher Mediclean forte) and treat with ultrasonic waves.
- Remove implant storage tray and rinse thoroughly with cold water.

Rev.: 04-12-2025

Cleaning: machine

For cleaning and disinfection in a CDD and for the protection against crosscontamination during automatic cleaning, the implant tray with the implants that are intended to be reprocessed, must always be separated from instruments which are contaminated with blood or secretion.

The cleaning and disinfection device (CDD) must meet DIN EN ISO 15883-1 requirements.

- Pre-rinse: 3 minutes with cold demineralised water, without additive;
- Emptying;
- Cleaning: with demineralised water, heat to 55°C and wash/clean for 5 minutes, add cleansing agent at 45°C, alkaline cleansing agent, strength 0.5%;
- Emptying.
- Neutralization: 3 minutes with warm water (>40°C) with addition of neutralizer, strength 1ml / l;
- Emptying;
- Final rinse: 2 minutes with warm deionised water (>40°C) (without any additive);
- Emptying

Disinfection

Thermal Disinfection A0 value 3000: Completely desalinated water, the thermal disinfection is carried out at temperatures > 80 °C and corresponding application time according to the A0 concept, DIN EN ISO 15883 and guideline DGKH, DGSV and AKI (e.g. A0 3000 = 90°C and 5 minutes application time). The operator is responsible for the implemented A0 value.

Drying

Complete drying must be ensured by the CDD, e.g. 60°C and 30 minutes. The instruments must be removed from the CDD promptly once the cleaning and disinfection program has finished. If necessary, compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Maintenance, inspection, and testing

After cleaning / disinfection the implants must be macroscopically clean, i.e., free from visible residues and soiling. Inspection is performed visually. Insufficiently cleaned implants must be cleaned again and then adequately rinsed and dried. Deformed or damaged implants must be removed and disposed of, as they can no longer be guaranteed to be safe to use.

Packing

The implant storage tray must be placed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868,
- DIN EN ISO 11607,
- suitable for steam sterilization (vapor permeable),
- adequate temperature resistance up to 138°C.

Sterilization equipment and sterilization wrapping must match both the wrap contents and the employed sterilization method.

Sterilization

Taking into account the respective national regulations, the following method must be employed for sterilization:

- Vacuum autoclave with triple vacuum and adequate drying of the products (Vacuum minimal 15 Minutes)
- Steam sterilizer complying with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes hold time at 134°C

It is essential to attain a Sterility Assurance Level of 10⁻⁶.

Storage

Reprocessed sterile implants must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semisterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments.

The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority.

Further information on reprocessing

Validated machine cleaning and disinfection is always preferred over manual cleaning because of the greater certainty of the method. Good cleaning helps to preserve value and is a precondition of successful sterilization. During machine processing, the following points should be noted:

- Correct loading of the trays for rinsing is a precondition for effective machine processing. Trays must not be overloaded.
- Rinsing shadows due to large instruments must be avoided.

The times and temperatures specified in these reprocessing instructions are minimum requirements and must not be less than those stated here. If they are to be reduced for technical reasons, this must be validated by the operating authority.

Exceeding the stated times and temperatures is always possible but leads to increased stress on the material, which may result in premature ageing of the instruments. The use of other sterilization methods is outside our responsibility.

Information on validation of the processing

Validation was performed with the following equipment, materials, and chemicals:
Cleaning and disinfection device: Type Miele Disinfector G 7735 CD and 7836 CD compartmented cart for surgical instruments
Cleaning neodisher® Fa. Dr. Weigert GmbH & Co. KG
Neutraliser: neodisher® Z, Dr. Weigert GmbH & Co. KG
Ultrasonic Sonorex DIGITEC DT 156 BH

bath: MMM Vaculab 969 S 3000,
Steriliser: MMM Selectomat S 3000,
Stiefenhofer KS 666-2ED, H+P Varioclav 400E

Note

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method.

If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device. Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

11. LIABILITY

In the event of discrepancies between the english and the german version of these instructions for use only the German version shall be applicable. Only the latest revision of the instructions for use applies. The current version date is printed on each edition of the instructions for use.
MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use, non-observance of post-operative instructions and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

CAUTION:
For USA, federal law restricts this device to purchase by or on behalf of physician or hospital.

12. SIGNS AND SYMBOLS

Symbol	Meaning
	Manufacturer
	Batch number
	Article number
	Non sterile
	Warning
	Follow instructions

	MR unsafe
	CE label accord. To directive 93/42/EEC
	Non reusable
	Prescription use only

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