Rev.: 10-12-2025

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 Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

MEDICON IMF screws are made of titanium. The material is biocompatible, corrosionresistant and non-toxic in the biological milieu. The surface is chemically passive and the material is antimtic. IMF screws are supplied unsterile.

IMF screws are available in the following designs:

	IMF screws with side-drilled holes and cross slotted head,		
	self-tapping	Cat. no.	Thread length
		68.94.30B	8 mm
		68.94.31B	10 mm
	Colour: Grey	68.94.32B	12 mm
	MF screws with side-drilled holes and cross slotted head,		
	self-tapping	Cat. no.	Thread length
		68.94.26B	10 mm
		68.94.27B	12 mm
	Colour: Green	68.94.28B	8 mm
	IMF-Schrauben mit Querbohrungen und Kreuzschlitzkopf, Low Profile		
	self-drilling	Cat. no.	Thread length
		68.94.38B	8 mm
	Colour: Green	68.94.41B	11 mm

Contents
1. General informations
2. Purpose and proper use/Indications
3. Contraindications
4. Possible side effects and complications
5. Single use product
6. MR instructions
7. Using and handling the IMF screws
8. Technical instructions
9. Warnings and precautions
10. Decontamination, cleaning and sterilization
11. Liability
12. Signs and symbols

1. GENERAL INFORMATIONS
The IMF screws are supplied unsterile and must be disinfected, cleaned and sterilised prior to use. Please observe the following instructions. This will guarantee that the screws function faultlessly and reliably.

2. PURPOSE AND PROPER USE/INDICATIONS
MEDICON IMF screws should only be used by surgeons with sufficient training and adequate experience in oral, orthodontic and maxillofacial surgery. IMF (intermaxillary fixation) screws are designed for temporary fixation of the mandible and maxilla and are used for perioperative stabilisation of simple fractures of the mandible and maxilla and for orthognathic procedures.
Indications include:

- Temporary fixation of simple, non-dislocated fractures of the mandible or maxilla
- Temporary fixation with orthognathic procedures
- Temporary support fixation for bone healing

The use of osteosynthesis plate systems in orthodontic surgery means that less time is generally required for intermaxillary fixation for patients with fractures of the mandible and maxilla. However, temporary jaw fixation during the operation is often required to provide support while repositioning and for establishing

normal occlusion. At the same time elastic tension immediately after the operation may be useful for correcting minor occlusal irregularities.

3. CONTRAINDICATIONS

- Combined fractures of the mandible and maxilla
- Dislocated fractures

- Multiple fractures and comminuted fractures
- Segmented dentition
- Dentition with reduced stability
- Paediatric surgery
- Patients, who are incapable of following instructions for postoperative care.
- Patients with: inadequate or reduced-quality bone tissue, acute and/or latent infection, circulatory disorders, metal allergy, reactions to foreign bodies and/or mental, physical or neurological disorders.

- 4. POSSIBLE SIDE EFFECTS AND COMPLICATIONS**
- Complaints, pain, abnormal sensations or palpability of the IMF screws.
 - Hypersensitivity of the patient due to the material of the foreign bodies in the form of allergic reactions.
 - Increased reaction of the connective tissue in the region of the fracture and/or the IMF screws.
 - Insufficient bone formation, osteolysis, osteoporosis, inhibited revascularization or infection that could cause loosening, bending or fracture of the IMF screws.
 - Fracture, bending or loosening of the IMF screw.

Apart from the above side effects and complications, the surgical procedure may cause other problems such as injury to nerves, infections and/or pain, which are not necessarily due to the IMF screws. If complications occur, they are often the result of incorrect patient screening, lack of practice or incorrect repositioning of the bone fragments rather than caused by the IMF screws themselves.

- 5. SINGLE USE PRODUCT**
-  Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

- 6. MR INSTRUCTIONS**
- Instruments**
-  The use of medical devices in the vicinity of an MR poses a hazard. Individual medical devices must not be located near the equipment during the application of these procedures.

Implants
The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MR environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MRI devices. the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device.
Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

7. USING AND HANDLING THE IMF SCREWS
Selection of the suitable IMF screws depends on the actual anatomical situation and is the responsibility of the surgeon. Before use the surgeon must make sure that the IMF screws are the correct design for the planned procedure. A lot number is printed on the package label.
We recommend recording the lot number in the patient's surgical record and in the patient's files, because the manufacturing process can be traced with this number. IMF screws are designed and constructed for single use only on a single patient and may not be reused. An explanted IMF screw must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist.

7.1 Self-tapping IMF screws
It is not necessary to tap a thread for self-tapping IMF screws. **Follow the instructions for drilling.**

7.2 Self-drilling IMF screws
It is usually not necessary to drill a hole for self-drilling IMF screws. It may however be necessary to drill a pilot hole in dense, hard bone. Follow the instructions for drilling for these cases.

7.3 Drilling
Use the specified drill bits with a red ring provided by MEDICON to drill the holes, to ensure that the hole is the correct diameter. Use only drill bits with sharp cutting edges. Blunt bits must not be used.

 **There is a danger during drilling of heat damaging the bone. Make sure that adequate cooling is provided at all times. Excessive force during drilling may cause the drill bit to break, causing danger to the patient, the surgeon or others.**

7.4 Screwdrivers
Use only the screwdrivers recommended for the specific IMF screws. The screwdriver must be applied to the IMF screw head with axial pressure to ensure that the blade is firmly seated in the head and that the IMF screw and screwdriver are correctly aligned.Signs of wear on the screwdriver blade adversely affect the firm connection between blade and IMF screw. Check the blade every time the screwdriver is used. Blades that cannot be seated firmly in the IMF screw head must not be used.

- 8. TECHNICAL INSTRUCTIONS**
- The junction between the attached and reflected mucosa has been found to be the most suitable position for the screw. IMF screws can be screwed in at any point on the dental arch, but preferably between the canine and the first premolar. One IMF screw per quadrant is normally sufficient, but more can be added if necessary.
 - When drilling make sure that the hole is between the roots of the teeth and that the roots are not damaged. If the drill touches a root this may make it more difficult to insert the IMF screw and may result in complications.
 - Screw in the IMF screw until the bottom of the screw head is in contact with the buccal mucosa. The IMF screw must screw in without resistance and without requiring force. If resistance is encountered, unscrew the IMF screw until the resistance ceases. If the resistance is not reduced, select a shorter IMF screw.
 - To prevent irritation of the soft tissue, the tip of the IMF screw must not penetrate the palatal or lingual mucosa.
 - The hole must be drilled at right angles to the mucosa and can be drilled through without requiring an incision of the gingiva.
 - Self-drilling IMF screws can also be screwed in without an incision.
 - Use stainless steel wires with a diameter of 0.4 mm to 0.6 mm for intermaxillary fixation. We recommend implant quality wires. Elastic bands can be used for tension.
 - After plating the fracture, the intermaxillary fixation is removed and the occlusion checked. If it can be seen that the fracture will heal correctly and the patient will retain normal occlusion, the IMF screws can be removed. This can be carried out as an outpatient procedure without requiring local anaesthetic.

- 9. WARNINGS AND PRECAUTIONS**
-  **The surgeon is responsible for correct selection of the IMF screws.**
- IMF screws must not be reused.
 - The surgeon must discuss in detail the planned operation with IMF screws and their possible risks with the patient, and inform the patient of what to do after the operation.
 - The patient must be instructed in postoperative oral hygiene and told to inform the surgeon of any changes at the surgical site.
 - The physician or surgeon in charge must then define the next stages and discuss them with the patient.

IMF screws are designed for use with the associated instruments. The use of IMF screws and instruments from other manufacturers with the MEDICON products involves unpredictable risks. Use or combine only the specified MEDICON products with one another.

- 10. DECONTAMINATION, CLEANING AND STERILIZATION**
-  **The following vCJD prion-specific protective measure is indicated when processing instruments.**

In the event of diagnosis of a definite or probable case of vCJD.
If it is not possible to use disposable products, the instrument used, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

If prion contamination is suspected:
In the event of suspected prion contamination, incineration of the instrument is recommended according to the vCJD task force final report.

If vCJD is excluded
Continue to use after instrument processing is completed. Otherwise, the instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

In the case of a non-identifiable vCJD illness.
Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions – e.g. mechanical alkaline cleaning combined with steam

sterilization.
If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (CNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134°C.

 **Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress-corrosion cracking. Only cleaned and disinfected implants may be sterilized.**

Reprocessing restrictions
Frequent reprocessing has a slight impact on implants

Location of use
Implants that have been contaminated by blood and/or secretions must not be placed back in the implant storage tray. They must be disposed of.

Preparation for decontamination
The implant storage tray must be placed in instrument holders suitable for machine cleaning. The instrument holders (e.g. wire trays) must allow subsequent cleaning in an ultrasonic bath or in the cleaning and disinfection device (CDD) where they are not impeded by acoustic or rinsing shadows.

- Pre-cleaning**
- Place implant storage tray for 15 minutes in ultrasonic bath at 40°C with 0.5% alkaline cleaner (Neodisher Mediclean forte) and treat with ultrasonic waves.
 - Remove implant storage tray and rinse thoroughly with cold water.

A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

- Cleaning: machine**
The cleaning and disinfection device (CDD) must meet DIN EN ISO 15883-1 requirements.
- Pre-rinse 1: 3 minutes with cold demineralised water, without additive;
 - Emptying;
 - Cleaning: with demineralised water, heat to 55°C and wash/clean for 5 minutes, add cleansing agent at 45°C, alkaline cleansing agent, strength 0.5%;
 - Emptying.
 - Neutralization: 3 minutes with warm water (>40°C) with addition of neutralizer, strength 1ml / l;
 - Emptying;
 - Final rinse: 2 minutes with warm deionised water (>40°C) (without any additive);
 - Emptying

Disinfection
"Thermal Disinfection A0 value 3000: Completely desalinated water, the thermal disinfection is carried out at temperatures > 80 °C and corresponding application time according to the A0 concept, DIN EN ISO 15883 and guideline DGKH, DGSV and AKI (e.g. A0 3000 = 90°C and 5 minutes application time). The operator is responsible for the implemented A0 value."

Drying
Complete drying must be ensured by the CDD, e.g. 60°C and 30 minutes. The instruments must be removed from the CDD promptly once the cleaning and disinfection program has finished. If necessary, compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Maintenance, inspection, and testing
After cleaning / disinfection the implants must be macroscopically clean, i.e., free from visible residues and soiling. Inspection is performed visually. Insufficiently cleaned implants must be cleaned again and then adequately rinsed and dried. Deformed or damaged implants must be removed and disposed of, as they can no longer be guaranteed to be safe to use.

Packing
The implant storage tray must be placed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868,
- DIN EN ISO 11607,
- suitable for steam sterilization (vapor permeable),

- 

12. SIGNS AND SYMBOLS

Symbol	Meaning
	Manufacturer
	Batch number
	Article number
	Warning
	Follow instructions
	MR unsafe
	CE label accord. To directive 93/42/EEC
	Non reusable
	Prescription use only

CON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines or instructions for use. The liability for defects by MEDICON eG is also not liable in case of changes or repairs to the product without prior written

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