Gebrauchsanweisung
Operating Instructions
Manual de instrucciones
Mode d'emploi
Istruzioni d'uso

2016-07-01

CE 0123

Rx only

Medicon eG · Gänsäcker 15
78532 Tuttlingen · Germany
Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50
E-Mail: sales@medicon.de
Internet: www.medicon.de

DE

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EN

MEDICON REFERENCE PIN ACCORDING TO KRETSCHMER

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

The MEDICON Kretschmer holders REF.: 68.84.04-06 are manufactured from stainless steel.

Table of contents

- General information
- Purpose and proper use/indications
- Contraindications
- Possible side effects and complications
- Single use product
- CT and MRI instructions
- Application/Handling of Pins
- Instrument
- Dekontamination, cleaning and sterilization
- Liability
- Signs and symbols

1. GENERAL INFORMATION

and must be disinfected, cleaned and sterilised prior to use. Please observe the following instructions. This will guarantee that the reference pins functions faultlessly and reliably.

2. PURPOSE AND PROPER USE/INDICATIONS

MEDICON reference pins acc. to Kretschmer REF.: 68.84.00 should only be used by surgeons with sufficient training and adequate experience in oral, orthodontic and maxillofacial surgery. MEDICON reference pins according to Kretschmer REF.: 68.84.00 are especially intended for insertion into the glabella during operations in the fields of oral and maxillofacial surgery. They were designed for monomaxillary and bimaxillary correction osteotomy for vertical reference measurement.

3. CONTRAINDICATIONS

Patients with insufficient or poor quality bone tissue in the area of the Glabella.

4. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

The possible complications are often the result of incorrect selection of the patient, lack of exercise or inadequate planning of the optimal implant position, rather than being caused by the implant itself. Apart from the side effects and complications mentioned above, the surgical intervention can also lead to problems such as nerve injuries, infections, pain, etc.

5. SINGLE USE PRODUCT

MEDICON reference-pins acc. to Kretschmer REF.: 68.84.00 are designed and constructed for single use only on a single patient and may not be reused. An explanted reference pin must never be reused. Even if the referencepin appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted reference-pin, reutilization incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

6. CT AND MRI INSTRUCTIONS

The reference pin acc. to Kretschmer poses a risk in the proximity of a CT or MRT and must not be close to such devices during their use.

7. APPLICATION/HANDLING OF PINS

The screwing of the reference-pin acc. to Kretschmer REF.: 68.84.00 can only be done with a special fixation REF.: 68.84.04-06. Any use of implants and instruments from other manufactures in combination with MEDICON reference pin according to Kretschmer REF.: 68.84.00 entails unpredictable risks, since the products are not adapted to each other. This may also cause danger to the patient, the surgeon or others. The implants have to be removed after surgery. We can therefore only accept liability for implants up to the end of surgery. We do not accept responsibility for any damage arising from Reference Pins that were not removed after healing. It cannot be excluded that implants will break, loosen, corrode, migrate through tissue and cause pain. The packaging label carries a "Lot No.". We recommend copying this "Lot No." to the patient records as this number allows the manufacturing process to be tracked.

8. INSTRUMENTS

See the operating instruction for MEDICON Fixation for the Reference-Pin according to Kretschmer REF.: 68.84.04-06 for more information about handling, maintenance, control etc.

9. DECONTAMINATION, CLEANING AND STERILIZATION

Matrix for the preparation procedure MPK 01				
Procedure			T	t
Place of use	01	Disposal	Wet	--
			Dry	--
Dekontamination	02	Pre-cleaning	manual	++
			Ultrasonic	40 °C 15 min. --

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able, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device. Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

10. LIABILITY

In the event of discrepancies between the English and the German version of these instructions for use only the German version shall be applicable. Only the latest revision of the instructions for use applies. Due to constant technical development the contents of these MEDICON instructions for use are updated regularly. Please use the MEDICON extranet www.medicon.de/extranet to ensure that you are using the current version. The version date of each edition of the instructions for use can be found in the footnote. MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the MEDICON Repair Service (MRS). Please contact MEDICON eG if you have further questions concerning the instruments.

CAUTION: For USA, federal law restricts this device to purchase by or on behalf of physician or hospital.

11. SIGNS AND SYMBOLS

- Follow instructions
- CE label accord. to directive 93/42/EEC
- Warning
- Single use product
- Not suitable for use with MRI
- Order No.
- Manufacturer
- Rx only** Prescription use only

Procedure			T	t	
04	Cleaning	Manual	45 °C	>10 min. ++	
		Mechanically	55 °C	10 min. ++	
		Alkaline cleaning agent > pH10	--	-- ++	
		Acidic cleaning agent	--	-- 0	
		Neutral cleaning agent	--	-- ++	
05	Rinsing	Desalinated water	>40 °C	2 min. ++	
06	Neutralisation	Citric acid	>40 °C	3 min. ++	
		Phosphoric acid	>40 °C	3 min. 0	
07	Disinfection	Chemical	20 °C	10 min. ++	
		Thermal	90°C	5 min. ++	
08	Thermal dry	Thermal	max. 120 °C	++	
		Medical compressed air	--	-- ++	
Maintenance and inspection	09	Functional test	Instrument oil REF 46.00.44	--	-- +
Packaging	10	Packaging	Sterile packaging	--	-- +
Sterilization	11	Steam	SAL 10 ⁻⁶	134 °C	5 min. ++
			SAL 10 ⁻⁶	134 °C	4 min. ++
			SAL 10 ⁻⁶	132 °C	4 min. ++
		Hot air	--	180 °C	<30 min. 0
			Gas	--	-- --
			Plasma	--	-- --

The up-to-date and comprehensive process description for reprocessing, information on prion protection measures (v) CJK as well as the validation references are available for download at www.medicon.de/extranet/

Note

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method. If the previously described equipment, materials, and chemicals are not avail-

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