DE

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Rx only

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EN

MEDICON HOLDER FOR REFERENCE PIN ACCORDING TO KRETSCHMER

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

The MEDICON Kretschmer holders REF.: 68.84.04-06 are manufactured from stainless steel.

Table of contents

1. General information

2. Purpose and proper use/Indications

3. Contraindications

4. Possible side effects and complications

5. CT and MRI instructions

6. Application/Handling of Pins

7. Instruments

8. Decontamination, cleaning and sterilization

9. Liability

10. Signs and symbols

1. GENERAL INFORMATION

The MEDICON Kretschmer holders REF.: 68.84.04-06 are supplied unsterile and must be disinfected, cleaned and sterilised prior to use. Please observe the following instructions. This will guarantee that the fixations will function faultlessly and reliably.

2. PURPOSE AND PROPER USE/INDICATIONS

The MEDICON Kretschmer holders REF.: 68.84.04-06 should only be used by surgeons with sufficient training and adequate experience in oral, orthodontic and maxillofacial surgery. The MEDICON Kretschmer holders REF.: 68.84.04-06 are especially used in the fields of oral and maxillofacial surgery for insertion of the reference pins according to Kretschmer REF.: 68.84.00 into the glabella. They also serve to determine the vertical dimension of the maxilla.

3. CONTRAINDICATIONS

Patients with insufficient or poor quality bone tissue in the area of the Glabella.

4. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

Potential complications are often a result of a lack of practice on the part of the physician and are usually not caused by the holders themselves. Incorrect insertion of the holders and the resulting incorrect insertion of the reference pins REF.: 68.84.00 can cause problems such as nerve damage, infections, pain etc.

5. CT AND MRI INSTRUCTIONS

The instruments intended for use with the MEDICON Kretschmer holders pose a risk in the surroundings of a CT or MRT. They must not be near these devices during the application of these procedures.

6. APPLICATION/HANDLING OF PINS

Exclusively the MEDICON Kretschmer holders REF.: 68.84.04-06 are to be clamped into the designated Medicon fittings. The combination of incorrect MEDICON fittings or products from other manufacturers entails unpredictable risks and/or holds the danger of contamination of the material as well as an insufficient adjustment of the fixation, fitting and implant. This may also cause danger to the patient, the surgeon or others.

7. INSTRUMENTS

The fixations are subject to wear and tear and mechanical strains even when used normally. To prevent failure of, or mechanical damage to the fixations during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use MEDICON Kretschmer holders REF.: 68.84.04-06 that were found to be damaged, deformed or not fully functional.

8. DECONTAMINATION, CLEANING AND STERILIZATION

Matrix for the preparation procedure MPK 01

Procedure			T	t	
Place of use	01	Disposal	Wet	–	–
			Dry	–	++
Decontamination	02	Pre-cleaning	manual		--
			Ultrasonic	40 °C	15 min. --
	03	Cleaning	Manual	45 °C	>10 min. ++
			Mechanically	55 °C	10 min. ++
			Alkaline cleaning agent > pH10	–	–
			Acidic cleaning agent	–	–
			Neutral cleaning agent	–	–
	05	Rinsing	Desalinated water	>40 °C	2 min. ++
			Citric acid	>40 °C	3 min. ++
	06	Neutralisation	Phosphoric acid	>40 °C	3 min. 0
			Chemical	20 °C	10 min. ++
	07	Disinfection	Thermal	90°C	5 min. ++

Procedure

08

Thermal dry

Thermal

max. 120 °C

++

Maintenance and inspection

09

Functional test

Instrument oil REF 46.00.44

–

–

+

Packaging

10

Packaging

Sterile packaging

–

–

++

Sterilization

11

Steam

SAL 10⁻⁶

134 °C

5 min.

++

SAL 10⁻⁶

134 °C

4 min.

++

SAL 10⁻⁶

132 °C

4 min.

++

Hot air

–

180 °C

<30 min.

0

Gas

–

–

–

–

Plasma

–

–

–

–

Required process steps, mechanically:

01-04-05-06-05-07-08-09-10-11

Required process steps, manual:

01-04-05-07-05-08-09-10-11

Legend:

T Temperature

t Soaking time

min Minutes

°C Degrees Celsius

++ Procedure validated by the manufacturer

+ Required process step

0 Non-validated procedure, check material compatibility before initial use

00 Not applicable

– No information

-- Not required

The up-to-date and comprehensive process description for reprocessing, information on prion protection measures (v) CJK as well as the validation references are available for download at [www.medicon.de/extranet/](#)

Note

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method. If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device. Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

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