Deutsch

Rx only

Rev.: 08-12-2025

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Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

MEDICON implants are manufactured from pure titanium or Ti6Al4V al-loy. Both materials are biocompatible, corrosionproof and nontoxic in the biological environment. The surface is chemically passive and the material is antimagnetic. The implants are supplied in non-sterile condition and must be disinfected, cleaned and sterilised prior to use.

Table of contents:

1. General information
2. Purpose and proper use/Indication
3. Contraindications
4. Possible side effects and complications
5. Single use product
6. MRI instructions
7. Implant selection
8. Application / Handling of plates
9. Application / Handling of screws
10. Instruments
11. Sterility
12. Liability
13. Signs and Symbols

1. GENERAL INFORMATION

The MEDICON Micro Neuro Flap Sets sterile consist of 3 microplates and 6 self-tapping microscrews and are supplied in a sterile condition.

2. PURPOSE AND PROPER USE/INDICATION

The MEDICON Micro Neuro Flap Sets sterile may be used only by sufficiently trained and experienced surgeons in craniofacial neurosurgery. The MEDICON Micro Neuro Flap Sets sterile are intended for osteosynthesis on the osseous structures of the cranial skull.

The implants are intended to support the normal bone healing process, but not to replace normal body structures or body weight in the event of inadequate bone healing. The indications include:

- Cranial fractures
- Fixation of the bone flap

3. CONTRAINDICATIONS

- Patients with unstable physical and/or mental health conditions.
- Patients that are unable to follow the instructions for post-operative care. This may be caused e.g. by psychological/mental or neurological problems of the patient.
- Patients with inadequate bone tissue, or bone tissue of insufficient quality, and patients with circulatory disorders or latent infections
- Extreme material sensitivity, i. e. foreign-body reaction of the patient. For this case, appropriate tests must be carried out under all circumstances prior to implantation (even when there is only a suspicion of such condition).
- Acute infections.

4. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

- If the implants were not adapted properly, bone fusion can be delayed, or fusion is insufficient or does not occur at all.
- Complaints, pain, abnormal sensations or palpability of the implant.
- Extreme material sensitivity of the patient (foreign-body sensitivity) manifested by allergic reactions, etc.
- Stronger connective-tissue reaction in the fracture region or in the vicinity of the implant.
- Delayed or insufficient fracture union, which can result in implant breakage.
- Inadequate healing.
- Insufficient bone formation, osteolysis, osteomyelitis, osteo-porosis, inhibited revascularization or infections that can result in loosening, bending, rupture or breakage of the implants.
- Breakage, bending, migration or loosening of the implant.
- Reduced bone density due to stress shielding.
- Loosening of the implant as a result of insufficient screw tightening.
- Osseous necrosis, osteoporosis, inhibited revascularizat on, bone-resorption, and poor bone formation, which can cause premature loss of fixation leading to non-union.
- Extreme or repeated bending of the plates in the same place can result in plate breakage.
- Unstable comminuted fractures can lead to an increased connec-tive-tissue reaction in the fracture region.

Apart from the side effects and complications mentioned above, the surgical intervention can also lead to problems such as nerve injuries, infections, pain, etc., which are not directly attributable to the implants.

5. SINGLE USE PRODUCT

 Implants are designed and constructed for single use only on a single patient and may not be reused.

An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wear out or failure. The user is liable for non-compliance with the instructions for use.

6. MRI INSTRUCTIONS**Instruments**

The use of medical devices in the vicinity of an MRI poses a hazard. Individual medical devices must not be located near the equipment during the application of these procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MRI environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MRI devices.

7. IMPLANT SELECTION

The surgeon is responsible for the correct selection of implants.

The following aspect are crucial for the correct selection of implants:

- Bone defect to be treated
- Body weight of the patient
- general health and degree of activity of the patient.

Incorrect implant selection can lead to premature implant loss and loosening, deformations or fractures of the implants.

Among other factors, the success of the intervention depends on how the implants are handled. The plates must be bent carefully and cautiously, avoiding repeated and extensive deformation. Damage and scratches reduce the stability of the product and result in early fatigue of the implant material. In cases of failed, delayed or insufficient fracture union, the incidence of implant bending, implant failure or implant breakage cannot be excluded. Therefore, immobilization of the fracture site must be ensured until there is solid bone fusion. Incessant load changes on the implant can result in fatigue fractures.

The implants and instruments are adapted to each other. Any use of implants and instruments from other manufactures in combination with MEDICON products entails unpredictable risks, since the products are not adapted to each other. To avoid risks, only such MEDICON products that are designated for combination with each other may be combined with each other.

Prior to applying the product, the surgeon must thoroughly discuss the intended surgical result with the patient. In such discussions, special attention must be given to the postoperative aspects and aftercare, which may become necessary. The patient must be informed that, due to the limited stability of the fixation, the fixation site must not be put under the full body weight. The patient must also be informed the non-observance of this advice can put at risk any safe bone union. The patient must be instructed to inform, without delay, the operating surgeon of unusual changes at the operating site. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing.

The implants can be removed after bone union is completed. An implant left in the body after complete healing may act as a load-baring element and can contribute to an increased risk of re-fracturing in active patients. For this reason, too, we accept liability for the implants only up to complete healing. We will not accept liability for any damage or injuries caused by implants that were not removed after complete healing of the treated fracture. The surgeon is responsible for deciding whether the implants shall remain in the body after complete healing of the fracture or whether they shall be removed. For his decision, the surgeon has to weigh risks and advantages of leaving the implants in the body. It cannot be excluded that implants break, loosen, migrate in tissue and cause pain.

The Micro Neuro Flap Sets sterile are provided with a packaging label. The packaging label carries a "Lot Number". We recommend transferring this "Lot Number" to the patient records as this number allows tracing the production history of the implant back to the raw material.

8. APPLICATION/HANDLING OF PLATES

The plates have to be adapted to the natural or planned bone contours as precisely as possible. For adapting the plates use the instruments provided for this purpose by the manufacturer. During the plate bending procedure, there occur cold stresses, under which the titanium hardens while, at the same time, its ductility is decreased. Therefore it is absolutely necessary to achieve the required plate shape in as few bending steps as possible. Excessive bending can lead to postoperative plate breakage. Overly aggressive application of the bending instruments can result in visible damage to the implant. When such damage has occurred, a new plate must be selected and adapted more cautiously. Deformed plate holes for the screws entail not only an increased risk of breakage at the

location of the holes; they also compromise the precise seating of the screw head in the plate. Therefore the plates must be bent carefully under all circumstances. Deformed plates must be inspected for notches, deformed screw holes and other mechanical damage prior to application.

9. APPLICATION/HANDLING OF SCREWS

Self-tapping screws do not require pre-cutting of threads. Still, in dense and hard bones, it may be necessary to prepare drill holes. Self-tapping screws are unsuitable for small bone fragments as the axial pressure on them can lead to fragment shifts or screw breakage.

To apply the drill holes, only use the respective drill designated by MEDICON for the specific screw so that the correct drill hole diameter for the screw is obtained. Only use drill bits with sharp cutting edges. During drilling, there is a risk of heat damage to the bone. Therefore, only work at low drill speed and always ensure sufficient cooling during the drilling procedure.

Excessive application of force during the drilling procedure can lead to drill breakage, endangering the patient, the surgeon and third parties. If, due to a failed drilling effort or a slack drill hole, the screw cannot be anchored firmly in the bone, the emergency screws of the MEDICON Micro System must be used.

Only use the designated drill handles and blades of the MEDICON Micro system. To ensure that the screwdriver blade is firmly seated in the screw head, axial force must be applied when inserting the screwdriver into the screw head. This ensures the correct axial alignment of screw and screwdriver, and prevents slipping of, or damage to the screw head or the screwdriver blade.

Application of excessive force in tightening the screws can lead to screw breakage during the intervention.

Wear of the screwdriver blade will impair the firm connection between blade and screw. Whenever such wear is detected, the screwdriver blade must be replaced by a new blade.

Once the implantation has been completed, the secure connection between every screw and plate must be checked. If necessary re-tighten the screws. Prior to implant removal, all screw heads must be cleaned carefully, using a scalpel or another suitable instrument, so that the screwdriver can be optimally seated in the screw head.

10. INSTRUMENTS

Only designated MEDICON instruments for the particular purpose shall be used. Before using the MEDICON Micro Neuro Flap Set sterile, ensure that all the instruments required for the implantation are available. These specifically include bending and adaptation tools, twist drills, screwdriver handles and bits. The instruments intended for application of the system are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force. To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional.

11. STERILITY

The implants are sterilized by radiation and must be stored in their original packaging until directly before use. The implants can then be removed from their original and protective packaging. Check the expiration date and the integrity of the sterile packaging prior to use. Implants that have exceeded their expiration date must not be used. Implants with damaged packaging must not be reused. Risk to patients and/or user infection and reduction in function as a result of reusing products. Contamination and/or impaired function of the products may lead to injury, illness or death.

12. LIABILITY

In the event of any discrepancies between the German version of these instructions for use and any other translated versions, only the German version shall be applicable. Medicon eG assumes no liability for damages caused by improper use, care or maintenance, non-observance of post-operative instructions or non-compliance with the restrictions for use and guidelines in the instructions for use. The liability for defects by Medicon eG is also not applicable in case of changes or repairs to the product without prior written consent from Medicon eG as well as in case of repairs which are not carried out by an authorized workshop of Medicon eG or the Medicon Repair Service (MRS).

CAUTION:

Federal law and regulations in the United States restrict the sale of this device to a physician or hospital or on prescription.

13. SIGNS AND SYMBOLS

	Instruction for use
	Caution
	CE marking according to Directive 93/42/EEC
	Manufacturer
	Date of manufacture
	Single use
	Do not re-sterilize
	Sterilized using irradiation
	MR unsafe
	Humidity limitation
	Protect from sunlight
	Prescription only (U.S. law)
	Batch code
	Article number
	Use by date
	Do not use if package is damaged
	Keep dry
	Temperature limitation

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