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Rx only

Carefully read the operating instructions and the "Distraction Devices" brochure prior to clinical application, and keep them safe and at hand. The instructions and notes contained therein must be followed.

TABLE OF CONTENTS

- General notes
- Intended use
- Indications
- Contraindications
- Possible adverse effects and complications
- Single-use product
- MR notes
- Application and handling
- Pre-operative and post-operative behaviour
- Performing the distraction
- Decontamination, cleaning and sterilisation
- Liability
- Description of symbols and icons

1. GENERAL NOTES

MEDICON intraoral mandibular distractors are made of titanium. The material is biocompatible, corrosionresistant and non-toxic in a biological environment. It permits almost artefact-free imaging with:

- Conventional X-ray imaging
- Computer tomography

The surface is chemically passive; the material is antimagnetic.

The implants and instruments are delivered in non-sterile condition. The implants must be cleaned, disinfected and sterilized before the single use. The instruments must be cleaned, disinfected and sterilized before the first application and before every further use. Please observe the following notes. These will ensure flawless and reliable functioning of the product.

2. INTENDED USE

MEDICON intraoral mandibular distractors are used for the correction of malpositions of the mandible by doctors and surgeons with sufficient experience in the field of oral and maxillofacial surgery in hospitals and clinical practices. The intraoral application of distractors enables a growthretarded mandible to be lengthened in an inconspicuous manner, whether unilaterally or bilaterally. The implants are single-use products and may be used only in conjunction with the specified MEDICON application instruments.

3. INDICATIONS

This results in the following indications:

Horizontal distractors:

- Hemifacial microsomia, possibly with airway obstruction
- Hemifacial microsomia, possibly with facial asymmetry
- Pierre-Robin syndrome
- Treacher-Collins syndrome
- Crouzon syndrome
- Extreme facial asymmetry
- Mandibular retrognathia

Vertical distractors:

- Hemifacial microsomia
- Short mandibular ramus
- Open bite with short mandibular ramus
- Pierre-Robin syndrome
- Craniofacial microsomia

Medial mandibular distractor:

- Transversal crowding of the dental arch
- Front teeth crowding

Mandibular distractor for alveolar ridge augmentation:

For this, the basal bone segment must meet at least the following requirements for adaptation of the distractor:

- Horizontal dimension: ≈ 15 mm
- Vertical dimension: ≈ 8 mm
- Transversal dimension: ≈ 4 to 5 mm

This can be the case for the following clinical pictures:

- Dentoalveolar deficit: Following resection of segments of the mandible in case of benign or malignant neoplasia or following trauma
- Mandibular atrophy of the anterior or posterior sector or of the entire mandible

4. CONTRAINDICATIONS

MEDICON intraoral mandibular distractors may not be used in case of:

- Patients who are unable to comply with the instructions for postoperative care. For example, due to psychological, mental or neurological problems

- Patients in unstable physical and/or mental condition
- Patients with inadequate or low-quality bone tissue, with perfusion problems or latent infections
- Proven allergy to titanium
- Material hypersensitivity, i.e. reaction of the patient to foreign bodies. Here appropriate tests are mandatory before implantation (even in case of mere suspicion)
- Acute infections

Additional contraindications for parallel mandibular distractors include:

- Open bite, facial asymmetry, laterognathia

5. POSSIBLE ADVERSE EFFECTS AND COMPLICATIONS

- In case of insufficient adaptation of the implants, bone healing may be delayed, or the bone may heal insufficiently or not at all
- Discomfort, pain, abnormal sensations due to the implant
- Material hypersensitivity of the patient due to foreign objects in the form of allergic reactions
- Increased reaction of the connective tissue in the area of the osteotomy and/or the implant
- Delayed or insufficient healing, which can result in breakage of the implant
- Insufficient bone formation, osteolysis, osteomyelitis, osteoporosis, inhibited revascularisation or infection that may lead to loosening, bending, tearing or breaking of the implants and endanger the fixation or the healing of the bone
- Breaking, bending, migration or loosening of the implant
- Decrease in bone density due to stress shielding
- Implant loosening caused by insufficient tightening of the screws
- Bone necrosis, osteoporosis and bone resorption can lead to non-union formation
- Extreme and/or multiple bending of the implants at the same point may cause breakages
- Increased connective tissue reactions in the area of the fracture which may present as deep or superficial early or late infections

In addition to the undesirable effects and complications already mentioned, the surgical procedure can lead to other problems such as nerve damage, infections, pain etc. that are not necessarily attributable to the implant.

6. SINGLE-USE PRODUCT

MEDICON implants are intended and designed for single use on one patient only and may not be reused. An explanted implant must never be reused. Even if the implant appears to be intact or fully functional, there may be wear, minor defects and attrition due to overuse that are not visible to the naked eye. As the impact of the forces and conditions inside the body on the stability, function and material quality of an explanted implant cannot be estimated, reimplantation would entail the risk of premature wear or failure and is therefore unacceptable. The user will be held liable if the operating instructions and the MEDICON brochure "Distraction Devices" are not observed.

7. MR NOTES

Instruments

The use of medical devices in the vicinity of an MR poses a hazard. Individual medical devices must not be located near the equipment during the application of these procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MR environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MR devices.

8. APPLICATION AND HANDLING

Implant selection:

The surgeon is responsible for correct selection and use of the implants.

The following aspects are critical for correct implant selection:

- Bone defect/deformity to be treated
- Patient's body weight
- Physical condition and level of activity of the patient.

Incorrect implant selection can result in premature loosening, bending or implant breakage and therefore loss of function. Only use of the correct components will result in optimally stable fixation, while a wrong decision may lead to loosening, bending or breakage of the implant and/or bone. Always choose the mandibular distraction device that best suits the specific situation and indication. All implants must be checked for function, deformation and mechanical integrity prior to implantation. Damaged implants must be replaced, as they can lead to partial or even total loss of function. Damage leads to reduction in the product's strength and premature fatigue of the implant. If no bone healing sets in, or if healing is delayed or inadequate, the possibility that the implant will bend, break or loosen cannot be ruled out. The constant load change the implants are subject to may lead to fatigue breakage. It cannot be ruled out that the implants may break, loosen, bend, migrate in the tissue and cause pain. Implants and instruments are developed and made to be compatible with each other. Using products by other manufacturers with MEDICON products may lead

to unforeseeable risks and/or danger of the material becoming contaminated, as well as the instrument and implant not being sufficiently compatible. Patients, users and/or third parties may be endangered as a result. The decision on whether the implants should remain in the body or be removed after complete healing is the responsibility of the surgeon and should be made after considering the risks and benefits of each approach.

Use and handling of plates:

The plates must be adapted to the natural or planned bone structure as precisely as possible. Use only the instruments recommended by MEDICON for adapting the plates. During the bending of the plate, cold straining occurs, under which the titanium increases in hardness but simultaneously loses deformability. It is therefore essential to use as few bending steps as possible to achieve the desired shape of the implant. Excessive bending can lead to post-operative breaking of the plates. The success of an operation depends, among other things, on how the products are handled.

The plates must be bent carefully, avoiding excessive deformation. The plates must not be bent to and for more than 2 to 3 times. Damage leads to reduction in strength and premature fatigue.

Excessively forceful use of the bending instruments can lead to visible implant damage. In such cases, a new plate must be used and adapted more carefully. Deformed plate holes for screw insertion not only lead to increased risk of breakage in these locations but also compromise the precise fixation of the screw head in the plate. The plates must therefore be bent carefully. Deformed plates must always be checked for notches, deformed screw holes or other mechanical damage before being used on the patient.

Damaged plates must be replaced. Always choose the plates that best suit the specific situation and indication. Should it be necessary to shorten a plate, use only the instruments provided for this purpose. The plate must be cut between the screw holes. Ensure that the cut-off part of the plate is not catapulted towards the patient, the surgeons or others, to prevent putting them at risk. The cutting edges of a shortened plate must be deburred in order to prevent tissue injury or irritation.

Use and handling of screws:

Self-tapping screws do not require drilling a threaded hole. Use only the drills specified by MEDICON to drill holes so that the correct hole diameter required for the screws can be achieved. The suitable screws can be found in the "Distraction Devices" brochure.

Use only drills with sharp cutting edges.

There is a risk of heat damage to the bone during drilling. Therefore, always ensure sufficient cooling during the drilling process, and use low drilling speeds. Excessive force during the drilling process can result in fracture of the tool, which can be dangerous for the surgeon, patient and third parties. If the screw cannot be securely anchored in the bone due to an incorrect or worn-out drill hole, the emergency screws intended for such circumstances are to be used.

The screwdriver must be inserted into the screw head with axial pressure to ensure that the blade sits tightly on the head. This guarantees correct longitudinal alignment of screw and screwdriver and prevents slipping of or damage to the screw head or the blade. If excessive force is used when fixing the screws, the screws may break during surgery, or the screw heads may tear off. Wear and tear effects to the screwdriver blade impair the firm connection between the blade and screw. In this case, replace the blade with a new one. Upon completion of the implantation, the secure connection between all screws and plates must be checked. If necessary, retighten the screws.

Before removing the implant, all screw heads must be thoroughly cleaned with a scalpel or other suitable instrument so that the screwdriver sits perfectly on the screw head.

Self-tapping screws:

When using self-tapping screws, you are usually not required to drill a hole. If the bone is dense and hard, it may still be necessary to drill holes, e.g. in the calvarium or mandible. Self-tapping screws are not suitable for small bone fragments, as the axial pressure may lead to displacement of the fragments or breakage of the screw.

Use and handling of instruments:

The instruments intended for application of the system are subject to wear and mechanical stresses even when used normally, but especially if used too forcefully. In order to prevent failure or mechanical damage to the instruments during surgery, they must be checked before each use to ensure that they are mechanically intact, that there are no deformations, and that the parts are fully functioning. Damaged instruments must not be used and are to be replaced. Use only accessories and instruments specifically approved by MEDICON, in order to avoid risks in connection with the compatibility of the products.

9. PRE-OPERATIVE AND POST-OPERATIVE CONDUCT

Pre-operatively:

Before using the products, the surgeon must thoroughly discuss the desired operation result with the patient. Particular attention must be paid to post-operative behaviour and potentially required follow-up care.

Post-operatively:

The patient must be instructed to immediately inform the surgeon about any unusual changes at the operation site. The surgeon must consider possible implications such as implant failure, and discuss required measures for further healing with the patient. The

decision on whether the implants should remain in the body or be removed after complete healing is the responsibility of the surgeon and should be made after considering the risks and benefits of each approach. If healing fails to set in, or is delayed or inadequate, the possibility that the implant will bend, fail or break cannot be ruled out. Therefore adequate immobilisation of the fracture site must be ensured until the bone has grown together again firmly. The constant load change the implants are subject to may lead to fatigue breakage. The treating physician decides on the type, duration and intensity of physical activities after the surgery. The patient must be informed that failure to comply with the physician's instructions may lead to non-foreseeable complications. This applies in particular if the patient manipulates the implant or the apparatus. The patient must be regularly monitored until removal of the distractors and until the final completion of the treatment. After explantation, the surgeon is responsible for disposal of the implants according to the standard. It cannot be ruled out that the implants may break, loosen, bend, corrode, migrate in the tissue or cause pain. A "LOT" number (batch number) is indicated on the label of the packaging. To ensure 100% traceability, we recommend noting down the "LOT" number in the patient records.

Disposal:

Observe the applicable national laws when disposing of MEDICON products and accessories!

10. PERFORMING THE DISTRACTION

PLEASE NOTE: Anti-clockwise = opening / Clockwise = closing

The distraction is performed by anti-clockwise rotation. Use the screwdriver provided to perform the distraction (Figure A). The arrow on the screwdriver handle shows the direction in which to turn in order to activate the distraction (Figure A).

One complete turn corresponds to a distraction distance of 0.4 mm, except for the mandibular distractor for alveolar ridge augmentation. In this case, a complete rotation corresponds to a distraction distance of 0.3 mm.

The patients or their family members, respectively, are to be informed comprehensively about the distraction process and the handling of the distractors, especially if the distraction is to be continued outside the hospital or clinic. The patient is to be provided with a schedule for activation of the distractors, indicating the distractions to be performed.

11. DECONTAMINATION, CLEANING AND STERILISATION

Opening the distractor:

The distraction devices must be opened for processing and sterilisation (see chapter 10 "Performing the distraction"). Avoid overwinding the distraction device when opening and closing it. For this reason, open the distraction device only until you reach the STOP mark (Figure B). Opening the distraction device beyond this mark can lead to damage.

The following prion-specific protective measure against (v)CJD for instrument preparation is indicated.

In case of diagnosis of definitive or probable (v)CJD:

- Where use of disposable items is not possible, the instruments used, where contamination has taken place or cannot be excluded, must be disposed of as waste for incineration.

In suspected cases:

- If prion contamination is suspected, it is recommended that any medical instruments be incinerated in accordance with the final report by the vCJD task force.

In case of exclusion of (v)CJD:

- Reuse after final processing. Otherwise all instruments where contamination has taken place or cannot be excluded must be disposed of as waste for incineration.

EN

CE₀₁₂₃

Rx only

Rev.: 10-12-2025

Instruction for use

If (v)CJD cannot be detected:

- Even if there is nothing to suggest presence of a prion disease, two procedures with at least partial effectiveness against prion contamination should be used for processing, e.g. alkaline mechanical cleaning combined with steam sterilisation.

If no alkaline mechanical cleaning process or other cleaning process with proven effectiveness against prion contamination is used, and the medical devices come into contact with risk tissue (CNS, eyes, lymphatic tissue), the Robert Koch Institute (RKI) recommends an extended sterilisation time of 18 minutes at 134 °C (273 °F).

Cleaning solutions to which hydrogen peroxide has been added and/or very alkaline washing solutions can cause colour changes. This may cause the coding function to be lost. Only cleaned and disinfected implants may be sterilised.

Limitations on reprocessing:

Frequent reprocessing has little effect on the implants.

Place of use:

Implants contaminated with blood and/or secretions must not be placed back into the implant storage tray. They must be disposed of.

Preparation for decontamination:

The implant storage tray must be appropriately placed on a washer-proof instrument carrier for rinsing. The instrument carriers (for example wire mesh trays) must be such that subsequent cleaning in the washer-disinfector (WD) is not impaired by "rinsing shadows" or "sound shadows" that may render any areas inaccessible to ultrasound or rinsing.

Pre-cleaning:

Always use a fresh and previously unused cleaning solution for manual pre-cleaning in the ultrasonic bath in order to prevent cross-contamination.

- Place the implant storage tray in the ultrasonic bath and treat with ultrasound for 15 minutes at 40 °C (104 °F) using 0.5 % alkaline cleaner (neodisher® MediClean forte).
- Remove the implant storage tray and thoroughly rinse it with cold water.

Automated cleaning:

To prevent cross-contamination during machine processing, always clean and disinfect the implant cassette containing the implants to be reprocessed in the washer-disinfector (WD) in a separate batch from instruments contaminated with blood or secretion. The washer-disinfector (WD) must meet the requirement of DIN EN ISO 15883-1.

- Prewashing: 3 minutes with softened cold water, without additives
- Emptying
- Cleaning: with softened water, heating to 55 °C (131 °F) and 5 minutes washing/ cleaning, addition of the cleaning agent at 45 °C (113 °F), alkaline cleaning agent (neodisher® MediClean forte), dosage 0.5 %.
- Emptying
- Neutralisation: 3 minutes with warm water (>40 °C (104 °F)) with addition of neutraliser, dosage 1 ml/l
- Emptying
- Final rinse: 2 minutes with warm, fully demineralised water (>40 °C (104 °F)) (with no other additives)
- Emptying

Disinfection:

Thermal disinfection A0 value 3000:

- Fully demineralised water, thermal disinfection is carried out at temperatures >80 °C (176 °F) with a corresponding application time according to the A0 concept of standard DIN EN ISO 15883 and DGHK, DGSV and AKI guidelines (for example A0 3000 = 90 °C (194 °F) and 5 minutes application time).
- The operator is responsible for the A0 value to be implemented.

Drying:

Sufficient drying must be ensured using the washer-disinfector device (WD), for example 60 °C (140 °F), 30 minutes. The instruments must be removed from the WD immediately after completion of the cleaning and disinfection programme.

If necessary, use of compressed air for drying is recommended due to its good and rapid effect (RKI recommendation).

Maintenance, control and inspection:

Following cleaning/disinfection, the implants must be macroscopically clean, i.e. free from visible contaminants and residues. The inspection is performed visually. Insufficiently cleaned implants must be cleaned again and then rinsed and dried sufficiently. Deformed or damaged implants must be sorted out and disposed of, as safe use cannot be guaranteed.

Packaging:

The implant storage tray is to be inserted into a suitable sterile barrier system. The sterile barrier system must meet the following requirements:

- DIN EN 868

- DIN EN ISO 11607
- suitable for steam sterilisation (steam permeability)
- sufficient temperature resistance up to 138 °C (281 °F)

Sterilisation accessories and sterilisation packaging must be adapted to the package content as well as to the sterilisation procedure to be applied.

Sterilisation:

For sterilisation, the following sterilisation procedure is to be used, taking into account all applicable national requirements:

- Fractionated vacuum procedure, fractionated three times and with sufficient product drying time
- Steam steriliser validated according to DIN EN 13060 or DIN EN 285 and according to DIN EN ISO 17665-1
- Sterilisation time and temperature: At least 5 minutes hold time at 134 °C (273 °F)

Reaching an SAL (Sterility Assurance Level) of 10⁻⁶ is indispensable.

Storage:

Reprocessed sterile implants are to be stored in dry, dust-protected, low-germ, dark and cool rooms, free of vermin, in an appropriate reusable sterilisation container. To prevent condensation, extreme temperature fluctuations should be avoided during storage.

No chemicals may be stored together with implants. Walls, floors and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must have a ground clearance of at least 30 cm (12"). The permissible on-site storage period depends on the type of sterile barrier system used and on the storage conditions. The permissible storage period is to be specified by the operator.

Additional information on reprocessing:

A validated automated cleaning and disinfection procedure is always to be preferred to manual cleaning because of the greater safety of the process. Thorough cleaning also helps to maintain value and is a prerequisite for successful sterilisation.

Please take note of the following points regarding automated processing:

- For effective automated processing, load the wire mesh baskets in such a way that they can easily be mechanically rinsed. Mesh baskets must not be overloaded
- Avoid "rinsing shadows" cast by large instruments

The times and temperatures stated in these reprocessing instructions are minimum requirements that must be adhered to. Should any reduction be needed for procedural reasons, this must be validated by the operator.

Exceeding the specified times and temperatures is generally possible. However, this will lead to increased stress on the material, which can cause premature wear. We accept no responsibility for use of other sterilisation procedures.

Information on validation of the processing procedure:

The validation was carried out with the following devices, materials and chemicals:

Washer-disinfector (WD):	Type Miele Disinfector G 7735 CD and 7836 CD Trolley for surgical instruments neodisher® FA Dr. Weigert GmbH & Co. KG neodisher® Z Dr. Weigert GmbH & Co. KG
Alkaline cleaning agent:	Sonorex DIGITEC DT 156 BH
Neutraliser:	MMM Vaculab 969 S 3000, MMM Selectomat S 3000 Stiefenhofer KS 666-2ED H+P Varioclav 400E
Ultrasonic bath:	
Sterilisers:	

Note:

The person processing the instrument is responsible for ensuring that the processing actually performed with the equipment, materials and staff at the reprocessing facility achieves the desired results. Validation and routine monitoring of the process are normally required. If the devices, materials and chemicals described above are not available, it is up to the persons processing the instrument to have their procedures appropriately validated. Please observe the notes and regulations of the relevant national laws and standards. MEDICON eG reserves the right to make changes to these instructions due to new findings.

12. LIABILITY

In the event of any discrepancies between the non-German and the German version of these operating instructions, only the German version is authoritative. Only the latest revision of the operating instructions applies. To ensure that you are using the current version. The version date of each edition of the operating instructions is included into the printout. MEDICON eG assumes no liability for damage caused by improper use, incorrect postoperative conduct, care or maintenance, or non-compliance with the restrictions on use and other guidelines in the operating instructions.

MEDICON eG assumes no liability for defects in case of changes or repairs to the product without prior written consent from MEDICON eG either, nor in case of repairs not carried out by workshops authorised by MEDICON eG or by the MEDICON Repair Service.

13. DESCRIPTION OF SYMBOLS AND ICONS

	Observe the operating instructions
	Please note
	CE label acc. to Directive 93/42/EEC
	Manufacturer
	Do not reuse
	MR Unsafe
Rx only	Subject to medical prescription (US federal law)
	Production lot number, batch
	Item number
	Non-sterile

Rx only

PLEASE NOTE: According to US federal law, in the USA this product may be bought only by a physician or hospital or upon prescription!

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