
TRANSTRANTRALE OBERKIEFERDISTRAKTOREN

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
TRANSTRANTRAL MAXILLARY DISTRACTORS

Carefully read the operating instructions and processing instructions prior to clinical application and keep them safe and at hand. The instructions contained therein must be followed.

Table of contents:

1. General notes
2. Intended use
3. Indications
4. Contraindications
5. Possible adverse effects and complications
6. Single use product
7. MR Notes
8. Application and handling
9. Pre-operative and post-operative behaviour
10. Performing the distraction
11. Decontamination, cleaning and sterilisation
12. Liability
13. Signs and Symbols

1 GENERAL NOTES

The transcranial maxillary distractors are made of titanium. The material is biocompatible, corrosion resistant and non-toxic in a biological environment. The surface is chemically passive; the material is antimagnetic.

MEDICON transcranial maxillary distractors may not be reused and are supplied in non-sterile condition. They must therefore be cleaned, disinfected and sterilised before the single use. Please observe the following notes and the processing instructions. These will ensure flawless and reliable functioning of the product.

2 INTENDED USE

MEDICON transcranial maxillary distractors are used for correction of malpositions of the maxilla by doctors and surgeons with sufficient training and experience in the field of oral and maxillofacial surgery in hospitals and clinical practices. The implants are single-use products and may be used only in combination with the specified MEDICON application instruments.

3 INDICATIONS

Indications for transcranial distraction are congenital and acquired maxillary retrognathia, such as maxillary retrognathia in conditions of orofacial cleft, maxillary hypoplasia in syndromes or conditions following severe midfacial comminuted fractures. A requirement for this procedure is the osseous continuity of the anterior nasal aperture, the canine fossa and the posterior wall of the maxillary sinus, as the distractors are fixed in place at these points. For Le Fort I osteotomy, mobilisation of the maxilla must therefore take place before the maxillary tuberosity. In other cases, osseous consistency must first be established with bone transplants.

4 CONTRAINDICATIONS

Transcranial maxillary distractors may not be used in case of:

- Patients who are unable to comply with the instructions for postoperative care. For example, due to psychological, mental or neurological problems
- Patients in unstable physical and/or mental condition
- Patients with inadequate or low-quality bone tissue, with perfusion problems or latent infections
- Proven allergy to titanium
- Material hypersensitivity, i.e. reaction of the patient to foreign bodies. Here appropriate tests are mandatory before implantation (even in case of mere suspicion!)
- Acute infections

5 POSSIBLE ADVERSE EFFECTS AND COMPLICATIONS

- In case of insufficient adaptation of the implants, bone healing can be delayed or the bone may heal insufficiently or not at all
- Discomfort, pain or abnormal sensations due to the implant
- Material hypersensitivity of the patient due to foreign objects in the form of allergic reactions
- Elevated reaction of the connective tissue in the area of the fracture and/or the implant
- Delayed or insufficient healing which can result in breakage of the implant
- Insufficient healing
- Insufficient bone formation, osteolysis, osteomyelitis, osteoporosis, inhibited revascularisation or infection that may lead to a loosening, bending, tearing or breaking of the implants and endanger the fixation or the healing of the bone

- Braking, bending, migration or loosening of the implants
- Decrease in bone density due to "stress shielding"
- Implant loosening caused by insufficient tightening of the screws
- Bone necrosis, osteoporosis and bone resorption can lead to pseudoarthrosis
- Extreme and/or multiple bending of the plates at the same point may cause the plates to break
- Increased connective tissue reactions in the area of the fracture which may present as deep or superficial early or late infections

In addition to the undesirable effects and complications already mentioned, the surgical procedure can lead to other problems, such as nerve damage, infections, pain etc. that may not necessarily be attributable to the implants.

6 SINGLE USE PRODUCT

The MEDICON implants are intended and designed for single use on one single patient only and may not be reused.

An explanted implant must never be reused. Even if the implant appears to be intact or fully functional, there may be wear, minor defects and attrition due to overuse which are not visible to the naked eye. As the impact of the forces and conditions inside the body on the stability, function and material quality of an explanted implant cannot be anticipated, reimplantation would entail the risk of premature wear or failure and is therefore unacceptable. The user will be held liable if the operating instructions and the MEDICON brochure "Transcranial Distraction" are not observed.

7 MR NOTES

Instruments

The use of medical devices in an MR environment poses a hazard. Individual medical devices must not be placed in the immediate vicinity of the equipment during the performance of MR procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Use in an MR environment may therefore pose a risk. It is recommended that patients and medical personnel be informed that the implants are not approved for use with MR equipment.

8 APPLICATION AND HANDLING

Implant selection:

The surgeon is responsible for the correct selection of implants.

The following aspects are crucial for the correct selection of implants:

- Bone defect to be treated or intended bone repositioning
- Body weight of the patient
- general health and degree of activity of the patient

breakage and therefore lead to function loss. All implants must be checked for function, deformation and mechanical integrity prior to implantation. Damaged implants must be replaced, as they can lead to partial or even total function loss. Damage leads to a reduction in the product's strength and to premature fatigue of the implant. If bone healing fails to materialise, or if healing is delayed or inadequate, the possibility that the implant will bend, break or loosen cannot be ruled out. The constant load change the implants are subjected to may lead to fatigue breakage.

It cannot be ruled out that the implants may break, loosen, bend, migrate into tissue or cause pain. Implants and instruments are developed and made so that they are compatible with each other. Using products by other manufacturers with MEDICON products may lead to unforeseeable risks and/or the risk of the material becoming contaminated as well as the instrument and implant not being sufficiently compatible. This may cause risks to patients, users and/or third parties. Whether the implants should remain in the body or be removed after complete healing is to be decided by the surgeon who takes the respective decision after weighing up the risks and benefits.

Use and handling of plates:

The plates must be adapted to the natural or planned bone structure as precisely as possible. Only use the instruments recommended by MEDICON for adapting the plate. While bending the plate, cold straining occurs under which the titanium increases in hardness but simultaneously loses deformability.

It is therefore essential to use as few bending steps as possible to achieve the desired shape of the implant. Excessive bending can lead to post-operative breaking of the plates. The success of an operation depends among other things on how the products are handled.

The plate must be bent carefully while avoiding excessive deformation. The plates must not be bent to and from more than 2 to 3 times. Damage leads to a reduction in strength and to premature fatigue. Excessively forceful use of the bending instruments can lead to visible implant damage. In such cases, a new plate must be used and adapted more carefully. Deformed plate holes for screw insertion not only lead to an increased risk of breakage in these locations but also compromise

EN

CE 0123

Rx only

Rev.: 15-12-2025

TRANSANTRAL MAXILLARY DISTRACTORS

the precise fixation of the screw head in the plate. The plates must therefore be bent carefully.

Deformed plates must always be checked for notches, deformed screw holes or other mechanical damage before being used on the patient. Damaged plates must be replaced. Always choose the plates that best suit the specific situation and indication. Should it be necessary to shorten a plate, only use the instruments provided for this purpose.

The plate must be cut between the screw holes. Ensure that the cut-off part of the plate is not thrust towards the patient, the surgeons or others to prevent putting them at risk. The cutting edges of a shortened plate must be deburred to prevent tissue injury or irritation.

Use and handling of screws:

Self-tapping screws do not require drilling a threaded hole. Only use the drills specified by MEDICON to drill holes so that the correct hole diameter required for the screws can be achieved.

Only use a drill with sharp cutting edges.

There is a risk of heat damage on the bone during drilling. Therefore always ensure sufficient cooling during the drilling process and use low drilling speeds. Excessive force during the drilling process can result in fracture of the tool, which can be dangerous for the surgeon, patient and third parties. If the screw cannot be securely anchored in the bone due to an incorrect or worn-out drill hole, the emergency screws intended for such circumstances are to be used.

The screwdriver must be inserted into the screw head with axial pressure to ensure that the blade sits tightly on the head. This guarantees the correct longitudinal alignment of screw and screwdriver and prevents slipping of or damage to the screw head or the blade.

If excessive force is used when fixing the screws, the screws may break during surgery or the screw heads may tear off. Wear and tear effects to the screwdriver blade impair the solid connection between the blade and screw. In this case, replace the blade with a new one.

Upon completion of the implantation, the secure connection between all screws and plates must be checked. If necessary, retighten the screws.

Before removing the implant, all screw heads must be thoroughly cleaned with a scalpel or another suitable instrument so that the screwdriver sits perfectly on the screw head.

Self-tapping screws:

When using self-tapping screws, you are usually not required to drill a hole. If the bone is dense and hard, it may still be necessary to drill holes, e.g. in the calvarium or mandible. Self-tapping screws are not suitable for small bone fragments, as the axial pressure may lead to displacement of the fragments or breakage of the screw.

Use and handling of instruments:

The instruments intended for application of the system are subject to wear and mechanical stresses even when used normally, but especially if used too forcefully.

In order to prevent failure or mechanical damage to the instruments during surgery, these parts should be checked before use to ensure that they are mechanically intact, that there are no deformations, and that the parts are fully functioning. Damaged instruments must not be used and are to be replaced.

Only use accessories and instruments specifically approved by MEDICON in order to avoid risks in connection with the compatibility of the products.

Use and handling of the locking rings:

Immediately after insertion of the distractors, the extensions (shafts with hinges) must be secured in the frontal maxillary vestibule until the end of the distraction phase. For this, before insertion into the vestibule, the larger ring of the double ring is pushed over the extension and gently compressed with flat pliers, so that the ring can no longer slide over the ends of the extension (Figure C).

However, the extension must still be able to move freely, i.e. the ring must not be pressed firmly onto the extension. The small ring is then sutured into the vestibule with non-resorbable suture material. After completion of the distraction phase, the extension is removed.

9 PRE-OPERATIVE AND POST-OPERATIVE BEHAVIOUR**Pre-operative:**

Before using the product, the surgeon should thoroughly discuss the desired operation result with the patient. Particular attention must be given to post-operative behaviour and potentially required follow-up care.

Post-operative:

The patient must be instructed to immediately inform the surgeon about unusual in situ changes. The surgeon considers possible implications such as implant failure and discusses required measures for further healing with the patient.

Whether the implants should remain in the body or be removed after complete healing is to be decided by the surgeon who takes the respective decision after weighing up the risks and benefits.

If healing fails to materialise, or is delayed or inadequate, the possibility that the implant will bend, fail or break cannot be ruled out. Therefore adequate immobilisation of the fracture or osteotomy line must be ensured until the bone grows together again tightly. The constant load change the implants are subjected to may lead to fatigue breakage. The treating physician takes decisions regarding the type, duration and intensity of physical activities after the surgery.

A "LOT" number (batch number) is indicated on the label of the packaging. To ensure 100% traceability, we recommend noting down the "LOT" number in the patient records.

Disposal:

Observe the applicable national laws when disposing of the MEDICON products and accessories.

10 PERFORMING THE DISTRACTION (FIGURE A)**Adjusting the extension screws and distraction process:**

To adjust the extension and perform the distraction, the double ended screwdriver must be used. The extension screw is operated with the square tip of the screwdriver. For distraction, the Allen tip is used. The distraction is performed by anticlockwise rotation. One turn corresponds to a distance of 0.4 mm. Please note the marking on the screwdriver. The patients or their family members, respectively, are to be informed comprehensively about the distraction process and the handling of the distractors, especially if the distraction is to be continued outside the hospital or clinic. The patient is to be provided with a schedule for the activation of the distractors, indicating the distractions to be performed

11 DECONTAMINATION, CLEANING AND STERILISATION**Opening the distractor:**

The threaded distraction rod ④ and the extension screw ① must be fully unscrewed for processing (Figure B).

Cleaning solutions to which hydrogen peroxide has been added and/or washing solutions with high alkalinity can cause colour changes. This may cause the coding function to be lost. Only cleaned and disinfected instruments may be sterilised.

LIMITATIONS ON REPROCESSING

Frequent reprocessing has little effect on these instruments.

PLACE OF USE

Implants with contaminations caused by blood and/or secretions must not be placed back into the implant storage tray. They must be disposed of.

PREPARATION FOR DECONTAMINATION

The implant storage tray must be appropriately placed on a washer-proof instrument carrier for rinsing. The instrument carriers (e.g. wire sieve trays) must be such that subsequent cleaning in the washing and disinfection device is not impaired by any areas that are inaccessible to ultrasonics or rinsing.

PRECLEANING

Always use fresh and unused cleaning solution for manual pre-cleaning in the ultrasonic bath to prevent cross-contamination.

- Place the implant storage tray in the ultrasonic bath and treat with ultrasound for 15 minutes at 40° C using 0.5% alkaline cleaner (Neodisher MediClean fort).
- Remove the implant storage tray and thoroughly rinse with cold water.

AUTOMATED CLEANING

To prevent cross-contamination during machine processing, always clean and disinfect the implant cassette including the implants to be reprocessed in the washing and disinfecting device (WD) separate from instruments contaminated with blood or secretion. The WD must comply with the requirements according to DIN EN ISO 15883-1.

- Prewashing: 3 minutes with softened water, without additives
- Emptying
- Cleaning: with softened water, heating to 55° C and 5 minutes washing/cleaning, adding of the cleaning agent at 45° C, alkaline cleaning agent (Neodisher MediClean fort), dosage 0.5%
- Emptying
- Neutralisation: 3 minutes with warm water (>40° C) including the addition of neutraliser, dosage 1 ml/l
- Emptying
- Final rinse: 2 minutes with warm, fully demineralised water (>40° C) (without any other additives)
- Emptying

DISINFECTION

Thermal disinfection A0 value 3000:

- Fully demineralised water, thermal disinfection is carried out at temperatures >80° C with a corresponding application time according to the A0 concept of standard DIN EN ISO 15883-1 and guideline DGHK, DGSV and AKI (e.g. A0 3000 = 90° C and 5 minutes application time). The operator is responsible for the A0 value to be implemented..

DRYING

Sufficient drying must be ensured by the WD (e.g. 60° C, 30 minutes). The instruments must be removed from the WD immediately after completion of the cleaning and disinfection programme.

If required, use of pressurised air for drying is recommended due to its good and quick effect (RKI recommendation).

MANUAL CLEANING

- Rinse the instruments for at least one minute with cold tap water to remove visible contamination and move the movable parts while rinsing the instrument. Use a soft bristle brush and lumen brush to remove dirt during the rinsing process and move the movable parts. Use a syringe to rinse the lumens (interior spaces, threads and holes) with rinsing water
- Prepare enzyme cleaner Enzol® with 1 ounce per gallon of tap water and soak the instruments in it for at least one minute
- Brush off the instruments after the soaking process using a soft bristle brush and lumen brush for at least one minute to remove residual dirt
- Use a syringe to clean interior spaces, threads and holes using an enzyme cleaner solution
- Remove the instruments from the enzyme cleaner solution and rinse them under tap water to remove residues of the cleaning solution
- Prepare the neutral cleaner Valsure® Neutral in an ultrasound unit using ¼ ounce per gallon of tap water (37-40° C). Immerse the instruments in the ultrasound unit and move them
- Use a syringe to rinse the lumens (interior spaces, threads, holes) using a clean-rinsing solution and treat the instruments with ultrasound for at least 10 minutes.
- Remove the instruments from the ultrasound unit and rinse them under cold demineralised water to remove residues of the cleaning solution
- Fully dry the instruments using pressurised air and a lint-free disposable cloth.

MAINTENANCE, CONTROL AND INSPECTION

Following cleaning/disinfection, the implants must be macroscopically clean, i.e. free from visible dirt and residues. The inspection is performed visually. Insufficiently cleaned implants must be cleaned again and then rinsed and dried sufficiently. Deformed or damaged implants must be sorted out and disposed of as safe use cannot be guaranteed.

PACKING

The implant storage tray must be packed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868
- DIN EN ISO 11607
- suitable for steam sterilization (steam permeable)
- adequate temperature resistance to 138° C

Sterilization accessories and sterilization packaging must be appropriate for the content of the package and for the sterilization process.

STERILIZATION

Use the following sterilization procedures subject to the applicable national requirements:

- Fractionated vacuum method with triple fractionation and sufficient product drying.
- Steam sterilization in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes holding time at 134° C.

It is essential that a SAL (Sterility Assurance Level) of 10⁻⁶ is reached

STORAGE

Processed sterile implants must be stored in appropriate reusable sterilisation containers which are dry, protected from dust and low-germ, in a dark and cool place, free from vermin.

To prevent condensation, extreme temperature fluctuations should be avoided during storage. No chemicals may be stored together with implants.

Walls, floors and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must have a ground clearance of at least 30 cm.

The permissible on-site storage period depends on the type of sterile barrier system used and on the storage conditions. The permissible storage period is to be specified by the operator.

ADDITIONAL INFORMATION ON REPROCESSING

A validated machine cleaning and disinfection procedure is always to be preferred to manual cleaning because of the greater safety of the process. Thorough cleaning also helps to maintain value and is a prerequisite for successful sterilisation.

Please take note of the following points during automated processing:

- For effective automated processing, load the wire mesh baskets in such a way that they can easily be mechanically rinsed. Mesh baskets must not be overloaded
- Avoid "rinsing shadows" caused by large instruments

The times and temperatures stated in these reprocessing instructions are minimum requirements which must be adhered to. Should any downward deviation be needed for procedural reasons, this must be validated by the operator.

Exceeding the specified times and temperatures is generally possible. However, this will lead to increased stress on the material, which can cause premature wear. We accept no responsibility for the use of other sterilisation procedures.

INFORMATION ON VALIDATING THE PROCESSING PROCEDURE

Validation is carried out with the following devices, materials and chemicals: Cleaning and disinfection equipment: Type Miele Desinfektor G 7735 CD Trolley for surgical instruments

Alkaline detergent: neodisher® MediClean forte Dr. Weigert GmbH & Co. KG
Neutraliser: neodisher® Z Dr. Weigert GmbH & Co. KG
Ultrasonic bath: Sonorex DIGITEC DT 156 BH
Sterilisers: MMM Vaculab 969 S 3000 MMM Selectomat S 3000 Stiefenhofer KS 666-2ED H+P Varioclav 400E

NOTE

The person processing the instrument bears the responsibility that the processing actually performed with the equipment, materials and staff at the reprocessing facility achieves the desired results. Validation and routine monitoring of the process are normally required.

If the devices, materials and chemicals described above are not available, it is up to the persons processing the instrument to have their procedures appropriately validated. Observe the notes and regulations of the relevant national guidelines and standards. MEDICON eG reserves the right to make changes to these instructions due to new findings.

12 LIABILITY

MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

13 SIGNS AND SYMBOLS

	Manufacturer
	Date of manufacture
	Lot number
	Item number
	Non-sterile
	Warning
	Consult instructions for use
	MR unsafe
	CE marking according to Directive 93/42/EEC
	Medical device
	Prescription only

Federal law and regulations in the United states restrict the sale of this device to a physician or hospital or on prescription!

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Detergente alcalino: neodisher® MediClean forte Ditta Dr. Weigert GmbH & Co. KG
Neutralizzatore: neodisher® Z Ditta Dr. Weigert GmbH & Co. KG