Rx only

Deutsch

Rev.: 04-12-2025

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Medicon implants are manufactured from pure titanium or Titanium alloy. Both materials are biocompatible, corrosion-proof and nontoxic in the biological environment. The surface is chemically passive and the material is antimagnetic. The implants are supplied in unsterile condition.

Table of contents

1. Indications
2. Contraindications
3. Possible side effects and complications
4. Single Use Product
5. MR Notes
6. Implant selection
7. Application/Handling of plates
8. Application/Handling of screws
9. Instruments
10. Decontamination, cleaning and sterilization
11. Note

These operating instructions must be read carefully and completely prior to clinical application of the products.

MEDICON eG as the supplier of these products is not liable for damage caused by, or arising from improper use or handling, especially due to non-observance of the operating instructions or improper cleaning or maintenance of the products. The Medicon plate systems may be used only by sur-geons of proper training and experience.

1. Indications

The Medicon CMS plate system was developed for the fixation of very thin bones of the calvaria and the facial skull, especially in infant and pediatric surgery. The plates of 0.4 mm thickness are intended for the fixation of osseous structures of the orbita regions, the frontonasal transition and other bone structures under thin soft tissue cover. The 0.55 mm plates are specially designed for segmental osteotomies and forward repositioning of the maxilla. The implants are intended to support normal bone healing, but not to replace normal anatomical structures or to bear the body weight in cases of insufficient bone healing.

2. Contraindications

- Patients who are unable to follow the instructions for postoperative care. Such inability may be caused by psychological/mental or neurological conditions of the patient.
- Patients with unstable physical and/or mental health conditions.
- Patients with insufficient or poor-quality bone tissue, circulatory disorders or latent infections.
- Oversensitivity to materials and foreign-body reactions. Preimplantation tests must be carried out even if such oversensitivity or reactions are only suspected.
- Acute infections.

3. Possible side effects and Complications

In many cases, possible complications are attributable not to the implants as such, but to incorrect patient selection, lack of practice or inadequate repositioning of the bone fragments.

- Insufficient adaptation of the implants can lead to bone union delay or failure.
- Complaints, pains, abnormal sensations or palpability of implants.
- Material/foreign-body sensitivity of the patient with allergic reactions.
- Stronger connective-tissue reaction in the fracture region or in the vicinity of the implant.
- Delayed or insufficient fracture healing leading to possible implant breakage.
- Inadequate healing.
- Insufficient bone formation, osteolysis, osteoporosis, reduced revascularization or infection, which can lead to loosening, bending, breakage or rupture of the implant.
- Breakage, bending, migration or loosening of the implant.
- Reduced bone density due to "stress shielding".
- Loosening of the implant as a result of insufficient screw tightening.
- Osseous necrosis, osteoporosis, inhibited revascularization, bone-resorption, and poor bone formation, which can cause premature loss of fixation leading to non-union.
- Extreme or repeated bending of the plates in the same place can result in plate breakage.
- Unstable comminuted fractures can lead to an increased connective-tissue reaction in the fracture region.

This surgical procedure may cause not only the above-mentioned side effects and complications but also problems such as injuries to nerves, in-fectious, pain etc., which are not necessarily caused by the implant.

4. Single use product

Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

5. MR Notes

Instruments

The use of medical devices in an MR environment poses a hazard. Individual medical devices must not be placed in the immediate vicinity of the equipment during the performance of MR procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Use in an MR environment may therefore pose a risk. It is recommended that patients and medical personnel be informed that the implants are not approved for use with MR equipment.

6. Implant Selection

The surgeon is responsible for the correct selection of implants.

The following aspects are critical for the correct selection of implants:

- bone defect to be treated / intended bone repositioning,
- body weight of the patient,
- general health condition, degree of activity of the patient.

Incorrect implant selection can lead to premature implant loss and cause loosening, deformations and fractures of the implants. Failure to use the proper component to maintain adequate blood supply and provide rigid fixation may result in loosening, bending, or fracturing of the product and/or bone. The success of an operation also depends on how the implants are handled. The plates must be bent carefully and cautiously, avoiding repeated and excessive bending. Scratches and other damage can reduce the stability of the product and result in early fatigue of the implant. Medicon implants and instruments have been designed for being used together. Using implants or instruments from other manufacturers in conjunction with Medicon products entails unpredictable risks since the products will not be matched to each other. To avoid such risks, only Medicon products that are intended for combination with each other may be used together. Prior to applying the products, the surgeon must thoroughly discuss with the patient the intended outcome of the operation. In such discussions, special attention must be paid to postoperative care measures that may become necessary. The patient must be informed that, due to its limited stability, the fixation must not put under the load of the full body weight. The patient must also be informed that non-observance of this advice can put at risk any safe bone union. The patient must be instructed to inform the surgeon immediately if they notice unusual changes in the situs. The patient must be monitored carefully if there is a manifest change in the fixation region. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing. In cases of failed, delayed or insufficient fracture union, the incidence of implant bending, implant failure or implant breakage cannot be excluded. Therefore, immobilization of the fracture site must be ensured until there is solid bone fusion. The incessant load changes, to which the implant will be subjected, can lead to fatigue fractures. The implants can be removed once the fracture has healed completely. An implant left in the body after complete healing may act as a load-bearing element and can contribute to an increased risk of re-fracturing in active patients. For this reason, we accept liability for the implants only until healing is complete. We will not accept liability for any damage or injuries caused by implants that were not removed as soon as healing of the treated fracture was complete. The surgeon is responsible for deciding whether the implants shall remain in the body after complete healing of the fracture or whether they shall be removed. For his decision, the surgeon has to weigh risks and advantages of leaving the implants in the body. It cannot be excluded that implants will break, loosen, corrode, migrate through tissue and cause pain. The packaging label carries a "Lot No.". We recommend copying this "Lot No." to the patient records as this number allows tracing the production hi-story of the implant back to the raw material.

7. Application and Handling of plates

The plates need to be adapted to the natural or planned bone contours as precisely as possible. Only use the instruments provided by the manufacturer for the purpose of adapting the plates. During the plate bending procedure the plates are subjected to cold stresses, under which the titanium hardens and its ductility is decreased. Therefore it is absolutely necessary to achieve the required plate shape in as few bending steps as possible. Excessive bending can lead to postoperative plate breakage. The plates must not be bent back and forth more than 2 – 3 times. Overaggressive application of the bending instruments can result in visible damage to the implant. When such damage has occurred, a new plate must be selected and adapted more cautiously. Deformed plate holes for the screws entail not only an increased risk of breakage at the holes; they also compromise the precise seating of the screw head in the plate. Therefore the plates must be bent cautiously. Adapted plates must be inspected for notches, deformed screw holes and other mechanical damage before they are applied to the patient. Always select the plates that provide the closest fit to the specific situation and indication. If a plate needs to be shortened, only use the instruments supplied for this purpose. When shortening a plate, care must be taken that the portion of the plate that is cut off is not flung towards the patient or towards third parties during the cutting procedure. The cut edges of a shortened plate must be deburred to prevent tissue injuries or irritations. Bending templates must not be implanted under any circumstances.

8. Application and Handling of Screws

Self-tapping screws do not require pre-cutting of threads. To apply the drill holes, only use the drills designated by Medicon for the respective screw so that the correct drill hole diameter for the screw is obtained. Only use drill bits with sharp cutting edges! During drilling, there is the risk of heat damage to the bone. Therefore, only work at low drill speed and always ensure sufficient cooling during the drilling procedure. Excessive application of force during the drilling procedure can lead to drill breakage, endangering the patient, the

surgeon and third parties. Use the emergency screws intended for this purpose if, due to a failed drilling effort or a worn drill hole, the original screw cannot be anchored firmly in the bone. To ensure that the screwdriver blade is firmly seated in the screw head, axial pressure must be applied when the screwdriver is inserted into the screw head. This ensures the correct axial alignment of screw and screwdriver, and prevents slipping of the screwdriver blade and damage to the screw head. Application of excessive force when tightening the screws can lead to intraoperative screw breakage. Wear of the screwdriver blade will impair the firm connection between blade and screw. Whenever such wear is detected, the screwdriver blade must be replaced by a new blade. Once the implantation has been completed, every screw and plate connection must be checked for stability. If necessary re-tighten the screws. Prior to implant removal, all screw heads must be cleaned carefully, using a scalpel or another suitable instrument, so that the screwdriver can be optimally seated in the screw head.

9. Instruments

The instruments intended for application of the system are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force. To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional.

10. Decontamination, cleaning and Sterilization

Implants of titanium and titanium alloys must not be cleaned with a process that forms active oxygen by dosage of hydrogen peroxide (e.g. Miele oXI Vario). Only cleaned and disinfected implants may be sterilized.

10.1 Preparations for decontamination

The implants must be sorted into the implant tray for cleaning in the cleaning and disinfection device.

10.2 Cleaning

Suitable pH-neutral or alkaline cleaners can be used. Selection of the cleaning agent is based on the properties of the implants and national guidelines and recommendations. The cleaning and disinfection device must conform to the requirements of DIN EN ISO 15883-1.

- Pre-rinse 1: 1 minute with demineralized cold water with no additives;
- Drain;
- Pre-rinse 2: 3 minutes with demineralized cold water with no additives;
- Drain;
- Cleaning: use demineralized water, heat to 55°C and wash and clean for 5 minutes, add cleaning agent at 45°C, alkaline cleaning agent, dosage 0.5%;
- Drain;
- Neutralization: 3 minutes with warm water (>40°C) with neutralizer, dosage 1 ml/l;
- Drain;
- Final-rinse: 2 minutes with warm demineralized water (>40°C) (no additives).
- Drain

The manufacturer's recommended parameters for the cleaning agent concentration, temperature and soaking time must be followed and automatic dosage devices must be adjustable.

10.3 Thermal disinfection

Carry out thermal disinfection, using fully desalinated water, at a temperature of 80 – 95 °C for an exposure time acc. to A0 draft, EN ISO 15883-1.

10.4 Drying

Sufficient drying time must be provided by the cleaning and disinfection device. If necessary, the use of compressed air for drying is recommended because of its effective and quick effect (RKI recommendation).

10.5 Maintenance, control and inspection

The implants must be macroscopically clean, i.e. free of visible residues, after cleaning and disinfection. The implants must be visually inspected in this respect. Critical areas require particularly careful inspection. Inadequately cleaned parts must be cleaned and thoroughly rinsed again.

10.6 Sterilization

Use the following sterilization procedures subject to the applicable national requirements:

- Fractionated vacuum method with triple fractionation and sufficient product drying
- Steam sterilization in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes holding time at 134°C

10.7 packing

The implants must be packed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868 and DIN EN ISO 11607
- suitable for steam sterilization (steam permeable)
- adequate temperature resistance to 138°C

Sterilization accessories and sterilization packaging must be appropriate for the content of

the package and for the sterilization process.

10.8 Storage

The implants must be stored in the implant tray in a dry, dust-protected, low-microbe, dark and cool place protected from vermin. To prevent the formation of condensation major temperature fluctuations should be avoided. Chemicals must not be stored with implants. The approved storage period depends on the type of sterile barrier system and the storage conditions. The approved storage period must be defined by the owner.

11. note

The user is responsible for validating the recommended parameters to take into account any differences with regard to sterilization chambers, packaging methods and loading configurations. The same applies if a procedure that differs from our recommendations is used. It is essential that a SAL (Sterility Assurance Level) of 10⁻⁶ is reached. Please contact MEDICON eG if you have further questions concerning the implants or instruments.

Caution: for USA, federal law restricts this device to purchase by or on behalf of physician or hospital.

12. SIGNS AND SYMBOLS

	Manufacturer
	Date of manufacture
	Lot number
	Item number
	Non-sterile
	Warning
	Consult instructions for use
	MR unsafe
	CE marking according to Directive 93/42/EEC
	Medical device
	Prescription only

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