## Deutsch

C0123

Rx only

Rev.: 11-12-2025

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**Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.**

MEDICON implants are manufactured from titanium. This material is biocompatible, corrosion-resistant and non-toxic in the biological environment. The surface is chemically passive and the material is antimagnetic. The implants are supplied in unsterile condition.

#### Table of contents:

1. Indications
2. Contraindications
3. Possible side effects and complications
4. Single Use Product
5. MR Note
6. Implant selection
7. Application/Handling of plates
8. Application/Handling of screws
9. Instruments
10. Instructions on reprocessing
11. Signs and symbols

MEDICON eG, as the original distributor of these products, does not accept any liability for immediate or consequent damage caused by improper use or handling, especially due to non-observance of the operating instructions or inappropriate care or maintenance of the products. The MEDICON plate systems may be used only by sufficiently trained and experienced surgeons in cranial neurosurgery.

#### 1. Indications

The MEDICON Low Profile implants, consisting of plates and screws, are intended for use in surgical interventions on the osseous structures of the skull. Drill hole cover plates are to be used for covering trepanation holes and fixing the bone flap. Indications include:

- fractures
- fixation of the bone flap

The Low Profile system is designed for use in the cranium for craniotomy fixation specifically in the front area to limit palpability. The implants have been designed for supporting the normal bone healing process, not for replacing normal body structures or bearing the body weight in cases of insufficient bone healing.

**Attention:** The Low Profile screws and plates have been specially developed for thin applications on the bone. For this reason, the Low Profile screws are suitable only for applications at which low stresses occur.

#### 2. Contraindications

- Patients that are unable to follow the instructions for postoperative care. This may be caused e.g. by psychological/mental or neurological problems of the patient.
- Patients with inadequate bone tissue, or bone tissue of insufficient quality, and patients with circulatory disorders or latent infections.
- Extreme material sensitivity, i.e. foreign-body reaction of the patient. For this case, appropriate tests must be carried out under all circumstances prior to implantation (even when there is only a suspicion of such condition).
- Acute infections.

#### 3. Possible side effects and complications

- Complaints, pain, abnormal sensations or palpability of the implant.
- Extreme material sensitivity of the patient (foreign-body sensitivity) manifested by allergic reactions, etc
- Heightened reaction of connective tissue in the region of the fracture and/or implant.
- Insufficient bone formation, osteolysis, osteomyelitis, osteoporosis, inhibited revascularization or infections that can result in loosening, bending, rupture or breakage of the implants.
- Delayed or insufficient fracture union, which can result in implant breakage.
- Bone necrosis.
- Inadequate healing.
- If the implants were not adapted properly, bone fusion can be delayed, or fusion is insufficient or does not occur at all.
- Breakage, bending, migration or loosening of the implant.
- Reduced bone density due to stress shielding.

Apart from the side effects and complications mentioned above, the surgical intervention can also lead to problems such as nerve injuries, infections, pain, etc.,

which are not directly attributable to the implant. The possible complications are often the result of incorrect selection of the patient, lack of exercise or inadequate repositioning of the bone fragments, rather than being caused by the implant itself. Extreme bending of the plates, or repeated bending at the same place, can lead to plate fractures. The implants can come loose if the screws are not tightened sufficiently.

#### 4. Single use product

Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wear out or failure. The user is liable for non-compliance with the instructions for use.

#### 5. MR Note

##### Instruments

The use of this medical devices in the vicinity of an MR system poses a potential hazard. The individual medical devices must not be located in a close proximity to the equipment during the procedure.

##### Implantats

The implants are not MRI-safe and have not been tested for MRI safety. Therefore, use in an MRI environment may present a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MRI systems.

#### 6. Implant Selection

**The surgeon is responsible for the correct selection of implants.**

The following aspect are crucial for the correct selection of implants:

- Bone defect to be treated
- Body weight of the patient
- general health and degree of activity of the patient.

Incorrect implant selection can lead to premature implant loss and loosening, deformations or fractures of the implants. Among other factors, the success of the intervention depends on how the implants are handled. The plates must be bent carefully and cautiously, avoiding repeated and extensive deformation. Damage and scratches reduce the stability of the product and result in early fatigue of the implant material. In cases of failed, delayed or insufficient fracture union, the incidence of implant bending, implant failure or implant breakage cannot be excluded. Therefore, immobilization of the fracture site must be ensured until there is solid bone fusion. Incessant load changes on the implant can result in fatigue fractures. The implants and instruments are adapted to each other. Any use of implants and instruments from other manufacturers in combination with MEDICON products entails unpredictable risks, since the products are not adapted to each other. To avoid risks, only such MEDICON products that are designated for combination with each other may be combined with each other. Prior to applying the product, the surgeon must thoroughly discuss the intended surgical result with the patient. In such discussions, special attention must be given to the postoperative aspects and aftercare, which may become necessary. The patient must be informed that, due to the limited stability of the fixation, the fixation site must not be put under the full body weight. The patient must also be informed the non-observance of this advice can put at risk any safe bone union. The patient must be instructed to inform, without delay, the operating surgeon of unusual changes at the operating site. The patient must be monitored carefully if there is a manifest change in the fixation region. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing. The implants can be removed after bone union is completed. An implant left in the body after complete healing may act as a load-bearing element and can contribute to an increased risk of refracturing in active patients. For this reason, too, we accept liability for the implants only up to complete healing. We will not accept liability for any damage or injuries caused by implants that were not removed after complete healing of the treated fracture. The surgeon is responsible for deciding whether the implants shall remain in the body after complete healing of the fracture or whether they shall be removed. For his decision, the surgeon has to weigh risks and advantages of leaving the implants in the body. It cannot be excluded that implants break, loosen, corrode, migrate in tissue and cause pain. The packaging label carries a lot number. We recommend transferring this lot number to the patient records as this number allows tracing the production history of the implant back to the raw material.

#### 7. Application/Handling of plates

The plates have to be adapted to the natural or planned bone contours as precisely as possible. For adapting the plates use the instruments provided for this purpose by the manufacturer. During the plate bending procedure, there occur cold stresses, under which the titanium hardens while, at the same time, its ductility is decreased. Therefore it is absolutely necessary to achieve the required plate shape in as few bending steps as possible. Excessive bending can lead to postoperative plate breakage. Overly aggressive application of the bending instruments can result in visible damage to the implant. When such damage has occurred, a new plate must be selected and adapted more cautiously. Deformed plate holes for the screws entail not only an increased risk of breakage at the location of the holes; they also compromise the precise seating of the screw head in the plate. Therefore the plates must be bent carefully under all circumstances. Deformed plates must be inspected for notches, deformed screw holes and other mechanical damage prior to application. Always select the plates that provide the closest fit to the specific situation and indication. If a plate needs to be shortened, only use the instruments supplied for this purpose. The plate must be cut between the holes. Care must be taken that the portion of the plate that is cut off is not flung towards the patient or towards third parties during the cutting procedure. The edges of a shortened plate must be deburred to prevent tissue injuries or irritations.

#### 8. Application/Handling of Screws

Self-tapping screws do not require pre-cutting of threads. To apply the drill holes, only use the respective drill designated by MEDICON for the specific screw so that the correct drill hole diameter for the screw is obtained. Only use drill bits with sharp cutting edges. During drilling, there is a risk of heat damage to the bone. Therefore, only work at low drill speed and always ensure sufficient cooling during the drilling procedure. Excessive application of force during the drilling procedure can lead to drill breakage, endangering the patient, the surgeon and third parties. If, due to a failed drilling effort or a slack drill hole, the screw cannot be anchored firmly in the bone, the emergency screws intended for this purpose must be used. To ensure that the screwdriver blade is firmly seated in the screw head, axial force must be applied when inserting the screwdriver into the screw head. This ensures the correct axial alignment of screw and screwdriver, and prevents slipping of, or damage to the screw head or the screwdriver blade. Application of excessive force in tightening the screws can lead to screw breakage during the intervention. Wear of the screwdriver blade will impair the firm connection between blade and screw. Whenever such wear is detected, the screwdriver blade must be replaced by a new blade. Once the implantation has been completed, the secure connection between every screw and plate must be checked. If necessary retighten the screws. Prior to implant removal, all screw heads must be cleaned carefully, using a scalpel or another suitable instrument, so that the screwdriver can be optimally seated in the screw head.

#### 8.1 Self-tapping screws

It is normally not necessary to drill a pilot hole for self-tapping screws. Still, in dense and hard bones, e.g. at the calvarium, it may be necessary to prepare drill holes. Self-tapping screws are unsuitable for small bone fragments as the axial pressure on them can lead to fragment shifts or screw breakage.

#### 9. Instruments

The instruments intended for application of the system are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force. To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional.

#### 10. Instructions on reprocessing

(Mechanical cleaning, alkaline, pH >10; with pre-cleaning and ultrasonic treatment).

**The following prion-specific protective measure (vcJD) is indicated when preparing the instrument.**

In the event of a diagnosis of **definite or probable** vcJD:

- If the use of disposable products is not possible, the instruments used (which had become contaminated or whose contamination cannot be ruled out) must be disposed of as waste for incineration.

#### In suspected cases:

- After the validated cleaning and disinfection process and until the diagnosis has been clarified, keep the instruments physically isolated and in an alkali-resistant and steam-sterilizable container (e.g. V4A stainless steel, material number 1.4401) with a tightly closable lid, or in a suitable

disposable container. Store under dry conditions, safeguard with a lead seal and label permanently.

#### If vcJD infection is ruled out:

- Re-use after final treatment. Otherwise the instruments used (that had become contaminated or where contamination cannot be ruled out) must be disposed of as waste for incineration.

In the event that vcJD infection is **not apparent**:

- the instrument(s) that have had (or could have had) possible direct contact with the CNS, or with the ocular fundus or lymphatic tissue must be subjected to the two-stage decontamination/preparation procedure described below.

**Instruments made of stainless steel may not be stored in physiological saline (nacl) solution, as prolonged exposure will lead to corrosion problems such as pitting and stress corrosion. only cleaned and disinfected instruments may be sterilised.**

#### 10.1 limitations on reprocessing

Frequent reprocessing has little effect on these instruments. Their service life normally comes to an end as a result of wear and damage by use.

#### 10.2 place of use

Implants with contaminations caused by blood and/or secretions must not be placed back into the implant storage tray. They must be disposed of.

#### 10.3 preparation for decontamination

The instruments must be placed on appropriate instrument carriers so that they can be rinsed in the washing and disinfection device. The instrument carriers (e.g. sieve trays) must be such that subsequent cleaning in an ultrasonic device or in the washing and disinfection device is not impaired by any areas that are inaccessible to ultrasonics or rinsing.

#### 10.4 precleaning

- Place the instruments in an ultrasonic bath for 15 min at 40°C with a 0.5% alkaline cleaner (pH 11.4 - 11.9) and subject to ultrasound.
- Remove the instruments and rinse with cold, softened water.

A high level of contamination in the ultrasonic bath will impair the cleaning effect and promote corrosion. The cleaning solution must be regularly replaced, depending on the conditions of use. The decisive criterion is optically recognisable contamination. In any case, the bath must be changed frequently, at least once a day. Comply with the applicable regulations for processing medical products in the country concerned.

#### 10.5 cleaning

The cleaning and disinfection device must conform to the requirements of DIN EN ISO 15883-1.

- Pre-rinse 1: 1 minute with cold demineralised water with no additives;
- Drain;
- Drainin 2: 3 minutes with cold demineralised water with no additives;
- Drain;
- Cleaning: use demineralised water, heat to >55°C, wash and clean for 10 minutes, add cleaning agent at 45°C, alkaline cleaning agent, dosage 0.5% (pH 11.4 - 11.9);
- Drain;
- Neutralization: 3 minutes with warm deionised water (>40°C) with neutralizer, dosage 1 ml/l;
- Drain;
- Final rinse: 2 minutes with warm deionised water (>40°C) (no additives);
- Drain

Carry out thermal disinfection, using fully desalinated water, at a temperature of 80 – 95 °C for an exposure time acc. to A0 draft, EN ISO 15883-1.

#### 10.6 Drying

Sufficient drying time must be provided by the cleaning and disinfection device. The instruments must be removed from the cleaning-disinfection device immediately after the cleaning and disinfection program has finished. If necessary, the use of compressed air for drying is recommended because of its effective and quick effect (RKI recommendation).

#### 10.7 Maintenance, control and inspection

The instruments must be macroscopically clean, i.e. free of visible residues, after cleaning and disinfection. The instruments must be visually inspected in this respect.

Inadequately cleaned parts must be cleaned, thoroughly rinsed and dried again. Defective instruments must be replaced, as they cannot fulfil their function, or cannot fulfil it safely enough.

### 10.8 Packing

The instruments must be packed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868
- DIN EN ISO 11607
- suitable for steam sterilization (steam permeable)
- adequate temperature resistance to 138°C

Sterilization accessories and sterilization packaging must be appropriate for the content of the package and for the sterilization process.

### 10.9 Sterilization

Use the following sterilization procedures subject to the applicable national requirements:

- Fractionated vacuum method with triple fractionation and sufficient product drying.
- Steam sterilization in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes holding time at 134°C.

It is essential that a SAL (Sterility Assurance Level) of 10<sup>-6</sup> is reached

### 10.10 Storage

Reprocessed sterile instruments must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments. The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority

### 10.10 Further Information on Reprocessing

A validated machine cleaning and disinfection procedure is always to be preferred to manual cleaning because of the greater safety of the process. Thorough cleaning also helps to maintain value and is a prerequisite for successful sterilisation. Please take notice of the following points regarding mechanical processing:

- Load the wire mesh baskets in such a way that they can easily be mechanically rinsed. Mesh baskets must not be overloaded.
- Avoid „rinsing shadows“ caused by large instruments.

The times and temperatures stated in these reprocessing instructions are minimum requirements, they must be reached. Should any downward deviation be needed for procedural reasons, this must be validated by the operator. Exceeding the specified times and temperatures is basically possible. However, this will lead to increased stress on the material, which can cause the instruments to wear out prematurely. We accept no responsibility for the use of other sterilisation procedures.

### 10.11 Information on validating the reprocessing procedure

Sterilisers: MMM Vaculab 969 S 3000, MMM Selectomat S 3000, Stiefenhofer KS 666-2ED, H+P Varioclav 400E

### 10.13 note

The person preparing the instrument bears responsibility that the preparation actually performed with the equipment, materials and staff at the reprocessing facility achieves the desired results. Validation and routine monitoring of the process are normally required. If the abovedescribed apparatus, materials and chemicals are not available, it is up to the persons preparing the instrument to have their own procedures validated. Please refer to the instructions and requirements of the relevant national regulations and standards as well as the instructions for use included with the medical device, if applicable. Medicon reserves the right to make alterations due to recent findings.

 **Caution: for USA, federal law restricts this device to purchase by or on behalf of physican or hospital.**

Rx only

### 11. Signs and Symbols

|                                                                                   |                                             |
|-----------------------------------------------------------------------------------|---------------------------------------------|
|  | Manufacturer                                |
|  | Date of manufacture                         |
|  | Lot number                                  |
|  | Item number                                 |
|  | Non-sterile                                 |
|  | Warning                                     |
|  | Consult instructions for use                |
|  | MR unsafe                                   |
|  | CE marking according to Directive 93/42/EEC |
|  | Medical device                              |
|  | Prescription only                           |

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