DE

CE 0123

Rx only

Rev.: 15-12-2025

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Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

Table of content:

1. General information
2. Purpose and proper use/ Indication
3. Contraindications
4. Possible side effects and complications
5. MR Note
6. Single use product
7. Implant selection
8. Use/handling of fastening screws
9. Instruments
10. Combination products
11. Decontamination, cleaning and sterilization
12. Liability
13. Signs and symbols

1 GENERAL INFORMATION

The MEDICON condyles are supplied unsterile and must be disinfected, cleaned and sterilised prior to use. Please observe the following instructions. This will guarantee that the screws function faultlessly and reliably. MEDICON condyles may be used only by surgeons with adequate training and sufficient experience in the treatment of mandibular joint defects.

MEDICON condyles are available in two lengths, each for the left and the right side of the mandible.

Length	Left side	Right side
40 mm	REF. 68.95.46	REF. 68.95.47
48 mm	REF. 68.95.48	REF. 68.95.49

2 PURPOSE AND PROPER USE/INDICATION

MEDICON condyles serve as a temporary joint body replacement in reconstructions of the mandible using the MEDICON titanium reconstruction plate system 2.4 (REF.: 68.95.10 - 68.95.19). The implants are used to bridge the time until final reconstruction.

MEDICON condyles are designed for temporary use and are not intended as a permanent replacement.

3 KONTRAINDIKATIONEN

- Patients that are unable to follow the instructions for postoperative care. This may be caused e.g. by psychological/mental or neurological problems of the patient
- Extreme material sensitivity, i. e. foreign-body reaction of the patient. For this case, appropriate tests must be carried out under all circumstances prior to implantation (even when there is only a suspicion of such condition)
- Acute infections
- Patients with unstable physical and/or mental health conditions.
- Patients with insufficient or low-quality bone tissue of the socket.

Direct contact between the condyle and the bone of the mandibular fossa is to be avoided.

4 POSSIBLE SIDE EFFECTS AND COMPLICATIONS

In many cases, possible complications are attributable not to the implants as such, but to incorrect patient selection, lack of practice or inadequate repositioning of the bone fragments.

- Complaints, pains, abnormal sensations or palpability of implants.
- Material/foreign-body sensitivity of the patient with allergic reactions.
- Stronger connective-tissue reaction in the fracture region or in the vicinity of the implant.
- Breakage, bending, migration or loosening of the implant.
- Loosening of the implant as a result of insufficient screw tightening.
- Discomfort in the area of the implant by prolonged usage of this temporary joint body replacement.
- Attaching the condyle at an incorrect angle can limit the function of the jaw on both sides.
- Massive bending and breakage of the implant.
- Osteonecrosis, osteoporosis, limited revascularization, bone resorption and poor bone formation can lead to early loss of fixation and thus to nonunion (pseudarthrosis).
- Inadequate healing.

- Reduced bone density due to "stress shielding".

Apart from the above side effects and complications, the surgical procedure may cause other problems such as injury to nerves, infections and/or pain, which are not necessarily due to the implant.

Complete restoration of the jaw function is not provided by the implant. Normal masticatory forces can weaken the implant and result in failure of the implant.

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Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

6 MR Note**Instruments**

The use of medical devices in an MR environment poses a hazard. Individual medical devices must not be placed in the immediate vicinity of the equipment during the performance of MR procedures. **Implants**

The implants are not MR-safe and have not been tested for MR safety. Use in an MR environment may therefore pose a risk. It is recommended that patients and medical personnel be informed that the implants are not approved for use with MR equipment.

7 IMPLANT SELECTION

The following aspects are critical for the correct selection of implants:

- bone defect to be treated / intended bone repositioning,
- body weight of the patient,
- general health condition, degree of activity of the patient.

Medicon implants and instruments have been designed for being used together. Using implants or instruments from other manufacturers in conjunction with Medicon products entails unpredictable risks since the products will not be matched to each other. To avoid such risks, only Medicon products that are intended for combination with each other may be used together. This may also cause danger to the patient, the surgeon or third parties. The condyles may be used only with the MEDICON titanium reconstruction plates (REF.: 68.95.10 - 68.95.19). Use with other plates is not permissible.

Prior to applying the products, the surgeon must thoroughly discuss with the patient the intended outcome of the operation. In such discussions, special attention must be paid to postoperative care measures that may become necessary. The patient must be instructed to inform the surgeon immediately if they notice unusual changes in the situs.

The patient must be monitored carefully if there is a manifest change in the fixation region. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing. It cannot be excluded that implants will break, loosen, corrode, migrate through tissue and cause pain.

The packaging label carries a "Lot No.". We recommend copying this "Lot No." to the patient records as this number allows tracing the production history of the implant back to the raw material.

For the fixation of the condyle on the reconstruction plate, only fastening screws REF.: 68.95.35 are to be used. The plate must be secured on the condyle with at least three screws (see Figure 1).

**8 USE/HANDLING OF FASTENING SCREWS**

The fastening screws must always be inserted at a right angle to the condyle and to the plate to ensure tight connection between condyle, plate and screws and to avoid damage to the thread.

The screw driver must be inserted into the screw head with axial pressure to ensure that the blade is seated tightly in the screw head. This ensures correct longitudinal alignment of screw and screw driver and avoids slipping of or damage to the screw head or the blade.

If excessive force is used when fixing the screws, the screws may break during surgery, or the screw heads may tear off. Wear and tear effects to the screw driver blade impair the solid connection between the blade and screw. In this case, replace the blade with a new one. The screws are self-retaining on the screw driver blade.

Excessive tightening of the fastening screws may cause the threads to pull out. In such a case, the defective fastening screw and / or the condyle is to be replaced with a new one. The damaged fastening screw and / or condyle may not be used further.

Stability of the connection between all screws and implants must be checked. If necessary, retighten the screws.

9 INSTRUMENTS

Only the instruments designed for the MEDICON condyles and MEDICON titanium reconstruction plates may be used. The instruments intended for application of the system are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force.

To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional..

10 KOMBINATIONSPRODUKTE

Condyles for reconstruction plates:

REF.: 68.95.46 und 68.95.47

REF.: 68.95.48 und 68.95.49

Combination options:

Plates: REF.: 68.95.10 bis 68.95.19

Screw: REF.: 68.95.35

If items not listed are added through product line extensions, in case of doubt please consult your retailer or MEDICON eG.

11 DECONTAMINATION, CLEANING AND STERILIZATION

The following vCJD prion-specific protective measure is indicated when processing instruments.

In the event of diagnosis of a definite or probable case of vCJD.

- If it is not possible to use disposable products, the instrument used, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

If prion contamination is suspected:

- In the event of suspected prion contamination, incineration of the instrument is recommended according to the vCJD task force final report.

If vCJD is excluded:

- Continue to use after instrument processing is completed. Otherwise, the instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

In the case of a non-identifiable vCJD illness.

- Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions - e.g. mechanical alkaline cleaning combined with steam sterilization.

If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (CNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134° C.

Cleaning solutions to which hydrogen peroxide is added and/or highly-alkaline washing solutions can cause colour changes. This can cancel out colour-coding functions. Only cleaned and disinfected implants may be sterilized.

Reprocessing restrictions

Frequent reprocessing has little effect on these Implants.

Location of use

Implants that have been contaminated by blood and/or secretions must not be placed back in the implant storage tray. They must be disposed of.

Preparation for decontamination

The implant storage tray must be placed in instrument holders suitable for machine cleaning. The instrument holders (e.g. wire trays) must allow subsequent cleaning in an ultrasonic bath or in the cleaning and disinfection device (CDD) where they are not impeded by acoustic or rinsing shadows.

Pre-cleaning

- Place implant storage tray for 15 minutes in ultrasonic bath at 40° C with 0.5% alkaline cleaner (Neodisher Mediclean forte) and treat with ultrasonic waves.
- Remove implant storage tray and rinse thoroughly with cold water.

A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

Cleaning: machine

The cleaning and disinfection device (CDD) must meet DIN EN ISO 15883-1 requirements.

- Pre-rinse 1: 3 minutes with cold demineralised water, without additive;
- Emptying;
- Cleaning: with demineralised water, heat to 55° C and wash/clean for 5 minutes, add cleansing agent at 45° C, alkaline cleansing agent, strength 0.5%;
- Emptying.
- Neutralization: 3 minutes with warm water (>40° C) with addition of neutralizer, strength 1ml / l;
- Emptying;
- Final rinse: 2 minutes with warm deionised water (>40° C) (without any additive);
- Emptying

Disinfection

Thermal Disinfection A0 value 3000: Completely desalinated water, the thermal disinfection is carried out at temperatures > 80° C and corresponding application time according to the A0 concept, DIN EN ISO 15883 and guideline DGKH, DGVS and AKI (e.g. A0 3000 = 90° C and 5 minutes application time). The operator is responsible for the implemented A0 value.

Drying

Complete drying must be ensured by the CDD (e.g. 60° C 30 min.). The instruments must be removed from the CDD promptly once the cleaning and disinfection program has finished. If necessary, compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Maintenance, inspection, and testing

After cleaning/disinfection the implants must be macroscopically clean, i.e., free from visible residues and soiling. Inspection is performed visually. Insufficiently cleaned implants must be cleaned again and then adequately rinsed and dried. Deformed or damaged implants must be removed and disposed of, as they can no longer be guaranteed to be safe to use.

Packing

The implant storage tray must be packed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868
- DIN EN ISO 11607
- suitable for steam sterilization (steam permeable)
- adequate temperature resistance to 138° C

Sterilization accessories and sterilization packaging must be appropriate for the content of the package and for the sterilization process.

Sterilization

Use the following sterilization procedures subject to the applicable national requirements:

- Fractionated vacuum method with triple fractionation and sufficient product drying.
- Steam sterilization in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes holding time at 134° C.

It is essential that a SAL (Sterility Assurance Level) of 10⁻⁶ is reached

Storage

Processed sterile implants must be stored in appropriate reusable sterilisation containers which are dry, protected from dust and low-germ, in a dark and cool place, free from vermin.

To prevent condensation, extreme temperature fluctuations should be avoided during storage. No chemicals may be stored together with implants.

Walls, floors and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must have a ground clearance of at least 30 cm.

The permissible on-site storage period depends on the type of sterile barrier system used and on the storage conditions. The permissible storage period is to be specified by the operator.

ADDITIONAL INFORMATION ON REPROCESSING

A validated machine cleaning and disinfection procedure is always to be preferred to manual cleaning because of the greater safety of the process. Thorough cleaning also helps to maintain value and is a prerequisite for successful sterilisation.

Please take note of the following points during automated processing:

- For effective automated processing, load the wire mesh baskets in such a way that they can easily be mechanically rinsed. Mesh baskets must not be overloaded
- Avoid "rinsing shadows" caused by large instruments

The times and temperatures stated in these reprocessing instructions are minimum requirements which must be adhered to. Should any downward deviation be needed for procedural reasons, this must be validated by the operator. Exceeding the specified times and temperatures is generally possible. However, this will lead to increased stress on the material, which can cause premature wear. We accept no responsibility for the use of other sterilisation procedures.

INFORMATION ON VALIDATING THE PROCESSING PROCEDURE

Validation is carried out with the following devices, materials and chemicals:
Cleaning and disinfection equipment: Type Miele Desinfektor G 7735 CD
Trolley for surgical instruments

Alkaline detergent: neodisher® MediClean forte Dr. Weigert GmbH & Co. KG
Neutraliser: neodisher® Z Dr. Weigert GmbH & Co. KG
Ultrasonic bath: Sonorex DIGITEC DT 156 BH
Sterilisers: MMM Vaculab 969 S 3000 MMM Selectomat S 3000 Stiefenhofer KS 666-2ED H+P Varioclav 400E

NOTE

The person processing the instrument bears the responsibility that the processing actually performed with the equipment, materials and staff at the reprocessing facility achieves the desired results. Validation and routine monitoring of the process are normally required.
If the devices, materials and chemicals described above are not available, it is up to the persons processing the instrument to have their procedures appropriately validated. Observe the notes and regulations of the relevant national guidelines and standards. MEDICON eG reserves the right to make changes to these instructions due to new findings.

12 LIABILITY

MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

13 SIGNS AND SYMBOLS

	Manufacturer
	Date of manufacture
	Lot number
	Item number
	Non-sterile
	Warning
	Consult instructions for use
	MR unsafe
	CE marking according to Directive 93/42/EEC
	Medical device
	Prescription only

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Detergente alcalino: neodisher® MediClean forte Ditta Dr. Weigert GmbH & Co. KG
Neutralizzatore: neodisher® Z Ditta Dr. Weigert GmbH & Co. KG