C€0123

Rx only

Deutsch

Rev.: 08-12-2025

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Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

MEDICON implants are manufactured from titanium alloy. The material is bio-compatible, corrosion-proof and nontoxic in the biological environment. The surface is chemically passive and the material is antimagnetic. The implants are supplied in unsterile condition.

Contents:

1. General Information
2. Purpose and proper use/indications
3. Contraindications
4. Possible Side Effects and Complications
5. Single use product
6. CT and MR instruction
7. Selection of screws
8. Handling
9. Application Instruments
10. Postoperative care
11. Decontamination, cleaning and sterilization
12. Liability
13. Signs and symbols

1. General Information

The MEDICON Aarhus Mini-Implant system is supplied unsterile and must be cleaned, disinfected and sterilised prior to use. Please observe the following instructions. This will guarantee that the implant functions faultlessly and reliably.

2. PURPOSE AND PROPER USE/INDICATIONS

The MEDICON Aarhus Mini-Implant system may be used only by orthodontists, dentists, oral surgeons and surgeons specializing in maxillofacial surgery of proper training and experience. All standard hygienic precautions required for invasive procedures for example a sterile working environment, sterile gloves, surgical masks must be taken before and during insertion of MEDICON Aarhus Mini-Implant screws.

The MEDICON Aarhus Mini-Implant System is designed for orthodontic treatment (anchorage of orthodontic appliances with wires, springs, rubber rings) for correction of incorrectly positioned teeth and jaws. It is used for treatment of patients where anchorage is impossible in the posterior field and as required with:

- Retraction and/or intrusion of the front teeth
- Intrusion of molars and/or premolars without occluding teeth
- Mesial movements of molars where no front retraction is acceptable
- Proclination of lower incisors where no posterior anchoring units are available or where no reactive forces can be tolerated
- Space closure against absolute anchorage

3. CONTRAINDICATIONS

The MEDICON Aarhus Mini-Implant System must not be used if the following contraindications are existing:

- Patients that are unable to follow the instructions for postoperative care. This may be caused e.g. by psychological/mental or neurological problems of the patient.
 - Patients with inadequate bone tissue, or bone tissue of insufficient quality, and patients with circulatory disorders of acute infections.
 - Patients with unstable physical and/or mental health conditions.
 - Extreme material sensitivity, i.e. foreign-body reaction of the patient. For this case appropriate tests must be carried out under all circumstances prior to implantation (even when there is only a suspicion of such condition).
 - Acute Infections
 - Radiotherapy in the head region and patients with tissue damaged by radiation
- Recurring diseases of the oral mucosa and poor oral hygiene

4. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

- Complaints, pain, abnormal sensations or palpability of the implant.
- Material hypersensitivity of the patient due to the foreign bodies in the form of allergic reactions.
- The use of different materials may cause corrosion.
- Increased reaction of the connective tissue in the area of the implant.
- Inadequate bone formation, osteolysis, osteomyelitis, osteoporosis, inhibited revascularization or infection that may cause loosening, bending, cracking or breakage of the implants.
- Breakage, bending, migration or loosening of the implant.
- Inadequate osseointegration that could cause loosening or breakage of the implant and failure of the orthodontic treatment.
- Decrease of bone density as a result of stress shielding.

This surgical procedure may cause not only the above-mentioned side-effects and complications but also problems such as injuries to nerves, infections, pain etc., which are not necessarily caused by the implant.

5. SINGLE USE PRODUCT

Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

6. MR Notes**Instruments**

The use of medical devices in an MR environment poses a hazard. Individual medical devices must not be located in the immediate vicinity of the equipment during the performance of these procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MRI environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MRI devices.

7. SELECTION OF SCREWS

Selection of the MEDICON Aarhus Mini-Implant screws is the responsibility of the maxillofacial or oral surgeon, dentist or orthodontist and must match the patient's anatomy.

Before using the Aarhus Mini-Implant screws, the maxillofacial or oral surgeon, dentist or orthodontist must explain the desired course of treatment and the expected result in detail to the patient. In such discussions, special attention must be paid to postoperative care measures that may become necessary. Selecting the incorrect implants may result in premature implant loss along with loosening, bending or breakage of the implant due to which the implant may lose its function. Damage and scratches will reduce the strength of the product and cause premature fatigue of the implant. This will adversely affect the success of the treatment. MEDICON implants and instruments have been designed for being used together. Using implants or instruments from other manufacturers in conjunction with MEDICON products entails unpredictable risks since the products will not be matched to each other. To avoid such risks, only MEDICON products that are intended for combination with each other may be used together.

8. HANDLING

MEDICON Aarhus Mini-Implant screws can be inserted in the maxilla in a vestibular of palatal position or in the alveolar ridge. Insertion in the mandible should be vestibular only. Based on current knowledge, we advise against lingual insertion and insertion in the vicinity of extraction wounds, dental follicles and first teeth. The following positions offer the best conditions for an insertion of the screw:

Maxilla:	Mandibula:
– The infrazygomatic crest	– The retromolar area
– The anterior nasal spine	– The symphysis
– The palate	– The alveolar process

The alveolar process

MEDICON Aarhus Mini-Implant screws have a self-tapping thread, which makes pilot drilling unnecessary. However, it is up to the surgeon to decide whether it is better to drill a pilot hole into the cortical bone of the lower jaw.

Only use drill bits with sharp cutting edges! In this case select the drill and length to match the screw diameter. Drill a cavity of at least 4.0 mm with adequate cooling and at a speed of 800 – 1500 rpm. Excessive application of force during the drilling procedure can lead to drill breakage, endangering the patient, the surgeon and third parties.

Application of excessive force when tightening the screws can lead to intraoperative screw breakage. The screw must be inserted in bone of appropriate quality only. If immediate therapeutic loading is not possible, an unloaded healing phase of at least three to four weeks is required. To ensure a safe function of the MEDICON Aarhus Mini-Implant screw, it must be securely anchored in the bone (primary stability) and the screw head must be placed in the region of the alveolar gingiva. When using the Aarhus Mini-Implant screw as an anchoring element, make sure that the head and the surrounding soft tissue are not subject to any unfavourable mechanical stresses (such as movement of the mucosa, influence by bands and/or tongue, manipulations).

Use only 50cN springs for anchorage, connected directly to the screw or to the teeth that are anchored by the screws.

The line of force must always pass through the screw if the screw is used as direct anchorage. Forces that induce torque around the longitudinal axis of the screw (identical to the force that is required to insert or remove the screw) must be avoided. A thorough force system must be developed before insertion of the screw. Other important information is covered in the various manuals on surgical technique.

9. APPLICATION INSTRUMENTS

When inserting implants, use only those instruments specifically designed and sold

for the MEDICON Aarhus Mini-Implant System (see the product brochure). The instruments intended for use of the system are subject to wear and mechanical stress even during normal use, and particularly if excessive force is used. In order to prevent failure or mechanical damage of the application instruments during surgery, they must always be checked before each use for mechanical integrity, full functional range, and possible deformation. Instruments with any damage must be discarded.

10. POSTOPERATIVE CARE

Periodic postoperative care and observation and also reduced physical activity after the operation and during the healing phase is very important for the success of the treatment. Implants may come loose, be displaced, bend or break as a result of overload. The supervising physician is responsible for deciding the type, duration and intensity of physical activity after the operation. The patient must be advised that ignoring medical advice regarding the above and other areas may result in unforeseeable complications. This is particularly applicable if the patient manipulates the implant or appliance himself. The patient must be instructed to inform, without delay, the operating surgeon of unusual changes at the operating site. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing. It cannot be excluded that implants will break, loosen, migrate through tissue and cause pain. The implants must be removed after the treatment has been completed. We can therefore only accept liability for implants up to the completion of treatment. We do not accept any responsibility for any damage arising from implants which remain in the body after conclusion of treatment.

The packaging label carries a Lot no. We recommend transferring this Lot no. to the patient records as this number allows tracing the production history of the implant back to the raw material.

11. DECONTAMINATION, CLEANING AND STERILIZATION

The following prion-specific protective measure (vCJD) is indicated when preparing the instrument.

In the event of a diagnosis of definite or probable vCJD:

If the use of disposable products is not possible, the instruments used (which had become contaminated or whose contamination cannot be ruled out) must be disposed of as waste for incineration.

In suspected cases:

If prion contamination is suspected, it is recommended that any medical instruments be incinerated in accordance with the final report by the nvCJD task force.

If vCJD infection is ruled out:

Re-use after final treatment. Otherwise the instruments used (that had become contaminated or where contamination cannot be ruled out) must be disposed of as waste for incineration.

In the event that vCJD infection is not apparent.

Even if there is no known existence of a prion disease, we should, in preparation, set up two procedures with at least partial effectiveness against prion contamination, for example, an alkaline mechanical cleaning process combined with steam sterilisation.

If an alkaline mechanical cleaning process is not used, or another cleaning process with proven effectiveness against prion contamination is deployed, and medical devices that come into contact with risk tissue (CNS, eyes, lymphatic tissue) are used, the Robert Koch Institute (IRK) recommends an extended sterilisation time of 18 minutes at 134°C.

Cleaning solutions to which hydrogen peroxide has been added and/or washing solutions with high alkalinity can cause colour changes. This may cause the coding function to be lost. Only cleaned and disinfected instruments may be sterilised.

Limitations on reprocessing

Frequent reprocessing has little effect on the implants.

Place of use

Implants with contaminations caused by blood and/or secretions must not be placed back into the implant storage tray. They must be disposed of.

Preparation for decontamination

The instruments must be placed on appropriate instrument carriers so that they can be rinsed in the washing and disinfection device. The instrument carriers (e.g. sieve trays) must be such that subsequent cleaning in an ultrasonic device or in the washing and disinfection device is not impaired by any areas that are inaccessible to ultrasonics or rinsing.

Precleaning

- Place the instruments in an ultrasonic bath for 15 min at 40°C with a 0.5% alkaline cleaner (pH 11.4 -11.9) and subject to ultrasound.
- Remove the instruments and rinse with cold, softened water.

A high level of contamination in the ultrasonic bath will impair the cleaning effect and promote corrosion. The cleaning solution must be regularly replaced, depending on the conditions of use. The decisive criterion is optically recognisable contamination. In any case, the bath must be changed frequently, at least once a day. Comply with the applicable regulations for processing medical products in the country concerned.

Cleaning

The cleaning and disinfection device must conform to the requirements of DIN EN ISO 15883-1.

- Pre-rinse 1: 1 minute with cold demineralised water with no additives;
- Drain;
- Pre-rinse 2: 3 minutes with cold demineralised water with no additives;
- Drain;
- Cleaning: use demineralised water, heat to >55°C, wash and clean for 10 minutes, add cleaning agent at 45°C, alkaline cleaning agent, dosage 0.5% (pH 11.4 – 11.9);
- Drain;
- Neutralization: 3 minutes with warm deionised water (>40°C) with neutralizer, dosage 1 ml/l;
- Drain;
- Final-rinse: 2 minutes with warm deionised water (>40°C) (no additives); Drain

Carry out thermal disinfection, using fully desalinated water, at a temperature of >80°C for an exposure time acc. to A0 draft, EN ISO 15883-1 and DGHK, DGSV and AKI (e.g. A0 3000 = 90°C 5 Min. exposure time). The operator is responsible for realizing the A0 value.

Drying

Sufficient drying time must be provided by the cleaning and disinfection device. The instruments must be removed from the cleaning-disinfection device immediately after the cleaning and disinfection program has finished. If necessary, the use of compressed air for drying is recommended because of its effective and quick effect (RKI recommendation).

Maintenance, control and inspection

The instruments must be macroscopically clean, i.e. free of visible residues, after cleaning and disinfection. The instruments must be visually inspected in this respect. Inadequately cleaned parts must be cleaned, thoroughly rinsed and dried again. Defective instruments must be replaced, as they cannot fulfil their function, or cannot fulfil it safely enough.

Packing

The instruments must be packed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868
- DIN EN ISO 11607
- suitable for steam sterilization (steam permeable)
- adequate temperature resistance to 138°C

Sterilization accessories and sterilization packaging must be appropriate for the content of the package and for the sterilization process.

Sterilization

Use the following sterilization procedures subject to the applicable national requirements:

- Fractionated vacuum method with triple fractionation and sufficient pro-duct drying.
- Steam sterilization in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes holding time at 134°C.

It is essential that a SAL (Sterility Assurance Level) of 10⁻⁶ is reached.

Storage

Reprocessed sterile instruments must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments. The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed

EN

CE₀₁₂₃

Rx only

Rev.: 08-12-2025

Instruction for use

and the storage conditions. This storage period must be established by the operating authority.

Further informations on reprocessing

A validated machine cleaning and disinfection procedure is always to be preferred to manual cleaning because of the greater safety of the process. Thorough cleaning also helps to maintain value and is a prerequisite for successful sterilisation.

Please take notice of the following points regarding mechanical processing:

- Load the wire mesh baskets in such a way that they can easily be mechanically rinsed. Mesh baskets must not be overloaded.
- Avoid „rinsing shadows“ caused by large instruments.

The times and temperatures stated in these reprocessing instructions are minimum requirements, they must be reached. Should any downward deviation be needed for procedural reasons, this must be validated by the operator. Exceeding the specified times and temperatures is basically possible. However, this will lead to increased stress on the material, which can cause the instruments to wear out prematurely. We accept no responsibility for the use of other sterilisation procedures.

Information on validating the reprocessing procedure

Validation is carried out with the following apparatus, materials and chemicals:

Cleaning and disinfection equipment: Type Miele Desinfektor G 7735 CD
Trolley for surgical instruments
Alkaline detergent: neodisher® MediClean forte, Fa. Dr. Weigert GmbH & Co. KG
Neutralizing detergent: neodisher® Z, Dr. Weigert GmbH & Co. KG
Ultrasonic bath: Sonorex DIGITEC DT 156 BH
Sterilisers: MMM Vaculab 969 S 3000, MMM Selectomat S 3000, Stiefenhofer KS 666-2ED, H+P Varioclav 400E

Note

The person preparing the instrument bears responsibility that the preparation actually performed with the equipment, materials and staff at the reprocessing facility achieves the desired results. Validation and routine monitoring of the process are normally required. If the above-described apparatus, materials and chemicals are not available, it is up to the persons preparing the instrument to have their own procedures validated. Please refer to the instructions and requirements of the relevant national regulations and standards as well as the instructions for use included with the medical device, if applicable. Medicon reserves the right to make alterations due to recent findings.

12. LIABILITY

In the event of discrepancies between the english and the german version of these instructions for use only the german version shall be applicable. The information in these instructions for use applies only to the product with which the instructions for use are supplied. Only the latest revision of the instructions for use applies. Due to constant technical development the contents of these MEDI-CON instructions for use are updated regularly.

MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

 **Federal law and regulations in the United States restrict the sale of this device to a physician or hospital or on prescription.**

13. SIGNS AND SYMBOLS

	Observe instructions for use
	Caution
CE ₀₁₂₃	CE marking according to Directive 93/42/EEC
	Manufacturer
	Single use
	MR unsafe
Rx only	Subject to medical prescription (US federal law)
	Batch code
	Article number

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