Mainzer Distraktor

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

Mainz Distractor

The Mainz Distractor was developed for augmentation of the alveolar process of the lower jaw where it is atrophied or lost due to trauma or tumour resection.

Indications

Reconstruction of the bony alveolar process in:

- Atrophy of the toothless or partially toothed lower jaw
- Alveolar process defects, e.g. due to tumour resection or trauma
- Primary alveolar process hypoplasia with congenital absence

Please note that the "Mainz Distractor" is a titanium implant for single use

The distractors are implants which are intended for single use. Re-use of explanted distractors is not permissible. Even a distractor which looks

undamaged after explantation may, as a result of preceding stresses, have signs of wear in the material and the drive which may lead to failure and endanger the success of an operation.
The distractors are non-sterile at delivery.

The label includes a "Lot-No." (batch number). We recommend inclusion of this Lot-No. in the surgery protocol, since this number can be used to trace the production of the distractor.

Medicon eG, as the distributor of the distractors, assumes no liability for indirect damage or consequential damage caused by improper use or handling, particularly including non-compliance with the operating instructions or improper care or maintenance.

The Medicon Distractors may only be used by surgeons with sufficient training and experience.

Carefully read the operating manual and user information in the brochure prior to use.

It is absolutely necessary to observe the information for use in our brochure "Mainz Distractor".

If you do not have this brochure, you will receive it upon request.

Treatment process

Distractor placement,
Commencement of distraction after 4-5 days,
Distraction end after another 5-18 days,
Implantation is possible immediately after completion of distraction,
Healing phase,
Clearing and placing load on implants after 4-5 months.

Possible complications and side effects

Complaints, pain, abnormal sensations or palpability due to the implants. Material hypersensitivity of the patient due to the foreign bodies in the form of allergic reactions.

Increased connective tissue reaction in the regions of the osteotomy or the implants.

Insufficient bone formation, osteolysis, osteoporosis, reduced revascularisation or infection which may lead to loosening, bending, breaking or failure of the implant.

Delayed or insufficient healing, which may cause implant breakage. Bone necrosis.

Insufficient healing.

Fracture of the distraction segment

The distraction cannot be continued. However, alternative bone augmentation, e.g. by means of interposition of a bone transplant, can be done. (Sandwich osteoplasty)

Fracture of remaining tissue

In such a situation, bridging osteosynthesis of the remaining tissue must first be done. If the osteosynthesis provides sufficient stability, it is biologically possible to continue distraction, however there is no experience in this regard. Experiences with secondary distraction of healed bone transplants show that after the fracture has healed, renewed distraction can be done.

Loosened distractor

If the distractor loosens in the late phase of distraction, it is not possible to further raise the distractor. If sufficient height has already been obtained, the distractor should be replaced immediately by an implant.

If the distractor loosens in the early phase, the distraction must be discontinued.

Wound dehiscence

If wound dehiscence occurs after segmental osteotomy, therapy does not have to be discontinued. Even though it has been shown that wound dehiscence can heal spontaneously even under the distraction, it is recommended to pause activation for 2 days.

Warnings

The intervention must be done solely with the instruments intended for this purpose.

See accessories.

The minimum height of the remaining tissue (basal segment) must not be less than 6 mm as a measure to prevent fractures.

The distraction segment should be 5 mm high if possible. Minimum height of 3 mm and minimum width of 5 mm are required.

Carrying out the distraction

The distraction is carried out by turning left. When activating the distraction screws, use the provided screwdriver, Ref. 69.00.10.

The arrow at the end of the screwdriver indicates the turning direction for activating the distraction screws.

One complete screw rotation corresponds to a distraction distance of 0.3 mm. Daily distraction is between 2 and 3 rotations, that is = 0.6 to 0.9 mm.

Preparation for decontamination

For cleaning and sterilization purposes the implants are supplied in a non closed position. In any case the implants must be activated up to the elevation top, using the screwdriver (69.00.10) prior to application. The implants must be sorted into the implant tray which is intended for the system.

Possible surface contamination from use on the surface of the implant tray must be removed with a single use cloth or paper towel.

Cleaning

Suitable pH-neutral or alkaline cleaning agents are used as cleaners. The selection of the cleaner depends on the material characteristics of the implant tray, as well as national guidelines and recommendations.

Thermal disinfection

Completely desalinated water; thermal disinfections is done at temperatures of 80-95°C with a corresponding acting time as per the A0 concept, EN ISO 15883-1.

Drying

Sufficient drying must be ensured by means of the cleaning and disinfection unit.

Cleaning agent

The cleaning agent manufacturer's stated parameters for concentration, temperature and acting time must be observed, and automatic dosing devices must be controllable.

Maintenance, control and verification

After cleaning / disinfection, the implants and trays must be macroscopically clean, that is, free of visible residues. The verification is visual. Critical areas require particularly careful control. Insufficiently cleaned parts must be re-cleaned and subsequently sufficiently rinsed.

Sterilization

Only cleaned and disinfected trays and implants are allowed to be sterilized.

Only the listed sterilization methods may be used for sterilization.

- Fractioned vacuum process with sufficient product drying
- Steam sterilizer validated as per DIN EN 13060 and/or DIN EN 13485 and as per DIN EN 554.
- Sterilization time and temperature:
 Minimum 5 minutes hold time at 134°C

Storage

In order to avoid the formation of condensation, larger temperature variations must be avoided. Chemicals may not be stored with instruments.

The duration of storage depends on the type of packaging.

The user is responsible for validating the recommended parameters in order to take possible differences between sterilisation chambers, packaging methods and load configurations into account. The same applies if procedures other than the recommendations are used. Reaching a SAL (Sterility Assurance Level) of 10⁻⁶ is mandatory.


Titan "Mainz Distraction Screws" and accessories

*Single use product



Art. No.	Description
69.00.05	Distractor "Model Mainz" lift height 5.0 mm, Ø 3.6mm*
69.00.07	Distractor "Model Mainz" lift height 7.0 mm, Ø 3.6mm*
69.00.09	Distractor "Model Mainz" lift height 9.0 mm, Ø 3.6mm*
69.00.11	Distractor "Model Mainz" lift height 11.0 mm, Ø 3.6mm*
69.00.03	Counter-nut for distractor "Model Mainz"
69.00.04	Flat wrench for counter-nut
69.00.10	Activation screwdriver
63.68.16	Drill with dental shaft 22 x 1.6 mm / Package of 10 items
69.00.20	Spiral drill with dental shaft 32 x 2.0 mm
69.00.22	Stage drill with dental shaft 26 x 2-3.0 mm
69.00.24	Spiral drill with dental shaft 32 x 3.0 mm
69.00.30	Directional indicators / pair
69.00.27	Thread cutter for ratchet
69.00.28	Thread cutter, long, 34.5 mm, for dental angle tool
69.00.29	Thread cutter, short, 26.0 mm, for dental angle tool
69.00.12	Placement instrument for ratchet
69.00.14	Placement instrument for dental angle tool
69.00.19	Placement instrument, manual
69.00.16	Ratchet
69.00.18	Ratchet handle, single
69.00.17	Ratchet insert, single
69.00.33	Gauge
69.00.35	Parallel drilling aid
69.00.40	Teflon sterilisation and storage cassette
49.40.50	C-Container 300 x 135 x 62 mm, cover with holes, for storage and sterilization. For carrying out the osteotomy, we recommend the micro-saw osseoscalpel "Model Sachse"
68.24.05	(special brochure upon request)

Surgical Procedure:


1. Mark position of distraction devices with a "Rosen" burr cat.-no. 63.68.16.


2. + 3. First drilling using a 2.0 mm twist drill cat.-no. 69.00.20 down to the basal segment. Assure for parallel alignment of the second bore hole by using the direction indicator cat.-no. 69.00.30.


4. Parallel enlargement of the drill holes by using graduated drill diameters 2/3 mm, cat.-no. 69.00.22.


5. Second drilling to the basal segment with twist drill diameters 3,0 mm, cat.-no. 69.00.24.


6. Cutting of thread with screw tap cat.-no. 69.00.27 alt. 69.00.28 or 69.00.29.

7. Realisation of a box shaped osteotomy with an appropriate distance to neighbouring structures, preserving a lingual soft tissues pedicle. The Osseoscalpel cat.-no. 68.24.05 or a small oscillating saw is recommended for this procedure.



8. + 9. Insertion of the distraction devices using the positioning instrument cat.-no. 69.00.12 alt. 69.00.14 or 69.00.19 and -if necessary- tighten the lock nut of the distraction devices to compensate for a possible distance between distraction device and bone. Test run of the distractor.

10. + 11. Distraction: Activation of the threaded spindle by clockwise turning the central screw clockwise. One turn leads to a distraction of 0.33 mm.

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