Deutsch

Rx only

Rev.: 12-12-2025

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Rev.: 12-12-2025

Instruction for use

Carefully read the operating instructions and the “Surgical technique and instruments” brochure prior to clinical application and keep them safe and at hand. The “Surgical technique and instruments” brochure can be requested from MEDICON. The instructions and notes contained therein must be followed. MEDICON eG will assume no liability for the application, unless appropriate training has been provided.

Table of Contents

- General notes
- Intended use and indication
- Contraindications
- Possible adverse effects and complications
- Single-use product
- MR Notes
- Application and handling
- Pre-operative and post-operative behaviour
- Sterility
- Liability
- Description of symbols and icons

1. General Notes

The MEDICON “MediRod® pedicle screw-rod system sterile” is made of titanium alloy. The material is biocompatible, corrosion-resistant and non-toxic in a biological environment. The surface is chemically passive; the material is antimagnetic. The implants are not reusable and are supplied sterile. Please observe the following notes. These will ensure flawless and reliable functioning of the product. The reusable surgical instruments are made of stainless materials. The material is corrosion-resistant and has excellent properties in a biological environment.

2. Intended Use and Indication

The MEDICON “MediRod® pedicle screw-rod system sterile” may be used only by surgeons with sufficient experience in spinal surgery (neurosurgeon, orthopaedist, trauma surgeon, etc.) who have been appropriately introduced to the system.

The MEDICON “MediRod® pedicle screw-rod system sterile” is used in hospitals and clinical practices to stiffen segments of the spine. The MEDICON “MediRod® pedicle screw-rod system sterile” is used to stabilize two or more vertebrae in conditions of:

- Degenerative diseases
- Fractures
- Correction of post-traumatic malpositions
- Correction of deformities, such as scolioses or kyphoses
- Spondylitis
- Removal of spinal tumours or metastases

The MEDICON implants in the “MediRod® pedicle screw-rod system sterile” are used for internal fixation of the spine until osseous vertebral fusion occurs. The MEDICON implants are single-use products and may be used only in combination with the specified MEDICON application instruments.

3. Contraindications

The MEDICON “MediRod® pedicle screw-rod system sterile” must not be used with:

- Severe osteoporosis
- Highly tumorous, destroyed regions of vertebrae
- Anatomically unsuitable parts of the spine
- Other skeletal sections for which it is not intended
- Patients who are unable to comply with the instructions for postoperative care, e.g. due to psychological, mental or neurological problems
- Patients with inadequate or low-quality bone tissue, with perfusion problems or latent infections
- Proven allergy to titanium
- Material hypersensitivity, i.e. reaction of the patient to foreign bodies. Here appropriate tests are mandatory before implantation (even in case of mere suspicion!)
- Acute infections

4. Possible adverse effects and Complications

Possible undesired effects and complications associated with the MEDICON “MediRod® pedicle screw-rod system sterile” include:

- Absence of bone fusion
- Failure of the fixation due to loosening in the bone
- Failure of the implant due to implant breakage at the screws and/or rods
- Failure of the implant due to insufficient bone stiffening of the area operated on and thus absence of mechanically stable conditions
- Mobility in the operated area due to remaining discs or delayed/insufficient healing of bone grafts/cages persists for several months, leading to a risk of follow-up operations or fatigue breakage

- In the event of a lack of cooperation on the part of the patient, or where cooperation cannot possibly be obtained, pseudarthrosis formation and/or implant failure may occur
- Loosened implants can lead to secondary irritations or injury to adjacent anatomical structures
- Incorrect implantation of the screws due to imprecise placement in the pedicle may lead to nerve root damage or even paraplegia
- Incorrect screw size selection can lead to breaking of the pedicle
- Incorrect screw length selection can injure the internal organs in the thorax and abdomen or, in the medium-term, lead to erosion bleeding when close to blood vessels
- Material hypersensitivity of the patient due to foreign objects in the form of allergic reactions
- Failure of the rod due to repeated or extreme bending at the same point
- Discomfort, pain, abnormal sensations due to the implant
- Breaking, bending, migration, loosening of the implant and/or further surgical intervention to remove the system
- Early or late infections
- Nerve damage resulting from surgical trauma or the presence of the implant. Neurological symptoms, including bowel and/or bladder dysfunctions, impotence, retrograde ejaculation and paraesthesia
- Bursitis
- Fissures in the dura mater occurring during the operation may require further surgical intervention to restore the dura, and may lead to continued leakage of cerebrospinal fluid or a fistula, and under certain circumstances even to meningitis
- Vascular damage resulting from surgical trauma or internal fixator vessel damage may lead to life-threatening or fatal haemorrhages. Incorrectly positioned implants in the vicinity of large blood vessels may erode these vessels and lead to life-threatening haemorrhaging in the late postoperative period
- Damage to lymph vessels and/or lymph exudation
- Impairment or damage to the spinal cord
- Bone fracture
- Degenerative changes or instability in segments adjacent to a stiffened spine

The implants should be removed after complete healing. A “MediRod® pedicle screw-rod system sterile” that remains in the body after complete healing can lead to the following complications, individually or in combination:

- Corrosion with localised tissue reaction or pain
- Changes in the position of the implant with consequent injuries
- Risk of additional injuries due to postoperative trauma
- Bending, loosening and/or fracture, making removal more difficult or impossible
- Pain, discomfort or non-physiological sensations due to the presence of the product
- Potentially increased risk of infection
- Bone loss caused by stress shielding

In the selection of patients, the following factors should be observed that promote possible adverse events and complications:

Patient weight:

- Excess weight or obesity of a patient can increase the load on the implant in such a way that failure becomes likely.

Occupation or activity of the patient:

- Heavy lifting, muscle strain, body rotation, repeated bending, stooping, running or manual labour should be avoided in private or occupational activities until the bones have healed completely. Even after complete healing, the patient may not be able to successfully perform the activities listed above any longer.

Senility, mental illness, alcoholism or drug abuse:

- These circumstances may contribute to the patient ignoring certain restrictions or precautions that the implant calls for. This can lead to failure of the implant or to other complications.

Certain degenerative diseases:

- In some cases, a degenerative disease can be so advanced at the time of implantation that the expected service life of

the MEDICON “MediRod® pedicle screw-rod system sterile” system is clearly reduced. In such cases, orthopaedic tools can only delay or temporarily halt the degeneration.

Sensitivity to foreign bodies:

- Preoperative testing cannot entirely exclude hypersensitivity or allergy. Such complications may occur even if the implant has already been in the body for some time.

Smoking:

- In smokers, a higher rate of pseudarthroses was found after surgical procedures where bone implants were used. Furthermore, diffuse disc degeneration has been observed in smokers. Smoking-induced progressive degeneration of adjacent segments can lead to later clinical failure (intermittent pain), even if initially successful stiffening of the bone has occurred and clinical improvement has been demonstrated.

In addition to the undesirable effects and complications already mentioned, the surgical procedure can lead to other problems such as nerve damage, infections, pain etc. that are not necessarily attributable to the implant.

5. Single Use Product

Implants are intended and designed for single use on one patient only and may not be reused.

An explanted implant must never be reused. Even if the implant appears to be intact or fully functional, there may be wear, minor defects and attrition due to overuse that are not visible to the naked eye. As the impact of the forces and conditions inside the body on the stability, function and material quality of an explanted implant cannot be estimated, reimplantation would entail the risk of premature wear or failure and is therefore unacceptable. The user will be held liable if the operating instructions are not observed.

6. MR Notes

Instruments

The use of medical devices in the vicinity of an MR poses a hazard. Individual medical devices must not be located near the equipment during the application of these procedures.

Implants

The implants are not MR-safe and have not been tested for MR safety. Therefore, use in an MR environment may pose a risk. It is recommended to inform patients and medical staff that the implants are not approved for use in MR devices.

7. Application and handling

Implant selection: The surgeon is responsible for correct patient-specific selection and use of the implants.

The following aspects are critical for correct implant selection:

- Indication
- Patient's body weight
- Patient's physical condition, age and level of activity

Incorrect selection of the implants with regard to the aforementioned selection criteria can result in overloading and associated premature implant failure. Only use of the correct components will result in optimally stable fixation, whereas an incorrect decision may lead to loosening, bending or breakage of the implant and/or bone. The success of an operation depends, among other things, on how the implants are handled. The rods must be bent carefully, and after being bent once, the rod must not be bent back at the same point. Damage leads to reduction in product strength and premature fatigue of the implant.

Implants and instruments are developed and made to be compatible with each other. Use of implants and instruments by other manufacturers in conjunction with the MEDICON products is associated with unforeseeable risks, as such products are not designed to be used together. For this reason, only the specially designed MEDICON products should be used together, in order to prevent risks to patients, users and/or third parties. If complete healing fails to occur, or is delayed or inadequate, the possibility that the implant will bend, fail or break cannot be ruled out. Therefore, adequate immobilisation of the fracture site must be ensured until the bone has grown together again firmly. The constant load change the implants are subject to may lead to fatigue breakage. It cannot be ruled out that implants may break, loosen, bend, corrode, migrate in the tissue or cause pain. The implants temporarily stiffen segments of the spine; due to the anatomical limitations and despite modern surgical materials, they cannot last for an unlimited period.

We can therefore assume liability for the implants only until complete healing. We cannot be held liable for potential damage by implants that have not been removed after healing.

Using and handling the “MediRod® pedicle screw-rod system sterile”:

The rods must be adapted to the natural or planned bone structure as precisely as possible. Use only the instruments recommended by MEDICON for bending the rods. During the bending of the rod, cold straining occurs, under which the titanium increases in hardness but simultaneously loses deformability. It is therefore essential to use as few

bending steps as possible to achieve the desired shape of the implant. Excessive bending and incorrect setting of the bending instrument can lead to breakage of the rod. The rod may not be bent back at the same point after being bent once. Excessively forceful use of the bending instruments can lead to visible rod damage. For these reasons, the rods must be checked for any damage after bending, and damaged rods must be replaced. All implants must be checked for function, deformation and mechanical integrity prior to implantation. Damaged implants must be replaced, as they can lead to partial or even total loss of function. Always choose the rod and screw lengths that suit the specific situation and indication best. Should it be necessary to shorten a rod, use only the instruments provided for this purpose. Ensure that the cut-off part of the rod is not catapulted towards the patient, the surgeons or others, to prevent putting them at risk.

Pedicle screws do not require drilling a threaded hole. For drilling, only the cortex drills provided by MEDICON that are intended for the respective screws and match

them in colour may be used, so that the appropriate diameter for the screw and corresponding hole depth are achieved. Make sure that the screws are positioned correctly with the aid of imaging procedures. Use only drills with sharp cutting edges. If subjected to great force or used at an unfavourable angle in an application, screws, rods, drills and other parts of the system may become overloaded and break. The tight and correct connection between screws, heads, rods and the inner locking screw must be checked and ensured during surgery. Rods may loosen as a result of insufficient tightening of the inner locking screw.

Use and handling of instruments:

The instruments intended for Application of the “MediRod® pedicle screw-rod system sterile” are subject to wear and mechanical stresses even when used normally, but especially if used too forcefully.

In order to prevent failure or mechanical damage to the instruments during surgery, they must be checked before each use to ensure that they are mechanically intact, that there are no deformations, and that the parts are fully functional. Damaged instruments must not be used and are to be replaced. Use only accessories and instruments specifically approved by MEDICON, in order to avoid risks in connection with the compatibility of the products.

Detailed information about the handling of instruments and safety information can be found in the MEDICON brochure “Surgical technique and instruments”.

The correct application and handling of the instruments is described in the brochure “Surgical technique and instrument preparation”. The instructions and notes in this document must be followed.

8. Pre-Operative and Post-Operative conduct

Pre-operatively:

Before using the product, the surgeon should thoroughly discuss the desired operation result with the patient. Particular attention must be paid to post-operative behaviour and potentially required follow-up care.

Post-operatively:

Post-operative care and the patient's ability to follow instructions are the most important aspects of successful bone healing. The patient must be instructed to immediately inform the surgeon about any unusual changes at the operation site. The patient must be carefully monitored if there are actually changes in the fixation area. The surgeon must consider possible consequences such as implant failure and discuss required measures for further healing with the patient. The patient must be aware of the limitations enforced by the implant and be instructed to avoid or limit any physical activities, especially lifting and turning movements, as well as participation in sporting activities. The patient must be aware that a metal implant is not as strong as normal healthy bone and that excess loading – particularly in case of incomplete bone healing – can lead to loosening, bending, and/or breakage.

Displaced or damaged implants may migrate and damage nerves or blood vessels. An active, weakened, or demented patient who cannot properly use disburdening walking aids is at particular danger in the post-operative rehabilitation phase. The decision on whether the implants should remain in the body or be removed after complete healing is the responsibility of the surgeon and should be made after considering the risks and benefits of each approach. After explantation, risks should be minimised with adequate post-operative follow-up care. In active patients, the risk may be increased by implants remaining in the body. After explantation, the surgeon is responsible for disposal of the implants according to the standard. A “LOT” number (batch number) is indicated on the label of the packaging. To ensure 100% traceability, we recommend noting down the “LOT” number in the patient records.

Disposal:

Observe the applicable national laws when disposing of the MEDICON system and accessories.

9. Sterility

The implants are sterilised by irradiation and must be stored in their original packaging until immediately before use. The implants may then be removed from the original and protective packaging. The expiry date and intactness of the sterile packaging must be checked before use. Implants must not be used after the expiry date. The implant must be free from packaging residues prior to use. If products are reused, there is a risk of infection of patients and/or users, and impairment of proper functioning. Contamination and/or impaired function of the product can lead to injury, illness or death.

The “MediRod® pedicle screw-rod system sterile” is a sterile product for single use. Products from damaged or broken sterile packaging may not be used, as sterility cannot be guaranteed. The product may not be re-sterilised.

The “MediRod® pedicle screw-rod system sterile” may not be re-sterilised.

10. Liability

In the event of any discrepancies between the non-German and the German version of these operating instructions, only the German version is authoritative. Only the latest revision of the operating instructions applies. Due to constant technical development, the content of these MEDICON operating instructions is updated regularly. The version date of each edition of the operating instructions is included into the printout. MEDICON eG assumes no liability for any damage caused by improper use, incorrect postoperative handling, care or maintenance,

ENG

CE 0123

mediRod® Tri

Rx only

STERILE R

MD

medicon®
 MEDICON eG
 Gänsäcker 15, 78532 Tuttlingen
 Germany
 Tel.: (49) 7462 / 2009-0
 Fax: (49) 7462 / 2009-50
 E-Mail: sales@medicon.de
 Internet: http://www.medicon.de

Rev.: 12-12-2025

Instruction for use

or non-compliance with the restrictions on use and other guidelines in the operating instructions. MEDICON eG assumes no liability for defects in case of changes or repairs to the product without prior written consent from MEDICON eG either, nor in case of repairs not carried out by workshops authorised by MEDICON eG or by the MEDICON Repair Service.

11. Liability Description of symbols and icons

	Observe the operating instructions
	Caution
	CE label acc. to Directive 93/42/EEC
	Manufacturer
	Date of manufacture
	Do not reuse
	Do not re-sterilise
	Sterilised by irradiation
	MR Unsafe
	Humidity limitation
	Protect from sunlight
	Subject to medical prescription (US federal law)
	Production lot number, batch
	Item number
	Use by
	Do not use if packaging is damaged
	Store in a dry place
	Limit temperature

Please contact MEDICON eG if you have any further questions regarding the implants or instruments. If you do not have the brochure "Surgical technique and instruments", MEDICON eG can provide an up-to-date copy upon request.

PLEASE NOTE: According to US federal law, in the USA this product may be bought only by a physician or hospital or upon prescription.