Gebruiksaanweisung
Operating Instructions
Manual de instrucciones
Mode d'emploi
Istruzioni d'uso
Rev.: 2017-01-26

Rx only

Medicon eG · Gänsäcker 15
78532 Tuttlingen · Germany
Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50
E-Mail: sales@medicon.de
Internet: www.medicon.de

DE **MEDICON DREHMOMENTHANDGRIFF**

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EN

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

The MEDICON torque handle may only be used by surgeons who have sufficient experience with this product or have received an appropriate introduction.

- Table of contents:**
- General information
 - Intended use/Indication
 - Possible side effects and complications
 - Technical information
 - Handling and safety instructions
 - Dekontamination, cleaning and sterilization
 - Liability
 - Signs and Symbols

1. GENERAL INFORMATION
 The MEDICON torque handle is not supplied sterile and must be cleaned, disinfected and sterilized before every use.

Please observe the following instructions. These will guarantee accurate and reliable functioning of the product.

2. INTENDED USE / INDICATION
Torque limiter for the locking screw REF 714.11.50 of the MediRod pedicle screw and rod system. Manual operation, reusable. Non-adjustable torque limiting instrument for protection against torque overload.

3. POSSIBLE SIDE EFFECTS AND COMPLICATIONS
Risks and side effects have to be considered in connection with the overall MediRod pedicle screw system. Please consult the instructions for use **G7141210**. These instructions for use can be downloaded on the Medicon extranet at www.medicon.de/extranet or requested directly from MEDICON eG.

4. TECHNICAL INFORMATION
The torque for the torque handle has been pre-set to 10 Nm by the manufacturer. This instrument cannot be disassembled and is supplied non-sterile. It is intended for use with a limited time period and can be reused within its working life.
The instrument has to be used manually, i.e. it may not be used with an active device such as a power drill.
This instrument is not a measuring device!

- Product precision: ±10%
- Max. rec. axial force: 100 N
- Radial torque: 0.3 Nm
- Working area clockwise with click at nominal torque 10 Nm
- Working life: Either 3 years or 6000 clicks or 300 sterilizations.

5. INSTRUCTIONS AND REPROCESSING
Appropriate use and a functioning torque limiter are the prerequisite for safe use of the instrument. Any use other than the intended use may damage the instrument.
Safe use cannot be guaranteed with damaged instruments. Over-tightening the locking screw on a non-functioning torque handle may compromise the stability of the implanted pedicle screw system and cause the implant to fail.

The torque handle can be used for tightening and loosening the locking screws. Use the screw driver handle for this, REF 714.13.47 in combination with the Torx insert, REF 714.13.29. As soon as the torque limit of 10 Nm has been reached, the torque limiter disengages and produces an audible "click" to prevent application of any torque values higher than the preset value.
Discontinue use of the instrument if the torque limiter does not trigger despite application of appropriate force.

Warning and precautions

Failure to comply with these warnings and precautions can result in injuries, malfunctions or other unexpected incidents:

- The instrument has to be completely cleaned, disinfected, sterilized and tested for proper functioning before first use and before each subsequent use.
- It is very important to examine the surgical instrument for visible damage and wear, e.g. cracks, fractures or defects, before every use.
- Never use damaged instruments.
- The instrument may only be used, maintained and prepared by professional and trained medical personnel.

6. DEKONTAMINATION, CLEANING AND STERILIZATION

The following vCJD prion-specific protective measure is indicated when processing instruments.

In the event of diagnosis of a definite or probable case of vCJD:

If it is not possible to use disposable products, the instrument used, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

If prion contamination is suspected:
In the event of suspected prion contamination, incineration of the instrument is recommended according to the vCJD task force final report.

If vCJD is excluded:
Continue to use after instrument processing is completed. Other-wise, the instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.
In the case of a non-identifiable vCJD illness.
Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions – e.g. mechanical alkaline cleaning combined with steam sterilization.

If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (CNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134°C.

Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress-corrosion cracking.
Only cleaned and disinfected instruments may be sterilized.

Reprocessing restrictions
Reconditioning is limited to 300 sterilizations.

MEDICON TORQUE HANDLE

Location of use
– Directly after use, gross contamination should be removed from the instrument with a disposable cloth/paper.
– Fixing agents or hot water (>40°C) must not be used as these leads to fixation of residues and can affect cleaning outcome.
– The instruments should be conveyed promptly to processing.
– Dry disposal preferred.

Preparation for decontamination
To disassemble the MEDICON PEN sction cannulae, follow steps 1 to 5 in reverse order, ensuring that the screwing movements in steps 3 and 5 are also carried out in the opposite direction.
Articulated instruments must be processed in open position. For rinsing, the instruments must be placed in instrument holders suitable for machine cleaning.
The instrument holders (e.g. wire trays) must allow subsequent cleaning in an ultrasonic bath or in the cleaning and disinfection device (CDD) where they are not impeded by acoustic or rinsing shadows.

Pre-cleaning
– Place the instruments in cold water for 5 minutes.
– Brush the instruments (plastic brushes) under cold water until all visible soiling is removed.
– Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again.

Cleaning: machine
The cleaning and disinfection device (CDD) must meet DIN EN ISO 15883-1 requirements.

- Pre-rinse 1: 1 minute with cold demineralised water, without additive,
- Emptying,
- Pre-rinse 2: 3 minutes with cold demineralised water, without additive,
- Emptying,
- Cleaning: with demineralised water, heat to 55°C and wash/ clean for 10 minutes, add cleansing agent at 45°C, alkaline cleansing agent, strength 0.5%,
- Emptying,
- Neutralization: 3 minutes with warm water (>40°C) with add-tion of neutralizer, strength 1ml/l
- Emptying,
- Final rinse: 2 minutes with warm deionised water (>40°C) (without any additive),
- Emptying

Disinfection
Thermal Disinfection A0 value 3000: Completely desalinated water, the thermal disinfection is carried out at temperatures > 80 °C and corresponding application time according to the A0 concept, DIN EN ISO 15883 (e.g. A0 3000 = 90°C and 5 minutes application time). The operator is responsible for the implemented A0 value.

Drying
Complete drying must be ensured by the CDD. The instruments must be removed from the CDD promptly once the cleaning and disinfection program has finished.
If necessary, compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Cleaning/disinfection: manual
The cleansing agents and disinfectants employed must be suitable for manual cleaning and disinfection of instruments and must be compatible with each other. The disinfectant must be of tested efficacy. When selecting the disinfectant agent and method, the relevant lists and recommendations of the Robert Koch Institute (RKI) and the German Society for Hygiene and Microbiology (DGHM) should be noted.

Pre-cleaning
– Place the instruments in cold water for 5 minutes;
– Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
– Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again;
– Thorough rinsing with demineralised water

Cleaning/disinfection
– Place the instruments in a bath with a tested cleansing and disinfectant agent;
– The instruments must be completely covered with the solution;
– The application times, temperatures, and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed;
– Remove the instruments and rinse for 2 minutes with cold deionised water;
– Repeat the cleansing process if visible contamination is still present on the instrument;
– If necessary, clean and disinfect the instruments in an ultrasonic bath.

Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed sooner. A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

Drying
Dry manually with compressed air and a lint free cloth. Compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

Maintenance, inspection and testing
After cleaning/disinfection the instruments must be macroscopically clean, i.e., free from visible residues and soiling. Inspection is performed visually. Insufficiently clean instruments must be cleaned again and then adequately rinsed and dried. Before functional testing, instruments with movable parts must have cooled down and been lubricated with instrument maintenance oil (Medicon REF 46.00.40). Instruments such as clamps and needle holders, which have a detent, should only be closed up to the first detent (risk of stress cracking). Defective instruments (hairline cracks, deformation or wear) must be replaced as they no longer fulfil their function or do so without adequate safety. Corroded instruments must also be replaced as they may corrode intact instruments by extraneous rust.

Packing
The instruments must be placed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:
– DIN EN 868,

- DIN EN ISO 11607,
 - suitable for steam sterilization (vapor permeable),
 - adequate temperature resistance up to 138°C
- Sterilization equipment and sterilization wrapping must match both the wrap contents and the employed sterilization method.

Sterilization
Taking into account the respective national regulations, the following method must be employed for sterilization:
– Vacuum autoclave with triple vacuum and adequate drying of the products (Vacuum minimal 15 Minutes)
– Steam sterilizing complying with DIN EN 13060 or DIN EN 285 and validated in accord-ance with DIN EN ISO 17665-1.
– Sterilization time and temperature: at least 5 minutes hold time at 134°C

It is essential to attain a Sterility Assurance Level of 10⁻⁶.

Storage
Reprocessed sterile instruments must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments. The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority.

Further information on reprocessing
Validated machine cleaning and disinfection is always preferred over manual cleaning because of the greater certainty of the method. Good cleaning helps to preserve value and is a precondition of successful sterilization. During machine processing, the following points should be noted:
– Correct loading of the trays for rinsing is a precondition for effective machine processing. Trays must not be overloaded.
– Rinsing shadows due to large instruments must be avoided.
– The instruments must be placed or stored based on their susceptibility to mechanical damage in order to prevent them from becoming damaged.

The times and temperatures specified in these reprocessing instructions are minimum requirements and must not be less than those stated here. If they are to be reduced for technical reasons, this must be validated by the operating authority. Exceeding the stated times and temperatures is always possible but leads to increased stress on the material, which may result in premature ageing of the instruments. The use of other sterilization methods is outside our responsibility.

Information on validation of the processing
Validation was performed with the following equipment, materials, and chemicals:
Cleaning and disinfection device: Type Miele Disinfector G 7735 CD and 7836 CD compartmented cart for surgical instruments
Cleaning agent alkaline: neodisher® FA, Dr. Weigert GmbH & Co. KG
Cleaning agent enzymatic: Endozime, Ruhof (enzymatic)
Decondex 23 Neutrazym, Borer, Switzerland neodisher® FA, Dr. Weigert GmbH & Co. KG

Water pistol: Selecta
Cleaning brush: Plastic/nylon brushes
Ultrasonic bath: Sonorex
Steriliser: MMM Vaculab 969 S 3000, MMM Selectomat S 3000, Siefenhofer KS 666-2ED, H+P Varioclav 400E

Note
The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method. If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device.

7. LIABILITY
In the event of discrepancies between the english and the german version of these instructions for use only the German version shall be applicable. Only the latest revision of the instructions for use applies. Due to constant technical development the contents of these MEDICON instructions for use are updated regularly.
To ensure that you

Gebruiksaanweisung Operating Instructions Manual de instrucciones Mode d'emploi Istruzioni d'uso

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