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Operating Instructions

Manual de instrucciones

Mode d'emploi

Istruzioni d'uso

DE

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EN

Please read these instructions carefully before clinical application and keep them safe and close to hand.


Table of contents

- General notes
- Indications
- Possible side effects and contraindications
- Warning and precautions
- Secure handling
- Assembly
- Cleaning and disinfection
- Sterilization
- Storage
- Liability
- Signs and Symbols

1. GENERAL NOTES

The endoscopes are used for illumination and visualisation for diagnostic purposes and/or surgical therapeutic treatment. The endoscopes have been developed for different, specific areas of application. These may only ever be used in the area of application intended for each specific model. When selecting a model, users have to ensure that the chosen model is suitable for the specific task. The endoscopes may only be used by appropriately trained users with the relevant experience. The instructions for use do not provide or explain any clinical applications. The instructions for use are not intended as instructions for carrying out endoscopic examinations and surgeries.

The endoscopes are supplied non-sterile and have to undergo a validated processing cycle before each use (see chapters 7 and 8).

As a rule, only endoscopes labelled "AUTOCLAVABLE" may be treated in an autoclave. Endoscopes without this label cannot be processed in an autoclave.

2. INDICATIONS

Clinical syndromes where minimally invasive procedures are appropriate for diagnostics or therapy and where the examination or surgery is carried out visually inside the patient. In our catalogues and brochures, the endoscopes are allocated to the appropriate field of application, e.g. laparoscopy, arthroscopy, otoscopy, sinuscopy, etc.

The main objective of endoscopic diagnostics or endoscopic surgery is to protect the tissue and therefore improve retention of function to a limited degree.

3. POSSIBLE SIDE EFFECTS AND CONTRAINDICATIONS

Infection risk

Infection risks are of great importance during endoscopic examinations. The risk factors for increased infection risk are divided into two groups:

Process-related risks

- type and extent of tissue damage during therapeutic procedures
- surrounding circumstances of the endoscopic procedure (emergency procedure or elective surgery)
- competence and experience of the person carrying out the examination/ surgery
- correct cleaning and disinfection of endoscopes and accessories

Patient-related risks

- reduced immune status or immune suppression in the patient (HIV, leukaemia, lymphoma, immune suppression therapy, advanced liver or kidney disease, old age)
- presence of special sources of infection or special anatomic circumstances
- situations which favour the adhesion of bacteria in the organ (valvular defect, heart valve replacement, endoprosthesis, permanent intravenous catheter)

Endoscopic examinations can lead to endogenic transmission of the body's own microorganisms with subsequent bacteremia. The national and international recommendations for the preventive administration of antibiotics before certain procedures (ESGE guidelines 1998) have to be noted here.

4. WARNING AND PRECAUTIONS

Endoscopes have to be processed according to these instructions before each use. This also applies to brand-new devices. A function test has to be carried out on the endoscopes before each use and they have to be checked for possible damage. Discontinue use of damaged endoscopes immediately. This also applies to the products used in combination with the endoscopes, e.g. trocar sleeves, adapters, fiberoptic cables, etc.

The endoscopes may only be used with the combination products intended and offered by the manufacturer for this purpose. Before using the endoscopes, the user has to ensure that the combination products are suitable for the intended purpose and pose no specific risk in case of combined use.

The high-energy light emitted from the fiberoptic cable at the distal end can reach a temperature of up to 41 °C up to a distance of 8 mm. Consequently, direct contact of the tip of the endoscope with body tissue or flammable substances has to be avoided. Reduce the light intensity when working in the direct vicinity of the object. When adapting with a camera and the camera adapter, the

proximal end of the endoscope has to be completely dry before adapting to prevent the endoscope from misting up during the procedure.

5. SECURE HANDLING

Process the brand-new endoscope before use (supplied non-sterile). Check the endoscope for sharp edges and bent, loose or broken parts before each use. Check the endoscope for functioning and damage after each processing cycle. Discard any damaged endoscopes immediately. When using the endoscope in a trocar, avoid bending stress during insertion and extraction. Endoscopes may only be used by the physician or by medical personnel under supervision of a physician. Use of the endoscope requires appropriate training, knowledge and experience with clinical application of endoscopic techniques.

6. ASSEMBLY

The following section describes the mechanical assembly of the endoscopes. The adapter required in each case has to be screwed onto the fiberoptic cable connector of the endoscope and onto the endoscope end of the fiberoptic cable. The threads can be greased if required. Remove any excess grease. Ensure that the surfaces of the optical fibres are clean.

Light Post Adapter Assembly

Adapters can be used to fit most light guides. Simply screw the appropriate adapter on or off the fibre optic light post to prepare the Scope for connection to the light guide (see illustration).

7. CLEANING AND DISINFECTION

Cleaning
The instruments must be placed on appropriate instrument carriers so that they can be rinsed in the washing and disinfection device. The instrument carriers (e.g. sieve trays) must be such that subsequent cleaning in an ultrasonic device or in the washing and disinfection device is not impaired by any areas that are inaccessible to ultrasonics or rinsing. Please do only use suitable pH-neutral or alkaline cleaning agents.

The choice of the cleaning agent depends on the properties of the instruments and on the national guidelines and recommendations. In case of elevated chloride concentration of the water piping or stress corrosion may occur. By using fully desalinated water or alkaline cleaning agents the danger of such corrosions will be minimized.

Thermal disinfection

Carry out thermal disinfection, using fully desalinated water, at a temperature of 80-95 °C for an exposure time acc. to AD draft, EN ISO 15883.

Drying

Sufficient drying time must be provided by the cleaning and disinfection device.

Cleaning agents

Please observe all parameters off the manufacturer concerning concentration, temperature and reacting time. Automatic doser must be controllable.

ENDOSCOPES

Control and Checking

After cleaning/disinfection the instruments must be macroscopically clean, i.e. free of visible residues. Visual checking. Instruments that are not sufficiently cleaned must be cleaned once again and rinsed sufficiently afterwards. Exchange defective instruments (fissures, deformation or abrasion) since they do not fulfil their function any more or not safely enough.

The Endoscope has three optical surfaces that must be cleaned and checked routinely to ensure both maximum transmission of illumination and a high quality image.

Cleaning

The Endoscope has three optical surfaces (see illustration):

- Proximal eye funnel (eye-piece)
- Distal tip
- Fibre optic connector

All should be cleaned with warm water and mild soap. If stains are still present, a mixture (1:1) of methyl alcohol and acetone may be used. A fine woven cloth or lens tissue should be used for cleaning. Dry with a soft woven cloth. By reflecting light off all optical surfaces, you can detect any foreign matter or scratches which would degrade image quality.

Note: Any foreign matter present on the fibre surface after cleaning may tend to burn and discolour the surface when the high intensity lamp is in use.

Removing deposits

Due to not sufficient cleaning or foreign matter in steam, deposits may develop on the three optical surfaces during the autoclaving process.

Any foreign matter remaining on the optical surfaces will reduce the optical performance of the endoscopes and may lead to a cloudy picture. These deposits can easily be removed using the bio-compatible polishing paste which is enclosed to each endoscope. To remove deposits, dab some polishing paste onto a clean cotton tipped swab. Press gently the swab onto the optical surface to be cleaned and scrub the window in circles. Rinse the optical surfaces with tap water to remove any remaining polishing paste and clean the surfaces as instructed in chapter 2.

Note: Cleaning with polishing paste should not be part of your routine cleaning procedures and should only be performed when the image - as viewed through the scope - is cloudy.

Chemical disinfection and steris process

The Scope is immesible and should be disinfected using the process obtained from the individual chemical manufacturer. HSW Scopes are materially compatible with the Steris process. The Endoscopes should be rinsed after soaking with sterile water and dried with a sterile cloth.

8. STERILIZATION

Important: Be sure that the Scope is cleaned thoroughly before sterilizing.

The Endoscopes should be sterilized or disinfected in a container which secures the instrument in place. Be sure the needle portion does not experience any force or stress which can destroy the delicate internal optics.

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