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
Arthroscope Shafts and Irrigation Cannulas

Contents

- General information
- Indication
- Application
- Instructions for reprocessing
- Disassembly and assembly

GENERAL INFORMATION

These instruments are not supplied in sterile condition and must be cleaned, disinfected and sterilised before being used for the first time and before every subsequent use.

Please observe the following instructions. They will ensure that the instruments will function correctly and reliably.

INDICATION

Arthroscope shafts are used to insert endoscopes and to irrigate the joint during arthroscopy.

Irrigation cannulas are used to irrigate the joint.

APPLICATION

Insert only the specified MEDICON trocars and endoscopes in correct diameters and lengths into the arthroscope shafts. Use only the specified trocars for the irrigation cannulas.

INSTRUCTIONS FOR REPROCESSING

WARNINGS

- Stainless-steel instruments must not be immersed in physiological saline solution (NaCl solution), because extended contact with the solution may cause corrosion such as pitting corrosion and stress corrosion.
- Only cleaned and disinfected instruments may be sterilised.

Restrictions for reprocessing

Frequent reprocessing has little effect on these instruments. The end of service life is usually obvious because of wear and damage caused by use.

Place of use.

- Immediately after use, remove coarse residue from the instruments with a disposable cloth or paper towel.
- Do not use fixing agents or hot water (>40° C), because this may fix residues and affect the cleaning results.
- The instruments must be processed immediately after use.
- Dry heat disposal is preferred.

Preparations for decontamination

Arthroscope shafts and irrigation cannulas must be disassembled, because otherwise areas that could have been in indirect contact with the patient may not be cleaned. This may pose a risk of infection. For manual processing, system-specific brushes and cleaning guns are required. The arthroscope shafts and irrigation cannulas should also be treated by ultrasonic cleaning if available.

The arthroscope shafts and irrigation cannulas must be disassembled and deposited ready for rinsing on trays that are suitable for machine washing. The instrument trays (such as sieve trays) must be constructed in such a way that subsequent cleaning by ultrasonics or cleaning and disinfecting machine will not be hampered by acoustic or cleaning shadows. The instruments must be rinsed by a connection to the tray, attaching them to special

nozzles etc. Internal lumina that are not connected will not be reached by the cleaning agent and will therefore not be cleaned and disinfected.

Cleaning

Please use only suitable pH-neutral or alkaline cleaning agents. The choice of the cleaning agent depends on the properties of the instruments and on the national guidelines and recommendations. If the water has an elevated chloride concentration, pitting or stress corrosion of the instruments may occur. The use of fully demineralised water or alkaline cleaning agents will minimise the danger of such corrosion.

The cleaning and disinfection machine must comply with the requirements of DIN EN ISO 15883-1.

1. Prerinse 1: 1 minute with softened cold water without additives
2. Drain
3. Prerinse 2: 3 minutes with softened cold water without additives
4. Drain
5. Cleaning: with softened water, heat to 55° C and wash and clean for 5 minutes, add cleaning agent at 45° C, alkaline cleaning agent 0.5% strength
6. Drain
7. Neutralisation: 3 minutes with warm water (>40° C) with neutraliser added at 1 ml/L
8. Drain
9. Final rinse: 2 minutes with warm demineralised water (>40° C) (without additives)
10. Drain

Please observe all parameters specified by the manufacturer for concentration, temperature and application time. Automatic dosers must be adjustable.

Thermal disinfection

Use demineralised water. The thermal disinfection process takes place at 80–95° C with an application time as specified by the AO concept, EN ISO 15883-1.

Drying

Please ensure sufficient drying time with the cleaning and disinfection machine. The temperature must not exceed 95° C. If necessary, compressed air is recommended for drying because of its fast and effective action (recommended by the RKL).

Cleaning agents

Please observe all parameters specified by the manufacturer for concentration, temperature and application time. Automatic dosers must be adjustable.

Maintenance, inspection and testing

After cleaning and disinfection the instruments must be macroscopically clean, i.e. free of visible contamination and residues. They are inspected visually. Critical areas such as joints and serrations require particularly careful inspection. Instruments that are not sufficiently or totally cleaned must be cleaned once again and thoroughly rinsed afterwards. Instruments with moving parts must be cool and oiled with instrument oil before testing their function. Instruments such as clamps and needle holders that have a locking function must only be closed to the first click to prevent tension cracks. Defective instruments (hairline cracks, deformation or wear) must be replaced, because they will no longer function correctly or will not function at all.

Sterilization

Only cleaned and disinfected instruments may be sterilized.

The following sterilization procedures must be used, subject the national requirements:

- Triple fractionated vacuum process with sufficient product drying
- Steam sterilizer conforming to EN 13060 or DIN EN 285 and validated in accordance with DIN EN 17665-1.
- Sterilization time and temperature:
at least 5 minutes hold time at 134° C

For sterilization use only the following sterilization procedures:

- Fractionated vacuum procedure with sufficient product drying
- Steam sterilizer conforming to EN 13060 or DIN EN 13485 and validated in accordance with DIN EN 554.
- Sterilization time and temperature:
at least 5 minutes hold time at 134° C

Packaging

The instruments must be stored in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868 and DIN EN ISO 11607
- suitable for steam sterilization (steam-permeable)
- sufficient temperature resistance to 138° C

Sterilization accessories and sterilization packaging must be designed for the package contents and the sterilization process.

Storage

Reprocessed sterile instruments must be stored in suitable reusable sterilization containers in dry, dark and cool rooms protected from dust, microbes and vermin. They should be protected from large temperature fluctuations to prevent the formation of condensate. Do not store chemicals together with instruments. The storage life depends on the type of packaging.

The approved storage time at the location depends on the type of sterile barrier system and the storage conditions. The owner must specify the approved storage period.

Note

The user is responsible for validating the recommended parameters to take account of variations in sterilization chambers, packaging methods and loading configurations.

The same applies if procedures are adopted that differ from those recommended. It is essential to reach a SAL (Sterility Assurance Level) of 10⁻⁶.

MEDICON eG, as supplier of these products, is not liable for direct or consequential defects caused by or arising from improper use or handling, especially through non-observance of the enclosed operating instructions or through improper care or maintenance.


Disassembly

Arthroscope shafts:

1. Unscrew knurled nut **A**.
2. Remove cover ring **B**.

Arthroscope shafts and irrigation cannulas:

3. Unscrew spring cap **C** and remove cock **D** from the tap. Remove both cocks from the model with two taps.

Assembly

Before assembly, spray all moving parts with spray oil (Medicon part no. 46.00.40). Carry out Assembly in reverse order as described under Disassembly.

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