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Please read these instructions carefully before clinical application and keep them safe and close to hand.

MEDICON implants are made of pure titanium (in accordance with ASTM F67) or Ti6Al4V alloy (ASTM F136). Both materials are biocompatible, corrosion-resistant and non-toxic in a biological milieu.

**Contents**

1. General notes
2. Purpose and proper use/indications
3. Contraindications
4. Possible side effects and complications
5. Single Use Product
6. MR Notes
7. Selection of implants
8. Use/handling of plates
9. Use/handling of screws
10. Instruments
11. Decontamination, cleaning and sterilization
12. Liability
13. Signs and symbols

**1 GENERAL NOTES**

The MEDICON Manus Modular Osteosynthesis System for hand surgery is supplied unsterile and must be disinfected, cleaned and sterilized prior to use.

Please observe the following instructions. This will guarantee that the implants function faultlessly and reliably.

**2 PURPOSE AND PROPER USE/INDICATIONS**

The Manus Modular Osteosynthesis System for hand surgery is used for osteosynthesis of bone fractures and the stabilisation of small and medium bones as part of hand, foot, paediatric, plastic and trauma surgery and orthopaedic procedures: fractures and luxations, osteotomies, bone implantations pseudoarthrosis. The implants are designed to support normal bone healing, but not to replace normal structures of the body or to bear body weight in the event of inadequate bone healing. MEDICON plate systems should only be used by surgeons with sufficient training and experience at doctor 's surgeries, hospitals and clinics.

**3 CONTRAINDICATIONS**

- The use of plates and screws in body parts and bones other than those for which they are intended.
- Use despite acute infections.
- Patients that are unable to follow the instructions for postoperative care. This may be caused e.g. by psychological/mental or neurological problems of the patient.
- Patients with inadequate bone tissue, or bone tissue of insufficient quality, and patients with circulatory disorders or latent infections.
- Extreme material sensitivity, i. e. foreign-body reaction of the patient. For this case, appropriate tests must be carried out under all circumstances prior to implantation (even when there is only a suspicion of such condition).
- Acute infections.
- Documented titanium intolerance.

**4 POSSIBLE SIDE EFFECTS AND COMPLICATIONS**

- Complaints, pain, abnormal sensations or palpability of the implant.
- Extreme material sensitivity of the patient (foreign-body sensitivity) manifested by allergic reactions, etc.
- Heightened reaction of connective tissue in the region of the fracture and/or implant.
- Insufficient bone formation, osteolysis, osteomyelitis, osteoporosis, inhibited revascularization or infections that can result in loosening, bending, rupture or breakage of the implants.
- Delayed or insufficient bone union, which can result in implant breakage.
- Bone necrosis.
- Inadequate healing.
- If the implants were not adapted properly, bone fusion can be delayed, or fusion is insufficient or does not occur at all.
- Breakage, bending, migration or loosening of the implant.
- Reduced bone density due to stress shielding.
- Tendon rupture due to implant erosion.

Apart from the side effects and complications mentioned above, the surgical intervention can also lead to problems such as nerve injuries, infections, pain, etc., which are not directly attributable to the implant. Extreme bending of the plates, or

repeated bending at the same place, can lead to plate fractures. The implants can come loose if the screws are not tightened sufficiently.

**5 SINGLE USE PRODUCT**

Implants are designed and constructed for single use only on a patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

**6 MR NOTES**

**Instruments**

The use of medical devices in an MR environment poses a hazard. Individual medical devices must not be placed in the immediate vicinity of the equipment during the performance of MR procedures.

**Implants**

The implants are not MR-safe and have not been tested for MR safety. Use in an MR environment may therefore pose a risk. It is recommended that patients and medical personnel be informed that the implants are not approved for use with MR equipment.

**7 SELECTION OF IMPLANTS**

The surgeon is responsible for the correct selection of implants.

The following aspect are crucial for the correct selection of implants:

- Bone defect to be treated or intended bone repositioning
- Body weight of the patient
- general health and degree of activity of the patient.

Incorrect implant selection can lead to premature implant loss and loosening, deformations or fractures of the implants. Among other factors, the success of the intervention depends on how the implants are handled. The plates must be bent carefully and cautiously, avoiding repeated and extensive deformation. Damage and scratches reduce the stability of the product and result in early fatigue of the implant material. In cases of failed, delayed or insufficient bone union, the incidence of implant bending, implant failure or implant breakage cannot be excluded. The implants and instruments are adapted to each other. Any use of implants and instruments from other manufacturers in combination with MEDICON products entails unpredictable risks, since the products are not adapted to each other. To avoid risks, only such MEDICON products that are designated for combination with each other may be combined with each other. Prior to applying the product, the surgeon must thoroughly discuss the intended surgical result with the patient. In such discussions, special attention must be given to the postoperative aspects and aftercare, which may become necessary. The patient must be informed that, due to the limited stability of the fixation, the fixation site must not be put under the full body weight. The patient must also be informed the non-observance of this advice can put at risk any safe bone union. The patient must be instructed to inform, without delay, the operating surgeon of unusual changes at the operating site. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing. The implants can be removed once the fracture has healed completely. An implant left in the body after complete healing may act as a loadbearing element and can contribute to an increased risk of refracturing in active patients. For this reason, too, we accept liability for the implants only up to complete healing. We will not accept liability for any damage or injuries caused by implants that were not removed after complete healing of the treated fracture. The surgeon is responsible for deciding whether the implants shall remain in the body after complete healing of the fracture or whether they shall be removed. For his decision, the surgeon has to weigh risks and advantages of leaving the implants in the body. It cannot be excluded that implants will break, loosen, migrate through tissue and cause pain. The packaging label carries a LOT no. We recommend transferring this lot no. to the patient records as this number allows tracing the production history of the implant back to the raw material.

**8 USE/HANDLING OF PLATES**

The plates have to be adapted to the natural or planned bone contours as precisely as possible. For adapting the plates use the instruments provided for this purpose by the manufacturer. During the plate bending procedure, there occur cold stresses, under which the titanium hardens while, at the same time, its ductility is decreased. Therefore it is absolutely necessary to achieve the re-quired plate shape in as few bending steps as possible. Excessive bending can lead to postoperative plate breakage. Overly aggressive application of the bending instruments can result in visible damage to the implant. When such damage has occurred, a new plate must be selected and adapted more cautiously. Deformed plate holes for the screws entail not only an increased risk of breakage at the location of the holes; they also compromise the precise seating of the screw head in the plate. Therefore the plates must be bent carefully under all circumstances. Deformed plates must be inspected for notches, deformed screw holes and other mechanical damage prior to application. Always select the plates that provide the closest fit to the specific situation and indication. If a plate needs to be

shortened, only use the instruments supplied for this purpose. The plate must be cut between the screw holes. Care must be taken that the portion of the plate that is cut off is not flung towards the patient or towards third parties during the cutting procedure. The edges of a shortened plate must be deburred to prevent tissue injuries or irritations.

**9 USE/HANDLING OF SCREWS**

Self-tapping screws do not require precutting of threads. To apply the drill holes, only use the respective drill designated by MEDICON for the specific screw so that the correct drill hole diameter for the screw is obtained. Only use drill bits with sharp cutting edges. During drilling, there is a risk of heat damage to the bone. Therefore, only work at low drill speed and always ensure sufficient cooling during the drilling procedure. Excessive application of force during the drilling procedure can lead to drill breakage, endangering the patient, the surgeon and third parties. If, due to a failed drilling effort or a slack drill hole, the screw cannot be anchored firmly in the bone, the next larger screw diameter intended for this purpose must be used. To ensure that the screwdriver blade is firmly seated in the screw head, axial force must be applied when inserting the screwdriver into the screw head. This ensures the correct axial alignment of screw and screwdriver, and prevents slipping of, or damage to the screw head or the screwdriver blade. Application of excessive force in tightening the screws can lead to screw breakage during the intervention.

Wear of the screwdriver blade will impair the firm connection between blade and screw. Whenever such wear is detected, the screwdriver blade must be replaced by a new blade. Once the implantation has been completed, the secure connection between every screw and plate must be checked. If necessary retighten the screws. Prior to implant removal, all screw heads must be cleaned carefully, using a scalpel or another suitable instrument, so that the screw-driver can be optimally seated in the screw head.

**10 INSTRUMENTS**

The MEDICON instruments intended for application of the system are subject to wear and tear and mechanical strains even when used normally, but even more so when applied with excessive force. To prevent failure of, or mechanical damage to the instruments during the operation, prior to each use the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional.

**11 DECONTAMINATION, CLEANING AND STERILIZATION**

The following vCJD (Creutzfeld-Jakob-Disease) prion specific protective measure is indicated when processing instruments.

**In the event of diagnosis of a definite or probable case of vCJD:**

- If it is not possible to use disposable products, the implant and instrument used, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

**If prion contamination is suspected:**

- In the event of suspected prion contamination, incineration of the instrument is recommended according to the vCJD task force final report.

**If vCJD is excluded:**

- Continue to use after instrument processing is completed. Otherwise, the implant and instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

**In the case of a non-identifiable vCJD illness:**

- Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions - e.g. mechanical alkaline cleaning combined with steam sterilization.

If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (CNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134° C.

Note: Please observe the relevant national regulations concerning the reprocessing of the products for patients with CJD or suspected CJD or possible variants of CJD.

Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress corrosion cracking. Only cleaned and disinfected implants may be sterilized.

Reprocessing restrictions  
Frequent reprocessing has a slight impact on implants and instruments.

**Location of use**

Implants that have been contaminated by blood and/or secretions or are suspected to be contaminated must not be placed back in the storage tray and reprocessed. They have to be disposed in a suitable container. STERILIZE Instruments that have been contaminated by blood and/or secretions or are suspected to be contaminated must not be cleaned in the implant/instrument storage tray. They have to be cleaned separately and after proper cleaning they can be placed in the implant/instrument storage tray.

**Preparation for decontamination**

Before manual / machine cleaning the instruments have to be taken out of the implant/instrument storage tray. The screwdriver blade has to be removed from the screwdriver handle. Place the lid back onto the implant/instrument storage tray. The implant/instrument storage tray must be placed in instrument holders suitable for machine cleaning. The instrument holders (e.g. wire basket) must allow subsequent cleaning in an ultrasonic bath or in the cleaning and disinfection device (CDD) where they are not impeded by acoustic or rinsing shadows.

**Machine cleaning - alkaline and disinfection**

- Place the implants in the implant storage tray for 15 minutes in ultrasonic bath at 40° C with 0.5% alkaline cleaner (Neodisher Mediclean forte) and treat with ultrasonic waves.
- Remove implant storage tray and rinse thoroughly with cold water.

A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

The cleaning and disinfection device (CDD) must meet DIN EN ISO 15883-1 requirements.

**Pre-rinse 1: 3 minutes with cold demineralised water, without additive;**

- Emptying;
- Cleaning: with demineralised water, heat to 55° C and wash / clean for 5 minutes, add cleansing agent at 45° C, alkaline cleansing agent, strength 0.5%;
- Emptying.
- Neutralization: 3 minutes with warm water (>40° C) with addition of neutralizer, strength 1ml / l;
- Emptying;
- Final rinse: 2 minutes with warm deionised water (>40° C) (without any additive);
- Emptying

Thermal Disinfection A0 value 3000: Completely desalinated water, the thermal disinfection is carried out at temperatures > 80° C and corresponding application time according to the A0 concept, DIN EN ISO 15883 and guideline DGKH, DGSV and AKI (e.g. A0 3000 = 90° C and 5 minutes application time). The operator is responsible for the implemented A0 value."

Complete drying must be ensured by the CDD, e.g. 60° C and 30 minutes. The instruments and implant tray must be removed from the CDD promptly once the cleaning and disinfection program has finished. If necessary, compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

**Manual Cleaning: Implant tray**

- Rinse the implants under cold running tap water to remove visible soil.
- Use a soft bristled brush and lumen brush to aid in removal of visible soil while rinsing.
- Using a syringe, flush hard to reach areas while rinsing
- Prepare EnzoZ® enzymatic detergent at 1oz. per gallon warm tap water.
- Immerse the tray to soak for 1 minute into EnzoZ® enzymatic detergent
- Following the soak time, brush the implant tray with a soft bristled brush and lumen brush for a minimum of 1 minute to remove soil residuals.
- Using a syringe, flush hard to reach areas with the detergent solution.
- Remove the implant tray from the enzymatic detergent solution and rinse under running tap water to remove detergent residuals.
- Prepare a neutral detergent Valsure® neutral in a sonication unit at ¼ oz. per gallon of warm (37 - 40° C) tap water. Immerse in the sonicator and actuate.
- Allow the implant tray to sonicate for 10 minutes

EN

CE<sub>0123</sub>

MEDICON  
implantArt  
OSTEOSYNTHESIS

Rx only

MD

## MANUS MODULAR OSTEOSYNTHESIS SYSTEM

Rev.: 15-12-2025

- Remove the implant tray from the sonicator and rinse them under cold distilled water to remove detergent residuals.
- Dry the implant tray completely using medical compressed air and disposable, lint-free cloths.

**Manual cleaning: Instruments**

- Rinse instruments under cold running tap water for at least one minute to remove visible soil and actuate movable parts while rinsing. Use a soft bristled brush and lumen brush to aid in removal of soil while rinsing and actuate movable parts while brushing. Use a syringe to flush the lumen (internal spaces, threads, and holes) with rinse water.
- Prepare Enzo!® enzymatic detergent at 8 ml/l of warm tap water and immerse instruments to soak for at least one minute.
- Following the soak time, brush the instruments with a soft bristled brush and lumen brush for a minimum of one minute to remove soil residuals.
- Using a syringe, flush internal spaces, threads, and holes with the enzymatic detergent solution
- Remove the instruments from the enzymatic detergent solution and rinse under running tap water to remove detergent residuals
- Prepare a neutral detergent Valsure® neutral in a sonication unit at 2 ml/l of warm (37 - 40° C) tap water. Immerse in the sonicator and actuate.
- Using a syringe, flush the lumen (internal spaces, threads, and holes) with the detergent solution and then allow the instruments to sonicate for at least 10 minutes
- Remove the instruments from the sonicator and rinse them under cold distilled water to remove detergent residuals.
- Dry the instruments completely using medical compressed air and disposable, lint-free cloths.

**Maintenance, Inspection, and Testing**

After cleaning the implants and instruments must be macroscopically clean; free from visible residues and soiling. Inspection is performed visually. Insufficiently cleaned implants and instruments must be cleaned again and then adequately rinsed and dried. Deformed or damaged implants and instruments must be removed and disposed of, as they can no longer be guaranteed to be safe to use.

**Packaging**

The implant- and instrument storage tray must be placed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868,
- DIN EN ISO 11607,
- suitable for steam sterilization (vapor permeable),
- adequate temperature resistance up to 138° C.

Sterilization equipment and sterilization wrapping must match both the wrap contents and the employed sterilization method.

**Sterilization**

Taking into account the respective national regulations, one of the following methods must be employed for sterilization:

- Vacuum autoclave with triple vacuum and adequate drying of the products (Vacuum minimal 15 Minutes)
- Steam sterilizer complying with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes hold time at 134° C/273° F

It is essential to attain a Sterility Assurance Level of 10<sup>-6</sup>.

**Storage**

After sterilization, the medical devices must be stored in the sterilization packaging in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin. To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with implants and instruments. The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority. Sterile packaging should be carefully examined prior to opening to ensure that package integrity has not been compromised.

**Disposal**

National regulations have to be considered when disposing the products, its components and its packaging.

**Further information on reprocessing**

Validated machine cleaning and disinfection is always preferred over manual cleaning because of the greater certainty of the method. Good cleaning helps to

preserve value and is a precondition of successful sterilization. During machine processing, the following points should be noted:

- Correct loading of the storage trays for rinsing is a precondition for effective machine processing. Storage trays must not be overloaded.
- Rinsing shadows due to large instruments must be avoided.

**The times and temperatures specified in these reprocessing instructions are minimum requirements and must not be less than those stated here. If they are to be reduced for technical reasons, this must be validated by the operating authority. Exceeding the stated times and temperatures is always possible but leads to increased stress on the material, which may result in premature ageing of the implants and instruments. The use of other sterilization methods is out-side our responsibility.**

**Machine cleaning - alkaline and sterilization:**

Cleaning and Type Miele Disinfector G 7735 CD and 7836 CD disinfection device:  
compartmented cart for surgical instruments Cleaning agent  
alkaline: neodisher® Fa. Dr. Weigert GmbH & Co. KG Neutraliser:  
neodisher® Z, Dr. Weigert GmbH & Co. KG  
Ultrasonic bath: Sonorex  
Steriliser: MMM Vaculab 969 S 3000, MMM Selectomat S 3000,  
Stiefenhofer KS 666-2ED, H+P Varioclav 400E

**Manual cleaning:**

Cleaning agent: Enzo!® enzymatic  
Neutralizer: Valsure® Neutral detergent  
Ultrasonic bath: Bronson sonicator  
Additional equipment: Lumen brushes, soft-bristled brushes, syringes

**Note**

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method. If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device. Please note that all instruments sent to the MEDICON Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

**12 LIABILITY**

MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use and other guidelines in the instructions for use. The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

**13 SIGNS AND SYMBOLS**

	Manufacturer
	Date of manufacture
	Lot number
	Item number
	Non-sterile
	Warning
	Consult instructions for use
	MR unsafe
	CE marking according to Directive 93/42/EEC
	Medical device
	Prescription only

Federal law and regulations in the United states restrict the sale of this device to a physician or hospital or on prescription!

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Detergenti alcalini: neodisher® Fa. Dr. Weigert GmbH & Co. KG
Neutralizzatore: neodisher® Fa. Dr. Weigert GmbH & Co. KG