



steam with holding time of at least 5 minutes to a maximum of 20 minutes and vacuum drying for at least 10 minutes. Sterilizer (Class B) must be in accordance with applicable national standards and regulations (e.g. DIN EN 13060 or DIN EN 285). The sterilization standards DIN EN ISO 17665 (Sterilization of health care products - Moist heat) must be observed. Electrosurgical electrodes shall not be sterilized with hot air, EtO gas, gamma radiation or plasma.

NOTE: Before use, the electrosurgical electrodes must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the electrosurgical electrodes after the sterilization process. In case the aforementioned described and recommended chemicals and machines for manual pre-cleaning and/or automated cleaning and disinfection are not available, it is the responsibility of the user to validate their procedure. Also, if another sterilization procedure is chosen, the procedure deviant from the procedure described in paragraph 5 has to be validated by the user accordingly.

Limitation of reconditioning:
The lifetime of the product will depend on wear, handling, and duration of use, damage and frequency of reconditioning.

6 Visual inspection



Prior to each use, a visual inspection must be carried out on the electrosurgical electrodes. In particular the insulation of the electrosurgical electrodes must be inspected for pressure points or damage. Electrosurgical electrodes with damage or pressure points shall not be used and must be replaced by new ones.

7 Exclusion of repair and modification

Defective electrosurgical electrodes shall not be repaired. They must be replaced by new electrosurgical electrodes. Furthermore, unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

For service and repair please contact your national representative of MEDICON eG

8 Packaging, storage, transportation, handling

The electrosurgical electrodes must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The electrosurgical electrodes must always be handled with the utmost care during transportation, cleaning, disinfection, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

9 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping. Use our Return Form for returns.

10 Liability

CAUTION: According to US federal law, in the USA this product may be bought only by a physician or hospital or upon prescription!
In the event of any discrepancies between the non-German and the German version of these operating instructions, only the German version shall be authoritative. Only the latest revision of the operating instructions applies.
Due to constant technical development, the content of these MEDICON operating instructions is updated regularly. Please use the IFU-Website ifu.medicon.de to ensure that you are using the current version. The version date of the respective edition of the operating instructions is included in the printout. MEDICON eG assumes no liability for damage caused by improper use, incorrect postoperative conduct, care or maintenance, or non-compliance with the restrictions on use and other guidelines in the operating instructions. MEDICON eG will not accept liability either for changes or repairs made to the product without the prior written authorisation of MEDICON eG, nor for repairs that were not carried out by workshops authorised by MEDICON eG or the MEDICON Repair Service (MRS).

11 Disposal

Disposal of damaged instruments or of instruments that are not functionally reliable: To avoid a risk of infection to third parties, the instruments must be cleaned and sterilised before disposal.
In addition, the instruments must be disposed of in specific and appropriately labelled containers, in order to protect third parties from cutting injury. Observe the applicable national laws when disposing of reusable surgical instruments!
The product, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

12 Description of symbols and icons

Icon	Description
	Manufacturer
	Production lot number, batch
	Item number
	Non-sterile
	Caution
	Observe the operating instructions
	CE label acc. to Directive 93/42/EEC
	Medical device
	Store in dry conditions
	Protect from sunlight