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Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

Medicon eG assumes no liability as the distributor of these products for direct or consequential damages caused by improper usage or handling, particularly due to non-compliance with the enclosed usage instructions or due to improper care or maintenance. The reusable surgical instruments are manufactured from stainless materials. The corrosion resistant material features excellent qualities in a biological environment.

Table of contents:

1. General Information
2. Purpose and proper use/indications
3. Contraindications
4. Possible side effects and complications
5. CT and MRI instructions
6. Application/Handling of instruments
7. Decontamination, cleaning and sterilization
8. Liability
9. Signs and symbols

1 GENERAL INFORMATION

The bipolar HF instruments are delivered in non-sterile state and must therefore be disinfected and sterilised before single application and before every sterile provision in the surgical area. Please observe the following instructions. These will guarantee accurate and reliable functioning of the product.

2 PURPOSE AND PROPER USE/INDICATION

The bipolar HF instruments are used for invasive surgical procedures and are applied in combination with a HF surgery device. The functions include cutting, separating, suction and preparing with simultaneous coagulation of biological tissue. In bipolar surgery, an active coagulation electrode is used. The high-frequency current flows through this bipolar active electrode into the tissue to be coagulated, and through a so-called neutral electrode out of the tissue again. The decisive factor is that the contact area of the active electrode on the tissue is significantly smaller than the contact area of the neutral electrode. Following appropriate procedures, the bipolar HF instruments can be reused. The bipolar HF instruments by MEDICON may be used only by surgeons who have sufficient experience with this product, or who have received an appropriate introduction.

3 CONTRAINDICATIONS

- Applications contrary to the intended use of the product.
- Reuse without reprocessing (disinfection, cleaning and sterilization)
- Use without prior function test by the user / operator
- If, due to the patients indications, use of the bipolar HF instruments would pose a potential risk, the instruments must not be used. This decision is to be made by the physician.

4 POSSIBLE SIDE EFFECTS AND COMPLICATIONS

- Incidents that have been reported in connection with the use of bipolar systems:
- Unintentional activation resulting in tissue damage at an improper point and/or damage to the equipment.
 - Fire in connection with drapes and other flammable materials
 - Alternating current paths leading to burns in locations where the patient or user comes into contact with uninsulated instruments.
 - Forceps tips can cause injuries!
 - After application, forceps tips may be so hot that they can cause injuries.
 - Explosions caused by sparking in the presence of flammable gases.
 - Endogenous burns are caused by excessive current density in the patient's tissue.
 - Perforation of organs

In addition to the undesirable effects and complications already mentioned, the surgical procedure may lead to other problems such as nerve damage, infections, pain etc.

5 CT & MRT HINWEIS

Bipolar HFI-instruments pose a risk in the proximity of a CT or MRT and must not be close to such devices during their use.

6 APPLICATION/HANDLING OF INSTRUMENTS


Even when used normally, bipolar HF instruments are subject to wear and mechanical stresses, but even more so if used too forcefully. In order to prevent failure or mechanical damage to the instruments during surgery, these parts must be checked before each use to ensure that they are mechanically intact, that there

are no deformations, and that the parts are fully functional. Damaged instruments must not be used, and are to be replaced. The bipolar HF instruments by MEDICON must be connected to a HF surgical device via a suitable bipolar cable and may be operated only with bipolar coagulation current. Observe the application and safety information by the manufacturer of the HF surgery device. The surgeon is responsible for correct patient-specific selection and use of the instruments.

Failure to comply with these warnings and precautions may result in injuries, malfunctions or other unexpected incidents.

- These instruments may be used, maintained and processed only by qualified and trained medical personnel.
- All instruments must be completely cleaned, disinfected, sterilized and tested for proper functioning before the first and each subsequent use.
- It is very important to examine each surgical instrument for visible damage and wear, e.g. cracks, fractures or defects, on the insulation before each use. Especially areas such as blades, tips, closures, catches and latches, as well as all moving parts, insulations and ceramic elements must be inspected carefully.
- Never use damaged instruments.
- The electric connectors of the bipolar HF instruments must be fully inserted into the cable connectors. Only push/pull the flat connector to connect and disconnect the connection cable.

- Do not use in the presence of flammable or explosive substances.
- The instrument must not be deposited on the patient.
- Coagulate only when the contact areas are in the field of vision and have a good contact with the tissue to be coagulated. Do not touch any other metallic instruments, trocar sleeves, optical elements, etc. Activation of the HF voltage can lead to capacitive coupling when the working end does not contact the tissue to be coagulated or is not properly positioned for the emission of energy.
- HF voltage and laser must never be activated simultaneously. The working end must be pulled back from the laser fibre while the laser is in use so that the laser cannot be accidentally focussed on the working end or the shaft insulation. Conversely, to avoid arcing, the laser fibre must be pulled back when the working end is activated, especially when the laser fibre is surrounded by metal. Observe the operating instructions of the laser system for correct application of the laser.
- Observe the application and safety information by the manufacturer of the HF surgical device.
- The plug connections between the bipolar MEDICON HF instruments, the HF connection cable and the HF surgical device may be established and removed only when the HF surgical device is switched off. Non-compliance may lead to electric shocks or burns.
- The output power of the HF surgical device must always be set to the lowest value that produces the desired surgical effect. If despite use of the default settings the usual coagulation effect cannot be achieved, the output power of the device must never be increased without checking the system components first. The maximum peak voltage of the bipolar HF instruments is 250 Vp.
- When using flushing instruments at the same time, non-conductive flushing liquid should be preferred.
- Pacemakers or other active implants can be damaged by HF currents. For patients with pacemakers, consult a cardiologist before the procedure

7 DECONTAMINATION, CLEANING AND STERILIZATION

The following vCJD prion-specific protective measure is indicated when processing instruments.

In the event of diagnosis of a definite or probable case of vCJD:

If it is not possible to use disposable products, the instrument used, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

If prion contamination is suspected:

In the event of suspected prion contamination, incineration of the instrument is recommended according to the vCJD task force final report.

If vCJD is excluded:

Continue to use after instrument processing is completed. Otherwise, the instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

In the case of a non-identifiable vCJD illness:

Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions - e.g. mechanical alkaline cleaning combined with steam sterilization.

If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (CNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134 ° C.

Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stresscorrosion cracking. Only cleaned and disinfected instruments may be sterilized.

Reprocessing restrictions

Frequent reprocessing has few effects on these instruments. The end of their working life is normally determined by wear and damage through use. If there are any reprocessing restrictions, this is noted in the instructions for use of the specific product.

LOCATION OF USE

- Directly after use, gross contamination should be removed from the instruments with a disposable cloth/paper.
- Fixing agents or hot water (>40 ° C) must not be used as this leads to fixation of residues and can affect cleaning outcome.
- The instruments should be conveyed promptly to processing.
- Dry disposal preferred: Direkt nach der Anwendung sind grobe Verschmutzungen mit einem Einmaltuch/-papier von den Instrumenten zu entfernen.

PREPARATION FOR DECONTAMINATION

Articulated instruments must be processed in open position.

For rinsing, the instruments must be placed in instrument holders suitable for machine cleaning.

The instrument holders (e.g. wire trays) must allow subsequent cleaning in an ultrasonic bath or in the cleaning and disinfection device (CDD) where they are not impeded by acoustic or rinsing shadows.

PRE-CLEANING

- Place the instruments in cold water for 5 minutes;
- Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
- Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again: Einlegen der Instrumente für 5 Minuten in kaltem Wasser;

CLEANING: MACHINE

The cleaning and disinfection device (CDD) must meet DIN EN ISO 15883-1 requirements.

- Pre-rinse 1: 1 minute with cold demineralised water, without additive;
- Emptying;
- Pre-rinse 2: 3 minutes with cold demineralised water, without additive;
- Emptying;
- Cleaning: with demineralised water, heat to 55 ° C and wash / clean for 10 minutes, add cleansing agent at 45 ° C, alkaline cleansing agent, strength 0.5%;
- Emptying.
- Neutralization: 3 minutes with warm water (>40 ° C) with addition of neutralizer, strength 1ml / l;
- Emptying;
- Final rinse: 2 minutes with warm deionised water (>40 ° C) (without any additive);
- Emptying

DISINFECTION

Thermal Disinfection A0 value 3000: Completely desalinated water, the thermal disinfection is carried out at temperatures > 80 ° C and corresponding application time according to the A0 concept, DIN EN ISO 15883 (e.g. A0 3000 = 90 ° C and 5 minutes application time). The operator is responsible for the implemented A0 value.*

CLEANING/DISINFECTION: MANUAL

The cleansing agents and disinfectants employed must be suitable for manual cleaning and disinfection of instruments and must be compatible with each other. The disinfectant must be of tested efficacy. When selecting the disinfectant agent and method, the relevant lists and recommendations of the Robert Koch Institute (RKI) and the German Society for Hygiene and Microbiology (DGHM) should be noted.

PRE-CLEANING

- Place the instruments in cold water for 5 minutes;
- Brush the instruments (plastic brushes) under cold water until all visible soiling is removed;
- Internal spaces, threads, and holes are each rinsed with the water pistol for 10 seconds and brushed again;
- Thorough rinsing with demineralised water

CLEANING/DISINFECTION

- Place the instruments in a bath with a tested cleansing and disinfectant agent;
- The instruments must be completely covered with the solution;
- The application times, temperatures, and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed;
- Remove the instruments and rinse for 2 minutes with cold deionised water;
- Repeat the cleansing process if visible contamination is still present on the instrument;
- If necessary, clean and disinfect the instruments in an ultrasonic bath.

Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed sooner. A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed. Instrumente in ein Bad mit einem geprüften Reinigungs- und Desinfektionsmittel einlegen.

DRYING

Dry manually with compressed air and a lint free cloth. Compressed-air drying is recommended because of its good and rapid action (recommendation by the German RKI).

MAINTENANCE, INSPECTION, AND TESTING

After cleaning/disinfection the instruments must be macroscopically clean, i.e., free from visible residues and soiling. Inspection is performed visually. Insufficiently clean instruments must be cleaned again and then adequately rinsed and dried. Before functional testing, instruments with movable parts must have cooled down and been lubricated with instrument maintenance oil (Medicon REF. 46.00.40). Instruments such as clamps and needle holders, which have a detent, should only be closed up to the first detent (risk of stress cracking). Defective instruments (hairline cracks, deformation or wear) must be replaced as they no longer fulfil their function or do so without adequate safety. Corroded instruments must also be replaced as they may corrode intact instruments by extraneous rust.

PACKING

The instruments must be placed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- DIN EN 868,
- DIN EN ISO 11607,
- suitable for steam sterilization (vapor permeable),
- adequate temperature resistance up to 138 ° C.

Sterilization equipment and sterilization wrapping must match both the wrap contents and the employed sterilization method.

STERILIZATION

Taking into account the respective national regulations, the following method must be employed for sterilization:

- Vacuum autoclave with triple vacuum and adequate drying of the products (Vacuum minimal 15 Minutes)
- Steam sterilizer complying with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN ISO 17665-1.
- Sterilization time and temperature: at least 5 minutes hold time at 134 ° C

It is essential to attain a Sterility Assurance Level of 10-6.

EN

Rx only

Rev.: 14-01-2026

HIGH FREQUENCY INSTRUMENTS – BIPOLAR

STORAGE

Reprocessed sterile instruments must be stored in a suitable reusable sterilization container in a dry, dark, cool, and semi-sterile place, protected from dust and free from vermin.

To avoid the development of condensation, major temperature fluctuations should be avoided during storage. Chemicals must not be stored together with instruments.

The walls, floors, and ceilings of the storage room should be smooth and easy to clean and disinfect. Shelves must be at least 30 cm off the floor. The duration of storage permitted depends on the type of sterile barrier system employed and the storage conditions. This storage period must be established by the operating authority.

FURTHER INFORMATION ON REPROCESSING

Validated machine cleaning and disinfection is always preferred over manual cleaning because of the greater certainty of the method.

Good cleaning helps to preserve value and is a precondition of successful sterilization. During machine processing, the following points should be noted:

- Correct loading of the trays for rinsing is a precondition for effective machine processing. Trays must not be overloaded.
- Rinsing shadows due to large instruments must be avoided.
- The instruments must be placed or stored based on their susceptibility to mechanical damage in order to prevent them from becoming damaged.

The times and temperatures specified in these reprocessing instructions are minimum requirements and must not be less than those stated here. If they are to be reduced for technical reasons, this must be validated by the operating authority. Exceeding the stated times and temperatures is always possible but leads to increased stress on the material, which may result in premature ageing of the instruments. The use of other sterilization methods is outside our responsibility.

INFORMATION ON VALIDATION OF THE PROCESSING

Validation was performed with the following equipment, materials, and chemicals:
Cleaning and disinfection device:

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Cleaning agent enzymatic:
Endozime, Ruhof (enzymatic)
Decondex 23 Neutrazym, Borer,
Switzerland

Neutraliser:
neodisher® Z
Fa. Dr. Weigert GmbH & Co. KG

Ultrasonic bath:
Sonorex DIGITEC DT 156 BH

Steriliser:
MMM Vaculab 969 S 3000
MMM Selectomat S 3000
Stiefenhofer KS 666-2ED H+P
Varioclav 400E

NOTE

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method. If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device. Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

8 LIABILITY

MEDICON eG assumes no liability for damages resulting from improper use, care, maintenance, or failure to comply with the usage restrictions and other instructions specified in the instructions for use.

Furthermore, MEDICON eG's liability for defects shall be excluded in the event of modifications or repairs to the product without prior written consent from MEDICON eG, as well as for repairs carried out by workshops not authorized by MEDICON eG or by the Medicon Repair Service (MRS).

9 SIGNS AND SYMBOLS

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Batch code, lot numvber
	Catalogue number / Article number
	Caution
	Consult instruction for use
	MR unsafe
	CE marking in aacordance with Directive 93/42/EEC
	Medical Device
	Prescription only

Federal law and regulations in the United States restrict the sale of this device to a physician or hospital or on prescription!

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