Gebruiksaanweisung
Operating Instructions
Instrucciones de uso
Manuel d'instructions
Istruzioni d'uso

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BIPOLAR NON-STICK FORCEPS

 **Please read these instructions carefully before clinical application and keep them safe and close to hand.**

MEDICON eG as the supplier of these products is not liable for damage caused by, or arising from improper use or handling, especially due to non-observance of the operating instructions or improper cleaning or maintenance of the products. The MEDICON bipolar non-stick forceps sterile may be used only by sufficiently trained and experienced surgeons.

- Table of contents:**
- General notes
 - Specific function
 - Possible side effects and complications
 - Technical Informationen
 - Application and handling
 - Reprocessing
 - Signs and symbols

1. GENERAL NOTES
MEDICON bipolar non-stick forceps are not supplied sterile and must be cleaned, disinfected and sterilized before every use. Please observe the following instructions. This will guarantee that the instrument functions faultlessly and reliably.

2. PURPOSE AND PROPER USE
MEDICON bipolar forceps are used to grasp, prepare and coagulate tissue.

3. RISKS AND UNDESIRABLE EFFECTS
Incidents that have been reported in relation to using bipolar systems:

- unintentional activation resulting in tissue damage in the wrong place and/or damage to the equipment;
- medical drapes or other flammable materials being set alight;
- alternating current paths that can lead to burns on sites where the patient or user comes into contact with uninsulated parts;
- explosions, caused by sparks in an environment with flammable gases;
- perforated organs;
- suddenly heavy bleeding;
- bipolar tweezers have proven to be ineffective for tubal sterilisation or tubal coagulation during the sterilisation process and should not be used for this purpose.

4. TECHNICAL INFORMATION
Maximum electrical operating voltage, Umax: 0,5 kVp
Appropriate connection cable:
MEDICON bipolar cable with flat plug (e.g. [REF](#) 88.00.21)

5. SAFETY NOTES
MEDICON bipolar forceps must be connected with an appropriate bipolar cable to the bipolar output of an HF generator and must only be operated with bipolar coagulation current.

 **Warnings and precautionary measures**

If the warnings and precautions for use are not heeded, injuries, malfunction or other unexpected incidents may occur.

- These instruments may only be used, maintained and set up by proficient, trained medical staff.

- All instruments must be fully cleaned, disinfected, sterilised and checked to make sure they are fit for operation before being used for the first time, and before being used each time thereafter.
- It is very important to examine every surgical instrument before it is used for visible damage, and wear and tear, such as fissures or cracks, or insulation defects. Particular attention must be paid to checking areas such as blades, points, ends and notches as well as all moveable parts, insulations and ceramic elements.
- Never use damaged instruments.
- The electric flat pins of the bipolar forceps must be inserted fully into the cable socket. Hold the flat plug only to connect and remove the connection cable.
- Do not use in the presence of flammable or explosive materials.
- The instrument must not be put down on the patient.
- Only coagulate if the contact surfaces are in the field of vision and are in good contact with the tissue selected for coagulation. When doing so, do not touch any other metal instruments, trocar sleeves, optical devices or such like.
- Please pay attention to the instructions for use and safety indications by the manufacturer of the HF surgical device.
- The plugs between the MEDICON bipolar forceps, the HF connection cable and the HF surgical device may only be put in and taken out when the HF surgical device is switched off. If this measure is not respected, there is a risk of an electric shock or burns.
- The output power of the HF surgical device is always set to the lowest value required to obtain the required surgical effect. If coagulation does not occur as usual at the standard setting, the output power must never be increased without checking the system components beforehand.
- Pacemakers or other active implants may be damaged by HF current. A cardiologist is to be brought in before any surgical procedures are carried out on patients fitted with a pacemaker.

6. CONDITIONING
Reconditioning
The product design, the materials used, and the purpose of the device are such, that it is impossible to determine an upper limit to the number of preparation cycles that can be carried out. The lifespan of the instruments depends on how they are used and how well they are looked after. Instruments for electrosurgery are naturally subject to a higher level of wear and tear, which in turn depends on the type and duration of use.

Preparation and Transport
Immediately after the use, remove any gross contamination from the instruments with a disposable cloth or paper. Do not use any fixing agent or hot water (>40 °C), as this results in the fixation of residues and may affect the success of the cleaning. The instruments must be stored and taken to the preparation site in a closed container to prevent damage to the instruments and environmental contamination.

6.1 Reprocessing: machine
Place the instruments on a perforated tray on the slide-in trolley, or on the insets of the MIC cart, and start the cleaning procedure.

- Prerinse 1 minute with cold water
- Emptying
- Prerinse 3 minutes with cold water
- Emptying
- Wash for 5 minutes, heat to 55°C add cleansing agent at 45°C, alkaline cleansing agent, strength 0.5%.
- Emptying

- Neutralization: 3 minutes with warm water (>40°C) with addition of neutralizer
- Emptying.
- Final rinse: 2 minutes with warm water (>40°C)
- Emptying

Disinfection
Carry out thermal disinfection acc. to A0 draft (EN ISO 15883-1). Comply with the applicable regulations for processing medical products in the country concerned.

Drying
The outside of the instruments should be dried using the drying cycle of the cleaning/disinfection device. If necessary, the instruments can be dried off manually using a lint-free cloth. Instrument cavities are to be dried with sterile compressed air.

6.2 Reprocessing: manual
Cleaning
Prepare the cleaning bath in accordance with the manufacturer's instructions.

- Rinse the products under cold tap water (< 40°C) until all visible soiling has been removed. Stubborn soiling should be removed with a soft brush.
- Products are to be placed fully in the prepared cleaning bath. Leave to act for the length of time specified by the manufacturer.
- Clean the instrument that has been placed in the bath manually with a soft brush. All surfaces must be brushed off several times.
- Rinse the products very thoroughly under running water to remove any remaining cleaning product.

Disinfection
Prepare the disinfection bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfection bath and respect the necessary application time. Rinse the products very thoroughly under demineralised water to remove any remaining disinfectant.

Drying
Dry manually with compressed air and a lint free cloth. Compressed-air drying is recommended because of its good and rapid action.

6.3 Control and packing
Check the products to make sure they are clean; and if applicable, carry out an assembly and functional test in accordance with the user manual. If necessary, repeat the preparation process until the instrument is clean to the eye. The instruments are to be packed for sterilisation according to ISO 11607 and EN 868

MEDICON Bipolar non-stick forceps
The polished precious metal tips of the MEDICON non-stick forceps can tarnish like silver. This does not in any way harm the performance of the product. The original shine can be brought back through rubbing with a cloth for cleaning silver.

6.4 Sterilization
Sterilisation with fractionated vacuum method (in accordance with ISO 13060/ ISO 17665) subject to the applicable national requirements:

- 3 pre-vacuum phases with at least 60 mbar of pressure
- Heat up to a sterilisation temperature of at least 132°C; max. 137°C
- Shortest holding time: 3 minutes
- Drying time: at least 10 minutes.

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6.5 Storage
The sterilised instruments are to be stored in a dry, clean, dust-free environment in a moderate temperature of 5°C to 40°C.

6.6 Repairs
Do not carry out any repairs yourself. Servicing and repairs are to be carried out by trained, qualified personnel only. Contact the manufacturer or your technical department if you have any questions about repairs.
Caution! Defective products must be fully reprocessed before being sent back for repairs.

6.7 Information on validating the preparation procedure
Validation is carried out with the following apparatus, materials and chemicals:

Cleaning agent (automatically):	Neodisher® Fa. Dr. Weigert (alkaline) Endozime, Fa. Ruhof (encymatic) Enzol Enzym. Detergent, Johnson&Johnson neodisher® Z; Dr. Weigert
Cleaning agent (manually):	Miele G 7735 CD
Neutraliser:	Trolley for surgical instruments E 327-06
Cleaning and disinfection equipment:	Miele MIC-Wagen E 450

Details see record
SMP GmbH # 01707011901-2 (Cleaning automatically)
Northview Laboratories #P8H066 (Cleaning manually, sterilization)
Nelson Labs # 200432706-02 (Sterilization)
MDS GmbH Testbericht 084183-10 (Sterilization)

If the above-mentioned chemicals and machines are not available, it is incumbent upon the user to validate the process accordingly.

6.8 Handling
All surgical instruments should always be handled with great care when being moved, cleaned, maintained, sterilised and stored. This applies in particular to blades, fine tips, and other sensitive areas.

6.9 Warranty
MEDICON eG, as the original distributor of these products, does not accept any liability for immediate or consequent damage caused by improper use or handling, especially due to non-observance of the operating instructions or inappropriate care or maintenance of the products. Do not make changes to the product under any circumstances. This is not permitted and will result in the exclusion of liability by MEDICON eG.

7. SIGNS AND SYMBOLS

 Follow instructions

 CE label

 Manufacturer

 Warning

 Order-No.

Please contact MEDICON eG if you have further questions concerning the instruments.

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SMP GmbH # 01707011901-2 (Pulizia meccanica)
MDS GmbH Testbericht 084183-10 (sterilizzazione)