R only

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Rev.: V3 23-02-2026

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Rx only

Rev.: V2 12-02-2026

INVASIVE AND NON-INVASIVE MEDICAL DEVICES

The instructions for use and the reprocessing instructions must be read carefully prior to clinical use and kept available in a safe and easily accessible location at all times.

Table of contents

- 1 General Safety Information
- 2 Possible Complications
- 3 MRI Information
- 4 Application and Handling
- 5 Cleaning, Disinfection and Sterilization
- 6 Service
- 7 Liability
- 8 Disposal
- 9 Symbols and Pictogram Explanation
- 10 Basic UDI-DI / Description / Intended Use / Groups / Indications / Contraindications

1 General Safety Information

- The products are reusable and supplied non-sterile. They must be reprocessed before first use and after each subsequent use in accordance with the section "Cleaning, Disinfection and Sterilization."
- The medical devices are made from corrosion-resistant, biocompatible materials suitable for surgical applications.
- The product may only be used by trained and qualified personnel.
- Use the product exclusively for its intended purpose (see Intended Use) and familiarize yourself with its handling prior to first use.
- New or unused products must be cleaned after removal from the transport packaging prior to the first sterilization cycle and then stored in a dry, clean and protected environment.
- Before each use, the products must be visually and functionally inspected (e.g. for damage, deformation, free movement, complete markings).
- Damaged or non-functional products must not be used and must be removed from service or repaired.
- Products designed to be disassembled must be fully disassembled for reprocessing.
- For patients with confirmed or suspected CJD/vCJD, national reprocessing guidelines must be strictly followed.
- Serious incidents related to the product must be reported to the manufacturer and the competent authority.

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3 MRI information

The medical devices must not be used outside their intended purpose or defined scope of application. Complications may arise from improper use, use outside the intended purpose, or from inadequate reprocessing.

4 Application and Handling

The surgeon is responsible for selecting the appropriate medical devices for the intended procedure. Reusable medical devices are subject to wear and mechanical stress during normal use, particularly when excessive force is applied. Use only accessories and instrument components specifically designated by MEDICON eG to avoid risks related to product compatibility.

5 Cleaning, Disinfection and Sterilization

Note on Prions (CJD/vCJD):

For reprocessing in cases of confirmed or suspected prion disease (CJD/vCJD), the applicable national guidelines must be followed, in particular those issued by the Robert Koch Institute (RKI).

The reprocessing procedures described in this document apply exclusively to standard reprocessing processes. For patients with confirmed or suspected CJD/vCJD, the operator must initiate reprocessing in accordance with the RKI recommendations or treat and/or dispose of the instrument set as specified in the respective guidelines.

5.1 Preparation for Decontamination

Before cleaning, the instruments must be removed from the instrument tray and placed in an arrangement suitable for rinsing on appropriate instrument holders (e.g., wire mesh trays). The positioning must ensure that cleaning in the washer-disinfector (WD) is not impaired by shadowing or covered surfaces. Gross contamination must be removed immediately after use to prevent organic material from drying.

This includes:

- Rinsing under cold water,
- Actuating all movable parts
- Flushing lumens, threads and bores using a syringe or suitable flushing device

Products with multiple components must be fully disassembled prior to cleaning (see the respective disassembly instructions). During transport and interim storage, the instruments must be kept moist (e.g., with damp cloths or suitable pre-treatment solutions) to prevent drying of contaminants.

Damaged or non-functional products must not be reprocessed and must be removed from use.

5.2 Cleaning and Disinfection

5.2.1 Automated Cleaning and Disinfection

The washer-disinfector (WD) must comply with the requirements of DIN EN ISO 15883-1.

- Pre-rinse 1: 1 minute with softened cold water (no additives)
- Drain
- Pre-rinse 2: 3 minutes with softened cold water (no additives)
- Drain
- Cleaning:
 - With softened water, heating to 55 °C
 - Dosage of the cleaning agent at 45 °C
 - 10 minutes washing/cleaning
 - Cleaning agent: neodisher MediClean forte, concentration 0.5%
- Drain
- Neutralization:
 - 3 minutes with warm water (>40 °C)
 - Neutralizing agent: neodisher Z, concentration 1 ml/l
- Drain
- Final rinse:
 - 2 minutes with warm demineralized water (>40 °C), without additives
- Drain

Water Quality According to Validation:

For the final rinse, low-microbial (<10 CFU/ml) and low-endotoxin water (<0.25 EU/ml) was used during validation.

In routine practice, demineralized water is typically used, provided it meets the technical requirements for WD reprocessing. Ensuring the appropriate water quality is the responsibility of the operator in accordance with the specifications of the WD and the facility.

Thermische Desinfektion:

Thermal disinfection is performed in accordance with the A₀ concept defined in DIN EN ISO 15883-1.

Minimum requirement:

– A₀ = 3000, e.g. 90 °C for 5 minutes, using demineralized water

The responsibility for achieving the required A₀ value lies with the operator.

Drying

Adequate drying must be ensured by the washer-disinfector (e.g., 60 °C / 30 minutes). The instruments must be removed from the WD immediately after completion of the cleaning and disinfection cycle. If necessary, the use of compressed air is recommended for drying due to its efficient and rapid effect (RKI recommendation).

5.2.2 Manual Cleaning

Note:

Manual cleaning was assessed as part of the manufacturer's validation. Its effectiveness is confirmed only within the validated parameters described below. Deviations (e.g., different concentrations, chemicals, times, or temperatures) must be validated by the operator in accordance with DIN EN ISO 17664-1.

1. Pre-cleaning

- Remove gross contamination under cold tap water.
- Actuate all movable parts completely.
- Flush lumens, threads, and bores at least 5 times using a syringe or flushing device.
- Keep the instruments moist until cleaning to prevent drying of organic residues.

2. Enzymatic Cleaning (Soaking + Brushing + Ultrasonic Cleaning)

- Prepare enzymatic cleaner (e.g., Cidezyme®/Enzol® at 0.8%).
- Soak the instruments for 1 minute without ultrasound.
- Brush all surfaces and lumens (at least 1 minute).
- Ultrasonic cleaning: 5 minutes with the bath activated, temperature 37–40 °C.

3. Rinsing

- Rinse thoroughly under running water.
- Flush lumens at least 5 times.
- Ensure that all cleaning agent residues are completely removed from the instruments.

4. Manual Disinfection

- Immerse the instruments in an OPA solution (e.g., Cidex® OPA).
- Observe the exposure time specified by the manufacturer.
- Fill lumens with the disinfectant solution.

5. Final Rinse

- Rinse the instruments at least 5 times with tap water.
- Perform a final rinse with demineralized or purified water to avoid residues.
- Flush lumens at least 3 times.

6. Drying

- Completely dry the instruments using compressed air.
- Actively blow through lumens.
- If necessary, finish drying with a lint-free cloth.

5.4 Inspection / Visual and Functional Check

After cleaning, disinfection, and drying, all products must be visually and functionally inspected. The inspection ensures complete cleanliness, integrity, and proper functionality of the medical device.

The following criteria must be checked:

- Cleanliness: no visible contamination, no residues of cleaning or disinfecting agents, no moisture in lumens or joints
- Surfaces: no corrosion, no discoloration, no chipping, no excessive scratches
- Mechanics: free movement of all parts, full functionality, correct alignment
- Lumens: fully patent, free of residues or blockages
- Markings: article numbers, laser markings, and codes must be fully legible

Actions in case of deviations:

- If visible residues are present: re-clean and re-disinfect the product
- In case of functional defects or damage: remove the product from service, repair, or replace it
- Products with illegible markings or missing identification must not be reused

Inspection must be carried out by qualified personnel prior to each sterilization cycle and before every clinical use.

5.5 Maintenance / Lubrication

After inspection, products with movable parts must be maintained according to their functional requirements.

Lubrication:

- Joints and movable parts may only be lubricated if required for proper function.
- Only biocompatible instrument oils suitable for steam sterilization may be used (e.g., white oils without additional additives).
- Lubricants must be applied sparingly and only to the designated functional areas.
- Lubricants must not be introduced into lumens, threads, or bores unless explicitly required.

Maintenance:

- Damaged, worn, or improperly functioning instruments must be removed from service and repaired if necessary.
- Products may only be sterilized in a fully assembled and functional condition.

The operator is responsible for the regular functional inspection, maintenance, and documentation.

5.6 Packaging

Before sterilization, the medical devices must be packaged in suitable sterile barrier systems that comply with the requirements of DIN EN ISO 11607-1/2.

Packaging requirements:

- Use a validated sterile barrier system (e.g., film/nonwoven pouches, Tyvek® pouches, container systems)
- Packaging materials must be suitable for the intended sterilization method (e.g., steam sterilization)
- Packaging must be undamaged, clean, and functional
- The packaging must enable complete sterilization and allow aseptic presentation

Loading Instructions (Operator Responsibility):

- Instruments must be completely dry, clean, and free of residues.
- The arrangement of instruments in trays, containers, or pouches must follow the validated loading instructions of the CSSD/AEMP.
- Instruments must not damage the packaging material.
- The load must ensure that steam can reach all surfaces to be sterilized.

Operator Responsibility:

- Selection of the sterile barrier system and validation of the packaging process lie with the operator.
- Packaging must be inspected for integrity before each sterilization cycle.

Only fully dried and correctly packaged products may be sterilized.

5.7 Sterilization

Country	Fractionated Vacuum Method
Deutschland	at least 5 min at 134 °C (273 °F)
USA mind.	at least 4 min at 132 °C (270 °F), minimum drying time 20 min
Frankreich	at least 5 min at 134 °C (273 °F)
Other countries	At least 5 min at 132 °C (270 °F) / 134 °C (273 °F)

Note:

The sterilization parameters listed above correspond to the conditions evaluated during the manufacturer's validation. The validation of the sterilization process and the respective load configurations is the responsibility of the operator.

Note on unsuitable sterilization methods:

Only fractionated pre-vacuum steam sterilization processes according to the parameters listed above are suitable for this product. Other sterilization methods such as ethylene oxide (EO), formaldehyde, gas plasma / hydrogen peroxide (H₂O₂), hot air, or other low-temperature processes are not permitted.

5.8 Information on Reprocessing Validation

The validation of the reprocessing procedure was carried out using the devices, chemicals, and process parameters listed below. Effectiveness has been demonstrated only within these conditions. Deviations must be validated by the operator in accordance with DIN EN ISO 17664-1.

Washer-Disinfectors (WD):

- Miele G 7735 CD
- Loading rack for surgical instruments

Automated Cleaning Agent:

- neodisher® MediClean forte (alkaline)
- Concentration: 0,5 %
- pH (working solution): 10,4 – 10,8

Neutralization:

- neodisher® Z (acidic)
- Dosage: 1 ml/l
- pH (working solution): 2 – 3

Manual Cleaning (Enzymatic):

- Cidezyme® / Enzol®
- Dosierung: 0,8 %
- pH (working solution): 7,5 – 8,5

Note on alternative procedures:

Manual cleaning procedures and the use of alternative cleaning agents (e.g., enzymatic cleaners, ultrasonic baths, alternative neutralizing agents) have not been validated. The use of such procedures lies within the responsibility of the operator and must be validated on site in accordance with DIN EN ISO 17664-1.

Note:

The use of other cleaning agents, WD models, or differing process parameters lies within the responsibility of the operator and must be evaluated as part of an operator validation in accordance with DIN EN ISO 17664-1.

5.9 Limitation of Reusability

The end of the service life of reusable medical devices is determined by wear, material fatigue, and damage that may occur both during intended use and through repeated reprocessing. Even under normal handling conditions, the products are subject to mechanical stress, particularly when excessive force is applied.

English
Operating instructions

Rev.: V2 12-02-2026

Before each use, the medical device must be carefully inspected by qualified personnel for mechanical integrity, deformation, and full functionality. A thorough visual and functional inspection is the most effective method for determining the end of the service life and preventing device failure during the procedure.

Typical signs of damage or wear include:

- Corrosion (e.g., rust, pitting)
- Discoloration or material changes
- Excessive scratches, abrasion, or chipping
- Cracks or fractures
- Bent or deformed instrument components
- Restricted movement of joints or mechanical parts
- Missing, damaged, or illegible markings (e.g., article numbers, laser engravings)Korrosion (z. B. Rost, Lochfraß)

Sorting out and replacement:
Medical devices exhibiting one or more of the above characteristics must not be reused. They must be replaced, repaired, or disposed of appropriately.

A fixed maximum number of reprocessing cycles is not specified, as actual service life depends on usage, handling, and reprocessing conditions.

6 Service

For service and repairs, please contact the national representative of MEDICON eG in your country.

7 Liability

CAUTION: Under U.S. law, this product may be purchased in the United States only by a physician, hospital, or on the order of a licensed healthcare professional.

Validity of the Instructions for Use:

Only the most recent revision is valid. Due to ongoing technical developments, these Instructions for Use are updated regularly. The version date and revision number are indicated on the document. Please ensure that you are using the current version.

Disclaimer:

MEDICON eG assumes no liability for damage resulting from:

- improper application,
- use outside the intended purpose,
- incorrect handling, care, or maintenance,
- failure to comply with these Instructions for Use,
- modifications or repairs without prior written approval from MEDICON eG,
- or repairs carried out by unauthorized service centers.

Warranty claims also lapse in the event of modifications or repairs performed without prior written approval by MEDICON eG, or if maintenance or repairs are carried out by unauthorized workshops.

8 Disposal

Medical devices must be cleaned and sterilized prior to disposal to prevent infection and injury risks to third parties. They must be disposed of in designated, appropriately labeled containers. National and local regulations must be observed.

9 Symbols and Pictogram Explanation

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Batch code / LOT number
	Catalogue number
	Non-sterile
	Caution
	Consult instructions for use
	MR unsafe
	CE marking
	Medical device
	Prescription only

INVASIVE AND NON-INVASIVE MEDICAL DEVICES

10 Basic UDI-DI / Description / Intended Use / Groups / Indications / Contraindications

Basic UDI-DI	4046826;N001;01;WC 4046826;N001;02;WF 4046826;N001;03;WJ 4046826;N001;04;WM 4046826;N001;05;WQ 4046826;N001;06;WT	4046826;N002;01;WP 4046826;N002;02;WS 4046826;N002;03;WV 4046826;N002;04;WY 4046826;N002;05;X3 4046826;N002;06;X6 4046826;N002;07;X9	4046826;N003;01;X2 4046826;N003;02;X5 4046826;N003;03;X8 4046826;N003;04;XB 4046826;N003;05;XE	4046826;N004;01;XD 4046826;N004;02;XG 4046826;N004;03;XK 4046826;N004;04;XN 4046826;N004;05;XR 4046826;N004;06;XU
Description	Reusable, non-invasive medical devices intended for temporary use in cutting, holding and/or grasping medical devices or tissue, shaping tissue or medical devices, or for autopsy purposes.	Reusable, non-invasive medical devices not used in combination with active medical devices, intended for temporary diagnostic use, holding or pressing intact tissue, or marking.	Reusable, non-invasive medical devices not used in combination with active medical devices, intended for storing, positioning and/or presenting tissue, fluids, other medical devices, or medical auxiliaries.	Reusable, invasive medical devices intended for temporary use in spreading, pressing, retracting, grasping, extracting, or dilating tissue or medical devices, as well as for diagnostic purposes.
Groups	Knives, scissors / saws, mills, files / forceps / tweezers / clamps / handles, holders, magnets, application instruments	Clamps / stethoscope / tuning fork / percussion hammer and bulbs / compressors / springs, markers / compasses	Instrument baskets and rinsing trays / silicone mats / inlays, trays, cassettes, holding and separating elements / containers / filters	Mirrors, optics / probes, catheters, cannulas / loops and ligatures / handles and holders / retractors, spreaders, spatulas / forceps, dissection instruments (dental / podiatry)
Intended Use	These medical devices are specialized tools intended for temporary use during medical procedures and can be reused multiple times after proper cleaning and sterilization. They are designed to enable precise and safe medical interventions and may be used invasively where applicable. (They may only be connected to active devices of Class I.)			
Indication	Cutting, holding, grasping medical devices or intact tissue; shaping tissue or medical devices.	Holding or pressing intact tissue; diagnostic use; marking.	Storage, positioning, and presentation of tissue, fluids, medical devices or medical auxiliaries.	Spreading, pressing, retracting, grasping, extracting, dilating intact tissue or medical devices; diagnostic use.
Contraindication	Invasive or surgically invasive use.		Surgically invasive use.	
	These products are not intended for direct contact with the heart and must not be used on the central nervous system (CNS) or the central circulatory system.			

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