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Carefully read the operating instructions and processing instructions prior to clinical application and keep them safe and at hand. The instructions and notes contained therein must be followed. On the product label you will find the processing instructions (MPK) valid for your product in the section 'Special Instructions'.

TABLE OF CONTENTS

- 1 General instructions for safe handling
- 2 Possible complications
- 3 MRI advice
- 4 Use and handling
- 5 Limitation of reusability (shelf life)
- 6 Service
- 7 Processing (cleaning, disinfection and sterilisation) of products
- 8 Liability
- 9 Disposal
- 10 Description of symbols and icons
- 11 Basic UDI-DI / Description / Intended purpose / Indication / Contra-indications

1 GENERAL INSTRUCTIONS FOR SAFE HANDLING

- The medical devices are reusable, are supplied non-sterile and, for this reason, must be cleaned, disinfected and sterilised according to the instructions in section 'Processing' prior to their first use.
 - Only allow the product and accessories to be operated and used by persons who have the necessary training, knowledge or experience with regard to application, function check and cleaning/sterilisation.
- The user as well as the relevant qualified personnel undertake to familiarise themselves with the products before their use.
- Read and observe the operating instructions.
- Only use the product as intended (see "Intended purpose").
- Clean brand-new product after removing the transport packaging and before the first sterilisation.
- Store brand-new or unused product in a dry, clean and protected location.
- Before each use, qualified personnel must:
 - visually check the product for loose, loose, bent, broken, cracked, worn and broken off parts.
 - check the product for proper function.
- Do not use a damaged or defective product. Immediately sort out damaged products or send them to the service centre specified in these operating instructions.
- Replace damaged parts immediately with original spare parts.
- All medical devices that can be disassembled must be disassembled for processing and sterilisation.
- All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is based.

2 POSSIBLE COMPLICATIONS

The medical devices must only be used for the intended use and purpose. Complications can be caused by non-functional or incorrectly processed medical devices.

3 MRI ADVICE

Use of the medical devices in the proximity of MRI systems is dangerous. The individual medical devices must not be in the vicinity of the devices while MRI procedures are being carried out.

4 APPLICATION AND HANDLING

The surgeon is responsible for selecting the correct medical devices to be used. Even when used normally, the reusable medical devices are subject to wear and mechanical stresses, but even more so if used too forcefully. In order to avoid risks in connection with the compatibility of the products, use only accessories and instruments specifically approved by MEDICON eG.

5 LIMITATION OF REUSABILITY

The service life of reusable medical devices normally comes to an end as a result of wear and damage from use and processing. Even when used normally, the reusable medical devices are subject to wear and

mechanical stresses, but even more so if used too forcefully. Careful inspection and functional check of the medical device before use is the best way to determine the end of service life of the medical device. In order to prevent failure or mechanical damage to the medical devices during surgery, they must be checked before each use by qualified personnel to ensure that they are mechanically intact, that there are no deformations, and that the parts are fully functioning. Evidence of damage and wear to a reusable medical device may include, but is not limited to, corrosion (i.e. rust, pitting), discolouration, excessive scratching, chipping, wear and cracking. Medical devices that are not functioning properly, medical devices with markings that are not readable by humans or machines, missing or removed (worn) part numbers, damaged and excessively worn parts must not be used and must be replaced, repaired or disposed of.

6 SERVICE

For service and repair, please contact your national MEDICON eG representative.

7 PROCESSING (CLEANING, DISINFECTION AND STERILISATION) OF MEDICAL DEVICES

For the detailed and product-specific processing (cleaning, disinfection and sterilisation) of products, please note the valid MPK in the section 'Special Instructions'. Here you will find one of the following processing instructions, please observe the specific processing instructions for your product:
MPK01-1 / MPK03-1 / MPK04-1 / MPK05-1 / MPK07-1 / MPK08-1
MPK11-1 / MPK12-1 / MPK13-1 / MPK14-1 / MPK16-1 / MPK20-1
MPK20-2 / MPK21-1 / MPK22-1 / MPK23-1 / MPK24-2

8 LIABILITY

CAUTION: According to US federal law, in the USA this product may be bought only by a physician or hospital or upon prescription. In the event of any discrepancies between the non-German and the German version of these operating instructions, only the German version shall be authoritative. Only the latest revision of the operating instructions applies. Due to constant technical development, the content of these operating instructions is updated regularly. Please ensure that you are using the current version. The version date of the respective edition of the operating instructions is included in the printout. MEDICON eG assumes no liability for damage caused by improper use, incorrect postoperative conduct, care or maintenance, or non-compliance with the restrictions on use and other guidelines in the operating instructions. MEDICON eG will accept liability neither for changes or repairs made to the product without the prior written authorisation of MEDICON eG, nor for repairs that were not carried out by workshops authorised by MEDICON eG or the MEDICON Repair Service (MRS).

9 DISPOSAL

To avoid a risk of infection to third parties, the medical devices must be cleaned and sterilised before disposal. In addition, the medical devices must be disposed of in specific and appropriately labelled containers, in order to protect third parties from cutting injury. Observe the applicable national laws when disposing of reusable medical devices!

10 DESCRIPTION OF SYMBOLS AND ICONS

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Production lot number, batch
	Item number
	Non-sterile
	Caution
	Observe the operating instructions

11 BASIC UDI-DI / DESCRIPTION / INTENDED PURPOSE / INDICATION / CONTRA-INDICATIONS

Basic UDI-DI	4046826;B;R001;01;aj;2U 4046826;B;R001;02;aj;37 4046826;B;R001;03;aj;3J 4046826;B;R001;04;aj;3V 4046826;B;R001;05;aj;48	4046826;B;R002;01;ah;3J 4046826;B;R002;02;ah;3V 4046826;B;R002;03;ah;48 4046826;B;R002;04;ah;4K 4046826;B;R002;05;ah;4W 4046826;B;R002;06;ah;59	4046826;B;R003;01;aj;48 4046826;B;R003;02;aj;4K 4046826;B;R003;03;aj;4W 4046826;B;R003;04;aj;59 4046826;B;R003;05;aj;5L 4046826;B;R003;06;aj;5X	4046826;B;R004;01;aj;4W 4046826;B;R004;02;aj;59 4046826;B;R004;03;aj;5L 4046826;B;R004;04;aj;5X 4046826;B;R004;05;aj;6A 4046826;B;R004;06;aj;6M	4046826;B;R005;01;ak;5L 4046826;B;R005;02;ak;5X 4046826;B;R005;03;ak;6A 4046826;B;R005;04;ak;6M
Description	All reusable medical devices for dissecting tissue and bone.	All reusable medical devices for holding tissue and bone.	All reusable medical devices for the abrasive removal of tissue and bone.	All reusable medical devices for the holding away, retracting, and dilating of tissue and bone.	All reusable medical devices for application.
Intended purpose	An instrument (scissors, knives, needles/awls, cutting forceps and trocars), not used in combination with an active medical device, intended for a surgical procedure, whose temporary use is the cutting of tissue and bone during surgical procedures. The medical devices may be used only by surgeons with sufficient experience regarding the purpose and function of the instruments during surgical procedures in hospitals and medical practices. They are used in patients regardless of age, weight, health condition or nationality. The surgical instrument can be reused after suitable cleaning and sterilisation procedures have been carried out.	An instrument (forceps, tweezers, clamps, clips, snares), not used in combination with an active medical device, intended for a surgical procedure, whose temporary use is to hold tissue and bone, in the case of needle holders to hold needles, during surgical procedures. The medical devices may be used only by surgeons with sufficient experience regarding the purpose and function of the instruments during surgical procedures in hospitals and medical practices. They are used in patients regardless of age, weight, health condition or nationality. The surgical instrument can be reused after suitable cleaning and sterilisation procedures have been carried out.	An instrument (chisels, osteotomes, spoons, curettes, files, raspatories, saws, drills, cutters, taps, punches, rongeurs), not used in combination with an active medical device, intended for a surgical procedure, whose temporary use is the abrasive removal of tissue and bone in surgical procedures. The medical devices may be used only by surgeons with sufficient experience regarding the purpose and function of the instruments during surgical procedures in hospitals and medical practices. They are used in patients regardless of age, weight, health condition or nationality. The surgical instrument can be reused after suitable cleaning and sterilisation procedures have been carried out.	An instrument (dilators, spreaders, gags, contractors, retractors, dissectors, elevators, probes) not used in combination with an active medical device, intended for surgical procedures, whose temporary use is to hold away, spread and dilate tissue and bone during surgical procedures. The medical devices may be used only by surgeons with sufficient experience regarding the purpose and function of the instruments during surgical procedures in hospitals and medical practices. They are used in patients regardless of age, weight, health condition or nationality. The surgical instrument can be reused after suitable cleaning and sterilisation procedures have been carried out.	An application instrument (forceps, templates, gauges, insertion devices, tampers, pestles and holders), not used in combination with an active medical device, intended for a surgical procedure, whose temporary use is to bend, hold, transfer shapes and sizes and insert implants and medical devices during surgical procedures. The medical devices may be used only by surgeons with sufficient experience regarding the purpose and function of the application instruments during surgical procedures in hospitals and medical practices. They are used in patients regardless of age, weight, health condition or nationality. The application instruments can be reused after suitable cleaning and sterilisation procedures have been carried out.
Indication	Use during surgical procedures to dissect or remove tissue.	Use during surgical procedures to hold tissue.	Use during surgical procedures for the abrasive removal of tissue and bone.	Use in surgical procedures to hold away, retract or dilate tissue and bone.	Use in surgical procedures to apply other instruments and implants.
Contra-indications	The medical devices are not intended for direct contact with the heart, the central nervous system (CNS) or the central circulatory system. Applications contrary to the intended purpose.	The medical devices are not intended for direct contact with the heart, the central nervous system (CNS) or the central circulatory system. Applications contrary to the intended purpose.	The medical devices are not intended for direct contact with the heart, the central nervous system (CNS) or the central circulatory system. Applications contrary to the intended purpose.	The medical devices are not intended for direct contact with the heart, the central nervous system (CNS) or the central circulatory system. Applications contrary to the intended purpose.	The medical devices are not intended for direct contact with the heart, the central nervous system (CNS) or the central circulatory system. Applications contrary to the intended purpose.

Rx only

CE 0297

medicon®
MEDICON eG · Gänsäcker 15
D-78532 Tuttlingen · Germany
Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50
Email: sales@medicon.de
Internet: www.medicon.de

English

Rev.: 2025-05-13

PROCESSING (CLEANING, DISINFECTION AND STERILISATION) OF MEDICAL DEVICES

Carefully read the operating instructions and processing instructions prior to clinical application and keep them safe and at hand.

TABLE OF CONTENTS

7 Processing (cleaning, disinfection and sterilisation) of products

7.1 General principles

7.2 Cleaning and disinfection

7.3 Machine cleaning/disinfection [washer-disinfector (WD)]

7.4 Manual cleaning and disinfection

7.5 Inspection

7.6 Maintenance

7.7 Packaging

7.8 Sterilisation

7.9 Storage

7.10 Material compatibility

7.11 Reusability

7.12 Safety information and warnings

7 PROCESSING (CLEANING, DISINFECTION AND STERILISATION) OF THE PRODUCTS

7.1 General principles

All products must be cleaned, disinfected and sterilised before every use; this applies in particular to the first use after delivery because all products are supplied non-sterile (cleaning and disinfection after removing the protective transport packaging; sterilisation after packaging). Effective cleaning and disinfection is an essential prerequisite for effective sterilisation. As part of your responsibility for the sterility of the products, please ensure during use:

- that as a matter of principle only validated methods that are sufficiently specific for equipment and product are used for cleaning / disinfection and sterilisation, that the devices used (washer/disinfector (WD), steriliser, etc.) are regularly maintained, checked and calibrated, and that the validated parameters are adhered to for each cycle.

Please also ensure during use that you collect soiled medical devices separately and do not place them back in the sterilisation tray to prevent heavier contamination of the loaded sterilisation tray.

Clean/disinfect the soiled medical devices, sort them onto the sterilisation tray and then sterilise the fully loaded sterilisation tray.

Please also follow the legal provisions that apply in your country as well as the hygiene regulations that apply in the medical practice or hospital.

Note:

The products are to be used by trained personnel only.

Processing may be carried out only by trained personnel in the central sterile supply department (CSSD) / processing unit for medical devices of the clinic or in the processing room of the medical practice. The clinic or medical practice is also responsible for the selection and use of the necessary protective equipment and hygiene measures. Consider deviating and/or additional specifications for some products in the table 'Special instructions'.

7.2 Cleaning and disinfection

Principles

For the cleaning and disinfection, an automated process [washer-disinfector (WD)] should be used where possible. A manual process, including one in which an ultrasonic bath is used (frequency range 35 kHz), should be used only if an automated process is not available or according to countryspecific requirements (e.g. in Germany an automated process is mandatory for critical B products) due to the considerably lower effectiveness and reproducibility. Pre-treatment must be carried out in both cases.

The concentrations, temperatures and application times specified by the manufacturer of the cleaning agent or detergent and disinfectant as well as the instructions for final rinsing must be strictly observed. Only use freshly prepared solutions, only use sterile or low-germ (maximum 10 germs/ml) and low-endotoxin (maximum 0.25 endotoxin units/ml) water (e.g. purified water / highly purified water) and only use a soft, clean and lint-free cloth (caution: care should be taken with products with rough surfaces, threads, sharp edges or similar, on which particles from the cloth can get stuck) and/or filtered air.

Pre-treatment

Immediately after use (within 2 h at most) coarse contamination must be removed from the products. If, due to the duration of the application or as a result of organisational aspects, this time cannot be observed, the user must determine and validate measures on their own responsibility in order to prevent the soiling from complete drying:

Process:

1. Disassemble the products as much as possible (see specific disassembly/assembly instructions).
2. Rinse the products for at least 1 min under running water (temperature < 35 °C/95 °F). Move movable parts back and forth at least three times during pre-rinsing. If applicable (see table 'Special instructions'): Rinse all lumens of the products at least three times (tools and minimum volume depend on the cavity to be rinsed).
3. Place the disassembled products in a sufficiently large pre-cleaning bath (in an ultrasonic bath that has not yet been activated) for the specified application time so that the products are completely covered. Ensure that the products do not touch during cleaning. Support the pre-cleaning by completely brushing all inner and outer surfaces (at the beginning of the application time, for tools see table 'Special instructions'). The brushes for the channels must be slightly larger than the inner diameter of the respective channel; the shaft length of the brush must be at least as long as the channel. Move all movable parts back and forth at least three times during the pre-cleaning. If applicable (see table 'Special instructions'): Rinse all lumens of the products at least three times at the end of the application time (tools and minimum volume depend on the cavity to be rinsed).
4. Activate the ultrasonic device for another minimum application time (but not less than 5 min).
5. Then remove the products from the pre-cleaning bath and rinse them thoroughly with water at least three times (at least 1 min). Move movable parts back and forth at least three times during the final rinse.

If applicable (see table 'Special instructions'):

Rinse all lumens of the products at least three times (tools and minimum volume depend on the cavity to be rinsed).

When selecting the cleaning agent to be used, 1 observe that

- it is generally suitable for cleaning invasive medical devices made of metals and plastics,
- the cleaning agent is suitable for ultrasonic cleaning (no foam formation),
- the cleaning agent is compatible with the products (see table "Material compatibility").

7.3 Machine cleaning/disinfection [washer-disinfector (WD)]

– When selecting the WD, ensure

- that the washer-disinfector has as a matter of principle proven effectiveness (e.g. DGHM or FDA approval/clearance/ registration or CE marking in accordance with DIN EN ISO 15883),
- that where possible, a tested programme for thermal disinfection (AO value > 3,000 or, for older devices, at least 5 min. at 90 °C/194 °F) is used (in the case of chemical disinfection, there is a risk of disinfectant residues on the products),
- that the programme used is suitable for the products and includes a sufficient number of rinsing cycles (at least three depleting steps after cleaning (or neutralisation, if used) or conductivity value control is recommended to effectively prevent detergent residues),
- that for the final rinse only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water) is used,
- that the air used for drying is filtered (oil-free, low in germs and particles) and that the washerdisinfector is regularly

maintained, checked and calibrated.

- When selecting the cleaning agent system, ensure that
- if thermal disinfection is not used, a suitable disinfectant with proven effectiveness (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is also used and that
- this is compatible with the cleaning agent used, and
- the chemicals used are compatible with the products (see section 'Material compatibility').

The concentrations, temperatures and application times specified by the manufacturer of the cleaning agent and, if applicable, disinfectant as well as the instructions for final rinsing must be strictly observed.

Process:

1. Disassemble the products as much as possible (see specific disassembly/assembly instructions).
2. Place the disassembled products in the washer-disinfector. Ensure that the products do not touch. If applicable (see table 'Special instructions'): Enable active flushing by establishing a connection to the flushing connection of the WD.
3. Start the program.

Step	Description	Medium	Temperature [°C]	Duration [min.]
1	Pre-rinse	Water	< 30	1
2	Emptying	-	-	-
3	Cleaning	Alkaline cleaning agent: Neodisher Med/Clean forte (Dr. Weigert GmbH & Co. KG, Hamburg) Concentration: 0.2 to 1 % (according to detergent manufacturer's specifications)	55	10
4	Emptying	-	-	-
5	Rinsing	Deionised water	< 30	1
6	Emptying	-	-	-
7	Disinfection (thermal)	-	95	5
8	Drying	Hot air	100	25

4. Disconnect the WD and remove the products at the end of the programme.

5. Check and pack the products as soon as possible after removal (see 'Inspection', 'Maintenance' and 'Packaging'); where applicable, after additional drying in a clean location. Proof of general suitability of the products for effective automated cleaning and disinfection was provided by an independent, officially accredited and recognised (§ 15 (5) MPG) test laboratory using a G 7836 CD washer-disinfector (thermal disinfection, Miele & Cie GmbH & Co., Gütersloh) and the Neodisher Med/Clean forte pre-cleaning and cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg). The procedure described above was taken into account for this.

7.4 Manual cleaning and disinfection

When selecting the cleaning agent and disinfectant, ensure that

- that these are generally suitable for cleaning or disinfecting medical devices made of metals and plastics,
- that the cleaning agent is suitable for ultrasonic cleaning (no foam formation),
- that a disinfectant with proven effectiveness (e.g. VAH/DGHM or FDA/EPA
- approval/clearance/registration or CE marking) is used and that this is compatible with the cleaning agent used, and
- that the chemicals used are compatible with the products (see section 'Material compatibility').

During manual cleaning, disinfection with a possible risk of injury and infection, further occupational health and safety measures (e.g. protective clothing, goggles, gloves; room air filtration) must be observed according to national regulations (e.g. TRBA 250 in Germany).

The concentrations, temperatures and application times specified by the manufacturer of the cleaning agent and disinfectant as well as the instructions for final rinsing must be strictly observed. Only use freshly prepared solutions, only use sterile or low-germ (maximum 10 germs/ml) and lowendotoxin (maximum 0.25 endotoxin units/ml) water (e.g. purified water / highly purified water) and only use a soft, clean and lint-free cloth (caution: care should be taken with products with rough surfaces, threads, sharp edges or similar, on which particles from the cloth can get stuck) and/or filtered air.

Process: Cleaning

1. Disassemble the products as much as possible (see specific disassembly/assembly instructions).
2. Place the disassembled products in a sufficiently large cleaning bath (in an ultrasonic bath that has not yet been activated) for the specified application time so that the products are completely covered. Ensure that the products do not touch during cleaning. Support cleaning by completely brushing all inner and outer surfaces with a soft brush. The brushes for the channels must be slightly larger than the inner diameter of the respective channel; the shaft length of the brush must be at least as long as the channel. Move all movable parts back and forth several times during cleaning. If applicable (see table 'Special instructions'): Rinse all lumens of the products at least five times at the beginning and end of the application time (tools and minimum volume depend on the cavity to be rinsed).
3. Activate the ultrasonic device for another minimum application time (but not less than 5 min).
4. Then remove the products from the cleaning bath and rinse them thoroughly with water at least three times (at least 1 min). Move movable parts back and forth several times during the final rinse. If applicable (see table 'Special instructions'): Rinse all lumens of the products at least five times (tools and minimum volume depend on the cavity to be rinsed).

5. Check the products (see 'Inspection' and 'Maintenance' chapters).

Disinfection

6. Place the disassembled, cleaned and inspected products in the disinfection bath for the specified application time so that the products are fully covered. Ensure that the products do not touch during cleaning. Move all movable parts back and forth several times during disinfection. If applicable (see table 'Special instructions'): Rinse all lumens of the products at least five times at the beginning and end of the application time (tools and minimum volume depend on the cavity to be rinsed).

7. Then remove the products from the disinfection bath and rinse them thoroughly with water at least five times (at least 1 min). Move movable parts back and forth several times during the final rinse. If applicable (see chapter 'Special instructions'): Rinse all lumens of the products at least five times (tools and minimum volume depend on the cavity to be rinsed).

8. Blow the products dry using filtered compressed air.

9. Pack the products as soon as possible after removal (see 'Packaging' chapter); where applicable, after additional drying in a clean location.

Proof of general suitability of the products for effective manual cleaning and disinfection was provided by an independent, officially accredited and recognised (§ 15 (5) MPG) test laboratory using Cidezyme/Enzol pre-cleaning and cleaning agent and CideX OPA disinfectant (Johnson & Johnson GmbH, Norderstedt). The procedure described above was taken into account for this.

7.5 Inspection

Inspect the products after cleaning or after cleaning/disinfection for corrosion, damaged surfaces, flaking, soiling and discolouration and separate any damaged products (for the restriction on the number of reprocessing cycles, see 'Re-usability' chapter) Products that are still soiled must be cleaned and disinfected again.

7.6 Maintenance

Reassemble the disassembled products (see specific disassembly/assembly instructions). Instrument oils or grease must not be used. Exception (only for special medical devices, see table 'Special instructions', not for implants):

In the case of lubricating joints, care should be taken to use only instrument oils (white oil, without further additives) which - taking into account the maximum sterilisation temperature applied - are approved for steam sterilisation and have a proven biocompatibility, and to apply only a small amount to the joints.

7.7 Packaging

Sort the cleaned and disinfected products into the corresponding sterilisation tray. Please pack the products or the sterilisation trays in sterilisation containers or very large products in single-use sterilisation packs (single or double packaging) that meet the following requirements (material/process):

- DIN EN ISO / ANSI AAMI ISO 11607 (for the USA: FDA clearance)
- Suitable for steam sterilisation (temperature resistance of at least 138 °C (280 °F) and adequate steam permeability)
- Sufficiently protect the products or sterilisation packages against mechanical damage
- Regular maintenance in accordance with the manufacturer's specifications (sterilisation container)
- A maximum weight of 10 kg per package/content of the sterilisation container must not be exceeded.

7.8 Sterilization

For the sterilisation, only the sterilisation methods listed below may be used; other sterilisation methods are not permitted. Steam sterilisation

- Fractionated vacuum method (with adequate product drying) Steam steriliser according to DIN EN 13060 / DIN EN 285 or ANSI AAMI ST79 (for the USA: FDA clearance)
- Validated in accordance with DIN EN ISO 17665 (current IQ/OQ (commissioning) and productspecific performance qualification (PO))
- Maximum sterilisation temperature 134 °C (273 °F plus tolerance according to DIN EN ISO 17665)
- Sterilisation time (exposure time at the sterilisation temperature):

Country	Fractional vacuum process	Gravitation process
Germany	at least 5 min. at 134 °C (273 °F)	not permitted ^d
USA at least	4 min. at 132 °C (270 °F), drying time at least 20 min.	not recommended ^d
France	at least 5 min. at 134 °C (273 °F)	not permitted ^d
Other countries	at least 5 min. at 132 °C (270 °F) / 134 °C (273 °F)	not permitted ^d

Proof of general suitability of the products for effective steam sterilisation was provided by an independent, officially accredited and recognised (§ 15 (5) MPG) test laboratory using the steam steriliser HST 6dx6 (Zirbus technology GmbH, Bad Grund) and using the fractionated vacuum process as well as the LAWTON MEDOL. Typical conditions in hospitals and medical practices as well as the method described above were taken into account for this.

The flash sterilisation method is not permitted as a matter of principle. Also do not use hot air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation or plasma sterilisation.

7.9 Storage

After sterilisation the products must be stored in the sterilisation packaging in a dry and dust-free environment.

7.10 Material compatibility

When selecting cleaning agents and disinfectants, ensure that these do not contain any of the following ingredients:

- organic, mineral or oxidising acids (minimum permissible pH value 5.5)
- alkaline solutions / strong alkaline solutions (neutral/enzymatic (max. permissible pH 8.5, mandatory for products made of aluminium or other alkali-sensitive materials, see table 'Special instructions') or alkaline cleaning agent (max. permissible pH value 11) recommended.
- organic solvents (e.g. alcohols, ethers, ketones, benzines)
- oxidising agents (e.g. hydrogen peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons

Never clean any products, sterilisation trays and sterilisation containers with metal brushes or steel wool. The products, sterilisation trays and sterilisation containers must not be exposed to temperatures above 138 °C (280 °F).

7.11 Reusability

The products can be reused with due care and provided they are not damaged or soiled; any further use or use of damaged and/or soiled products is the responsibility of the user. In case of nonobservance, any liability is excluded.

7.12 Safety information and warnings

- Ensure that all equipment involved in the processing procedure is qualified and maintained according to the manufacturer's specifications.
- Risk of injury and infection during transport packaging for the processing of sharp-edged or pointed medical devices, as well as in case of carelessness of the personnel during manual cleaning or disassembly processes, when loading or unloading the loading trolleys of a washerdisinfector.
- Repair work may be carried out only by the manufacturer. (Caution: Otherwise, the person who carries out the repair automatically becomes the responsible manufacturer)
- If unsuitable process parameters or media are used or if washer-disinfectors (including ultrasonic devices) are loaded inappropriately, residues from the processing procedures may remain on the medical device.
- When cleaning in the ultrasonic bath, the basic rules applicable for this purpose (no excessive loading, degassing, etc.) must be observed.
- The media used in the functional test must not cause any re-contamination or crosscontamination. Testing facilities are to be included in regular hygiene measures.
- Instrument surfaces optimised for ultrasound imaging (e.g. sound-reflecting surfaces on biopsy cannulas) may have their reflectance reduced by deposits, surface roughening or deformation and contrast less well in the ultrasound image.
- Electrical radiation: Systems compatible with magnetic resonance can produce artefacts in magnetic resonance imaging by changing their properties in the form of magnetisation or changes in internal resistances, or can damage tissue by inducing eddy currents.
- Suitable, hygienically specified and regularly tested media must be used for drying.

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Production lot number, batch
	MRT unsafe
	CE marking
	Medical device
	Prescription only
	Item number
	Non-sterile
	Caution
	Observe the operating instructions

Legend

1. If you use a cleaning and disinfecting agent for this purpose (e.g. for reasons of occupational safety), please take into account that it should be aldehyde-free (otherwise blood contamination would be fixed), have a proven effectiveness (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for the disinfection of the products and be compatible with the products (see 'Material compatibility' chapter). Please note that the disinfectant used during the pre-treatment serves only to protect personnel and cannot replace the later disinfection step that must be carried out after cleaning is completed.
2. If you consider a lower water quality to be sufficient on the basis of national recommendations (e.g. KRINKO/RKI/BfArM recommendation on processing in Germany), this is your sole responsibility.
3. At least three vacuum steps
4. The use of the gravitational method is not permitted in the European Union. The use of the less effective gravitational method outside the European Union is only permitted if the fractionated vacuum method is not available. It requires significantly longer sterilisation times and must be validated by the user specifically for the product, device, process and parameters.
5. The actual drying time required depends directly on parameters that are the sole responsibility of the user (load configuration and density, steriliser condition, ...) and must therefore be determined by the user. Nevertheless, the drying time should not fall below 20 minutes.

		Special instructions																
Worst case (code)		MPK01-1	MPK03-1	MPK04-1	MPK05-1	MPK07-1	MPK08-1	MPK11-1	MPK12-1	MPK13-1	MPK14-1	MPK16-1	MPK20-1	MPK20-2	MPK21-1	MPK22-1	MPK23-1	MPK24-1
Geometric features		Products with longer/ lower ring lumen (standard) Disassembly for cleaning/disinfection possible, direct connection not possible	Products with longer/lower ring lumen (standard), direct connection not possible	Products with longer/lower ring lumen (standard) Disassembly for cleaning/disinfect on possible, direct connection not possible	Products with longer/ lower ring lumen (standard, not flushable)	Products with longer/ narrower ring lumen With flushing connection (LuerLock), with protective cap, no disassembly possible	Tube, shaft Products with LuerLock, no disassembly possible	Tube, shaft Products with LuerLock, no disassembly possible	Products with longer/narrower single lumen Luer/LuerLock	Products with very short/ wide ring lumen without LuerLock Disassembly for cleaning/disinfection possible, direct connection not possible	Segmented products with longer/narrow ring-shaped cannulas Disassembly for cleaning/disinfection possible, direct connection not possible	Products with longer/ lower single lumen, with flushing hose olive	Products with very narrow cavities without LuerLock, minimally open	Segmented products with longer/extremely narrow ring-shaped cannulas and complex mechanism on the inside Disassembly for cleaning/disinfection possible, direct connection not possible	Self-opening joint instruments with sliding surfaces, can be disassembled	Products with joint, can be disassembled	Products with joint / thread	Self-closing joint instruments
Flushing volume		50 ml (disposable syringe) with attached cannula (for backflushing the blind lumen)	50 ml (disposable syringe), flushing gun	5 ml (disposable syringe)	n/a	10 ml (disposable syringe)	10 ml (disposable syringe)	5 ml (disposable syringe)	50 ml (disposable syringe) / flushing gun	n/a	50 ml (disposable syringe) / flushing gun	10 ml (disposable syringe)	50 ml (disposable syringe) / flushing gun	n/a	n/a	n/a	n/a	n/a
Brush		Standard brushes	Standard brushes Long brush (length >320 mm, diameter approx. 6 mm)	Standard brushes Long brush (length >320 mm, diameter approx. 4 mm)	Standard brushes	Standard brushes	Standard brushes	Standard brushes	Standard brushes Long brush (length >510 mm; diameter approx. 4 mm)	Standard brushes	Standard brushes Flexible long brush (length >700 mm, diameter approx. 4 mm)	Standard brushes	Standard brushes Long brush (length >510 mm, diameter approx. 4 and 5 mm)	Standard brushes	Standard brushes Long brush (length >400 mm, diameter approx. 12 mm)	Standard brushes	Standard brushes	Standard brushes
Special/additional procedure for	Disassembled, flush five times on the inside, brush on the outside and inside	Standard brushes	Disassemble, open and close five times, insert in open position for ultrasound treatment	Disassemble, flush five times on the inside, on the outside, move the wire eyelet back and forth five times.	Brush on the outside, open and close five times, insert in open position for ultrasound treatment	Open and close five times, insert in open position for ultrasound treatment	Brush on the outside, flush on the inside at least 5 times, move the joint five times during soaking and flushing	Flush five times on the inside, connection via LuerLock, brush on the outside	Disassemble, flush five times on the inside, brush on the inside and outside	Disassemble brush on the inside and outside flush on the inside and outside at least 5 times	Brush on the inside and outside flush on the inside at least 5 times	Flush five times on the inside, brush on the outside (do not brush on the inside!)	Disassemble, brush on the inside and outside, flush on the inside and outside at least 5 times	Disassemble, brush on the inside and outside, flush on the inside and outside at least 5 times	Disassembled, flush handle tubes with running water, open and close joint five times, for ultrasound treatment insert joint in open position	Brush outside and in the gaps, open and close five times, for ultrasound treatment insert joint in open position	Brush on the inside and outside, open and close five times, insert joint in open position for ultrasonic treatment (if no lock provided, keep joint open by rubber band/stainless steel wire loop around handles).	Disassembled, brush on the inside and outside, move crank at least five times during soaking and flushing
	Disassembled, flush five times on the inside, brush on the outside and inside	Disassembled, flush five times on the inside, brush on the outside and inside	Disassembled, open and close five times, insert in open position for ultrasound treatment	Flush five times on the inside, on the outside, move the wire eyelet back and forth five times.	Open and close five times, insert in open position for ultrasound treatment	Open and close five times, insert in open position for ultrasound treatment	Brush on the outside, flush on the inside at least 5 times, move the joint five times during soaking and flushing	Flush five times on the inside, connection via LuerLock, brush on the outside	Disassembled, flush five times on the inside, brush on the outside	Disassembled brush on the inside and outside flush on the inside and outside at least 5 times	Brush on the inside and outside flush on the inside at least 5 times	Flush five times on the inside, brush on the outside	Disassembled, brush on the inside and outside, flush on the inside and outside at least 5 times	Brush on the inside and outside, open and close five times, insert in disassembled state during treatment in the ultrasonic bath	Disassembled, joint, open and close five times, insert joint in open position for ultrasound treatment	Brush outside and in the gaps, open and close five times, for ultrasound treatment insert joint in open position	Brush on the inside and outside, open and close five times, insert joint in open position for ultrasonic treatment (if no lock provided, keep joint open by rubber band/ stainless steel wire loop around handles).	Disassembled, brush on the inside and outside, move crank at least five times during soaking and flushing
	Disassembled, attach outer tube to flushing mandrel, place small parts in small parts basket, standard basket for other parts	Disassembled, flush five times on the inside, brush on the outside and inside	Standard, disassembled	Standard, disassembled	In standard basket	Connection to flushing connection	Connection to flushing connection, mouth part in open position	Connection via LuerLock	Disassembled, standard basket, sleeve of coupling in small parts basket	Disassembled small parts basket	Flushing hose, attach to olive	Attach flushing hose (hose connection)	Disassembled: Use the flushing lance for cannulated parts, standard basket for other parts	Standard basket, place punch in disassembled state on the flushing basket.	Standard basket, disassembled, insert joint in open position	Standard basket, insert joint in open position	Standard basket, insert joint in open position (if no lock provided, hold joint open by rubber band/stainless steel wire loop around handles)	Standard basket, disassembled
Packaging		Attach loosely, lubricate inner shaft and thread	Reassemble, no lubrication permitted	Re-attach (flushing wire only, does not apply to protective cap for flushing connection), no lubrication permitted	Loop shifted backwards no lubrication permitted	Standard procedure: Lubricate joint, open cap on flushing hose, attach protective cap	Lubricate joints, close open Luer-Lock mouth part	Reassemble, no lubrication permitted	Reassemble, lubricate threads and inner shafts	Almost closed but still assembled loosely, lubricate threads	No lubrication permitted	No lubrication permitted	Reassemble, lubricate all movable parts	Lubricate the joints and sliding surfaces according to GA572800.	Install (handle tubes) lubricate joints	Not completely closed, lubricate the joints and the thread of the spindle	Lubricate the joints in slightly open position	Lubrication not permitted
Sterilisation		Assembled loosely, lubricated	Assembled	Assembled (flushing wire only, does not apply to protective cap for flushing connection)	Standard	Standard procedure: Protective cap attached, open cap on flushing hose, mouth parts closed	Lubricated, protective cap closed, mouth part closed	Standard	Assembled, lubricated	Almost closed but still assembled loosely, lubricated	Standard	Standard	Assembled (does not apply to loading and conical part), lubricated	Sterilisation in assembled state	Assembled lubricated	Not fully closed	Joints lubricated in slightly open position	Assembled
Permitted maximum number of cycles		As long as the functional check has been carried out by the user without problems and the legibility of the labelling including the coding is ensured, the medical device can be processed in accordance with these validated processing instructions. A limitation of the number of processing cycles is not necessary due to the investigations carried out by us and their results.																
Classification recommendation according to KRINKO/RKI/ BfArM recommendation (only Germany, if used as intended)		Critical B																

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