MEDICON eG GERMANY · 47° 57' 31.4" N · 8° 46' 58.6" E

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mini mized

Minimal-invasive Zugangstechniken
Minimally Invasive Approaches

medicon[®]


Not all products illustrated in this documentation are available for sale in all countries.
Not for all declared products in this documentation Medicon is the legal manufacturer.
Please contact your local distributor for more information concerning availability and
regulatory information (e.g. classification) of the mentioned products.

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
The anterior cervical approach is a well accepted procedure in the midst of the spine surgery community. Perhaps more than in other procedures, in anterior cervical surgery, a painstaking attention to the detail is one of the key factors to the therapeutic success.

This is demonstrated by the broad range of incidence of complications in the different series, in contrast to the quite homogenous incidence of recurrent herniation in lumbar disc surgery. [1]

The "cervicalino" spinal spreading system is suited for anterior cervical approach techniques.

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


[1] Cheung KM, Mak K, Luk KD (2011) Anterior Approach to Cervical Spine. Spine (Oct 24)

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

- ▶ Miniaturized spreading system with "side load" feature for minimal invasive approach to cervical spine
- ▶ Miniaturized model improves freedom of movement for the surgeon in the surgical area
- ▶ Possible attachment to a holding arm allows for a secure and stable fixation
- ▶ Simplified handling when positioning the valves and the retractor

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- ▶ Miniaturized spreading system for minimal invasive approach to cervical spine
- ▶ Miniaturized model improves freedom of movement for the surgeon in the surgical area
- ▶ Possible attachment to a holding arm allows for a secure and stable fixation
- ▶ Artifact-free visualisation by PEEK-blades

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- ▶ Artifact-free visualisation by PEEK-blades and blade supports
- ▶ Miniaturized spreading system for minimal invasive approach to cervical spine
- ▶ Miniaturized model improves freedom of movement for the surgeon in the surgical area

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Approach techniques in spinal surgery are focusing more and more on the postoperative importance of preserving the innervation of the paravertebral muscles as much as possible. [1,2]

This concept is realized by the percutaneous insertion of pedicle screws and the use of tubular retractors in lumbar and posterior cervical discectomy. The self-retaining speculum-counter-retractor developed by Wolfhard Caspar in the late seventies has proven itself, with minor modifications, over the last decades. [3]

Today, the length of the skin incision and the extent of paravertebral muscle dissection are reduced due to the dimensions of the endoscopic working channels (6 mm) and the use of tubular retractors (11-15 mm). [4]


This decrease of muscle tissue traumatization leads to less need for postoperative analgesics. [5,6]

The minimally invasive transmuscular approaches, particularly the endoscopic variants, necessitate a learning curve and an expenditure of resources which have hampered the widespread use of these techniques.

- [1] Taylor, McGregor, Medhizadeh et al (2002) - The impact of self-retaining retractors on the para-spinal muscles during posterior spinal surgery. Spine; 27:2758-62
- [2] Weber, Grob, Dvorak et al (1997) - Posterior surgical approach to the lumbar spine and its effect on the multifidus muscle. Spine; 22:1765-72
- [3] Caspar (1977) - A new surgical procedure for lumbar disc herniation causing less tissue damage through a microsurgical approach. In: Advances in Neurosurgery; Wüllenweber, Brock, Hamer (eds); Vol 4: 74-77; Springer-Verlag, Berlin
- [4] Rütten, Komp, Godolias (2006) - A new full-endoscopic technique for the interlaminar operation of lumbar disc herniations using 6-mm endoscopes: prospective 2-year results of 331 patients. Minim Invasive Neurosurg 49(2): 80-7
- [5] Greiner-Perth, Boehm, Ezzati et al (2006) - Comparison of soft tissue trauma in microsurgical nucleotomy to a new technique: Microscopically Assisted Percutaneous Nucleotomy. An MRI-Study. Eur Spine J 15 (Suppl.4): 475
- [6] Brock, Kunkel, Papavero (2008) - Lumbar microdiscectomy: subperiosteal versus transmuscular approach and influence on the postoperative analgesic consumption. Eur Spine J, Vol 17,4: 518-522

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The "piccolino" spinal spreading system is well suited for miniaturized subperiosteal paravertebral as well as transmuscular approaches. When compared with conventional techniques, the reduced dimensions of the retractor do not limit access to the surgical field since the speculum flares out distally.

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







Aluminium-Valven
Aluminium blades

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The Interlaminar Approach

Positioning: According to surgeon's preference (Wilson frame; knee-elbow position).

X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the disc herniation (DH) and projecting on the intervertebral space (IVS). If necessary, an additional marker can be placed at the level of the extruded disc fragment.


Skin incision: Under the microscope, approx. 18 mm long, 5 mm off the midline, centred on the IVS or DH.

Fascia incision: Semilunar incision, pedicled at the midline, a few millimeters longer than the skin incision, 3 stay sutures.

Muscle dissection: Transect the tendinous insertion of the erector spinae muscle, if necessary. Whenever possible the subperiosteal dissection should be blunt. Retract the muscles laterally with a Langenbeck retractor up to the medial part of the joint complex.

Insertion of the "piccolino" speculum: Insert a speculum of sufficient length flush with the skin incision. The length of the lateral PEEK blade is determined by trial insertion, then fixed to the counter-retractor and placed in its definite position.

Fluoroscopy: The radiolucent PEEK blades ease the fluoroscopic check of the correct level.

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Interlaminar fenestration: Drill off or punch out the lower border of the upper hemilamina and/or the medial part of the joint. Narrow burr-handpieces facilitate visual control. The yellow ligament is opened. The disc fragment is removed, and, if necessary, discectomy is performed by standard microsurgical technique. The nerve root retractor can be replaced by the flattened sucker in right or left version (3 sizes).

Decompression of the spinal stenosis in cross-over technique: Tilt the operating table about 10° away from you. After drilling off the base of the cranial spinous process and a few millimeters of the inferior rim of the cranial hemilamina the seam between the two yellow ligaments is exposed.

Resect the ipsilateral yellow ligament first. Then drill off the medial aspect of the normally hypertrophic joint complex. At this point the table is tilted further (e.g., 30°) away from the surgeon. Resect the deep part of the interspinous ligament. Resect also the contralateral yellow ligament. Drill off the contralateral medial joint facet, if needed. Protect the dura with the flattened suction cannulas tip to the left and right respectively.

Wound closure: Remove the counter-retractor first and then the "piccolino" speculum. Perform the hemostasis carefully. The access tunnel is filled with warm Ringer solution. 2-3 fascial sutures. 2 subcutaneous sutures. Intracutaneous suture or adhesive strips.

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The Translaminar Approach

Positioning: According to surgeon's preference (Wilson frame; knee-elbow position). The lumbar hemilamina (HL) dips downwards obliquely in the caudo-cranial direction. By tilting the operating table about 15° head upwards, the HL becomes horizontal. That facilitates drilling of the translaminar hole (TH).


X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the disc herniation (DH) and projecting on the upper rim of the intervertebral space (IVS). The second marker is placed at the level of the extruded disc fragment, usually between the lower border of the cranial pedicle and the upper rim of the IVS.

Skin incision: Under the microscope, approx. 18 mm long, 5 mm off the midline, centered on the DH.

Fascia incision: Semilunar incision, pedicled at the midline, a few millimeters longer than the skin incision, 3 stay sutures.

Muscle dissection: The tendinous insertion of the erector spinae muscle can be transected, if necessary. Whenever possible the subperiosteal dissection should be blunt. Retract the muscles laterally with a Langenbeck retractor up to the medial aspect of the joint complex.

Insertion of the "piccolino" speculum: Insert a speculum of sufficient length flush with the skin incision. The length of the lateral PEEK blade is determined by trial insertion, then fixed to the counter-retractor and placed in its definite position. Its lateral and inferior border the posterior arch and its median junction with the spinous process can be seen.

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Fluoroscopy: Insert a dissector into the open "piccolino" and pack with a drawn out moist pad. The radiolucent PEEK blade helps when checking with lateral fluoroscopy. The dissector should point to the disc fragment dislocated superiorly.

Translaminar hole: The second skin marking, between the lower border of the cranial pedicle and the superior rim of the IVS, can be marked on the HL using a small diamond burr.

Centered on this, a TH about 10mm in diameter is drilled. Anatomical remarks: (a) the HLs become narrower in the caudocranial direction, i.e. the TH turns from round to oval; the IVS "ascends" increasingly under the HL and its superior rim may be inspected through the TH. (b) Keep a distance of at least 3 mm between the lateral border of the TH and the lateral rim of the HL in order to avoid a fracture of the pars interarticularis. (c) When drilling through the HL, there are 3 successive layers from the surface down: white (lamina externa), red (spongy bone), and white again (lamina interna).

The latter should be drilled with a diamond burr as there is no protective yellow ligament beneath the lamina interna. The TL should ideally be centered on the lateral border of the dural sac and the axilla of the exiting root. However, the root is often displaced cranial by disc fragments. The DH is removed, and, if necessary, the IVS is cleared using standard microsurgical technique. Due to the small surgical field, the nerve root retractor can be replaced by the flattened suction cannulas in right and left versions (3 sizes).

Wound closure: Remove the counter-retractor first and then the "piccolino" speculum. Any minor bleeding from muscle is coagulated. The access tunnel is filled with warm Ringer solution. 2-3 fascial sutures. 2 subcutaneous sutures. Intracutaneous suture or adhesive strips.

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The extraforaminal approach

Positioning: According to surgeon's preference (Wilson frame; knee-elbow position).


X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the disc herniation (DH) and projecting on the intervertebral space (IVS). AP-view: A rectangle is formed by 2 horizontal lines: (1) upper rim of the IVS (2) line projected on the middle third of the cranial transverse processes and by 2 vertical lines: (3) paramedian to the spinous process (4) lateral parapedicular line. The lateral line of the rectangle corresponds to the skin incision.

Skin incision: Under the microscope, approx. 25 mm long and approx. 40 mm off the midline. The surgeon's microscope oculars are tilted 5°-10° towards the midline.

Fascia incision: Straight, a few millimetres longer than the skin incision, 2-4 stay sutures.

Muscle dilation: Blunt muscle-splitting with the forefinger until the medial third of the transverse processes is reached. Expose the target field with small Langenbeck retractors.

Insertion of the "piccolino" speculum: Insert a speculum of appropriate length until both transverse processes are exposed. Measure the length of the lateral PEEK blade (skin → transverse processes) and of the medial PEEK blade (skin → facet joint). Usually the medial blade is 15 mm shorter than the lateral blade. If needed, bolster the handle of the "piccolino" speculum with pads in order to achieve additional medial convergence.

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Fluoroscopy: The radiolucent PEEK blades enable the fluoroscopic check of the correct level. After this check, tilt the operating table about 20° away from you in order to obtain a better view of the extraforaminal compartment.

Extraforaminal fenestration: Depending on the segment, particularly at L5/S1, partial burring of the lateral facet joint is indicated. Bipolar coagulation and incision of the medial insertion of the intertransverse muscle. Identification of the DH in the extraforaminal fat of the nerve. The latter is usually displaced by the DH in a cranial and lateral direction as well as towards the surface. The extruded disc fragment is removed. If necessary, the IVS is partially cleared using standard microsurgical technique. The nerve root retractor can be replaced by the flattened suction cannulas in right and left versions (3 sizes).

Wound closure: First remove the PEEK blades and then the "piccolino" speculum. Hemostasis is performed carefully. The access tunnel is filled with warm Ringer solution. 2-3 fascial sutures. 2 subcutaneous sutures. Intracutaneous suture or adhesive strips.

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The dorsolateral foraminotomy

Positioning: According to surgeon's preference: sitting, prone or so-called Concorde position (our preference).





X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the disc herniation (DH) and projecting on the intervertebral space (IVS). At the lower levels (C6-T1), the arms may be pulled with adhesive strapping or weights.

Skin incision: Under the microscope, approx. 20 mm long, 5 mm off the midline, centered on the IVS.

Fascia incision: Semilunar incision, pedicled at the midline, a few millimeters longer than the skin incision, 3 stay sutures.

Muscle dissection: Dissect the various muscle layers and fascias off the spinous process by blunt and sharp dissection. Retract the muscles laterally with a narrow Langenbeck retractor up to the laminofacet junction. Quite typically this target area is characterized by the yellow ligament tapering off laterally in V-shaped fashion.

Insertion of the "piccolino" speculum: Insert a speculum of sufficient length flush with the skin incision. Select an appropriate PEEK blade and insert the counter-retractor.

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Interlaminar fenestration: Drill off and/or punch out the lower rim of the superior hemilamina and/or the medial part of the joint complex (key-hole). Slim burr-handpieces facilitate visual control of the target structures. The yellow ligament is opened, the axilla of the nerve root is located, and the disc fragments are removed by standard microsurgical technique. The nerve root retractor can be replaced by the flattened suction cannulas in right and left versions (3 sizes).

Wound closure: Remove the counter-retractor first and then the speculum. Careful hemostasis is performed. The access tunnel is filled with warm Ringer solution. 2-3 fascial sutures. 2 subcutaneous sutures. Intracutaneous suture or adhesive strips.

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At the end of the 1970s, independently of one another, Williams, Caspar, and Yasargil introduced micro-surgery for the operative treatment of lumbar disc herniation. The advantages of minimal tissue trauma combined with early postoperative mobilization of the patients helped this method to become popular. As in conventional surgery, microdiscectomy requires dissection of the paravertebral muscles from the spinous process and hemilamina. Retraction due to scarring as well as denervation of the muscles, particularly rich in proprioceptors in the cervical region, may result in troublesome local syndromes. [1-3]

Endoscopic removal of disc fragments is performed through a cylindrical working channel (approx. Ø 7 mm) introduced by muscle-splitting. The technique, developed by Destandeau and Kambin and refined by Foley and Rütten, is very atraumatic to the muscles. However, its disadvantages are the two-dimensional view of the surgical field, the difficulty of using a burr and impairment of the optics by bleeding. [4-7]


- [1] Williams (1978) Microlumbar discectomy. A conservative surgical approach to the virgin herniated disc. Spine, 3(2): 175-182
- [2] Caspar (1977) A new surgical procedure for lumbar disc herniation causing less tissue damage through a microsurgical approach. In: Advances in Neurosurgery; Wüllenweber, Brock, Hamer (eds); Vol 4: 74-77; Springer-Verlag, Berlin
- [3] Yasargil (1977) Microsurgical operation of the herniated lumbar disc. In: Advances in Neurosurgery; Wüllenweber, Brock, Hamer (eds); Vol 4: 81-88; Springer-Verlag, Berlin
- [4] Destandeau (1999) A special device for endoscopic surgery of lumbar disc herniation. Neurol Res 21: 39-41
- [5] Kambin, Zhou (1997) Arthroscopic discectomy of the lumbar spine. Clin Orthop Relat Res 337: 49-57
- [6] Foley, Smith, Rampersaud (1999) Microendoscopic approach to farlateral lumbar disc herniation. Neurosurg Focus 7(5): e5
- [7] Rütten, Komp, Godolias (2006) A new full-endoscopic technique for the interlaminar operation of lumbar disc herniations using 6-mm endoscopes: prospective 2-year results of 331 patients. Minim Invasive Neurosurg 49(2): 80-7

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The Microdisc-XS technique combines the less traumatic "endoscopic" muscle-splitting approach with the proven "microsurgical" advantages in the surgical field.

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Skin incision: Under the microscope, approx. 18 mm long, 15 mm from the midline, centred on the IVS or DH. The surgeon's oculars are tilted 5°-10° towards the midline.

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Muscle splitting: With the index finger until contact is made with the lamino-facet junction (target area). Alternatively: sequential dilators. The target area is exposed with narrow Langenbeck retractors. A few insertions of the multifidus muscle are coagulated and divided.

Insertion of the spreading sleeve (SS): Replace the Langenbeck retractors with Love nerve root retractors, and advance a SS of suitable length with alternating rotating movements until contact is made with bone. Alternatively: after removing the Langenbeck retractors, advance the SS through a dilator as far as the target area.

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Interlaminar fenestration: Burr or punch out the lower border of the upper hemilamina and/or the medial part of the joint. Narrow burr-handpieces facilitate visual control. The yellow ligament is opened. The disc fragment is removed, and, if necessary, discectomy is performed by standard microsurgical technique. The nerve root retractor can be replaced by the flattened sucker in right or left version (3 sizes).

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

X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the DH and projecting on the upper rim of the IVS. The second marker is placed at the level of the extruded disc fragment, usually between the lower border of the cranial pedicle and the upper rim of the IVS.

Skin incision: Under the microscope, approx. 18 mm long, 10 mm from the midline, centered on the DH. The surgeon's oculars are tilted 5°-10° towards the midline.

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Muscle dilatation: With the index finger until contact is made with the HL and its inferior border. Alternatively: sequential dilators. The HL is exposed with narrow Langenbeck retractors. Its lateral and inferior border as well its junction with the spinous process should be identified.

Insertion of the spreading sleeve (SS): Replace the Langenbeck retractors with Love nerve root retractors and advance a SS of suitable length with alternating rotating movements until contact is made with the HL. Alternatively: after removing the Langenbeck retractors, the SS is advanced through a dilator as far as the HL.

Fixing the SS: The SS docks the contralaterally fixed snake arm and is locked in the desired position. (1) Spread the SS. (2) Introduce a dissector into the SS and tamponade with a moist unfolded sponge. (3) By means of lateral fluoroscopy adjust if necessary position and/or angle of the SS. The DH should project into the geometric centre of the oval opening of the SS.

Translaminar hole: The second skin marking, between the lower border of the cranial pedicle and the superior rim of the IVS, can be transferred to the HL using a small diamond burr. Centered on this, a TH about 10mm in diameter is drilled. Anatomical remarks: (a) The HLs become narrower in the caudocephalad direction, i.e. the TH turns from round to ovale; the IVS "ascends" increasingly under the HL and its superior rim may be inspected through the TH. (b) Keep a distance of at least 3 mm between the lateral border of the TH and the lateral rim of the HL in order to avoid a fracture of the pars interarticularis. (c) When drilling through the HL, there are 3 successive layers from the surface down: white (lamina externa), red (spongy bone), and white again (lamina interna). The latter should be drilled with a diamond burr as there is no protective yellow ligament beneath the lamina interna. The TL should ideally be centered on the lateral border of the dural sac and the axilla of the exiting root. However, the root is often displaced cephalad by disc fragments. The DH is removed, and, if necessary, the IVS is cleared using standard microsurgical technique. Due to the small surgical field, the nerve root retractor can be replaced by the flattened sucker in right and left versions (3 sizes).

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X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the DH and projecting on the IVS. AP-view: A rectangle is formed by 2 horizontal lines: (1) upper rim of the IVS (2) line projected on the middle third of the cranial transverse processes and by 2 vertical lines: (3) paramedian to the spinous process (4) lateral parapedicular line. The lateral line of the rectangle corresponds to the skin incision.

Skin incision: Under the microscope, approx. 25 mm long and approx. 40 mm from the midline. The surgeon's oculars are tilted 5°-10° towards the midline.

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Muscle splitting: With the index finger until contact is made with the junction of the medial third of the transverse process and the lateral joint facet. Alternatively: sequential dilators. The target area is exposed with narrow Langenbeck retractors.

Insertion of the spreading sleeve (SS): Advance a SS of suitable length between the Langenbeck retractors with alternating rotating movements until contact is made with the corner between the lateral joint facet (medially) and the superior transverse process. Alternatively: after removing the Langenbeck retractors, the SS is advanced through a dilator as far as the target area.

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Extraforaminal fenestration: Depending on the segment, particularly at L5/S1, partial burring of the lateral facet joint is indicated. Bipolar coagulation and incision of the medial insertion of the intertransverse muscle. Identification of the DH in the extraforaminal fat of the nerve. The latter is usually displaced by the DH in a cephalad and lateral direction as well as towards the surface. The extruded disc fragment is removed. If necessary, the IVS is partially cleared using standard microsurgical technique. The nerve root retractor can be replaced by the flattened sucker in right and left versions (3 sizes).

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
X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the DH and projecting on the IVS. At the lower levels (C6/T1), the arms may be pulled with adhesive strapping or weights.

Skin incision: Under the microscope, approx. 18 mm long, 15 mm from the midline, centered on the lamino-facet junction. The surgeon's oculars are tilted 5°-10° towards the midline.

Fascia incision: Straight, a few millimetres longer than the skin incision, 2 retaining sutures.

Muscle splitting: The fascia and the single muscle layers are split with scissors or with a dissector according to their fibre direction. The lamino-facet junction is aimed at between two narrow Langenbecks. In this target area, the yellow ligament typically runs laterally in a V shape.

Insertion of the spreading sleeve (SS): Replace the Langenbeck retractors with Love nerve root retractors and advance a SS of suitable length with alternating rotating movements until contact is made with bone. Alternatively: after removing the Langenbeck retractors, advance the SS through a dilator as far as the target area. With both techniques, the SS is centered on the V-shaped yellow ligament.

Fixing the SS: The SS docks the contralaterally fixed snake arm and is locked in the desired position. (1) Spread the SS. (2) Introduce a dissector into the SS and tampon with a moist unfolded sponge. (3) By means of lateral fluoroscopy adjust if necessary position and/or angle of the SS. The DH should project into the geometric centre of the oval opening of the SS.

Interlaminar fenestration: Burr and/or punch out the lower rim of the superior hemilamina and/or the medial part of the joint complex (key-hole). Slim burr-handpieces facilitate visual control of the target structures. The yellow ligament is opened, the axilla of the nerve root is located, and the disc fragments are removed by standard microsurgical technique. The nerve root retractor can be replaced by the flattened sucker in right and left versions (3 sizes).

Wound closure: The SS is closed and withdrawn slowly. Any minor bleeding from muscle is coagulated. The access tunnel is filled with warm Ringer solution. 2-3 fascial sutures. 2 subcutaneous sutures. Intracutaneous suture or adhesive strips.

MINIMIZED^{TO} MAXIMIZE

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The image shows a surgical instrument used for minimally invasive spine surgery. It features a long, thin, flexible shaft with a series of small, rounded protrusions along its length. The shaft is connected to a handle with a textured, knurled grip. The handle is black and has a circular opening in the center. The instrument is shown against a background of a human spine, with a specific vertebra highlighted in a light blue color.

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With an incidence rate of 87 per 100 000 inhabitants per year in Germany, lumbar intervertebral disc surgery is one of the most common procedures in the field of Neurosurgery and Orthopaedics. This high frequency rate requires procedures that, on the one hand, afford extensive, qualitative safety during surgery and, on the other hand, keep access-related trauma to a minimum, thus promoting rapid patient rehabilitation. Healthcare costs can be reduced by shortening the length of time spent in hospital and through the availability of an alternative, minimally invasive intervertebral disk surgical technique performed on an outpatient basis. [1]

The previous open intervertebral disc procedure was replaced in the 1970s by the microsurgical technique developed by Yasargil and Caspar, and which nowadays is considered as the "Gold standard". The advantage of this method compared with the "open" technique is the smaller access required and the excellent visualisation of the surgical field using a surgical microscope. Another important step towards reducing access trauma was achieved with the development of the endoscopically supported method by Destandeau and Foley for the surgical correction of lumbar herniated vertebral disks. The disadvantage of this method mainly lies in the two-dimensional depiction of the surgical field and the troublesome contamination of the optics due to bleeding. [2-5]

- [1] Kast, Antoniadis, Richter - Epidemiologie von Bandscheibenoperationen in der BRD. Zentralbl Neurochir 61 (2000): 22-25
- [2] Yasargil - Microsurgical operation of the herniated lumbar disc. In: Wüllerweber, Brock, Hamer (eds.) Advances in Neurosurgery, Springer-Verlag, Berlin 1977
- [3] Caspar - A new surgical procedure for lumbar disc herniation causing less tissue damage through a microsurgical approach. In: Wüllerweber, Brock, Hamer (eds.) Advances in Neurosurgery, Springer-Verlag, Berlin 1977
- [4] Destandeau - A special device for endoscopic surgery of lumbar disc herniation. Neurol. Res 21 (1999): 39-41
- [5] Foley, Smith - Microendoscopic discectomy. Tech Neurosurg 3 (1997): 301-307

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This technique, referred to as microscopically assisted percutaneous nucleotomy allows optimal, three-dimensional visualisation of the surgical field through the surgical microscope. Access to the spinal canal is provided by a working channel, introduced through transmuscular dilatation.

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
Pictures 1-5: The patient leans on his/her knees and elbows, thus ensuring that the lumbar spine is not curved.

A skin incision 15 mm long and approximately 2 cm from the midline is made at the level of the marker needle on the affected side.

The level is marked with the needle on the uninvolved side, monitoring the lateral C-arm image. The needle should be directed towards the lesion (e.g. cranially or caudally extruded sequestrum).

Insertion of the first dilator and confirmation by lateral C-arm image. This prevents intervention at the wrong level. The target point is the Ligamentum flavum in the interlaminar window.

Medial convergence is approximately 10°. The lower edge of the upper Lamina medial from the joint mass is a good anatomical landmark that is palpable with the dilator.

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
Pictures 6-10: Insertion of the second dilator with image intensifier control and removal of the first. Once again, the lower edge of the upper Lamina should be palpated as a safety precaution.

The external screw allows the working channel to be "gently screwed in" through the paraspinal musculature above the second dilator. Image intensifier control should then be repeated as a precautionary measure.

The handle is placed on the working channel, thus facilitating straightforward positioning of the working channel in the direction of the affected site. Soft tissue is removed from the Ligamentum flavum beneath the microscope.

Blunt opening of the Ligamentum flavum using a dissector.

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Pictures 11-15: Resection of parts of the flavum to extend access into the spinal canal to reveal the dural tube and/or nerve root.

Mobilisation of nerve root or dural sack. Protection of the neural structures using the nerve root retractor. Removal of the herniated disk.

Removal or "unscrewing" of the working channel under observation in order to stem any bleeding from the paraspinal musculature.

Examination beneath the surgical microscope shows how the muscle fibres overlaid each other on retracting the working channel. There is no need to insert a suction tube.

The wound is closed by one single suture of the fascia and one for the subcutis. Skin adaptation by a bandaid. The figure illustrates the cosmetic outcome 3 months after surgery.

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











Nerve Root Retractors and Alligator Type Scissors

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







Schmale Stanzplatte
Narrow footplate

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Medicon offers various versions of the dismantable rongeurs. The "CLASSIC" and "COMFORT" lines are each available with thick, narrow or extra-narrow footplate. The rongeurs of the "COMFORT" line are additionally equipped with a three-step control element, that permits selecting from among three different levels of grip width and jaw opening.

Extra schmale Stanzplatte
Extra narrow footplate

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© Copyright 12/2024, MEDICON eG, Tuttlingen, Germany