Lipo-Injektion Dosierhandgriff
REF 09.87.90/09.87.93/09.87.95

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Lipo-Injection Dosage-handle REF 09.87.90/09.87.93/09.87.9

Table of contents:

- Safety precautions
- Intendend use
- Putting into operation
- Handling
- Assembling- / Dismantling-instruction
- Instruction for preparation and sterilization

SAFETY PRECAUTIONS

Medicon instruments are manufactured of high-quality materials. They are developed, designed and produced with the greatest care. Constant further development and strictest quality monitoring are standard procedure at Medicon.

Only proper use of these quality instruments enables the best possible working results and a long service life of the instruments. To ensure this, the following operation instructions and safety precautions must always be observed and complied with.

It is the responsibility of the user to check the product on his/her own responsibility for

This operating instruction accompany the following instruments:

- 09.87.90
- 09.87.93
- 09.87.95

INTENDEND USE

The dosage-handle was developed specifically for re-injection of fat in Plastic and Aesthetic Surgery.

PUTTING INTO OPERATION

The Lipo-Injection dosage-handle is delivered as unsterile. Brand-new instruments must be cleaned, disinfected and sterilised before initial use and after each additional use

Before use it must be ensured that the tools employed are in faultless technical and sterile condition. Generally, tools which are soiled or damaged must not be used.

HANDLING



ILLUSTRATION:

- 1 Device to fix rubber stop
- 2 Bayonet lock
- 3 Unlocking device
- 4 Ratched rod
- 5 Pulling ring
- 6 Handle
- 7 Lever

ASSEMBLING- / DISMANTLING-INSTRUCTION

For preparation and sterilizing the bayonet lock (2) has to be screwed off the handle (6) in direction a.

After sterilizing resp. before usage fix the bayonet lock (2) onto the handle (6) again and turn in direction b.

PREPARATION OF LIPO-TRANSFER


WARNING AND PRECAUTIONS

Lipo-Transfer as well as its preparation may be effected only in sterile condition!


A: Dismantle disposable plastic syringe 09.87.92/ 09.87.94 / 09.87.97


B: Mount rubber stop on metal fixing device (1)


C: Push the barrel over the rubber stop up to the bayonet lock (2)

D: Lock barrel by turning in direction of arrow


E: Push rubber Stopp into the barrel by pressing the handle (7)


Lipo-Transfer (Example)


Fix connection between vacuum syringe 09.87.82, transfer-adaptor 09.87.83 and plastic syringe 09.87.92 / 09.87.94 / 09.87.97.

Keep pressed unlocking device (3) during Lipo-Transfer

Fix LL injection needle 09.88.01 -09.88.04 on disposable plastic syringe

INSTRUCTION FOR PREPARATION AND STERILIZATION

ADVICE

Please consider our general instruction for preparation MPK01 for thermostable medical devices. If it is not present at yours, please ask us for receiving it

WARNING AND PRECAUTIONS

- Instruments of stainless steels must not be deposit in physiologic saline solution (NaCl-solution), because a longish contact may cause corruptions like pitting and stress corrosion.
- Only cleaned and disinfected instruments may be sterilized.
- The cleaning, disinfection and sterilization has to be made in dismantled condition.
- Brand-new instruments have to be procured to a manual resp. automatic preparation before initial use.

RESTRICTION OF REPREPARATION:

Frequent re-preparation has slight effects on the instruments. The end of the life cycle is normally determined by wear and damages by usage.

- Directly after the usage gross contaminations as to be removed of the instrument with a disposable sheet / paper.
- No fixing mediums or hot water (>40°C) must be used, as this causes fixation of residuals and may have effect on the cleaning-result.
- Instruments have to be procured to the preparations-process immediately. Dry removal is preferred.

PREPARATION FOR DECONTAMINATION

Instruments have to be dismantled and placed in a suitable and machine-fitting instrument-tray.

The trays (e.g. wire baskets) must be so designed, that the following cleaning in ultrasonic bath or cleaning/disinfection device will not be interfered by sonic- or rinsing-shades.

CLEANING

pH-adjusted or alkaline cleaner are employed as cleaner. The choice of the cleaning agents depends on the characteristics of the instruments and national directives and recommendations.

If increased chloride-concentrations are existing in the water, pitting and stress corruptions may occur at the instruments. By using fully desalted water or alkaline cleaning agents such corruptions can be minimized.

The cleaning/disinfection device must meet the requirements of DIN EN ISO 15883-1.

- 1. Pre-rinsing 1:** 1 minute with softened cold water, without addition
- 2. Emptying**
- 3. Pre-rinsing 2:** 3 minutes with softened cold water, without addition
- 4. Emptying**
- 5. Cleaning:** with softened water, heating-up to 55°C and 5 minutes washing / cleaning, addition of cleaning agent at 45°C, alkaline cleaning agent, dosage 0,5%
- 6. Emptying**
- 7. Neutralisation:** 3 minutes with warm water (>40°C) with addition of neutralisator, dosage 1 ml/l
- 8. Emptying**
- 9. Final rinsing:** 2 minutes with fully desalted water (>40°C), without additions
- 10. Emptying**

The parameters for concentration, temperature and application time, given by the manufacturer of the cleaning agent, must be complied and automatic metering devices must be controllable.

THERMIC DISINFECTION

Use fully desalted water. The thermic disinfection has to be performed at temperatures of 80-95°C and corresponding application time acc. to A0-Konzept, DIN EN ISO 15883-1.

DRYING

An adequate drying has to be assured by the cleaning/disinfection device. The temperature must not exceed 95°C at this. If necessary, the usage of compressed air for drying is recommended (RKI-recommendation) because of its good and quick effect.

MAINTENANCE, CHECKING AND INSPECTION

After the cleaning / disinfection the instruments have to be macroscopic clean, i.e. free from visible contaminations and residuals. The checking is made visual. Critical areas like hinges, compression spring and ratched rod

have to be checked very carefully.

Insufficient cleaned instruments or instruments with blocked passages have to be recleaned and then rinsed off sufficiently. Instruments with movable parts have to be cooled down and lubricated with instrument maintenance oil prior to function test.

Defective instruments (hairline cracks, deformations or wears) has to be replaces, as the function is no more or no more sufficient fulfilled.

STERILIZATION

Only cleaned and disinfected instruments may be sterilized. For the sterilization following method of sterilization has to be used considering the respective national requirements:

- Fractional vacuum-procedure triple fractionated and with sufficient product-drying
- Steam sterilizer acc. to DIN EN 13060 resp. DIN EN 285 and validated acc. to DIN EN ISO 17665-1.
- Duration of sterilization and temperature:
Minimum dwell time 5 minute at 134°C

PACKAGING

The instruments have to be deposited in a suitable sterile barrier-system.

The sterile barrier-system has to meet following criterias:

- DIN EN 868 and DIN EN ISO 11607
- Suitable to steam sterilization (vapour permeability)
- Adequate temperature-stability up to 138°C
Accessory for sterilization and packaging for sterilization have to be adjusted to the content of packaging as well as to the used sterilization-procedure.

STORAGE

Re-prepared sterile instruments has to be stored in a suitable reusable sterilizing-case dry, dustproof, almost sterile, dark and in cool rooms, free from vermin. For avoiding the generation of condensate, greater variations in temperature should be avoided. Chemicals must not be stored together with instruments. The duration of storage acts on the style of packaging.

The admissible duration of storage acts on the style of the used sterile barrier-system and the storage conditions. The admissible durations of storage has to be defined by the operator.

ADVICE

The operator is accountable for validating the recommended parameters for accommodating potential differences regarding sterilization-cells, methods of packaging and loading-configurations.

The same applies, if a process is used which is not corresponding to the recommendations. It is essential to reach a SAL (Sterility Assurance Level) of 10⁻⁶.

MEDICON eG, being the distributor of these products, does not assume liability for direct damages or secondary damages, which arise out of improper usage or handling especially by disregarding the attached operating instruction or by improper care or maintenance.