
G517902

Otis Urethrotom

CE 0297


REF 51.79.02

medicon®

Medicon eG
Gänsäcker 15
78532 Tuttlingen · Germany
Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50
e-mail: sales@medicon.de
www.medicon.de

Deutsch	(DE)	3
English	(EN)	15
Español	(ES)	27
Français	(FR)	39
Italiano	(IT)	51

Otis Urethrotom

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DE

Produktbeschreibung

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Please read these instructions for use carefully before clinical application and store them in a safe place where they can be consulted easily.

Table of contents

1. General information	16
2. Scope of delivery	16
3. Intended use	16
4. Indications	17
5. Contraindications	17
6. Possible side effects and complications	17
7. Before use	17 – 18
8. Assembly	18
9. Disassembly	19
10. Use and handling	19
11. Instructions on reprocessing	20 – 25
12. Liability	26
13. Explanation of symbols and pictograms	26

1. General information

The Medicon Otis urethrotomes are not sterile when delivered and must be cleaned, disinfected and sterilised prior to their first and every subsequent use. Please observe the information below. This will guarantee the correct, reliable functioning of the devices.

2. Scope of delivery

Medicon Otis urethrotome, complete **REF 51.79.02**

comprising:

REF	51.79.05	Dilator
REF	51.79.10	Knife
REF	51.79.11	Attachment, conical
REF	51.79.12	Attachment, flat
REF	51.79.13	Attachment, ball

3. Intended use

The reusable Medicon Otis urethrotome has been developed for urethrotomy in medical facilities. In this procedure, the urethrotome is introduced into the urethra and an incision made at the constricted point. The Medicon Otis urethrotome may only be employed by qualified urologists with sufficient experience in this endoscopic treatment procedure.

4. Indications

- Short urethra strictures, constriction of the urethra, without pronounced fibrotic reaction
- Urolithiasis / urethral obstruction

5. Contraindications

Urinary tract infections (UTIs) or coagulation disorders

6. Possible side effects and complications

Incisions which are made too deeply can damage the corpus spongiosum (spongy body).

This results in haemorrhaging and longer periods of catheterisation to stem the bleeding. Erectile dysfunction and sphincter injuries with incontinence occur more rarely.

All surgical interventions are associated with certain complication risks. However, these are more often than not caused by insufficient preparation and not the instruments themselves.

The complications which may occur can often be traced back to incorrect patient selection or lack of practice on the part of the doctor.

Complications:




- Damage to nerves, tissue or vessels, infections or other accidental injuries during the surgical intervention.
- Using the instruments for purposes other than their intended use can damage them. Damaged instruments cannot always be used safely.

7. Before use

Warnings and safety measures:
The knife blade ⑦ poses a risk of injury.

The accessory articles intended for use with the instrument are subject to wear and tear and mechanical loads even when used normally. To prevent failures or mechanical damage to the instrument during the operation, all parts should be inspected for mechanical integrity, distortion, cracks, scratches, corrosion and complete functioning before every use. The knife blade must not be damaged. Instruments which are not in perfect condition must not be used.

Functional test:


- Smooth operation of the Set screw ⑤;
- When the set screw ⑤ is turned, the shaft ② must open and close (spread).
- The indicator on the scale ③ must change accordingly.

8. Assembly

Select and screw the distal attachment ① in place (Fig. 1). Insert the knife blade ⑦ in the guide channel ⑩ before pressing the knife shaft ⑧ gently against the instrument. Hold the knife ④ by the handle ⑨ and insert the shaft into the guide channel ⑩ (Fig. 2). The knife ④ must be inserted until the knife blade ⑦ disappears completely into the guide channel ⑩. When inserting the knife ④ into the guide channel, attention must be paid to the fact that the curvature of the handle ⑨ points away from the instrument (Fig. 3).

Only use the accessories specially approved by Medicon in order to avoid risks associated with the compatibility of the accessory articles with the urethrotome.

9. Disassembly

The knife ④ must be removed from the guide channel ⑩ before the urethrotome is prepared. This is done by holding the knife ④ by the handle ⑨ and pulling it out of the guide channel ⑩. The distal attachment must be unscrewed for cleaning. The instrument can be opened by turning the set screw ⑤ clockwise until the indicator shows 45 Fr.

10. Use and handling

The instrument must be covered with a medical lubricant, then the distal end of the instrument is inserted up to the required position in the ure-

thra. Now the spreading width can be set, which is done by turning the set screw ⑤ clockwise. The spreading width can be between 15 Fr and 45 Fr and can be read off the scale ③. Turning the set screw ⑤ opens the shaft ② of the instrument (Fig. 4).

The incision is performed when the knife is pulled out.

Once the knife ④ has been pulled out, the shaft ② is closed again by turning the set screw ⑤ anticlockwise. The indicator on the scale ③ should be at 15 Fr. The urethrotome can now be carefully removed from the urethra.

11. instructions on Reprocessing

(Machine cleaning, alkaline pH >10; manual cleaning, neutral enzymatic; without precleaning)

The following prion-specific safety measure (v)CJD for instrument preparation is indicated.

In cases of diagnosed definite or probable (v)CJD.

If it not possible to use disposable articles, the instrument set where contamination has been identified or cannot be excluded must be incinerated.

In suspected cases:

If prion contamination is suspected, the Task Force vCJD final report recommends incinerating the instrument set.

If (v)CJD is excluded:

Instrument can be used again after being prepared completely. In all other cases, the instrument set where contamination has been identified or cannot be excluded must be incinerated.

In the case of a non-detectable (v)CJD infection.

Even if there is no indication of the presence of a prion infection, two procedures with at least partial prion efficacy should be used

for preparation – e.g., alkaline machine cleaning in combination with steam sterilisation.

If no alkaline machine cleaning or other cleaning procedure with documented prion efficacy is used and the instruments in question are medical devices which come into contact with risk tissues (CNS, eyes, lymphatic tissue), the Robert Koch Institute (RKI) recommends an extended sterilisation time of 18 minutes at 134°C.

Instruments made of stainless steel must not be immersed in physiological saline solution (NaCl solution), as prolonged contact can result in forms of corrosion such as pitting and stress cracking corrosion. Only cleaned and disinfected instruments may be sterilised.

Preparation restrictions

Frequent preparation has little effect on the instruments. The end of the service life is normally determined by wear and tear and damage resulting from use.

Restrictions on reprocessing

Place of use

- Directly after use, heavy soiling should be removed with a disposable cloth/paper towel.
- Disassemble the instrument in accordance with the disassembly instructions
- No fixing agents or warm water ($>40^{\circ}\text{C}$) must be used as they result in the fixing of residues and remnants, which can impede the success of the cleaning.
- The instruments must be prepared **immediately**.
- Dry disposal is preferable.

Preparation for decontamination

Instruments with joints should be opened for preparation. The instruments should be placed in instrument trays suitable for use in machines in such a way that they are rinsed completely.

The instrument trays (e.g., wire baskets) must be structured in such a way that subsequent cleaning in an ultrasound bath or a washer-disinfector is not impeded by areas where the ultrasound or rinsing cannot reach.

Cleaning: machine

The washer-disinfector must comply with the requirements of DIN EN ISO 15883-1.

- rerinsing 1: 1 minute with softened, cold water, with no additives.
- Drainage.
- Prerinsing 2: 3 minutes with softened, cold water, with no additives.
- Drainage.
- Cleaning: with softened water, heating to 55°C and 10 minutes washing/cleaning, addition of cleaning agent at 45°C , alkaline cleaning agent, dose 0.5%.
- Drainage.
- Neutralisation: 3 minutes with warm water ($>40^{\circ}\text{C}$) adding neutralising agent, dose 1 ml /l.
- Drainage.
- Final rinse: 2 minutes with warm, completely demineralised water ($>40^{\circ}\text{C}$) (without any other additives).
- Drainage.

Disinfection

Thermal disinfection A0 value 3000: Completely demineralised water, thermal disinfection is performed at temperatures $>80^{\circ}\text{C}$ with corresponding exposure times in accordance with the A0 concept of the DIN

EN ISO 15883 standard and DGKH, DGSV and AKI guideline (e.g., A0 3000 = 90°C and 5 minutes exposure time). The operator bears the responsibility for the A0 value to be achieved.

Drying

It must be ensured that the washer-disinfector dries properly. The instruments should be removed from the washer-disinfector as soon as the cleaning and disinfection programme finishes. If necessary, compressed air is recommended for drying due to its great and rapid efficacy (RKI recommendation).

Cleaning/disinfection – manual

The cleaning and disinfectant agents must be thoroughly suitable for manual cleaning and disinfection of instruments and compatible with each other. The disinfectant must possess verified efficacy.

The pertinent lists and recommendations of the Robert Koch Institute (RKI) and the German Society for Hygiene and Microbiology (DGHM) must be taken into consideration when selecting the disinfectant and disinfection procedure.

Precleaning

- Immerse the instruments in cold water for 5 minutes.
- Brush (using plastic brushes) the instruments under cold water until all visible soiling has been removed.
- Internal cavities, threads and holes must each be rinsed with a water gun for 10 seconds and then brushed again.
- Remove the instruments and rinse them with cold, completely demineralised water.

Cleaning/disinfection

- Immerse the instruments in a bath with a tested cleaning and disinfectant agent.
- The instruments must be completely covered with the solution.
- The exposure times, temperatures and concentration specified by the cleaning/disinfectant agent manufacturer must always be observed.
- Remove the instruments and rinse them with cold, completely demineralised water for 2 minutes.
- Repeat the cleaning process if any soiling is still visible on the instrument.
- If necessary, perform the cleaning and disinfection ultrasonically in an ultrasound bath.

The solutions must be prepared fresh each day. The solutions must be

changed regularly if they become heavily soiled. Excessive soiling in the ultrasound bath affects the efficacy of the cleaning and increases the risk of corrosion. The cleaning solution must be renewed regularly depending on the application conditions. This should be done when soiling is visible. In any case, the liquid should be changed regularly – at least once a day. The national guidelines must be followed.

Drying

Manual drying with compressed air and a lint-free cloth. Compressed air is recommended for drying due to its great and rapid efficacy (RKI recommendation).

Maintenance, checks and inspections

Following cleaning/disinfection, the instruments must be macroscopically clean, i.e., free from visible soiling and residues/remnants. They are checked visually. Insufficiently cleaned instruments must be recleaned and then adequately rinsed and dried. Instruments with moving parts must be allowed to cool and must be lubricated with instrument care oil (e.g., Medicon **REF 46.00.40**) before they undergo a functional test. Instruments such as clamps and needle holders which have a notch must only be closed up to the first notch (risk of stress cracking). Defective instruments (hairline cracks, deformations or wear and tear) must be

replaced as they can no longer perform their function or cannot do so adequately. Corroded instruments must also be removed as they can trigger corrosion of intact instruments via the transfer of extraneous rust.

Packaging

The instruments must be inserted in a suitable sterile barrier system. The sterile barrier system must comply with the following criteria:

- DIN EN 868
- DIN EN ISO 11607
- Suitable for steam sterilisation (steam permeable)
- Heat resistant up to 138°C

Sterilisation accessories and sterilisation packaging must be compatible with both the contents of the packaging and the sterilisation procedure employed.

Sterilisation

The following sterilisation procedure should be applied for sterilisation taking the respective national requirements into consideration:

- Triple fractionation vacuum procedure with sufficient drying of the product
- Steam sterilizer in accordance with DIN EN 13060/DIN EN 285 and

validated in accordance with DIN EN ISO 17665-1.

- Sterilisation time and temperature: Minimum 5 minutes holding time at 134°C

Achieving an SAL (sterility assurance level) of 10^{-6} is indispensable.

Storage

Prepared, sterile instruments should be stored in a suitable, reusable sterilisation container dry, protected from dust, sterile and in dark, cool rooms free from vermin. To avoid the formation of condensation, considerable temperature fluctuations must be avoided during storage. No chemicals should be stored together with the instruments.

The walls, floors and ceilings of the storage room must be smooth and easily to clean and disinfect. Shelves must be a minimum of 30 cm above the floor. The permissible storage period on site depends on the type of sterile barrier system employed and the storage conditions. The permissible storage period must be specified by the operator.

Further information on preparation

A validated machine cleaning and disinfectant procedure should always be given preference over manual cleaning based on the higher degree of safety associated with the process. Good cleaning preserves the value

of the instruments and is a prerequisite for successful sterilisation.

The following points must be taken into account during machine preparation:

- The instruments must be loaded in the instrument trays in such a way as to allow proper rinsing – this is a prerequisite for efficient machine preparation. Instrument trays must not be overloaded.
- It must be ensured that large instruments do not cover areas preventing them from being properly rinsed.
- The instruments must be positioned/stored in accordance with their mechanical stability in order to exclude damage.

The times and temperatures specified in these instructions for preparation are minimum requirements and must be maintained at all times. If lower values need to be used for procedural reasons, the operator must validate them accordingly. It is generally possible to exceed the times and temperatures specified, but this exposes the materials to higher loads, which can cause the instruments to age prematurely. We are not responsible for the use of other sterilisation procedures.

Information on the validation of the preparation

The validation was performed with the following devices, materials and chemicals:

Washer-Disinfector:	Miele disinfector G 7735 CD and 7836 CD
Alkaline cleaning agent:	Container trolley for surgical instruments
Neutral cleaning agent:	neodisher® from Dr. Weigert GmbH & Co. KG
	Endozime, Ruhof (enzymatic)
	Decondex 23 Neutrazym from Borer Switzerland
Neutralising agent:	neodisher® Z from Dr. Weigert GmbH & Co. KG
Water gun:	Selecta
Cleaning brushes:	Brushes made of plastic/nylon
Ultrasound bath:	Sonorex
Sterilizers	MMM Vaculab 969 S 3000, MMM Selectomat S 3000 Stiefenhofer KS 666-2ED, H+P Varioclav 400E

Note

The person preparing the instruments is responsible for the fact that the preparation actually performed with the equipment, materials and personnel employed in the preparation facility achieves the required results. This usually requires validation and routine monitoring of the procedure. If the devices, materials and chemicals listed above are not available, the person preparing the instruments is responsible for validating his procedure accordingly. Please observe the information and specifications of the pertinent, national statutory regulations and standards. Please note that all instruments sent to the Medicon repair service (MRS) for repair must be cleaned and sterilised before being sent. Medicon eG reserves the right to make changes to these instructions to reflect new findings and discoveries.

12. Liability

The information in these instructions for use refers to the product with which the instructions for use are supplied. The instructions for use for the respectively current product can be downloaded from the Medicon Extranet at www.medicon.de/extranet or requested from MEDICON eG directly.

MEDICON eG provides no guarantees for damages resulting from improper use, care, maintenance or failure to observe the restrictions of use and other specifications in the instructions for use. MEDICON eG's liability for defects is also excluded if the product is modified or repaired without the advanced, written approval of MEDICON eG and in the case of repairs performed in workshops not authorised by MEDICON eG or not performed by the Medicon repair service (MRS).

13. Explanation of symbols and pictograms

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© Copyright 07/2025, MEDICON eG, Tuttlingen, Germany