G097770

Gebrauchsanweisung
Operating Instructions
Manual de instrucciones
Mode d'emploi
Istruzioni d'uso

CE 0123

Rx only

medicon®
Medicon eG · Gänsäcker 15
78532 Tuttlingen · Germany
Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50
E-Mail: sales@medicon.de
Internet: www.medicon.de

DE — VERESS KANÜLEN

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EN — VERESS CANNULAS

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

Contents

1. General notes
2. Intended use/Indications
3. Application/Handling
4. Disassembly
5. Assembly
6. Special requirements for preparation an pre-cleaning
7. Decontamination, cleaning and sterilization
8. Note
9. Liability
10. Signs and symbols

1. GENERAL NOTES

MEDICON Veress cannulas are supplied unsterile and must be disinfected, cleaned and sterilised prior to use.
Please observe the following instructions. This will guarantee that the cannulas will function faultlessly and reliably.

2. INTENDED USE/INDICATIONS

Medicon Veress cannulas are used mainly in laparoscopic, gynaecological, urological and arthroscopical surgeries. They are used for initial puncturing with subsequent gas insufflation.

3. APPLICATION/HANDLING

Observe the Instructions for use of the insufflation device to be used. Before puncturing, the snapping mechanism must be checked by pulling back the needle and letting it spring forth several times. During puncturing, hold the needle by the handle portion so that the inner cannula can move freely and return to its initial position after the puncture. Close the cock to stop the gas supply.

4. DISASSEMBLY

To unscrew the sharp-edged outer cannula, the grooved handle sleeve must be turned in counter-clockwise direction. Make sure that the sensitive tip of the outer cannula is particularly well-protected against damage during this process and during the subsequent cleaning and sterilization.
To remove the stopcock from the cock housing, unscrew the spring screw cap from the bottom of the cock housing.

5. ASSEMBLY

When re-inserting the stopcock into the cock housing, ensure that the lateral stop pin of the stopcock is positioned in the recess on the upper edge of the cock housing. Please grease the stopcock slightly. In order to mount the stopcock, screw the spring screw cap onto the thread on the bottom of the cock housing.
The sharp-edged outer cannula is now slid onto the blunt obturator and screwed on in clockwise direction, holding it by the handle sleeve. Functional testing is mandatory after the mounting process.

6. SPECIAL REQUIREMENTS FOR PREPARATION AND PRE-CLEANING

The insufflation needle should be flushed with cold mains water using a water gun for at least 10 seconds in pulse mode. Then dismantle into four parts and flush the tube separately likewise for at least 10 seconds using the water gun.
Connect the Luer-Lock connector of the cannula to a flushing connection of the trolley for MIC instruments and place the tube in a basket. Place the remaining parts into an instrument tray. Then start one of the following processes:

7. DECONTAMINATION, CLEANING AND STERILIZATION

Matrix for the preparation procedure MPK 02

Procedure				T	t	
Place of use	01	Disposal	Wet	–	–	--
			Dry	–	–	++
Decontamination	02	Pre-cleaning	manual	<30 °C	10 min.	++
	03		Ultrasonic	40 °C	15 min.	--
	04	Cleaning	Manual	45 °C	>10 min.	++
			Mechanically	55 °C	10 min.	++
			Alkaline cleaning agent > pH10	–	–	++
			Acidic cleaning agent	–	–	0
			Neutral cleaning agent	–	–	++
	05	Rinsing	Desalinated water	>40 °C	2 min.	++
	06	Neutralisation	Citric acid	>40 °C	3 min.	++
			Phosphoric acid	>40 °C	3 min.	0
	07	Disinfection	Chemical	20 °C	10 min.	++
			Thermal	90°C	5 min.	++
	08	Thermal dry	Thermal	max. 120 °C		++
			Medical compressed air	–	–	++
Maintenance and inspection	09	Functional test	Instrument oil REF 46.00.44	–	–	+
Packaging	10	Packaging	Sterile packaging	–	–	+
Sterilization	11	Steam	Sterile packaging	–	–	+
			SAL 10^{-6}	134 °C	5 min.	++
			SAL 10^{-6}	134 °C	4 min.	++
			SAL 10^{-6}	132 °C	4 min.	++
		Hot air	–	180 °C	<30 min.	0
		Gas	–	–	–	–
		Plasma	–	–	–	–

Required process steps, mechanically: 01-02-04-05-06-05-07-08-09-10-11
Required process steps, manual: 01-02-04-05-07-05-08-09-10-11

Legend:

T	Temperature	++	Procedure validated by the manufacturer
t	Soaking time	+	Required process step
min	Minutes	0	Non-validated procedure, check material compatibility before initial use
°C	Degrees Celsius	00	Not applicable
		–	No information
		--	Not required

The up-to-date and comprehensive process description for reprocessing, information on prion protection measures (v) CJK as well as the validation references are available for download at ***www.medicon.de/extranet/***

8. NOTE

The user is responsible for the actual processing achieving the desired results with the equipment, materials, and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method.

If the previously described equipment, materials, and chemicals are not available, it is the responsibility of the user to validate his method accordingly. Please note the instructions and regulations of the relevant national regulations and standards and any instructions for use accompanying the medical device.

Please note that all instruments sent to the Medicon Repair Service (MRS) for repair must be cleaned and sterilized prior to dispatch. Medicon eG reserves the right to modify these instructions whenever new information is obtained.

9. LIABILITY

In the event of discrepancies between the english and the german version of these instructions for use only the German version shall be applicable. Only the latest revision of the instructions for use applies.

Due to constant technical development the contents of these MEDICON instructions for use are updated regularly. Please use the MEDICON extranet ***www.medicon.de/extranet*** to ensure that you are using the current version. The current version date is printed on each edition of the instructions for use. MEDICON eG assumes no liability for damages caused by improper use, care or maintenance or non-compliance with the restrictions for use, non-observance of post-operative instructions and other guidelines in the instructions for use.

The liability for defects by MEDICON eG is also not applicable in case of changes or repairs to the product without prior written consent from MEDICON eG as well as in case of repairs which were not carried out by workshops authorised by MEDICON eG or the Medicon Repair Service (MRS).

Caution!
Federal law and regulations in the United States restrict the sale of this device to a physician or hospital or on prescription.

10. SIGNS AND SYMBOLS

	Follow instructions
	Warning
CE 0123	CE label accord. to directive 93/42/EEC
NON STERILE	Non-sterile
	Manufacturer
LOT	Batch number
REF	Order number
Rx only	Prescription use only (US federal law)

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