



Product data sheet REF 69.02.35

11.12.2025

| General data | | | |
|---|-----------------------------------|--|---|
| article number | 69.02.35 | implant tray Mandibular Reco Lock Syst 2.4 | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 12143 | GTIN | 04046826338763 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 12-142 |
| Classification | 1 (EU MDR) | Base UDI-DI | 4046826;N003;03;X8 |
| Net weight (kg) | 2,150 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|  | | | |
| size | | Coating | Anodizing is the process of enriching titanium with a hard anodized layer on aluminum components. |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.3547 (AlMg4,5Mn) | | |
| Screw | 1.4305 (X8CrNiS18-9) | | |
| Others | PPSU (Polyphenylsulfon) | | |