



Product data sheet REF 04.75.73

03.06.2026

General data			
article number	04.75.73	Gorney Freeman HM scrs cvd 23cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	38727	GTIN	04046826004835
Risk category	1r (EU MDR)	UMDNS-Code	13-500
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;01;ag;2U
Net weight (kg)	0,070	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.53 (catalogue cardio) 450.00.60 (catalogue plastic surgery) 451.03.23 (brochure medicut scrs)		
Complete picture		Detail images	
			
size	23cm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4117 (X38CrMoV15)		
Screw	1.4305 (X8CrNiS18-9)		
Others	WSG_20-GO-55_CRTZ (RCoCr-C)		