






## Product data sheet REF 09.21.16

11.12.2025

| General data  |  |   |                         |
|---|--|---|-------------------------|
| article number  | 09.21.16   | Lathbury (Farrell) cotton carrier 16cm serr |                         |
| Unit  | piece  | Customs code (HS code)                      | 90189084                |
| GMDN  | 32873  | GTIN  | 04046826010904          |
| Risk category   | 1r (EU MDR)  | UMDNS-Code                                  | 13-698                  |
| Classification  | 6 (EU MDR)   | Base UDI-DI                                 | 4046826;B;R002;01;ah;3J |
| Net weight (kg)   | 0,005  | Operating instructions                      | GATDR00X                |
| Latex free  | yes  | Preparation instructions                    |                         |
|   |  |   |                         |
| Medicon media   | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.40 (Katalog Nr 5 Dental )<br>450.00.46 (catalogue sinuscopy ) |   |                         |
| Complete picture  |  | Detail images                               |                         |
|   |  |   |                         |
| size  | 16cm   | Coating                                     |                         |
| Storage conditions  |  |   |                         |
| No special requirements   |  |   |                         |
| Materials   |  |   |                         |
|   | Material number (DIN designation)  |   |                         |
| Main part   | 1.4305 (X8CrNiS18-9)   |   |                         |