



Product data sheet REF 09.21.76

11.12.2025

General data			
article number	09.21.76	cotton carrier Ø 1.2mm 16cm triang serr	
Unit	piece	Customs code (HS code)	90189084
GMDN	32873	GTIN	04046826011031
Risk category	1r (EU MDR)	UMDNS-Code	13-698
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;01;ah;3J
Net weight (kg)	0,006	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT)		
Complete picture		Detail images	
			
size	Ø 1.2mm 16cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4305 (X8CrNiS18-9)		