






Product data sheet REF 09.73.14

11.12.2025

General data			
article number	09.73.14	Tuohy needle 1.50x80mm LL	
Unit	piece	Customs code (HS code)	90183900
GMDN	38749	GTIN	04046826012014
Risk category	1r (EU MDR)	UMDNS-Code	15-206
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;03;ag;3J
Net weight (kg)	0,026	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	1.50x80mm	Coating	Electroplated chrome layer on the surface
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	1.4310 (X10CrNi18-8)		
Others	2.0401 (CuZn39Pb3)		