



Product data sheet REF 09.73.93

11.12.2025

| General data | | | |
|---|--|----------------------------|---|
| article number | 09.73.93 | Bier needle 1.20x90mm Luer | |
| Unit | piece | Customs code (HS code) | 90183900 |
| GMDN | 38749 | GTIN | 04046826012045 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-206 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;03;ag;3J |
| Net weight (kg) | 0,032 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 1.20x90mm | Coating | Electroplated chrome layer on the surface |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |
| Screw | 1.4310 (X10CrNi18-8) | | |
| Others | 2.0401 (CuZn39Pb3) | | |