






Product data sheet REF 09.75.84

11.12.2025

General data			
article number	09.75.84	trocar 105mm f Ø 1.40mm	
Unit	piece	Customs code (HS code)	90183900
GMDN	46899	GTIN	04046826012267
Risk category	1r (EU MDR)	UMDNS-Code	14-154
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;05;ag;48
Net weight (kg)	0,058	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	105mm Ø 1.40mm	Coating	Electroplated chrome layer on the surface
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	1.4310 (X10CrNi18-8)		
Others	2.0401 (CuZn39Pb3)		