



## Product data sheet REF 09.77.88

11.12.2025

General data			
article number	09.77.88	3 way stopcock	
Unit	piece	Customs code (HS code)	90189084
GMDN	35384	GTIN	04046826012458
Risk category	2a (EU MDD)	UMDNS-Code	11-729
Classification	2 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,030	Operating instructions	
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size		Coating	Electroplated chrome layer on the surface
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	2.0401 (CuZn39Pb3)		