



## Product data sheet REF 09.83.38

11.12.2025

| General data   |                                   |                                   |                |
|--|-----------------------------------|-----------------------------------|----------------|
| article number   | 09.83.38                          | Frazier stylet only f 8 Fr 12.5cm |                |
| Unit   | piece                             | Customs code (HS code)            | 90189084       |
| GMDN   | 60909                             | GTIN                              | 04046826013332 |
| Risk category  | KEIN MP (EU MDD)                  | UMDNS-Code                        | 00-000         |
| Classification   | 6 (EU MDD)                        | Base UDI-DI                       |                |
| Net weight (kg)  | 0,035                             | Operating instructions            |                |
| Latex free   | yes                               | Preparation instructions          | MPK03          |
|   |                                   |                                   |                |
| Medicon media  |                                   |                                   |                |
| Complete picture   |                                   | Detail images                     |                |
|  |                                   |                                   |                |
| size   | 8 Charr 12.5cm                    | Coating                           |                |
| Storage conditions   |                                   |                                   |                |
| No special requirements  |                                   |                                   |                |
| Materials  |                                   |                                   |                |
|  | Material number (DIN designation) |                                   |                |
| Main part  | 1.4310 (X10CrNi18-8)              |                                   |                |
| Others   | 1.4305 (X8CrNiS18-9)              |                                   |                |