



## Product data sheet REF 09.84.08

25.08.2026

| General data                                                                                  |                                                                                |                                   |                |
|-----------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|-----------------------------------|----------------|
| article number                                                                                | 09.84.08                                                                       | Plester suct tube 8 Fr ang 13.5cm |                |
| Unit                                                                                          | piece                                                                          | Customs code (HS code)            | 90189084       |
| GMDN                                                                                          | 38749                                                                          | GTIN                              | 04046826013530 |
| Risk category                                                                                 | 2a (EU MDD)                                                                    | UMDNS-Code                        | 10-566         |
| Classification                                                                                | 6 (EU MDD)                                                                     | Base UDI-DI                       |                |
| Net weight (kg)                                                                               | 0,038                                                                          | Operating instructions            | G098395        |
| Latex free                                                                                    | yes                                                                            | Preparation instructions          | MPK03          |
| <div></div>  |                                                                                |                                   |                |
| Medicon media                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT ) |                                   |                |
| Complete picture                                                                              |                                                                                | Detail images                     |                |
| <div></div> |                                                                                |                                   |                |
| size                                                                                          | 8 Charr 13.5cm                                                                 | Coating                           |                |
| Storage conditions                                                                            |                                                                                |                                   |                |
| No special requirements                                                                       |                                                                                |                                   |                |
| Materials                                                                                     |                                                                                |                                   |                |
|                                                                                               | Material number (DIN designation)                                              |                                   |                |
| Main part                                                                                     | 1.4301 (X5CrNi18-10)                                                           |                                   |                |
| Others                                                                                        | 1.4305 (X8CrNiS18-9)                                                           |                                   |                |