



## Product data sheet REF 09.85.00

25.08.2026

| General data                                                                                  |                                                                                                                 |                          |                |
|-----------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|--------------------------|----------------|
| article number                                                                                | 09.85.00                                                                                                        | Yankauer suct tube 28cm  |                |
| Unit                                                                                          | piece                                                                                                           | Customs code (HS code)   | 90189084       |
| GMDN                                                                                          | 38749                                                                                                           | GTIN                     | 04046826013721 |
| Risk category                                                                                 | 2a (EU MDD)                                                                                                     | UMDNS-Code               | 10-566         |
| Classification                                                                                | 6 (EU MDD)                                                                                                      | Base UDI-DI              |                |
| Net weight (kg)                                                                               | 0,076                                                                                                           | Operating instructions   | G098395        |
| Latex free                                                                                    | yes                                                                                                             | Preparation instructions | MPK03          |
| <div></div>  |                                                                                                                 |                          |                |
| Medicon media                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.53 (catalogue cardio ) |                          |                |
| Complete picture                                                                              |                                                                                                                 | Detail images            |                |
| <div></div> |                                                                                                                 |                          |                |
| size                                                                                          | 28cm                                                                                                            | Coating                  |                |
| Storage conditions                                                                            |                                                                                                                 |                          |                |
| No special requirements                                                                       |                                                                                                                 |                          |                |
| Materials                                                                                     |                                                                                                                 |                          |                |
|                                                                                               | Material number (DIN designation)                                                                               |                          |                |
| Main part                                                                                     | 1.4301 (X5CrNi18-10)                                                                                            |                          |                |
| Handle                                                                                        | 1.4301 (X5CrNi18-10)                                                                                            |                          |                |
| Others                                                                                        | 1.4305 (X8CrNiS18-9)                                                                                            |                          |                |