



Product data sheet REF 09.85.01

25.08.2026

General data			
article number	09.85.01	Yankauer suct tube 31cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	38749	GTIN	04046826013738
Risk category	2a (EU MDD)	UMDNS-Code	10-566
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,065	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK03
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
			
size	31cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	1.4301 (X5CrNi18-10)		
Others	1.4305 (X8CrNiS18-9)		