



## Product data sheet REF 09.85.07

11.12.2025

| General data   |   |                          |                |
|--|---|--------------------------|----------------|
| article number   | 09.85.07  | DeBakey suct tube 29cm   |                |
| Unit   | piece   | Customs code (HS code)   | 90189084       |
| GMDN   | 38749   | GTIN                     | 04046826013752 |
| Risk category  | 2a (EU MDD)   | UMDNS-Code               | 10-566         |
| Classification   | 6 (EU MDD)  | Base UDI-DI              |                |
| Net weight (kg)  | 0,058   | Operating instructions   |                |
| Latex free   | yes   | Preparation instructions | MPK03          |
|   |   |                          |                |
| Medicon media  | 450.00.03 (general catalogue the art of surgery)<br>450.00.53 (catalogue cardio ) |                          |                |
| Complete picture   |   | Detail images            |                |
|  |   |                          |                |
| size   | 29cm  | Coating                  |                |
| Storage conditions   |   |                          |                |
| No special requirements  |   |                          |                |
| Materials  |   |                          |                |
|  | Material number (DIN designation)   |                          |                |
| Main part  | 1.4301 (X5CrNi18-10)  |                          |                |
| Others   | 1.4305 (X8CrNiS18-9)  |                          |                |