






## Product data sheet REF 09.87.10

03.06.2026

General data			
article number	09.87.10	liposuct cann 30cm Ø 5mm hdl conn	
Unit	piece	Customs code (HS code)	90189084
GMDN	45414	GTIN	04046826013998
Risk category	2a (EU MDD)	UMDNS-Code	17-640
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,075	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK11
			
Medicon media	450.00.60 (catalogue plastic surgery ) 451.03.82 (brochure liposuction )		
Complete picture		Detail images	
 			
size	30cm Ø 5mm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	1.4305 (X8CrNiS18-9)		