





## Product data sheet REF 09.87.24

03.06.2026

General data			
article number	09.87.24	Toledo liposuct cann 15cm Ø 3mm hdl conn	
Unit	piece	Customs code (HS code)	90189084
GMDN	45414	GTIN	04046826014124
Risk category	2a (EU MDD)	UMDNS-Code	17-640
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,056	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK11
			
Medicon media	450.00.60 (catalogue plastic surgery ) 451.03.82 (brochure liposuction )		
Complete picture		Detail images	
			
size	15cm Ø 3mm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	1.4305 (X8CrNiS18-9)		
Others	1.4021 (X20Cr13)		