






## Product data sheet REF 09.87.37

03.06.2026

| General data   |   |                                   |  |
|--|---|-----------------------------------|--|
| article number   | 09.87.37  | liposuct cann 20cm Ø 3mm hdl conn |  |
| Unit   | piece   | Customs code (HS code)            | 90189084   |
| GMDN   | 45414   | GTIN                              | 04046826014230   |
| Risk category  | 2a (EU MDD)   | UMDNS-Code                        | 17-640   |
| Classification   | 6 (EU MDD)  | Base UDI-DI                       |  |
| Net weight (kg)  | 0,059   | Operating instructions            | G098395  |
| Latex free   | yes   | Preparation instructions          | MPK11  |
|   |   |                                   |  |
| Medicon media  | 450.00.60 (catalogue plastic surgery )<br>451.03.82 (brochure liposuction ) |                                   |  |
| Complete picture   |   | Detail images                     |  |
|   |   |                                   |  |
| size   | 20cm Ø 3mm  | Coating                           | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions   |   |                                   |  |
| No special requirements  |   |                                   |  |
| Materials  |   |                                   |  |
|  | Material number (DIN designation)   |                                   |  |
| Main part  | 1.4301 (X5CrNi18-10)  |                                   |  |
| Handle   | 1.4305 (X8CrNiS18-9)  |                                   |  |